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Instructions for reprocessing according EN ISO 17664:2004 for suction cannulas

Article REF 0400 – 0469

Important information for reprocessing

1. Use suitable personal protective equipment when reprocessing the suction cannulas.
2. Incorrect reprocessing of the suction cannulas can lead to health hazards.
3. Cannulas are thermo-sensitive products.

Sterilize only in a mixture containing at least 1.5 kg metal instruments and do not place close to the wall of the boiler vessel to avoid projection of the heat onto the product that might melt it. Place the cannulas on a sheet of filter paper and not directly on the trays. If necessary, sterilize in sterilization boxes.

4. A specialized rack equipped to clean the hollow spaces is required for machine cleaning and disinfection of hollow items.
5. Prior to first use, cannula suction tips must be prepared according to these instructions.

Restrictions for reprocessing

Basically no restrictions. Dispose of any cannula that may be damaged.

Preparation for cleaning

Cleaning agent: enzymatic cleaning agent or aldehyde-free disinfectant with cleaning power. Scrub the suction tips using an appropriate cleaning brush and the liquid surface of the cleaning solution. Then visually inspect the lumen; further preparation is required until the lumen no longer contains foreign material.

Cleaning and disinfection (machine)

Equipment: Cleaning/disinfection machine in accordance with DIN EN ISO 15883

Cleaning agent: Mildly alkaline cleaning agent with pH > 10 at 55° C

Neutralization: Phosphoric acid. The suction cannulas must be adapted to rinse inside the lumen on the loading rack.

1. Check the adaptor is securely in place before starting the machine and once the process is over.

2. Validated program parameters:

5 min	Cold pre-rinse	> 30°C	Tap water
10 min	Cleaning phase	60°C	De-ionized water
2 min	Neutralization	60°C	Mixed water
3 min	Disinfection	93°C	De-ionized water

3. Inspect the products visually for cleanliness and residual moisture when removing them from the cleaning/disinfection machine. If necessary, perform a follow-up cleaning by hand including a subsequent rinsing using de-mineralized, low-bacteria water and/or subsequent drying with sterile compressed air.

Cleaning and disinfection (manual)

According to KRINKO/RKI/BfarM recommendations, hollow bodies of risk class "semi-critical B" should preferably be reprocessed using a machine process validated by the operator.

Inspection and function test

Inspect the suction cannulas for damage. Damaged suction cannulas must not be used on patients.

Packaging

Individually: Use a sterile barrier system (e.g. paper film packaging) in accordance with EN ISO 11607; this system must be intended by the manufacturer for steam sterilization. Suction cannulas may be packaged individually or in pairs. The sealed seam must not be under pressure. After heat sealing, inspect the seam for any problems. Re-package and re-seal the suction cannulas if there are any problems.

Sets: Sort the Suction cannulas into the provided sieves, so that they will not be damaged during transport. The weight of the sieves must not exceed 10 kg. Sterile barrier systems in accordance with EN ISO 11607 or sterile containers in accordance with DIN 58953-9 are suitable for packaging the sieves.

Sterilization

Equipment: Sterilizer in accordance with DIN EN 285 or small steam sterilizer in accordance with DIN EN 13060 with a Type B process.

Procedure: Steam sterilization with fractionated pre-vacuum, sterilization temperature 134°C, dwell time of at least 3 min. All validated fractionated steam sterilization processes for hollow medical instruments are suitable. The process parameters to be used in the process validation are listed below.

Evacuation phase:	1. Pre-vacuum	200 mbar	
	2. Pre-vacuum	200 mbar	
	3. Pre-vacuum	200 mbar	
Sterilization phase:	Dwell time	5 min	134°C
	Drying	16 min	200 mbar

After sterilization, check suction cannulas and their packaging for damage, contamination, and moisture. Dispose of damaged suction cannulas. Suction cannulas in damaged, dirty, or moist packaging must be reprocessed, packed, and sterilized again.

Storage

Store in a dry and dust-free location with no major fluctuations in pressure or changes in humidity. The storage period depends on the packaging system used and must be determined independently for each facility.

Peridenta care Dental Produkte e.K., Ernst-Grote-Str. 37a, 30916 Isernhagen, Tel. 0511-6960090, Fax 0511-6960091