Deutsch

English

Español

Français

Italiano

REF 90-324-58-40

Revision 5

Date of Release: 2018-04

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Contents

1	Product Liability and Warranty	24
1.1	General Information	24
1.2	Intended Purpose.....	24
1.3	Warranty	25
1.4	User's Inspection	25
1.5	Hotline.....	25
2	Notices Concerning this Document.....	26
2.1	Validity of this Document	26
2.2	Symbols Used in this Document.....	27
3	Application	28
3.1	Performance Data	28
3.2	Fixing the Neutral Electrode on the Connection Cable	29
4	Safety Notes.....	30
4.1	Checking for Operational Readiness.....	31
4.2	Risks of Mutual Disturbances	31
4.3	Before Use	32
4.4	Precautionary Measures During Application	36
4.5	After Application	37
5	Processing, Cleaning, Disinfection and Sterilization	38
5.1	Storage and Transport.....	38
6	Disposal	39

1 Product Liability and Warranty

1.1 General Information

We thank you for having decided to buy one of our products. This product bears the CE-marking, which means that it satisfies the essential requirements laid down in the EC Directive concerning medical devices.

We are the manufacturer of this product:

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · D-78532 Tuttlingen · Germany

Postfach 60 · D-78501 Tuttlingen · Germany

Tel. +49 7461 706-0 · Fax +49 7461 706-193

info@klsmartin.com · www.klsmartin.com

1.2 Intended Purpose

Neutral electrodes from Gebrüder Martin are intended for the use with electrosurgical devices from Gebrüder Martin. Neutral electrodes serve for the purpose to close the HF circuit during monopolar coagulation and cutting of biological tissue. Neutral electrodes are an electrosurgical accessory and assist in the treatment of a disease.

Disposable neutral electrodes consist of a carrier made of soft plastic, onto which the actual electrode surface is mounted as a metal foil. Upon this, an electrically conductive gel layer is applied, which establishes skin contact and attaches the electrode to the patient's skin.

Connection to the electrosurgery unit is established via a reusable cable with connector clip.

An electrosurgery unit that comprises a corresponding NE monitoring system (e.g. PCS, Patient Control System) can check the contact quality of the neutral electrode on the skin and sound an NE alarm if necessary. Observe the Instructions for Use of the respective electrosurgery unit.

Gebrüder Martin recommends use of a disposable neutral electrode with split contact surface. It provides better protection against unintentional burns under the neutral electrode.

The connection cables to be used (see chapter 2.1 "Validity of this Document", page 26) are not included in the scope of delivery.

No modifications on the product are allowed under any circumstances.

NOTICE

The disposable neutral electrode is a **non-sterile** single-use product and may be used only once!

⚠ WARNING

Danger of injuries from insufficient qualification of the operator!

Safe application of bipolar forceps requires the user to be familiar with the technology and application modes of electrosurgery.

Accessories for electrosurgery may be used only by persons specially trained or instructed in their use!

1.3 Warranty

Our Standard Terms and Conditions of Sale effective at the time shall apply. Agreements diverging from these Standard Terms and Conditions do not restrict the legal rights of the buyer.

Any warranty exceeding the above provisions shall require a contractual form.

1.4 User's Inspection

- Immediately upon receipt, the goods must be checked for completeness and potential damage in transit.
- Notice of any such damage must be given immediately.
- Disposable neutral electrodes are packed in a fully hygienic but non-sterile packaging.

1.5 Hotline

- Should you have any questions on how to handle the product or use it for clinical applications, please do not hesitate to contact the Product Management:
Tel: +49 7461 706-243
Fax: +49 7461 706-190
- Should you have any technical questions, please do not hesitate to turn to our Martin Service Center:
Tel: +49 7461 706-343
Fax: +49 7461 706-408
E-mail: service@klsmartin.com
- Should you have any questions concerning maintenance contracts or training courses, please contact our Technical Service Manager:
Tel: +49 7461 706-332
E-mail: service@klsmartin.com

NOTICE

Each packaging and in some cases also the product is marked with a specific lot number (LOT) and a reference number (REF). Please indicate LOT and REF in case of any complaint.

2 Notices Concerning this Document

Non-observance of this document can lead to serious or even lethal patient injury!

Improper handling and care as well as non-intended use can lead to premature wear and/or pose a risk to patients and users!

It is the operator's responsibility to ensure that all personnel handling the product have understood and do observe the notes and instructions in this document.

- Every user is required to read this document completely and follow it carefully.
- In particular, be sure to heed all cautions, warnings and danger notices.
- Keep this document accessible to users at all times.
- Please also observe the Instructions for Use of the electrosurgery unit. In particular the safety notes therein must always be observed!

2.1 Validity of this Document

This document applies to the following products:

Designation	REF
Disposable neutral electrode PCS-DuoSafe for adults and children, split contact surface, 107 cm ²	80-344-09-04
Connection cable for disposable neutral electrodes by Gebrüder Martin	80-294-40-04
Connection cable for disposable neutral electrodes maxium® e-version (Erbe)	80-294-43-04
Connection cable for disposable neutral electrodes maxium® i-version (Valleylab)	80-294-44-04

2.2 Symbols Used in this Document

Throughout this document, important information (such as general or safety-related notices) is marked with the following symbols and signal words:

WARNING

Life hazard or serious injury!

Indicates a situation which, if not avoided, **could** result in death or serious injury!

CAUTION

Risk of injury!

Indicates a situation which, if not avoided, could result in minor or moderate injury!

NOTICE

Risk of material damage!

Indicates a situation which, if not avoided, could lead to material damage (loss of time, data loss, device/machine failure, etc.)!

3 Application

3.1 Performance Data

Parameter	Range
HF power transfer	up to 400 W
Increase of skin temperature	< 6 °K
Permissible patient weight	> 5 kg (children and adults)


Designation	Rated voltage [kVp]	Monitoring of contact quality possible	Sterilizability	REF
Disposable neutral electrode PCS-DuoSafe	n. a.	Compatible: Minicutter ME 102 ME MB 1 ME MB 2 ME MB 3 ME 411 maxium® type ME 402 Not compatible: MD 62	No	80-344-09-04
Connection cable for disposable neutral electrodes (Martin)	0.5	n. a.	Yes	80-294-40-04
Connection cable for disposable neutral electrodes maxium® e-version (Erbe)	0.5	n. a.	Yes	80-294-43-04
Connection cable for disposable neutral electrodes maxium® i-version (Valleylab)	0.5	n. a.	Yes	80-294-44-04

3.2 Fixing the Neutral Electrode on the Connection Cable

- 1 Cable connection clip
- 2 Disposable neutral electrode

The disposable neutral electrode must be connected to the cable connection clip before use

- Open the latch (3) of the cable connection clip (1) and insert the tab (4) of the disposable neutral electrode into the cable connection clip.

- Close the latch and check secure fit of the cable connection clip.

The cable is tightly connected to the tab of the disposable neutral electrode.

4 Safety Notes

WARNING

Failure to observe the safety measures leads to risk of burns of the patient and user!

- A **one-part** neutral electrode cannot be monitored. In case of insufficient contact, there will be **no warning signal!** Generally we recommend use of a disposable neutral electrode with a split contact surface.
- The indicated maximum power (see chapter 3.1 "Performance Data", page 28) applies to usual application times of a few seconds. Interspersed cooling periods may be necessary.
- There is risk of electric shock if the output voltage of the electrosurgery unit is higher than the rated voltage of the accessory.
 - Make sure that the output voltage of the electrosurgery unit cannot exceed the rated voltage of the accessory!
- Monitoring of the contact quality is possible only when using a split contact surface. However, intra-operative monitoring is required even in this case!
 - Ensure that the contact quality monitor is compatible with the neutral electrode to be monitored.
- Regularly check connection cables and accessories for damage, when required also using a magnifying glass. Damaged accessories must not be used!
- There are risks that can be caused by neuromuscular stimulation. This may occur in particular when using modes that produce arcs between the active electrode and the tissue.
- Make sure that the disposable neutral electrode is safely attached to the skin over the entire contact surface, and that the cable connection clip is safely attached to the tab of the neutral electrode.
- Guide the supply line to the disposable neutral electrode in such a manner that neither the patient nor other supply lines are touched. Avoid cable loops!
- After use, remove the disposable neutral electrode from the skin by softly pulling on it. Begin on the corner and not on the connection cable or the tab in order to avoid injury to the skin or defects of the accessory.
- Always remove the connection cable from the electrosurgery unit by pulling on the connector. Never pull on the connection cable! This can lead to damage inside the connection cable that is difficult to notice and will result in premature failure of the product.
- To remove the connection cable from the disposable neutral electrode, open the cable connection clip and carefully pull off the connection cable.

- An audible alarm due to insufficient contact between neutral electrode and patient is sounded only if the disposable neutral electrode is used together with a compatible electrosurgery unit, see 3.1 "Performance Data", page 28.
- Do not bend the contact surface.
- Do not apply mechanical stress to the bend protector areas of the supply cable.
- Do not stress the product mechanically during storage either (e. g. no pressure from above), to avoid functional defects that can lead to burns.

4.1 Checking for Operational Readiness

When using the disposable neutral electrode with split contact surface, the quality of contact to the patient's skin must be checked. The color changes from red to green (status light next to the NE connection socket on the electrosurgery unit) only when the contact is flawless.

- Red = contact quality insufficient
- Green = contact quality good

The disposable neutral electrode is ready for operation only if the status light is green!

4.2 Risks of Mutual Disturbances

Failure to comply with the following instructions may result in death or serious injury!

- Pacemakers
can be damaged or disturbed by application of electrosurgery. Consult a cardiologist before the operation.
 - Metallic shunts in the HF current path
Implants can change the current density, leading to unnoticed endogenous burns.
 - Cables installed in parallel
Cables that are installed in parallel to the connection cable can inject noise and lead e.g. to disturbance of monitor images.
 - Earthed objects
Upon contact with an operation table with insufficient isolation of the patient, unintentional leakage currents can increase to such an extent that the device function can change or burns result.
 - Disinfection/cleaning agents
Insufficiently dried agents can lead to burns.
-

4.3 Before Use

- Check expiration date of the packaging. Do not use disposable neutral electrodes after expiry.
 - Open packaging only immediately before application.
 - Before use, perform a visual inspection of the disposable neutral electrode and its connections.
 - Damaged, modified or expired disposable neutral electrodes must not be used, as application safety may be impaired.
-

⚠ WARNING

Danger of burns!

Connection cables with defective or broken insulation can lead to burns of the user and patient or cause fires.

- If damaged, discard disposable neutral electrode immediately!
-

⚠ CAUTION

Danger of burns!

Failure to observe the following instructions can lead to injuries!

- Endogenous gases can be ignited by sparks that occur during the intended use of the electrosurgery unit.
 - Some materials saturated with oxygen, such as cotton or gauze, may catch fire from the sparks that occur during the intended use of the electrosurgical unit.
 - Use of flammable anesthetics or oxidizing gases such as nitrous oxide (N₂O) and oxygen should be avoided when an operation in the chest or head area is performed, unless these substances are removed by aspiration.
 - Whenever possible, non-flammable substances should be used for cleaning and disinfection. Combustible substances that are used for cleaning or disinfection or as solvents for adhesives should be evaporated before application of electrosurgery.
-

⚠ CAUTION

Failure to observe the following instructions can lead to injuries!

- To avoid image interference on the monitor, do not install the connection cable directly in parallel to camera cables.
 - At any rate, we recommend monitoring systems that feature devices for limiting HF currents.
 - Set the output power to the lowest setting possible for the respective purpose.
 - Connection cables must never come into direct contact with the patient's skin.
 - Do not lay out connection cables in loops.
 - Place patient on a dry, electrically insulating surface.
 - Earth operation table.
 - Skin-to-skin contact (e.g. between the arms and the body of the patient) must be avoided, for example by placing dry gauze in between.
 - In case of surgical procedures in which HF current could flow through body parts with a relatively small cross section, application of bipolar technology can be preferable to avoid unintentional tissue damage.
 - Do not reduce the conductive contact surface by cutting. The metal inlays that may be exposed during cutting could come into direct contact with skin and result in a shunt circuit.
 - If necessary, attach a smaller disposable neutral electrode.
 - Neutral electrodes whose gel layer is injured can cause 2nd or 3rd degree burns.
 - Disposable neutral electrodes whose gel layer is injured must not be used.
 - Ensure that the gel-free connection tabs are fully covered by the connection clip so that they do not come into contact with the patient's skin.
 - Do not use conductive gel to improve the contact. As the application thickness cannot be controlled and the gel may dry unevenly, there is the risk of distortion of the even contact surface, which can lead to so-called hotspots.
 - If necessary, attach a new disposable neutral electrode.
 - Attach disposable neutral electrode without overlapping!
-

⚠ CAUTION**Danger of injuries from unsuitable application area!**

- Select application area carefully and attach disposable neutral electrode over its full surface!
 - Select a convex skin area that is muscular or well supplied with blood in the vicinity of the operation field but at least 20 cm away from it. We recommend the upper arms or thighs as application site for adults and the torso as application site for children.
 - Select application so that no currents flow in the area of the heart.
 - Avoid skin areas with strong hair growth, scar tissue, bony protrusions, and inflamed or injured areas.
 - The disposable neutral electrode must be attached closer to the operation field than any ECG electrode.
 - Keep foreseeable current path in the body as short as possible. Never lead current paths through the entire body.
 - Under no circumstances lead current paths across the area of the heart.
 - Avoid application sites where fluid could accumulate. There is the risk of accumulations of combustible liquid under the patient or in depressions such as the navel or in body cavities such as the vagina. Remove any fluid that has accumulated in these locations before applying HF current.
 - Guide the supply cable to the disposable neutral electrode in such a way that it does not come into contact with other cables or the patient.
 - Avoid cable loops.
 - Only then connect neutral electrode, foot switch and connection cable to the electrosurgery unit.
-

NOTICE

- Prepare application site carefully:
 - Shave evenly
 - Clean entirely (ensure site is free from grease, and remove cosmetics)
 - Dry thoroughly
-

⚠ WARNING

Risk of burns from improper handling of the neutral electrode!

The risk of burns under of the neutral electrode is particularly high if monopolar cutting or contact coagulation currents with particularly high power and longer duration are used (e. g. TUR-P, endometrium ablation).

In order to minimize the risk of burns or other problems in the area of the neutral electrode, during application of the neutral electrode please mind the following:

- The neutral electrode must butt against the patient's body as close to the operation field as possible, reliably, and with its whole surface. For an operation field on the torso, the upper arms and thighs are good application sites.
 - Safe contact of the neutral electrode must be ensured for the total duration of the high frequency application.
 - Application of the neutral electrode to a limb must not impair perfusion. Particularly for longer operation times it must be made sure that the patient does not lie on the cable connection clip of the neutral electrode (risk of pressure necrosis).
 - The current paths in the body should be as short as possible and run in the longitudinal or diagonal direction of the body, not across it, the latter particularly not on the chest. Any metal parts in or on the body should be removed if possible, insulated, or paid special attention to.
-

⚠ WARNING

Risk of burns from high HF current concentration!

The hazard potential of HF surgery increases with the power applied. Therefore please note the following:

- The high-frequency power should be set to the lowest possible setting for the respective application.
 - An insufficient effect while using the standard setting can, for example, be caused by poorly applied neutral electrodes, insufficient contact in plug connections, cables that are damaged under the insulation or encrusted electrodes. This must be checked before the HF power setting is increased.
 - A defect in the electrosurgery unit may cause undesirable increases in output power.
-

4.4 Precautionary Measures During Application

WARNING

Danger of injury from damage to active implants!

In patients with active implants such as pacemakers or implanted electrodes, the application of the electrosurgical unit poses a risk. Application may result in irreparable damage to the active implant or impairment of its function.

Do not apply the neutral electrode above implants or other metal parts, nor above bone protrusions or scarred tissue. If necessary, prepare the application site by cleaning and degreasing it; strong hair-growth is to be removed. For removal do not use substances that desiccate the skin (e. g. alcohol).

CAUTION

Danger of burns from power modifications!

- Increasing the output of the electrosurgery unit during application can lead to severe burns!
-

If despite default setting of the electrosurgical device the usual coagulation or cutting power should not be achieved, or if the set power decreases during operation or the electrosurgery unit breaks down (with or without NE alarm), follow these instructions:

- Even if the coagulation or cutting effect is reduced during operation, never increase the output power of the device without prior examination of the surgical situation! Specifically, check the following:
 - Flawless contact of all HF plugs and connection cables
 - Functioning of the foot switch or the finger switch on the handle, respectively
 - Insulation of the HF components (external damage)
 - Cleanliness and wear of the distal end of the active electrode(s)
 - Correct application of the neutral electrode
- When repositioning the patient during a surgical procedure, ensure that the neutral electrode is still in full contact with the skin over the entire contact surface!

⚠ WARNING

Risk of infection from contamination!

- Do not reposition the disposable neutral electrode!
 - Repositioning of disposable neutral electrodes can lead to cross-contamination in the patient. If necessary, attach a new disposable neutral electrode.
- Reuse of disposable neutral electrodes can lead to cross-contamination between patients.
 - Do not reuse disposable neutral electrodes!

NOTICE

- On electrosurgery units with handles that feature separate finger switches for activation of the HF operating modes **Cutting** and **Coagulation**, different settings regarding current type and power output can be selected for the respective operating mode. By pressing the activation buttons **Cutting** or **Coagulation** during regular use, the electrosurgery unit can output different levels of power. This is not a safety issue and no reason to stop a required surgical procedure.

4.5 After Application

- After application, remove disposable neutral electrode from the attachment point on the connection cable. Do not pull on the connection cable.
- After use, remove the disposable neutral electrode from the skin by softly pulling on it. Begin on the corner (not on the connection cable or the tab) in order to avoid damage to the skin or defects of the accessory.
- Check the contact surface of the disposable neutral electrode pulled off.
 - If the contact surface is intact, dispose of the disposable neutral electrode correctly according the guidelines, see chapter 6 "Disposal", page 39.
 - If it has burned areas, keep disposable neutral electrode in an airtight bag for later purposes.

5 Processing, Cleaning, Disinfection and Sterilization

- Disposable neutral electrodes **must not** be reused! They **must not** be disinfected and **must not** be sterilized!
 - Connection cables can be reprocessed, see chapter 3.1 "Performance Data", page 28.
-

NOTICE

Potential damage from impermissible reprocessing of connection cables!

Failure to observe the reprocessing instructions can lead to destruction of the connection cables or to the following damage:

- Loss of all planned electrical or mechanical properties
 - Deformation of the plastics
 - Increasing porosity of the soldered joints in the connector
-

5.1 Storage and Transport

- Store products in a clean, cool and dry place.
- Protect against direct sunlight (min. 5°C (41°F), max. 30°C (86°F)).
- In case of return shipments send only clean and disinfected products in sterile packing.
- Observe the imprinted expiration date. Unopened packages can be stored for two years.
- Use up opened packages within seven days.
- Defective disposable neutral electrodes and connection cables must **not** be repaired. Replacements are to be provided in due time.

6 Disposal

WARNING

Risk of infection from contamination!

- Avoid skin contact during disposal, and use only clearly marked disposal containers.
 - After exceedance of the product lifetime and/or in case of defects, we recommend disinfecting the disposable neutral electrode and storing it in a closed container. Only thus can transmission of microbes or carryover of infection (e. g. skin diseases) be avoided.
 - Discard disposable neutral electrode into hazardous hospital waste.
-

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Martin Nederland/Marned B.V.
KLS Martin UK Ltd.