How to use

FILPIN
restoration retention system

1 Remove the required amount of pins and the drill from dispenser pack and autoclave at 134°C for 4 minutes before use. Always follow the manufacturers instructions for use. Latch the correct sized drill head (for your chosen size of FILPIN) onto a slow-speed handpiece and make a pilot hole up to the drill's self-limiting depth of 2mm. Repeat as necessary. NB: Ensure the pilot hole is drilled at a speed of 1500 rpm or less.

2 Latch a FILPIN onto the handpiece and insert into the hole - hardly any pressure is needed to encourage the pin to insert. Once the correct depth of 2mm is reached, FILPIN will self-shear at exactly the right point. Repeat as necessary. NB: Ensure the pin is inserted at a speed of 1000 rpm or less.


3 If required, bend the protruding pin/s with a Filpin Pin Bender or other suitable instrument to suit the restoration. Final build up of the restoration can now take place.

CE
0088

4 The final restoration.

FILRATCH Pin-Insertor

For inserting FILPIN when the dentine is too soft or the depth inadequate as it gives excellent tactile control.

FILPIN is from FILHOL DENTAL, specialising in the manufacture of high quality retention pins for 40 years.

Available in two sizes in unique dispensers each containing 35 x FILPINs and 1 x drill

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