
Fit Kofferdam

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A handwritten signature in black ink. The name 'Silvia' is written in a cursive script, with the 'i' and 'a' connected. Below it, the name 'Feiger' is written in a similar cursive style, with a long, sweeping tail on the 'g' that extends downwards.

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Table of Contents:	Page
Preface	32
Fit for the Rubberdam Technique	33
The whole purpose of a Rubberdam	33
Advantages for the Treating Team	33
• Protection	33
• Eases your Work	33
Disadvantages of the Rubberdam	33
Important Citations	34
Rubberdam Products and their Properties	34
Quality Inspection	34
Colours and Qualities (Thicknesses) of Rubberdam	34
The five Rubberdam Qualities (Thicknesses)	35
FIT NON-LATEX	35
Rubberdam Templates	36
Rubberdam Stamps	36
Rubberdam Punches	36
The FIT PUNCH	37
Punching	37
Rubberdam Clamps	37 – 40
Quick Clamp Orientation	41
Rubberdam Forceps	41
Rubberdam Frames	41
The FIT RUBBERDAM FRAME (acc. to Young)	42
The FIT FOLDABLE FRAME	42
A short Excursion into History	42
Important Accessories!	43
• Interdental Wedges	43
• Wedjets	43
• Rubberdam Napkins	43
• Dental Floss	43
• Heidemann Spatula	43
Examples of the Rubberdam Technique for various	
Areas of Indication	44 – 48

Isolating the Maxillary Anterior Teeth	46 – 48
The Isolation	48
• of an already cemented Bridge	48
• of children's first and second Molars for Fissure Sealing	48
• for the Removal of Amalgam	48
How to Punch an oblong Hole	49
Isolating the mandibular Teeth	50
Isolating Cervical Cavities (Class V)	50 – 51
Hatch cervical Clamp No. 214	51
Additional Tips and Tricks	52
Tray Equipment for the endodontic Treatment	53
A Tray for the conservative Treatment	53
How to position the Patient and yourself for a convenient	
Placement of the Rubberdam	53
How should I handle Clamps and Clamp Forceps?	54
How should the Clamp be positioned on the Tooth?	54
What should I observe while processing Composites under the Rubberdam?	54
How can I prevent the Rubberdam from Ripping at the Contact Areas?	54
How can I seal off a holey or ripped Rubberdam?	55
What can I do if the Rubberdam does not easily slide over the Teeth?	55
How can I fit a Rubberdam on Patients with sensitive Reactions	
to natural Rubber?	55
What should I observe during the Removal of the Rubberdam?	55
How can I fit a Rubberdam over Wisdom Teeth?	56
Closing Words	57

Introduction

Safety aspects and the endeavour to achieve optimal results in the field of adhesive dentistry and endodontics are reasons for using rubberdam in the dental practice. The rubberdam is a 'conditio sine qua non' during the placement of ceramic and composite inlays, condensed gold fillings, lateral composite fillings and adhesive bridges. Besides that, hygienic concerns for practice team and patients, protecting the patient's mucosa and gingiva against the irritating effects of chemical substances during root canal rinses and use of certain bleaching methods, for gingival retraction and the protection of cheeks and lips make the use of a rubberdam strongly recommendable. And on top of all this, an efficient use of rubberdam saves time and work and provides a clear work field.

Every dentist knows the rubberdam and was forced to deal with this subject during his studies. Yet in Germany only 5 – 15% of all dentists regularly use a rubberdam, while in the USA, Switzerland, and Scandinavia a rubberdam is fitted four to seven times as often. The reasons stated for the rejection in Germany are the time it takes to fit the rubberdam and its complicated use. J. K. Ingle paraphrased this argument as follows: "The most time – consuming about a rubberdam consists in convincing a colleague of its advantages and indispensability."

But the true reasons are mainly found in insufficient knowledge and lack of practicing the rubberdam techniques. The rubberdam is a simple, yet efficient tool that should help and not burden the user. Its advantages should more than compensate for the disadvantage of the additionally required effort. The problem can only be solved if the time needed to fit a rubberdam on the patient can be minimized. And that is only achievable by utilizing the simple rubberdam techniques to begin with, and with continuous practice. Course and information material could be a viable introduction to this subject.

With their courses, the video film, a CD and an information brochure, Hager & Werken render outstanding services to the subject and diversification of the absolute isolation as well. After fundamentally revising and supplementing the first and second edition, the third "FIT RUBBERDAM" brochure is now available. The historic aspects and both, advantages and disadvantages, are discussed introductorily, followed by a description of the constituent parts of the rubberdam set. The basic techniques and use of the rubberdam are thoroughly and systematically described and complemented by photographs. This information brochure contributes to an extensive and quick realization of the rubberdam techniques through comprehensible language and pictorial documentation.

Prof. Dr. med. dent. habil. Eckhard Beetke
Polyclinic for Conservative Dentistry University of Rostock/Germany
Department of Cariology and Endodontics

Fit for the Rubberdam Technique

The Inventor of the Rubberdam

Dr. La Roche (a Frenchman) declared in 1883 that he had applied the rubberdam technique already since 1857 and consequently claimed to be the inventor of this method. But it was actually Sanford Christie Barnum (1836 – 1885), a New York dentist, who on 15th March 1864 used a rubberdam for the first time, and later demonstrated its use during the Dental Society Congress in New York in June 1864. Initial reports were published in August 1864, and by 1867 the rubberdam technique was already widely spread.

The whole Purpose of a Rubberdam

Advantages for the treating team:

Protection

- Prevents swallowing of objects
- No aspiration
- No mucosal damage by rinses
- Protects the treating team from infections
- Retraction of soft tissue (gingiva, lips, cheek, tongue)

Eases your work

- Work field remains dry
- Clear visibility
- Sterile work process
- Gain of time approx. 20%
- Mouth always open
- Discussion interrupted
- No rinsing
- No change of cotton rolls

The biggest advantage for the patient is the comfort of feeling that the treatment happens outside of his mouth. Now he is no longer tormented by the otherwise intensive aspiration (which results in a very dry mucosa) and cotton rolls, or worse by the operator's fingers holding them down (which often triggers a gag reflex). By isolating the oral cavity, the patient's every-day oral conditions are maintained and he can easily swallow and breathe.

Disadvantages of the rubberdam:

- Loss of axial orientation during the preparation of the access cavity
- Possible traumatization of papillae
- Higher demands during X-ray taking
- Possible allergies (corrective: use of silicone squares; FIT NON-LATEX)

Important Citations

"The most time-consuming about a rubberdam is convincing a colleague of its advantages and indispensability." (J. I. Ingle)

"With a little practice it is possible to fit a rubberdam in a normal treatment situation within approximately one minute (two minutes maximally)." (J. I. Ingle)

Rubberdam – why and what for?

Nobody answered this question shorter and more precise than G. V. Black in 1908. He wrote:

"The rubberdam serves to keep the working area around the teeth clean, dry, and antiseptic, if necessary. The latter is especially desirable when it comes to root canal treatments." (G. V. Black, 1908)

Rubberdam Products and their Properties

Rubberdam – A natural latex sheet that is commercially available in rolls, or in 15 x 15 cm squares. Rubberdam is highly elastic, which is necessary for its purpose and expandability. Unfortunately, it does not retain these optimal properties indefinitely (approx. 1 year), then it becomes brittle, rips easily and does not isolate sufficiently any longer. If the latex sheet is always stored in a refrigerator or even freezer, it is still usable after a long time (far more than a year if stored as described).

Quality Inspection

The latex sheet is optimally serviceable if it can be stretched between the hands until it becomes absolutely transparent, regardless of its age.

Rubberdam – Colours and Qualities (Thicknesses)

Colours:

- Beige – light: preferred for endodontic treatments because of its transparency, but only for processing composites (poor display of outlines).
- Brown or even grey – dark: good colour contrast, reduction of light reflectance (problems may occur during the axial orientation).
- Green: pleasant, friendly, calming shade that provides a good display of outlines (light reflectance), non-glaring under OT lamps (reduction/no light reflectance), smells of mint for the patient's comfort.

Additional colours: blue, light blue, pink, lilac as alternative to the above mentioned light, dark, and green shades.

The five Rubberdam Qualities = Thicknesses

- Thin (0.13 – 0.18 mm) – The thinner the rubberdam the easier it is to apply, but the easier it also rips and it does not isolate as well as a heavy quality.
- Medium (0.18 – 0.23 mm) – The most commonly used quality, which is especially easy to handle, even for novices of the rubberdam technique, and of course for conservative treatments.
- Heavy (0.23 – 0.29 mm) – This quality allows a good gingival retraction and hardly ever rips.
- X-strong (0.29 – 0.34 mm) – It is extremely tear-proof and allows a maximal gingival retraction, but is hard to position in situ.
- Sp-heavy (0.34 – 0.39 mm) – This rubberdam is used only if special protection of the tissue is required.

Because the thicknesses of the rubberdam overlap, the highest measurement of one specification is at the same time always the lowest one of the next higher specification. But rubberdam users "feel" even these slight differences.

FIT NON-LATEX

- Absolutely free of powder and odour-free
- Does not stick
- Very good gingival retraction strength, also seals off well around the tooth
- Fibre-reinforced for extra tensile strength
- Especially pleasant touch
- Approx. 0.3 mm thick (x-strong), green
- Its reflective colour creates a clear work area
- Especially flexible, very close to latex material in its application
- Helps to eliminate powder and latex allergies, and is also recommendable for the use on Class V cavities

Summary:

If composites are processed under rubberdam, it is mandatory to select the colour prior to fitting the rubberdam. Afterwards the teeth desiccate rapidly and strongly, making a proper colour selection impossible. During the initial phase of "befriending" with rubberdam, a thin quality according to the material description should be chosen.

But practice perfects dexterity over the course of time, and it is recommendable to change over to a thicker quality later. Two rubberdam qualities are used in most practices. A thin one in the beginning, which is later exchanged for the medium to strong, and finally a very strong quality.

Worth knowing:

The stronger the rubberdam the better it isolates and the better its gingival retraction capacity.


Important:

One side of the rubberdam is smooth, the other one powdered. The smooth side should always point towards the oral cavity, and the powdered side towards the dentist. This has two reasons:

- The smooth side is easier to pull over the teeth that are to be separated (advantage for the dentist).
- The tongue will not be coated with rice or corn flour if it touches the rubberdam (comfort for the patient).

Rubberdam Template

It is usually made of white vinyl and enables an exact marking of the tooth position on the rubberdam. For this procedure the smooth side of the rubberdam is placed onto the template, the powdered side marked with a pen (to ensure the marking is not laterally reversed). The use of a template has the advantage that only the teeth to be separated are marked, which makes the orientation in situ easier, because no drawings are on the rubberdam. (Fig. 1)


Rubberdam Stamp

It is possible to transfer complete arches onto the rubberdam; the standard positions of the teeth are simply punched out where required.

Rubberdam Punches

Their 5 or 6 sizes of holes warrant the punch of perfect hole at any time. The following applies for the punch acc. to Ainsworth (Fig. 2):

- Hole No. 5 (largest)**
recommended for clamp teeth (at the end of the dental arch and for very large molars)
- Hole No. 4 (large)**
universally for molars
- Hole No. 3 (medium)**
for canines and premolars in the UJ and LW
- Hole No. 2 (small)**
for anterior teeth of the upper jaw
- Hole No. 1 (smallest)**
for the very slender anterior teeth of the lower jaw


An absolutely tight seal against moisture is ensured if the hole sizes are varied according to the diagram above.

The FIT PUNCH

The FIT PUNCH (Fig. 3) always punches perfectly neat holes into the rubberdam. The movable punch spike simultaneously catches the entire border and punches a tear-proof hole.

Fig. 3

Punching

In order to create a perfectly punched hole, the rubberdam square is stretched between thumb and ring finger as well as index and middle finger of the left hand (prevents wrinkles and subsequently undesired double perforation). The plate of the intact, sharp punch is now brought under the rubberdam and the punch handles are squeezed, pushing the spike into the chosen perforation of the plate to cut a perfectly round hole.

If the visual examination reveals holes with tears, cuts or jags, the rubberdam will rip at this imperfect hole once it gets stretched over the teeth. But if the perforation is perfectly round, the rubberdam can be stretched without ripping until it is almost transparent.

Rubberdam Clamps

Generally it is distinguished between two basic clamp shapes. They were developed to meet the requirements of the different rubberdam application methods. Differentiated is between wingless clamps and winged clamps.

Wingless Clamps (Fig. 4)


Their jaws are short and rounded. This type of clamp is secured prior to aspiration (loop dental floss distally around the clamp). The clamp is then placed onto the tooth and the rubberdam stretched over it. Finally the frame is fitted. The letter W in front of the number generally indicates a wingless clamp.

Fig. 4

Winged Clamps (Figs. 5/6)

Their jaws have "appendage wings". Their wings are pulled through the hole in the rubberdam, and both are placed on the tooth together. Corresponding to academic opinions and practice relevance some clamps emerged as being most commonly used.

Labial Clamps = double bow Clamps (Fig. 7)

In endodontics it is usually sufficient to isolate only the tooth that will be treated. Double bow clamps are especially suitable for single-rooted teeth.

- No. 9
Universal winged double bow clamp for labial cavities on anterior teeth and premolars (endodontics).
- No. 210
Mainly used for labial cavities and endodontic treatments in the anterior upper jaw. Often used for premolars as well, and sometimes on the 2nd molar of the lower jaw.
- No. 211
Predominantly used for labial cavities on anterior teeth of the lower jaw, and as endo-clamp for anterior teeth and premolars of the lower jaw.
- No. 00
Suitable for very small teeth mainly used in the lower jaw (34, 33, 43, 44). The clamp is also convenient for teeth that are irregularly positioned in the arch. Placed on lower incisors, they prevent (due to the closely situated lower lip) the rubberdam from shifting.
- No. 0
Its high distal bow facilitates placement on long-crowned teeth. Its narrowly shaped jaw is predestined for elongated teeth or narrow dental necks.

15 14 24 25
35 34 44 45

- No. 1
Standard clamp for upper premolars. Their distinctively rounded jaws enable a good fit on the dental arch on the gingival margin. It will not rotate even on very round-shaped premolars (15, 14, 24, 25).
- No. 2
Standard clamp for premolars, mainly in the lower jaw. Their flat jaws prevent gingival damage. Preferred for 35, 34, 44, 45.
- No. 2 A
For premolars with strong dental necks. Its wide jaws make it especially suitable for premolars of the upper and lower jaw.

15 14 24 25
35 34 44 45
- No. 206
Narrow jaws, for premolars with thin necks (closely resembles No. 1); (15, 14, 24, 25).
- No. 207
Closely resembles clamp 2A and therefore often used on

15 14 24 25
35 34 44 45
- No. 5
Universal clamp for molars in the upper jaw. The clamp with its deep jaws fits especially well at the gingival margin of round-shaped teeth and prevents any rotation (18, 17, 16, 26, 27, 28).
- No. 56
This clamp is shaped like No. 5, but it is bigger and therefore suitable mainly for large upper molars.
- No. 201
Large molar clamp for slightly bulbous teeth (mainly 18, 17, 16, 26, 27, 28).
- No. 7
An excellent standard clamp for molars in the lower jaw. The flat design of the jaws prevents gingival damage (38, 37, 36, 46, 47, 48).
- No. 202 – Large clamp for molars in the lower jaw (also see clamp No. 7) with an extended distal bow. This enables a practice relevant treatment of distal cavities on larger molars.

Useful Clamps to complete the Rubberdam Indication Range

- No. 8 A
Their strong point lies in clamping down on small upper and lower jaw molars that are not completely erupted or irregularly shaped. The four tips can be placed under the gingiva without damaging any tissue. Their jaw inclines downwards and thus enables the fitting of rubberdam below the gingival margin (deep gripping).
- No. 14 A
Especially suitable for incompletely erupted, underdeveloped or irregularly shaped molars. The tight fit of the clamp is warranted through a 4-point-contact and the inclination of its jaws. A narrow bow characterizes the clamp. Its shape corresponds to No. 14, but it is bigger (deep gripping).
- No. 212
Cervical clamp (= tissue retraction) for cervical cavities (class V fillings). This indication requires extra strong rubberdam to achieve the highest possible gingival retraction. It is imperative to punch the hole for the tooth to be treated approx. 2 mm further towards cervical compared to the arch of holes for the adjoining teeth. It is recommendable to punch the hole for this tooth one size larger than the others. KERR or Stent's composition is used to stabilize the clamp.

- No. 214
Hatch cervical clamp. The large, rigid double bow (= small jaw) and a more movable clamp arm (= large jaw) provide four contact points and form two anatomic jaws. The double bow with its small jaw is fitted vestibularly resp. labially, so that the gingiva is pushed back into the desired position. The movable clamp arm with its large jaw is fixated by tightening the knurled screw.

Caution:

Anaesthetics should be administered prior to fitting this clamp in order to push both rubberdam and gingiva back towards apical for the treatment of class V cavities.

Advantage compared to No. 212: Does not require KERR or Stent's composition for stabilization and fixation.

- No. 138 + 139
Left and right model, both wingless. Originally developed for wisdom teeth, these clamps are well proven for many areas of application, particularly for lower molars. Corresponding to the shape of lower jaw molars, the buccal jaw is more distinct than the lingual jaw. In addition, the lingual jaw sits deeper and the jaws are riffled. The arrangement of jaws on different levels allows the clamp to be fitted at almost any angle. The riffles provide a secure grip, which makes the clamps suitable for incompletely erupted teeth (e.g. fissure sealing on children's teeth).

FIT RUBBERDAM – TIGER CLAMPS

- 7 different clamps for the labial, premolar and molar area
- With special retention for better fixation on the tooth
- These rubberdam clamps correspond in shape and range of indication to the initially described article numbers. Their tiger-pattered riffles provide even better fitting to the tooth.
- The following FIT RUBBERDAM T-CLAMPS are available: 9T, 1T, 1AT, 2T, 2AT, 14T, 56T

No. 1 AT – A very interesting and useful clamp for anterior teeth and premolars that attaches the rubberdam to the root of the premolar. The downwardly inclined jaws grip the root tightly, even if subgingival cavities are present. Not suitable for healthy teeth in healthy gums (deep gripping).

FIT KOFFERDAM – CHILDREN'S CLAMPS

- Especially developed for the anatomic shape of children's molar
- Non-symmetric beak shape grips children's molars tightly and securely
- Flat distal bow designed for children's narrow and small mouths
- 8 different types in shiny finish
- The following clamps are available (the numbers indicate which quadrant and which tooth the clamp is meant for): 55 E, 54 D, 64 D, 65 E, 74 D, 75 E, 85 E, 84 D.

QUICK CLAMP ORIENTATION

Clamp:

No. 9, 214 (Hatch) (0, 00, 210, 211, 212, 1 AT)

No. 2, 1, (0, 00, 2A 206, 207)

No. 7, 8, 8A, 14A, (5, 201, 202, 56)

No. 7 (202)

No. 8 (201, 56, 5)

No. 8A

No. 14A

Tooth Category:

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Premolars

Molars

Universal lower jaw molars

Universal upper jaw molars

Small molars, incompletely erupted,
deeply destroyed or partially retained

Large molars, incompletely erupted,
deeply destroyed or partially retained

Rubberdam Forceps

(Fig. 8) Thanks to its excellent retention, rubberdam forceps grip any clamp bow securely and without tilt. They enable the adaptation of clamps to the tooth. Besides that, they are used to spread and fixate the clamps, and to remove them after the work is finished.

Rubberdam Frames

(Figs. 9/10) Stretching a rubberdam with rubberdam weights is not necessary any longer. Nowadays, a rubberdam frame with some sensibly placed, well-gripping retainer spikes allows for a quick and successful fixation of the rubberdam. Today's rubberdam frames are made of metal or synthetic material. Synthetic frames that fold in the middle ease the accessibility significantly; for that reason they are especially suitable for endodontic treatments.

Fig. 8

The FIT RUBBERDAM FRAME (according to Young)



(Fig. 9) This rust-proof and flexible U-shaped rubberdam frame is equipped with few, but sensibly arranged, well-gripping cylindrical retainer spikes to secure the rubberdam.

The FIT FOLDING FRAME

(Fig. 10) The Fit folding frame with its middle hinge is made of radiolucent, autoclavable synthetic material which makes it predestined for endodontics. It is used like a regular rubberdam frame, but one half can be folded away to obtain a good position for taking control x-rays. The enforced part of the frame provides good orientation, as it always points towards the chin.

A short Excursion into History ...

(Fig. 11) Stretching the rubberdam successfully with rubberdam weights requires a great amount of practice. The upper lappet must be stretched with an elastic spring behind the patient's head, and the weights must be attached to the lower lappets.

Important Accessories!

Interdental Wedges

Interdental wedges are used to fixate matrices, and also the rubberdam in the patient's mouth. The wedges are made of hickory (wood) or nylon (synthetic material). Their ideal anatomic shape ensures perfect adaptation to the tooth. Hickory wedges are placed approximately in order to push both, tissue and rubberdam, towards apical, and to prevent the rubberdam from ripping when it is wrapped around the drill for the treatment of extensive cervical defects.

Wedjets

A disposable rubberdam stabilization ligature made of natural latex. It is available in three strengths; each piece has a length of approx. 2 m.

- light blue = extra thin
- yellow = thin
- orange = thick

Various suppliers offer basic sets, which – besides the items that are available in the practice anyhow – contain the equipment that is needed to start. All items, clamps, and accessories can also be purchased individually.

Rubberdam Napkins

A must for patients with sensitive skin (they are practical for the operator and very comfortable for the patient) that should not get rationalized out of existence – especially during extensive treatments. The rubberdam is pulled through the hole in the napkin and fitted onto the frame. The napkins absorb saliva, water, and sweat.

Dental Floss

Dental floss is mainly used to control the accessibility of interdental spaces. But it makes an excellent tool to skirt a rubberdam around difficult contact areas as well. A dental floss ligature is often helpful in attaching a rubberdam to the tooth (especially to deciduous molars and canines). Ligatures should be fixated buccally with a surgical knot. A long ligature makes subsequent manipulation and removal easier.

Caution:

The use of non-waxed dental floss may result in cutting the dam. Solution: Use waxed, or even better PTFE-coated dental floss (e.g. MIRAFLOSS TAPE*).

*MIRAFLOSS TAPE slides easily through even the narrowest interdental spaces and unfolds optimally in there, allowing the rubberdam to be efficiently reworked. Furthermore it is extremely tear-proof and does not fray.

In case a clamp is ill fitting, disturbing or even dispensable, a ligature made of dental floss, or a double or multiple string of floss is threaded through the contact point and remains there as tie-down. A narrow rubberdam strip pulled into the interdental area also makes an excellent anchor. (ideal size: 1 cm x 3 cm).

Heidemann Spatula

This instrument is excellently suitable for a short separation of the contact points, and the inversion of the rubberdam around the necks of teeth. A continuous air stream aimed towards the sulcus makes the inversion an easy task.

Miradont Gel:

In addition to its original postoperative application, the gel also makes an excellent rubberdam lubricant. For that purpose it is thinly applied behind the rubberdam perforations to ease its application. The gel is absolutely water-proof and neutral in taste.

Examples of the Rubberdam Technique for various Areas of Indication

Endodontic isolation = Isolation of a single tooth

Rubberdam and rubberdam clamp are simultaneously brought onto the tooth to be isolated, meaning a suitable winged rubberdam clamp is clamped into the dam, for endodontic isolation in the anterior upper jaw e.g. clamp No. 210. For this procedure the first wing is placed behind the perforation from the powdered side.

The perforation is then stretched a little and the second wing can be pulled through. The clamp is now located on the powdered side of the dam, and the little wings on the smooth back side.

The wings keep the rubberdam perforation open and the tooth, on which the clamp will be placed, is well visible.

This places the clamp in to the field of vision, so to speak. In practice it means that the clamp, which is fixated into the rubberdam, is now brought into the clamp forceps.

Caution! Clamp forceps work like telescope forceps.

This means its jaws open if pressure is applied respectively when the lock is set. Due to that fact the clamp is easily removed from the clamp board, placed onto the tooth and removed again after the treatment.

In case of endodontic isolation in the anterior upper jaw it is recommendable to attach the clamp in the gingival margin from vestibular, to pull it over the strongest point of the tooth past the tuberculum and to carefully deposit it palatally (cp. Fig. 12).




The clamp forceps can now be put aside, rubberdam and clamps are in place.

(Fig. 12) The clamp forceps is attached and both clamp and rubberdam brought towards the mouth in one step. In order to make the pertinent tooth visible, the clamp is spread and then placed onto the tooth.

(Fig. 13) The rubberdam is now easily pushed from the clamp wings with a Heidemann spatula whose curved side is placed over the edge. The rubberdam slides under the clamp jaw and lies snugly against the tooth.


In order to stretch the rubberdam it is pulled through the opening of the FIT FOLDING FRAME and fixated on the well arranged retaining spikes.

In case that sulcal fluid still seeps on the dam, the dam is not slipped far enough under the clamp jaws. An approx. 30 cm long piece of waxed or PTFE-coated dental floss, initially placed from palatal behind the compact clamp (Fig. 14), will eliminate this problem.

Each end is then threaded under the clamp bow of the corresponding interdental space (Fig. 15). The floss is now intersected from vestibular under the clamp jaw and carefully tightened. (Fig. 16) The rubberdam fits snugly on the gingival margin and warrants the absolute isolation. The ends of the floss are either cut, or the floss is completely pulled out.

Now the rubberdam is fitted perfectly and the actual endodontic treatment can begin. During the opening of the canal (under water-cooling) it is recommendable to have an assistant aspirate the water mist with an aspiration canula that has a wide and faultless pelotte. (Caution! Burr formation on the pelotte may cause the dam to rip). Very effective is the use of the MIRASUC 3 P aspiration canula (aspiration canula with exchangeable pelottes, available in three shapes). The conical pelotte should be mounted for that purpose. (Fig. 17)

The patient can swallow saliva as always. Once the treatment has progressed to the point at which a control x-ray should be taken – thanks to the FIT FOLDING FRAME this is easily done. (Fig. 18)

Helpful items for this task are the EMMENIX FILM HOLDER (Fig. 18) or an artery clamp (the more gracile CROWN CLIPPER II – curved is very convenient as well).

In case of an endodontic treatment it usually suffices to isolate only the tooth that will be treated. This shows that endodontic cases are especially suitable for getting acquainted with the rubberdam.

Once the sitting is over the rubberdam clamp is returned to the clamp forceps. The fastest and easiest way of doing that is to place the forceps' retention pin into the visible perforation of the clamp first, and to spread the forceps a little afterwards while moving it over the tooth in order to place the second retention pin into the perforation of the second clamp.

The clamp forceps is then locked, whereby the clamp opens and can easily be removed from the tooth. After a single tooth isolation the patient can bear having rubberdam and frame removed in one go. The only way of performing that task gently and with little splatter is by placing the thumb from vestibular under both, rubberdam and frame, and placing the index finger of the right hand on the rubberdam under the cutting edge of the tooth, while the middle finger also lies on the dam and the ring finger and little finger are brought under the frame. The rubberdam can now be slightly lifted off the patient's face and splatter-freely pulled off towards the cutting edge.



Isolating the Maxillary Anterior Teeth

In order to provide a clear field of work and vision, it is recommendable to generally isolate teeth 13 - 24 with rubberdam when working on the anterior upper jaw. If this basic principle is followed, the assistant can always use the patient-free time to prepare and perforate the rubberdam squares for the pertaining indication and keep it readily available in the refrigerator. The single-stage technique enables the operator to achieve his goal, e.g. the above mentioned indication, very quickly, as rubberdam, rubberdam clamps, and rubberdam frame are applied at the same time.


The rubberdam, which is already perforated from 13 - 24, is taken by the left hand (powdered side visible on top), and the large rubberdam metal frame is placed approx. 1.5 cm from the left margin, during which the upper end of the frame should be flush with the upper edge of the rubberdam. Rubberdam and rubberdam frame are now in the left hand, the right hand is used to stretch the rubberdam until it can be fixated onto the cylindrical retainer spikes of the frame. (Fig. 19)

The right hand holds now both, rubberdam and frame, while the left hand is used to stretch the rubberdam far enough to be fixated onto those retainer spikes as well. The excess rubberdam left, right, and on the lower end of the frame can be used to fold a small collecting pocket, which is convenient to collect water and particles of old fillings during the preparation.

The attached sheet should now be folded lengthwise, so that this side of the frame is covered with rubberdam. The hand that is on the folded side now grabs the frame at the lower corner, so that the fold cannot turn back again. Rubberdam and frame are now in that hand. The free hand grips the folded-over end of the dam and pulls it up high enough to be fixated onto the upper edge of the frame on this side. (Figs. 20/21) These steps are now repeated on the opposite side. Due to the heavy stretching, the pocket may limit the field of vision. A perfect pocket can be created by gripping a lappet of the rubberdam from the middle of the pocket and to fixate this on the cylindrical retainer spike in the middle of the U-shaped crossbar. (Figs. 22/23)

A winged premolar clamp suitable for tooth 24 (e.g. No. 2, 2A, 209 or 206) is now brought onto the clamp forceps, and its mesial wings are inserted into the size 4 perforation hole in such a way, that its distal bow actually points towards distal. This also helps in marking tooth 25 in case of initial orientation problems. If the bow is located above this guide, everything is fine when fitting the dam on the patient.


The rubberdam clamp is now attached buccally on the gingival margin. Slight pressure of the work hand causes the clamp to spread open so it can be pulled over the occlusal plane and the strongest point of the tooth. The forceps' lock is released by decreasing pressure, and the clamp is gently deposited palatally. (Fig. 25)

The anchor tooth 24 is now clamped and located outside the rubberdam, and the clamp forceps are put aside. The rubberdam is now gripped with both hands in such a way that it can be steeply positioned into the interdental space (like a matrix band).

Now the perforation is pulled crosswise over the incisal edge into the next interdental space. Again, the latex sheet with the next perforation is inserted interdentally as steeply as possible, and the perforation is pulled crosswise over the incisal edge into the next interdental space. As each tooth, which is supposed to be isolated, appears outside of the latex sheet, the reverse side (distal of 13) is fixated with a scrap of rubberdam, an interdental wedge, a Wedjet, a string of dental floss, a ligature or – on 13 – with a clamp. (Fig. 26)

A ball-pointed plugger or a Heidemann spatula – its curved side over the plane – is brought between latex and clamp. The instrument is slightly turned between the fingers towards left, respectively right, causing the rubberdam to find its way under the clamp as if by magic. (Fig. 24)

In case sulcal fluid still seeps through and moistens the teeth, the procedure should be finished by another work step, the so-called inversion. This means, the edges of the rubberdam are wrapped around the necks of the teeth. This is the only way to warrant an absolute seal against moisture. If the operator inverts with the curved side of a Heidemann spatula placed over the edge, this side of the instrument is brought under the perforation and carefully deposited on the gingiva. During the entire procedure an assistant should aim an even air stream towards the sulcus.

If the operator now uses index and middle finger of the left hand to push the dam slightly onto the upper lip/gingival area, and removes the spatula while taking advantage of the inflating effect, the rubberdam can easily be turned over.

But if the operator prefers an inversion by means of dental floss, he threads a 30 – 40 cm long string of waxed dental floss into the interdental space on both sides in such a way, that it rests snugly against the gingival margin from palatal. Now the ends of the floss only have to be crossed from vestibular and slightly tightened, and the rubberdam will turn over. (Fig. 28)

The powdered side of the rubberdam lies on the pertaining tooth during the inversion procedure. And once the little wings of the clamp, which is positioned on tooth 24, are brought onto the rubberdam a perfect isolation is achieved.





Fig. 25 Application of rubberdam clamp and frame
Fig. 26 Mounting the rubberdam
Fig. 27 Moving the clamps onto the dam
Fig. 28 Inversion (possibly with floss only)

The Isolation


- of an already cemented bridge
- of 1st and 2nd molars for fissure sealing
- for amalgam removal

A third technique is also well suitable for this indication. But some facts should be taken into consideration before the start: the heavier the chosen rubberdam sheet, the better it isolates. Since there is no chance of pulling a rubberdam into the interdental space of a bridge, a so-called oblong hole is created. The oblong hole is a proven method for fissure sealing under rubberdam, and even for the removal of extremely close amalgam fillings.

The rubberdam sheet is placed with its smooth side onto the template; the area that is to be isolated is marked on the powdered side, and an oblong hole size 2 hole is punched in that area. (Figs. 29/30)

How do I punch an oblong Hole?

3–4 overlapping size 2 perforations are punched in order to prevent the rubberdam from being pulled in approximal.

A fourth technique is especially recommendable when working on children. The rubberdam sheet is first slipped over the teeth that are to be isolated. A scrap of rubberdam is then pulled distally from the last isolated tooth, and the dam is fixated on this side. In case this causes the dam to roll off one or two teeth, now is the time to reattach it. At that time the assistant should already have chosen a molar clamp (e.g. Nos. 7, 8, 200 or 18) and brought it into the clamp forceps.

If this method is not used to isolate each single tooth but a cemented bridge construction, or sealing fissures on tooth 26 or 27, the best procedure is to start the oblong hole perforation distally on tooth 27 and to extend the hole from tooth 25 to mesial, so that teeth 27, 26 and 25 are positioned within the oblong hole but already outside of the rubberdam. A strip of rubberdam or something similar is now threaded mesially from tooth 25 to fixate the rubberdam at this side. A molar clamp is placed onto tooth 27 for final fixation.

For better visibility in the molar area it is recommendable to attach the clamp first palatally on the gingival margin. The work hand applies a little pressure to the forceps to lead the spread clamp over the occlusal plane and the strongest point of the tooth onto the buccal gingival margin.

The rubberdam is now fixated, and the rubberdam sheet, that hangs loosely over the oral cavity, can be stretched over the frame. If no pocket is created it is up to the operator to place the frame under the rubberdam or on top of it for stretching.

Important: The frame should never be pushed into the patient's face during fitting. After the treatment the clamp is removed first, then the rubberdam strip, and finally both, rubberdam and frame.

Isolating mandibular Teeth

For work in the lower jaw it is recommendable to employ a technique that allows visibility of the clamps.

A suitable wingless molar clamp (e.g. No. W7, W8, 18, 28) should be selected. An approx. 40 cm long piece of dental floss is wrapped around the distal bow to prevent the clamp from being swallowed. (Fig. 31)

In order to bring the clamp into the clamp forceps, the ends of the floss are wrapped around the index finger of the left hand once or twice. The clamp can now be held between middle finger and thumb of the left hand. Holding the clamp forceps in the right hand the retention pin of the clamp forceps is first brought into the round perforation of the clamp arm that is located at the thumb. This way the forceps can rest on the thumb and spread widely enough to place the second retention pin into the round perforation of the opposite clamp arm at the middle finger. Now index finger, middle finger, and thumb of the right hand = work hand can lock the forceps.


Fig. 31

The secured clamp is now attached from lingual to the gingival margin of tooth 37. The right hand applies a little pressure to the forceps, causing the clamp to open widely enough to be lead over the occlusal plane and the strongest point of the tooth. Some pressure applied to the forceps releases the lock and places the clamp gently from buccally onto the gingival margin of the tooth. The fit of the clamp should be examined. For this purpose the index finger is used to apply a little pressure on the distal bow. If the clamp starts to swing during the stability test it is too big for the tooth and should be replaced by another one.

Isolating cervical Cavities (Class V)

If a tooth needs to be treated at its neck, it is best to isolate this tooth and in addition two distal ones on each side. The chosen rubberdam sheet should be of heavy or even better x-heavy quality to obtain the greatest possible gingival retraction.


Furthermore it is important to punch the hole for the tooth to be isolated approx. 2 mm further towards apical, compared to the perforations for the adjacent teeth.

The perforations for the adjacent teeth are punched as usual, according to the template and the given dental arch. (Figs. 32/33)

It is recommendable to punch the hole for the tooth to be isolated one size bigger than the others (Fig. 33). Both recommendations should be closely observed, especially if very deep subgingival damage is present.

Lets take a Look at the Isolation of a Class V Cavity on Tooth 11.

In this case, the x-heavy rubberdam sheet is best prepared to isolate from tooth 13 to tooth 24 (standard isolation of the anterior maxillar area). The smooth side of the latex is placed onto the template, and teeth 13, 12, 21, 22, 23, 24 are marked according to the template on the powdered side. Only tooth 11 is marked approx. 1–3 mm above the intended position (Fig. 32). Size 2 holes should be punched for teeth 13, 12, 21, 22, and 23, while a size 4 perforation is recommendable for the anchor tooth 24. A size 3 hole is suitable for tooth 11, which is to be treated. (Fig. 33).

Before the actual work field (tooth 11) is uncovered with the cervical clamp No. 214 (= Hatch Cervical Clamp), the rubberdam should be slipped over all teeth that are to be isolated. Tooth 24 is equipped with an adequate clamp. A suitable accessory (hickory wedge, rubberdam strip, Wedjet or a length of dental floss) is used to secure the dam distally of tooth 13.

Rubberdam and gingiva are now forced to retract (pushed towards apical) with the aid of clamp No. 214. An administration of anaesthetics is mandatory prior to this procedure. (Fig. 35)

Cervical Clamp No. 214 (Hatch)

A large, rigid double bow forms a small jaw that is attached vestibularly resp. labially. The movable part of the clamp forms a larger jaw that is attached lingually resp. palatinally (Fig. 34). Under slight pressure the four contact points slide gently under the gingiva without causing damage or injury, so that both, gingiva and rubberdam, are forced to retreat into the desired position and can then be fixated by tightening the screw. (Fig. 35)

Additional Tips and Tricks for Mounting a Rubberdam

Preparation of the work area


In order to work effectively it is recommendable to start the treatment with properly prepared and equipped trays.

A tray that is prepared for an endodontic treatment should be equipped with the following items:

- heavy rubberdam with universal perforation for endodontic treatments
- a universally usable premolar or molar clamp with wings, resp. a double bow clamp with wings that is suitable for anterior teeth and premolars (e.g. No. 210 for the upper jaw, No. 211 for the lower jaw, No. 9 for both upper and lower jaw, No. 6 for anterior teeth as well as premolars and the second lower molar)
- a radiolucent, FIT FOLDING FRAME
- rubberdam forceps
- a length of waxed dental floss (approx. 40 cm)
- a Heidemann spatula
- a complete set of basic dental instruments
- the EMMENIX FILM HOLDER
- dental film
- if desired, a slender artery clamp, resp. the CROWN CLIPPER II, curved

The EMMENIX FILM HOLDER provides useful services by holding films for control x-rays, during which the patient holds the film in situ by biting down on the holder. This prevents holder and film from slipping. Root canal instruments can be clamped in the CROWN CLIPPER. This is especially helpful for divergent root canals, since not even one's own fingers obscure the field of vision. The CROWN CLIPPER is also suitable for measurement exposures and for holding the film.

Fig. 36


A tray that is prepared for a conservative Treatment should be equipped slightly different.

The equipment consists of (Fig. 37):

- rubberdam of medium quality, already perforated for the row of teeth that will be isolated
- an adequate premolar or molar clamp with or without wings, depending on the area of application and preferred rubberdam application method
- a large U-shaped FIT RUBBERDAM STEEL FRAME
- rubberdam forceps
- a length of waxed dental floss (approx. 40 cm)
- a Heidemann spatula
- a complete set of basic dental instruments
- a wedge holder equipped with hickory wedges, or a scrap of rubberdam, maybe Wedjets

Fig. 37

How to position the Patient and yourself for the convenient Placement of the Rubberdam?

The operator should position himself and his patients in a way that enables him to pull the rubberdam towards himself. This means the fitting of rubberdam in the anterior maxillar area e.g. can best be accomplished while sitting by the reclined patient in the 12 o'clock position.

The fitting of rubberdam in the lateral anterior maxillar area is easier realized while sitting by the reclined patient in the 10 o'clock position.

Sitting by the reclined patient in the 2 o'clock position is recommendable for fittings in the left lateral anterior maxillar area.

Fitting in the entire mandibular area is easiest accomplished while sitting by the reclined patient in the 8 o'clock position. If those recommendations are observed, it is always possible to stretch the rubberdam towards oneself, so good results will be achieved quickly.

If one attempts to push the rubberdam away from oneself it takes much longer to accomplish the task successfully.

How should I handle Clamps and Clamp Forceps?

First of all it is advisable to get familiar with the way the clamp forceps works. It functions like a telescope forceps, so its "jaws" = work tips open when the handles are pressed together. Completely different compared to extraction pliers or technical forceps, whose way of functioning became second nature to the operator.


Now let us assume a right hander, who would best pick up the rubberdam clamp with thumb and index finger of the left hand while the distal bow points towards the palm.

Holding the clamp forceps in the right hand, the retention pin of the forceps is now inserted into the round perforation of the clamp arm that is grasped by the thumb. The clamp forceps can now rest on the left thumb while slight pressure is applied to the forceps with the right hand, causing the clamp to spread widely enough to insert the second retention pin into the round perforation of the opposite clamp arm. The forceps is now locked with index finger, middle finger, and thumb of the right hand.

How should the Clamp be positioned on the Tooth?

The clamp should grasp the tooth on 4 points in order to find a secure hold. If the clamp has only a 3-point contact, the tension of the rubberdam will cause the clamp to swing, rotate and wobble until it finally falls off the tooth. It is important to select a clamp that is not too big, or else the contact points lie on the papillae and probably cause injuries.

Important: Rubberdam clamps should generally correspond to the anatomy of the teeth to be clamped.


What should I observe while Processing Composites under the Rubberdam?

It is important to determine the colour PRIOR to fitting the rubberdam. As teeth desiccate very quickly during absolute isolation, a later decision would lead to the selection of significantly lighter colours.

How can I prevent the Rubberdam from Ripping at the Contact Areas?

Rubberdam often rips if one tries to force the rubber between the contact points of two teeth (in the occlusal area). In order to prevent this, it is recommendable to pull the rubber between the contact areas, starting with only one edge. This is done on all "easy" contact areas, while dental floss is used for "difficult" ones. First the edge of the rubber is pulled through the contact areas as described, and then the dental floss is pulled through to carry part of the rubber along. The dental floss is now returned from lingual, so that it is doubled up and both ends point towards buccal. (Fig. 39)

One end of the floss is carefully pulled out. During this process only a small amount of rubber should be carried along with the floss. The procedure should be repeated until the rubberdam has a perfect fit.

How can I seal off a holey or ripped Rubberdam?

In case the rubberdam is damaged by an instrument during the course of work, the hole can be closed with a cotton pellet drenched in cavity varnish. If larger rips exist it is better to attach a new rubberdam over the old one.

What can I do if the Rubberdam does not easily slide over the Teeth?

In most cases the patient's saliva suffices, and the rubberdam slips over the teeth easily. But sometimes the use of a lubricant is necessary to make the strips between the holes slide easier through the interdental spaces. The chosen lubricant must be water soluble and free of fat. Vaseline is not suitable. It would contaminate dental materials and agglutinate the gingiva. Besides that, the selected lubricant should not taste worse than the rinses that are used for endodontic treatments or the etching liquids or gels that are used when processing composites. A touch of MIRADONT GEL applied behind the perforations on the smooth side of the rubberdam helps the strips to slide easily through the approximal contacts.

How can I fit a Rubberdam on a Patient with sensitive Reactions to natural Rubber?

Rubberdam napkins – a must for treating patients with latex allergies. It is convenient for the operator and comfortable for the patient (especially during extensive treatments) if rubberdam napkins are not rationalized out of existence. The rubberdam is pulled through the hole in the napkin and fitted on the frame. The napkin absorbs saliva, water, and sweat. If necessary, the operator can also resort to Flexidam and FIT NON LATEX rubberdam.

What should I observe during the Removal of the Rubberdam?

In any case it should be done very carefully. If ligatures were applied, the threads should be cut with scissors. When removing the rubberdam it is important that no residue remains in the pocket region to prevent possible inflammations. Depending on the area in which it was attached, the rubberdam sheet is pulled out towards vestibular, labial resp. buccal and the interdental strips are cut with gum scissors (Fig. 40 Suture Scissors ES12).


Fig. 40

How can I fit a Rubberdam over Wisdom Teeth?

Due to the limited space and the close contact that wisdom teeth have to the soft tissue it is recommendable to attach the rubberdam with the distal bow inserted. The distal bow of the chosen clamp (suitable are winged and wingless clamps) is pulled through the perforation from the non-powdered back side of the latex sheet. This means the clamp is taken up with the forceps, both forceps and clamp are in the right hand (= work hand), while the latex sheet is held with the left hand. The clamp forceps with the integrated clamp is now brought under the sheet, and the sheet is stretched over the forceps with the right index finger. Now the left hand pulls at the rubberdam in such a way that the distal bow of the clamp slides through the latex. (Figs. 41/42)

The sheet is then gathered by the left hand and held laterally, so that clamp and rubberdam become clearly visible. (Fig. 43)

Once the clamp is positioned, the gathered sheet is spread again and pulled from the bow over tooth and clamp arms. Additional teeth are isolated in the usual way, and the dam is secured at the second end of the isolation. The dam may be fitted on a frame before the treatment is started. (Figs. 44/45)

Summary

Over the years that have passed since the first issue of this brochure, the interest in rubberdam techniques has steadily increased. The knowledge gained during numerous rubberdam workshops that were held so far lead to the desire to explain the rubberdam techniques even more detailed and graphic in this revised brochure.

On behalf of Hager & Werken, the author would like to take this opportunity to sincerely thank Prof. Dr. med. dent. habil. Beetke, who perused the manuscript and wrote this brochure's preface.

A handwritten signature in black ink. The name 'Silvia' is written in a cursive style with a large loop at the end of the 'i'. Below it, the name 'Peiger' is written in a similar cursive style, with a long, sweeping tail that extends downwards.

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HAGER & WERKEN GmbH & Co. KG