FR

COMPOSITION

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EN

COMPOSITION

Per sponge:
Native, non-denaturated, freeze-dried collagen
of bovine origin0.027 g

TSE/BSE

The risk assessment for virus transmission and animal transmissible spongiform encephalopathies has been performed in compliance with the provisions of Regulation 722/2012.

INDICATIONS

Local hemostasis for dental use when control of bleeding by ligaturing or other conventional means is ineffective or impractical.

CONTRA INDICATIONS

Patients with a genetic predisposition to allergy or a known history of allergy to collagen of bovine origin.

SIDES EFFETS

Allergic reaction,
Alveolalgia when used for packing of dental extraction sockets.

INSTRUCTIONS FOR USE

Before placing the hemostatic sponge, re-examine the socket thoroughly to remove all granulation tissue and bone debris.
Apply a dry sponge with or without pressure depending on the amount of bleeding, making sure that the socket does not dry out.
If necessary, the sponge can be held in place with a cross stitch over the socket.

The sponge will be completely and naturally resorbed.

WARNINGS AND PRECAUTIONS FOR USE

- Do not use in infected areas.
- Do not moisten the sponge.
- Any non-used piece of sponge should be disposed of since it is no longer sterile.
- It is recommended:
 - to use dry, smooth instruments,
 - to apply small pieces of sponge, just big enough to fill the socket.
- For single use only. Any reuse exposes other patients and the user to critical risks, including infection.
- Any excess not removed at the time of surgery may lead to a foreign body reaction that may present with or without clinical signs and symptoms as lesion or postoperative abscess formation upon imaging. Imaging may initially not be capable of distinguishing the difference.
- Check package integrity before use.

Do not use if package is open or damaged before first use.

- Do not resterilize.

STORAGE

So as to keep the sponges sterile and prevent them from regaining moisture, do not remove the heat-sealed film until use.

Keep away from sources of heat.

PRESENTATION

Box containing 3 blisters of 8 single use sponges.
Sponges radio-sterilised by β-ray irradiation.

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