ES

Estimado cliente

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EN

Dear customer

Thank you for choosing a quality product from Dentaaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current instructions for use supplied with the product and stored in the internet at www.dentaaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaaurum GmbH & Co. KG
Turnstr. 31 | 75228 Ispringen | Germany

2. General product description

The hyrax® click screw (Hygienic Rapid Expansion Screw) is a special expansion screw for rapid palatal expansion. It has a patented screw locking feature (Fig. 1) to prevent unwanted reversing. The retention legs of this expansion screw are welded onto preformed bands (Fig. 2) or onto so-called cap splints made of metal (Fig. 3). Alternatively they can also be polymerized into acrylic splints (Fig. 4). The bands or the splints are cemented into place as usual. Follow the instructions from the cement manufacturer.

3. Area of application

Used for rapid palatal expansion with a fixed appliance.

4. Fabrication of the hyrax® click appliance with preformed bands

Preformed bands are fitted on the maxilla left and right first bicuspid and left and right first molars. Using an impression, these bands are either transferred to or fitted on the hard plaster model. The retention legs of the hyrax® click screw are then bent to fit onto the four bands so that they can be welded into place (Fig. 2 + 5).

The hyrax® click screw should be “suspended” 1 – 2 mm above the palatal vault/gingiva (see Fig. 6) and its dorsal side should be in line with the mesio-palatal cusps of the first molars. In order to avoid leverage, it must be situated parallel to the occlusal plane. In order to make the task of inserting the safety key and activating the screw easier, the hyrax® click screw is fitted with the direction of rotation facing towards dorsal (Fig. 2 – 4).

5. Tips for working with hyrax® click screws

The bands can be connected by means of a remaloy® wire ø 1.0 mm (REF 528-100-00) either on the palatal or the buccal side as a means of reinforcement. These reinforcement wires are laser-welded above the retention legs. The palatal reinforcement wires should extend slightly over the bands and be supported by the adjacent teeth in order to distribute the expansion forces more effectively.

In order to bend the retention legs of the hyrax® click screw gently and easily, we recommend using the universal bending key (REF 611-122-00, Fig. 7 + 9) without the aid of pliers (lever principle).

If this key is not available, the retention legs must be stabilized by holding the expansion screw body with a pair of sturdy pliers, e.g. wire bending pliers Maxi (REF 013-522-00), in order to prevent the retention legs from breaking (Fig. 8).

When using the soldering technique please protect the body of the screw against excessive heat, e.g. with heat screen (REF 380-070-00).

6. Composition

The composition is included in the materials list; please refer to the catalog or www.dentaaurum.com.

7. Safety instructions

Always use only the blue safety key supplied to prevent swallowing (Fig. 10).

8. Information for single use products

The product is intended for single use only. Reconditioning of a product that has already been used (recycling) and its reuse on a patient is not permitted.

9. Product range

Description	straight legs	angulated legs	Quantity
hyrax® click Mini – 7 (Expansion 7 mm)	602-830-10 602-830-30 602-830-50	602-831-10 602-831-30 602-831-50	1 piece 10 pieces 50 pieces
hyrax® click Medium – 10 (Expansion 10 mm)	602-832-10 602-832-30 602-832-50	602-833-10 602-833-30 602-833-50	1 piece 10 pieces 50 pieces
hyrax® click Maxi – 12 (Expansion 12 mm)	602-834-10 602-834-30 602-834-50	602-835-10 602-835-30 602-835-50	1 piece 10 pieces 50 pieces

Packs (1 piece/ 50 pieces) with 1 safety key

Packs (10 pieces) with 10 safety keys

Safety key	611-120-01	10 pieces
Universal bending key	611-122-00	1 piece

10. Quality

Dentaaurum ensures faultless quality of its products. The recommendations in these instructions are based upon Dentaaurum Implant's own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as Dentaaurum has no influence on the processing on site.

11. Explanation of symbols used on the label

Please refer to the label. Additional information can be found at www.dentaaurum.com (Explanation of symbols REF 989-313-00).

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