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Privacy notice

EN

We process the personal data collected with this form (hereinafter „the data“) in order to ensure optimum support for our customers with regards to the device guarantee granted by us, and to be able to send them – as described in detail below – current information and offers on our products and services from time to time. Accordingly, METASYS Medizintechnik GmbH is responsible for compliance with the provisions of the data protection law.

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1. Purposes and legal bases of data processing

We shall use the given company master data (company, address and contact data, etc.) to register your company in our customer database. This shall in particular be used for the purpose of fulfilling contractual obligations in connection with the device guarantee granted by us. The legal basis for this is Art. 6 para. 1 item b) of the General Data Protection Regulation (GDPR). We shall process this data to be able to send our customers individual information and current offers on METASYS products and services from time to time. This data processing for the purpose of direct advertising shall also be in our legitimate interest. The legal bases for this are Art. 6 para. 1 items b) and f) GDPR.

2. Data security/passing on of data to third parties/duration of storage

The data you provide will of course be treated confidentially and protected by suitable technical and organisational precautions against loss and changes as well as unauthorised access by third parties. For internal administrative purposes, personal data of our customers are processed and stored in a central data management system of the METASYS Group, i.e. METASYS Medizintechnik GmbH and its EU affiliates. METASYS Medizintechnik GmbH is responsible for compliance with the provisions of this privacy policy and applicable data protection law. This data processing for administrative simplification and optimization is in the legitimate interest of METASYS Medizintechnik GmbH. The legal basis for this is Article 6 (1) lit. f) GDPR. We will not disclose your data to third parties outside the METASYS Group, unless we make use of specially selected and contractually obligated service providers (e.g. technicians, delivery services, printers) to fulfil the purposes stated in this privacy statement. The confidentiality and security of your data will be thereby be ensured in accordance with our privacy policy and applicable law. We shall store your data for as long as this is necessary for the aforementioned purposes and due to statutory documentation and retention obligations. If you no longer wish to receive information and offers from us, you can inform us at any time in accordance with section 3. In this case, we shall store your data for as long as this is legally required with regard to subsequent documentation or accountability obligations or with regard to the limitation of any mutual legal claims. The data will not be processed outside the European Union, unless the special data protection requirements of the GDPR are met.

3. Your rights

If and insofar as the respective legal requirements are met, you shall be entitled to request information from us about personal data or data processing concerning you (Art. 15 GDPR), rectification, erasure and restriction of personal data or data processing concerning you (Art. 16 to 18 GDPR) and transfer of personal data concerning you (Art. 20 GDPR). Furthermore, if the requirements of Art. 21 GDPR are met, you shall have a right of objection to data processing that relates to a „legitimate interest“ of the person responsible pursuant to Art. 6 para. 1 item f) GDPR or that is used for direct advertising. To exercise the aforementioned rights or if you wish to withdraw your consent to the processing of your personal data in whole or in part with effect for the future, you can contact us via the aforementioned contact details, in particular by using the e-mail address datenmanager@metasys.com. Art. 77 para. 1 GDPR also gives you the right of appeal to a supervisory authority, in particular in the Member State in which you are staying, in which you are working or where you suspect infringement of the GDPR to have taken place, if you believe that the processing of personal data concerning you is in contravention of the GDPR.

The data protection authority responsible for us in Austria is:

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Hygiene system

H1

Equipment Logbook
Assembly, operation and maintenance

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Information on guarantee claims**EN**

METASYS delivers exclusively on the basis of the present terms and conditions, insofar as contractual partners have not struck a deviating agreement in writing.

1. Guarantee claims on our products are valid for a period from 12 to 36 months (depending on the product).
2. This guarantee includes all claims regarding material failures that affect the function of the device. Damages caused by wrong or improper handling or normal abrasion and the exchange of the collection amalgam container are excluded from the guarantee claim. Compulsory maintenance rates have to be met precisely. The producer reserves to himself the right to request the corresponding product documents to ascertain service rates.
3. To assert the validity of the guarantee claim it is necessary to return the attached installation proof to the producer right after the installation. Installation and return of installation proof must be completed within 24 months after date of sale from METASYS.
4. Both the installation and the required inspections have to be carried out by a technician trained by METASYS and have to be registered in the METASYS equipment logbook.
5. We point out that guarantee claims for amalgam separators can only be asserted if the device has been treated with the suction system disinfectant METASYS GREEN&CLEAN M2 For air/water separations and central suction systems guarantee claims can only be asserted if the devices have been treated with the suction system disinfectant GREEN&CLEAN M2 or the cleaner GREEN&CLEAN CL according to the directions of use. In case of water decontamination system WEK/WEK Light, the unexceptional use of the water decontamination preparation GREEN&CLEAN WK is essential for the validity of a guarantee claim. This regulation does also apply to the hygiene system H1, which must exclusively be operated with GREEN&CLEAN H1. Implied warranty remains unaffected of this.

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


1. Index

The header on each page defines the user group particular information is aimed at.

Chapter	Page
1. Index	2
2. Explanation of the pictograms	2
3. General information	3
4. Product description	4
5. Explanation of the type plate	4
6. Technical data	4
7. Installation requirements	5
8. Modular assembly	6
9. Functional description	7
10. Cleaning cycles	9
11. Hygiene System H1 - maintenance	10
12. Hygiene System H1 - control panel	11
13. Electronic circuit board	12
14. Inspections	14
15. Cleaning and disinfection	14
16. Refilling the H1 cartridge	15
17. Disposal	15

2. Explanation of the pictograms

Information

Caution!

General warning sign

Refer to instruction manual/booklet

General information

Practice personnel, technicians


3. General information


The safety, reliability and performance of the appliance is only guaranteed by METASYS if the following instructions are observed:

- Assembly, alterations or repairs may exclusively be carried out by authorized service personnel in compliance with EN Standard 60601-1 (International Standard for Medical Electrical Apparatus, in particular Part 1: General Rules for Safety).
- The electrical installation must comply with the regulations of the IEC (International Commission for Electrical Engineering).
- The apparatus must exclusively be used in conformity with the instructions for installation, operation and maintenance.
- Only original parts may be used for repairs or replacements.
- All the guidelines provided by the manufacturer on the correct use of the Hygiene system H1 must be followed.
- All inspection and service work must be entered into the device documentation at the first pages of this manual
- When requested by an authorized engineer, METASYS agrees to make all documents available for the use of technically qualified service personnel.
- METASYS accepts no responsibility for damages caused due to external factors, such as wrong installation, improper use of the apparatus or unauthorized technical intervention.
- Users must study equipment and assure themselves of its good condition before every use.
- When the apparatus is dismantled at the end of its service life, it must be disposed of properly.

EN

1

2

4. Product description

The Hygiene System H1 **1** provides automatic and continuous cleaning, disinfection, defoaming and deodorisation of the entire suction system during operation, with integrated selective function of the suction hoses.

The system can be retrofitted with a tube holder for two suction hoses and two extra retainers with selectable configuration.

5. Explanation of the type plate

2 See illustration

- 2.1** Serial number
- 2.2** Address of the manufacturer
- 2.3** Mains supply data
- 2.4** Protection class II
- 2.5** Type BF sign

6. Technical data

Permissible water pressure	2–4 bar (working pressure: 2 bar)
Water connection	PU tube, nominal width 2 mm, external Ø 4 mm
Voltage	24 V AC
max. current consumption	450 mA
Frequency	50/60 Hz
max. power consumption	11 VA
Volume of the H1 cartridge	125 ml (pouch: 500 ml)
Suction system	wet or dry vacuum pump
Permissible vacuum	100 mbar – 250 mbar
Dimensions (H x W x D)	128 x 132 x 95 mm (excl. housing and suction tube holder)
Space requirements (H x W x D)	180 x 260 x 110 mm (excl. housing and suction tube holder)

Installation requirements

Practice personnel, technicians



7. Installation requirements

- A water filter **3.9** (40 µm) and a pressure regulator **3.8** (2 bar) must be attached to the mains water supply **3.7** (included in the scope of delivery).
- The underpressure may not exceed 250 mbar at the cannula.
- Voltage to terminal X2 must be 24 V AC.
- The manufacturer's regulations regarding the installation of particular units must be adhered to.

3 See illustration


- 3.1** Suction line
- 3.2** Suction hose 1
- 3.3** Suction hose 2
- 3.4** Vacuum switch
- 3.5** Connection cable (4- and 6-pole)
(control unit, support signal)
- 3.6** Voltage 24 V AC (red/black)
start signal for suction system (yellow)
- 3.7** Mains water supply
PU tube, nominal width 2 mm, external Ø 4 mm
- 3.8** Pressure regulator
- 3.9** Water filter

4 Fastening on the rear, see illustration

Modular assembly

Technicians


8. Modular assembly

To ensure maximum serviceability and to allow variations in assembly, the Hygiene System H1 consists of five modules:

5 Module 1: H1 cartridge & adapter plate or connection of pouch


The exchangeable H1 cartridge contains the disinfectant supply for an approx. 12-15 hours of suction time and is simply slid onto the adapter plate.

6 Module 2: H1 valve block, double membrane pump & DVGW (air gap) container

The valve block, in combination with the double membrane pump, is the heart of the Hygiene System H1. It determines the ratio of water to disinfectant, dispenses the cleaning agent to the respective handpiece and controls the elbow closure pieces.


In order to prevent any backflow of process water into the drinking water network, the fresh water is taken from a container with free outlet (per DVGW).

7 Module 3: Filter housing, elbow closure pieces & suction tubes


The filter housing contains a removable non-drip rake.

The elbow closure pieces (one per suction tube) serve as locking devices for the suction flow and in order to interrupt the supply of disinfectant if vacuum pressure drops. They are rotatably hung and form the connection of the suction tubes which are fitted with a thin inner tube that supplies the cleaning agent to the handpiece. The power of the suction flow can be adjusted with the regulator on the filter housing.

8 Module 4: Electronic circuit board & housing

The circuit board housing serves as a protection of the electronic circuit board. This forms the central control unit for the Hygiene System H1.

9 Module 5: Control panel

The control panel displays the main operating states and can be mounted where desired or as part of the METASYS suction tube holder.

For a more detailed description, refer to page 11.

Functional description

Technicians





9. Functional description

10 Function diagram, see illustration

- A** Suction flow
- B** Disinfectant GREEN&CLEAN H1
- C** Water
- 1** H1 cartridge or pouch
- 2** DVGW container (water)
- 3** Double membrane pump
- 4** Pump chamber
- 5** Suction flow closure
- 6** Chemical closure
- 7** Filter housing
- 8** Handpiece
- 9** Suction tube with inner tube
- 10** Amalgam separator or automatic separator
- 11** Suction system

Lifting the handpiece **8** of a suction hose results in a vacuum causing the respective suction flow sealing membrane **5** to be opened by an associated solenoid valve (L/R) **12a**.


The suction flow seal is designed as a pipe elbow and as such it is highly streamlined and is also insensitive to contaminants **12b**. The cleaning cycles start at the same time as the suction process.


EN

Functional description


Technicians

13

14

15

9. Functional description

13 + 14 See illustration

Water is drawn out of the DVGW container, or GREEN&CLEAN H1 is drawn out of the cartridge or pouch by means of a water pressure-driven double membrane pump.

15 See illustration

Once drawn off, the disinfectant is pumped through the knee sealing membrane and through lines in the suction hoses to the handpiece in use.


The cleaning mixture is then distributed evenly over the inner surface of the suction tube by the suction flow and thus ensures optimal cleaning and disinfection of the whole suction system.

16 Adjusting the suction flow rate, see illustration

The suction flow rate can be adjusted by a regulator above the filter cup.

A screwdriver or a coin may be used to adjust the regulator.

Please observe the markings on the regulator!

16

Cleaning cycles

Practice personnel, technicians

10. Cleaning cycles

A cleaning cycle only takes place with suction. As the vacuum decreases, the ongoing pumping of the cleaning medium is interrupted.

17 Initial cleaning, see illustration

- Initial cleaning at the start of each treatment stage. This program is activated by switching on the main device switch on the treatment unit. The increased disinfectant dosing creates the ideal starting prerequisites for the pending treatments with respect to defoaming and blood solution.
- The valve sequence (H1 disinfectant) with 10 x H₂O or 2 x H1, 10 x H₂O per suction hose is shown flashing on the external display.
- Sequence: 10 / 2 / 10 per hose
- The program duration is 1 min 5 sec per suction hose.

18 Normal operation


- takes place with a raised suction hose
- Sequence: 1 x H1 / 413 x H₂O / 1 x H1 per hose
- This sequence is halted with an interruption to the suction operation and is continued when the suction operation resumes or starts anew.

19 Special cleaning

- The special cleaning program is started by pressing the respective button on the operating interface.
- This program should be used at the end of each treatment day and after particularly intense loading of the suction system (e.g. after lengthy surgical procedures).
- Sequence: 40 x H₂O / 4 x H1 / 8 x H₂O per hose
- The program duration is 2 min 35 sec per suction hose.

This intensive cleaning rinses any residues or deposits out of suction system. Problems arising from periods of inactivity (such as evenings or holidays) are avoided. A final disinfection is guaranteed.









During these programs the Hygiene System H1 controls the suction system, which means that the suction tubes must remain in their holders.

EN

Hygiene system H1 - maintenance

Practice personnel, technicians

11. Hygiene system H1 - maintenance

Daily cleaning of the filter cup:

Please wear protective gloves when performing this task!

- 20** Turn down the container clip of the Hygiene System H1.
- 21** Pull out the filter cup horizontally.
- 22** Lift out filter rake.
- 23** Empty the amalgam particles into the designated collection container (METASYS ECOCENTER).
- 24** Clean or replace the filter rakes and grease the sealing lip, if necessary clean the underside of the housing.
- 25** Place filter cup with inserted filter rake horizontally back into the previously cleaned housing.
- 26** Turn container clip upwards and lock in place.

A change of the filter rake is recommended every 2 to 3 months. The filter rake must always be in place and the container clip closed when the suction system is in operation. If the performance of the suction system declines, check the filter rake for contamination.

Changing the H1 cartridge/pouch (see illustration 27)

Please also refer to the H1 instructions on the refill package and on page 11!

If the yellow monitor light no. 2 (empty tank warning) on the control panel flashes and a short alarm is audible every time the suction tube is lifted, please exchange the H1 cartridge or renew the pouch as follows:

Removing the H1 cartridge:

- Turn off the main switch of the unit.
- Press onto the front end of the cartridge briefly to release the locking mechanism (PUSH). The cartridge will slide into a position from which it can be pulled out on its slide rails horizontally.

Inserting the H1 cartridge:

- Fit the cartridge in the guide and slide in until it latches in the locking mechanism.
- Turn on the main switch of the unit.
- The initial cleaning program starts automatically (duration: 1 min 5 sec/suction hose).

Hygiene system H1 - control panel

Practice personnel, technicians

12. Hygiene system H1 - control panel

28 See illustration

28.1 Control light 1 (red): Malfunction

- flashing red and 5-times beeping in the event of a fault; additional short beep each time the suction hose is lifted, red flashing signal remains; hygiene system not functional, suction operation still possible:

Membrane break - suctioning functional but no cleaning and disinfection

- red light flashing without signal: Press the PUSH button once; if this does not work then an electrical fault exists:

Check electronics!

28.2 Control light 2: Empty tank warning

- yellow light and single beep whenever the suction tubes are lifted:


Change the H1 cartridge or H1 pouch!

28.3 Control light 3: Program in use



- green light on, suction hose in place:
ready for operation
- green light on, suction hose in use:
continuous cleaning in progress
- green light flashing slowly (approx. 3 sec):
initial cleaning in progress
- green light flashing quickly (approx. 1 sec):
special cleaning in progress

28.4 Start button for special cleaning program

- Press PUSH once to activate.

EN

13. Electronic circuit board

29 Wiring diagram, see illustration

- 29.1** Membrane rupture sensor (X1 - MB- SENS)
- 29.2** Voltage supply 24 V ~ (X2 - 24V ~)
- 29.3** Start of suction motor (X2 - SM)
- 29.4** Suction tube holder 1, 2 and 3 (X3 - ABL)
- 29.5** Sensor chemical tank (X4 - CH- SENS)
- 29.6** Valves V1 - V5 (X5 - VENTILE), 6-pole
- 29.7** Control panel (X6 - BED- TEIL)
- 29.8** Vacuum switch

30 Transition board, see illustration

- V1** H₂O
- V2** GREEN&CLEAN H1
- V3** Knee valve, right
- V4** Knee valve, left
- V5** Additional valve for third suction hose

13. Electronic circuit board

31.1 See illustration

The voltage supply must be provided by a safety transformer that complies with IEC 601-1 / VDE 0750 part 1 or DIN VDE 0551.

If the output contact X2 should switch electrical connections which are not of type BF, a cut-off relay that complies with IEC 601-1, section 20.2, insulation B-a must be used.

This cut-off relay is switched through the output contact X2, whereupon the operating voltage used must be the same as of the H1 system.



31.1 H1 mode setting

31.2 SIC1. fuse (750m A/T; wired in)

32.1 H1 mode setting, see illustration


Switch 1:

- 32.1** Standard
1 - two suction hoses
- 32.2** 3 suction hoses
0 - three suction hoses

Inspections • Cleaning and disinfection

Practice personnel, technicians


33

14. Inspections

We recommend an annual service check of the Hygiene System H1 by an authorised technician.

15. Cleaning and disinfection

The below mentioned parts can be cleaned and disinfected as follows:

34

33 Handpieces, see illustration

(This applies to METASYS handpieces only!)

The handpieces may be disassembled as follows: Switch off the unit and pull the handpieces from the connecting piece of the tube with a gentle twist.

Upon reassembly make sure that the O-rings are lubricated and in place, and slide back the handpieces with a twist.

Small suction tube:


33.1 Handpiece, small - order no: 50010047

Large suction tube:

33.2 Handpiece, large - order no: 50010048

33.3 Distributor, large - Order no: 50010001

Grease the o-rings with vaseline after cleaning!

35

34 Prefilter, see illustration

34.1 Filter rake - order no: 55010001

34.2 Filter cup - order no: 50010008

35 Lubricate the sealing lip with vaseline and place back properly!

The entire suction tubes can be disinfected with alcohol-free surface disinfectants, such as GREEN&CLEAN MK or GREEN&CLEAN SK (wipe disinfection). It is also possible to clean the suction tubes under running water.

36 Cleaning the cuspidor, see illustration

We recommend GREEN&CLEAN MB to clean the cuspidor.

Do not use household cleaners because these can lead to malfunctions of the suction system!

36

Refilling the H1 cartridge · Disposal

Practice personnel, technicians

16. Refilling the H1 cartridge

A refill set is included in the scope of delivery!

37 Open the closure of cartridge

- Attach key onto bayonet socket of cartridge.
- Turn the key counter-clockwise.

38 Remove cover lid of the cartridge

- Remove the cap.
- Pull out plunger with cartridge seal.

39 Attach the filling adapter




- Open the screw cap of the refill pouch.
- Attach the short end of the adapter to the refill pouch.
- Attach the long end of the adapter to the cartridge.

40 Refill the cartridge

- Fill up the cartridge with GREEN&CLEAN H1 until the maximum level is reached.
- Remove adapter after use and rinse with water.
- Close the refill pouch.
- Insert plunger with cartridge seal in filling opening.
- Attach the bayonet socket and turn the key clockwise.

17. Disposal

When the complete device is de-installed at the end of its service life it must be returned to the manufacturer or disposed of appropriately in accordance with the national regulations.

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