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Application and safety instructions for the MEISINGER

EN

Micro Screw System

Micro Screw System Basic

Screw System for safe Stabilization of Cortical Bone Graft Transplants
'Olsberger Konzept' by Prof. Dr. Fouad Khoury

Art.-No.: BMS00, BMSBA

The Bone Management® Micro Screw Systems are used for the safe fixation and stabilisation of cortical bone grafts. The special feature of these systems are the osteosynthesis screws, which have a diameter of 1,0 mm and 1,2 mm. These Micro Screws are made of surgical stainless steel, which gives them high stability despite their reduced diameter.

Indications for Use

The Bone Management® Micro Screw Systems are indicated for the treatment of bone fractures, especially for the fixation of transplanted bone blocks during the augmentation process in the oral cavity and maxilomandibular surgical field.

Note: Micro Screws are not intended to remain in the body permanently. After they have fulfilled their supportive function such as is the case after healing of a transplant, or healing of a fracture, for example, they need to be removed completely.

Utilization

The recommended speeds for the application and the maximum speeds are summarized in the instru- ment overview under „Content“. For proper, safe use, these must be observed.

Placement of the pilot hole through bone cylinder and recipient bone


For the application of the Micro Screws, a pilot hole through bone graft and recipient bone is necessary. This is done using a twist drill with a 0,2 mm smaller diameter:

- 1,0 mm Micro Screw $\triangleq$ twist drill MSPB1 or initial bur 202RF 008 ($\varnothing$ 0,8 mm)
- 1,2 mm Micro Screw $\triangleq$ twist drill MSPB2 or initial bur 202RF 010 ($\varnothing$ 1,0 mm)


The respective drilling should take place intermittently and under constant cooling with sterile, physiological saline solution. Since the Micro Screws have a self-tapping thread due to specially shaped thread flanks, no expansion drilling (gliding hole) to the outer diameter of the screws is necessary.

Selection of the screwdriver and picking up of the Micro Screws

Depending on the selected procedure (machine or manual screw insertion), the appropriate screw- driver must be selected. The screwdrivers MSSDE (screwdriver basic) or MSSDM (optionally with claw MSGSD) are used for manual insertion. For machine insertion with the angle piece, the MSSDW, screwdriver RA, is used.

Fixation of the bone graft

possible, fixation with at least 2 screws is recommended to prevent any rotation.

Removing the Micro Screws

Taking into account the patient-specific healing time and after clinically and radiographically confirmed complete bone healing as well as functional restoration, the Micro Screws are to be removed counter-clockwise using a screwdriver from the Micro Screw System. The removal of the Micro Screws is also necessary in the case of screw loosening, infection or weakening of the bone tissue due to the inserted bone screw.


Sterilisation of the Micro Screws

The instruments and Micro Screws are non-sterile when delivered and must be processed (cleaned / disinfected / sterilised) according to validated instructions for processing before use. For proper pro- cessing of the articles, please observe the instructions for processing (cleaning, disinfecting and sterilis- ing) of medical devices produced by Hager & Meisinger GmbH.

Instructions for disassembly and assembly of the screwdriver with claw

Screwdriver ①, Claw ②, Groove 1 ③, Groove 2 ④

1. Attach the claw ② to the screwdriver blade ①, as shown. Push the claw ② until it snaps into groove 1 ③, open position).

2. Now the screwdriver with claw is ready to use. To secure a screw, push the claw ② in the direction of the screwdriver blade over the screw until it snaps into groove 2 ④, closed position). For Disassembly push and remove the claw ② into the direction of the screwdriver blade.

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