Sebond Grip

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Sebond Grip

Product specification:

Sebond Grip is a light-curing universal priming opaquer for applying between Sebond Smart and light-curing powder or paste opaques currently available on the market.

Sebond Grip improves the chemical bonding between primers and opaques and contributes towards simple and secure adhesion.

Preparations and Adhesion:

The framework is prepared as usual after casting. The surface is sandblasted and blown off using oil-free compressed air.

Sebond Smart is used as a primer in accordance with the instructions. To achieve optimal adhesion, apply a thin coating of Sebond Grip using a brush with short, stiff bristles. If necessary, apply a second coat.

The flowability can be optimised by mixing on a mixing block. Care must be taken to ensure that a particularly thin layer is applied in the area of the retentions.

The material has very good light-curing properties that ensure optimal curing.

Application of opaquer:

Then apply the opaquer (powder/liquid or paste) according to the directions of the manufacturer and polymerise.

Sources of faults/avoidance of faults:

The surfaces of the framework must be fat-free and any remnants of polish must be cleaned off. Sebond Grip may only be applied in this layers. Do not use any intermediate polymerisation tools for polymerising.

Processing time of the mass:

1-3 min depending on the light.

Storage:

The material should be stored at 10-25°C. Close the syringe well immediately after use and protect from direct light.

Shelf life:

The maximum shelf life is printed on the label of each syringe. Do not use after the expiry date.

Guarantee:

Our technical instructions for use, whether provided orally, in writing or by practical instruction, are based on our own experience and can therefore be regarded only as guidelines. Our products undergo continuous development. We therefore reserve the right to make alterations.

Side effects:

Undesirable side effects of this medicinal product can be expected extremely rarely if it is handled and used correctly. However, immune reactions (e.g. allergy) or local discomfort cannot be completely excluded. If you become aware of undesirable side effects, even in doubtful cases, please notify us.

Contraindications / interactions:

This product must not be used or must be used only under the strict supervision of the treating doctor or dentist in the event of patient hypersensitivity to one of the constituents. In these cases, the composition of the medicinal product supplied by us is available on request. Known cross reactions or interactions between the medicinal product and other materials already present in the mouth must be taken into account by the dentist when using it.

Note:

Please give all the above information to the treating dentist if you are employing this medicinal product for a special model. When using it, please follow the existing safety data sheets.

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