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c/u de 10 ml

02018

**Precision impression material, C-silicone based,
medium consistency – medium bodied**

Fig. 1

Fig. 2

1. Preparation of the impression

Depending on the impression technique select an appropriate tray. For optimal adhesion of the impression to the tray apply a thin layer of **sili Adhesive Spray** or **sili Adhesive** liquid onto the tray and let dry for 60 sec. (Fig. 1).

2. Mixing

Extrude **silone®** and **cat^p universal FUTUR** from tubes in equal length onto mixing pad (Fig. 2). With dosage of liquid catalyst 1 cm **silone®** to 1 drop of catalyst. Mixing 30 sec. with spatula until homogeneous mix and uniform colour is achieved.

3. Disinfection

After removal rinse impression under lukewarm water. Disinfection for 15 min. in 2 % glutaraldehyde.

4. Model casting

Casting between 30 minutes and 24 hours after removal of the impression. Recommended materials: Dental plasters of class III and IV.

5. Galvanisation

Impressions can be electroplated with commercially available copper and silver plating baths.

6. Cleaning of trays

Remove set material with a blunt instrument. **sili Spray** or **sili liquid** film can be easily removed with **sili Spray Solvent** (free of FCHC). Work in well ventilated rooms. Clean and disinfect as usual.

Important working hints

- Increased temperatures accelerate, decreased temperatures retard the setting time.
- Do not intermix with addition curing silicones.
- Catalyst bottles or tubes should be closed tightly immediately after use.
- Don't leave any residual material in the patient's mouth.
- Avoid staining of clothing. Cured impression materials cannot be removed.
- Avoid contact of catalyst with eyes, may cause irritation. If necessary flush eyes with plenty of water and seek immediate medical care.
- Please follow the instructions of the safety data sheet!


Avoid contact with eyes and skin!
Wear suitable eye protection and gloves!

Further information:

Silicone based materials are proven a million times. On condition of a proper application, **undesired effects** are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to get tested the allergenic potential before the application of the material.

For use by trained specialists.

silone® / cat^p universal FUTUR:

Indications for use:

- Universal material for dentulous impression
- Partial and copping and inlay impressions

Technical Data:

DIN EN ISO 4823 – Type 2

- **Mixing ratio:**
silone® / cat^p universal FUTUR
1 cm silone® = 1.00 g /
1 cm cat^p = 0.05 g
silone® / sta-seal catⁱ
1 cm silone® = 1.10 g /
1 drop catⁱ = 0.02 g
- **Colour code:**
silone®: red-brown
cat^p universal FUTUR: blue
sta-seal catⁱ: light blue
- **Mixing time:**
approx. 30 sec.
- **Working time:**
approx. 1 min. 15 sec.*
- **Setting time:**
approx. 3 min. 45 sec.*
- **Time in mouth:**
approx. 2 min. 30 sec.
- **Strain in compression:**
9.50 %
- **Recovery from deformation:**
98.70 %
- **Linear dimensional change:**
0.45 %
- **Application:**
At 23 °C ± 2 °C /
73 °F ± 4 °F, 50 ± 5 % rel. humidity
- **Storage:**

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

Ordering information:

silone®	
Standard packing	02008
160 ml tube	
tube 3-pack	02009
3 x 160 ml tubes	
cat^p universal FUTUR	02015
5 x 35 ml paste	
sta-seal catⁱ	02018
5 x 10 ml bottles	

DETAX

GmbH & Co. KG
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
Telefon: 0 72 43/510-0 · Fax: 0 72 43/510-100
www.detax.de · post@detax.de

Made in Germany

04/2016

"C" silicone system

putty,
heavy bodied

kneadable,
very high
consistency

silaplast FUTUR

high final hardness
normal setting

Preliminary impression
material for corrections of
preparations, tray material,
overimpressions, wax-bite

Wash
material
light bodied

low viscosity
low con-
sistency

silasoft® N

normal setting

Correction impression
material for the double
impression and double
mix technique, relinings

Wash
material
light bodied

low viscosity
low con-
sistency

silasoft® S

normal setting

Correction impression
material for the double
impression and double
mix technique, relinings

Functional
heavy bodied

high con-
sistency and
plasticity

sta-seal f

medium sfinal hardness
plastic phase
extended setting time

Functional impressions,
functional border moulding

Catalyst
paste

Catalyst
liquid

cat^p / cat^f FUTUR

inert system
non-irritating

cat^p universal paste,
cat^f silaplast liquid,
silasoft® cat^f, liquid,

Special
instrument

Accessory

Deta-Cut

sili adhesive

sili solvent

Sterilizable, curved, for
cutting the preliminary
impression

adhesive for C-silicones
removal of the adhesive
and cleaning of trays

All DETAX condensation-
curing silicones are
portioned with dispensing
aids (tray, measuring bowl,
measuring scale etc.) and
kneaded by hand or mixed
homogeneously on a
mixing pad.

Ordering information:

silaplast FUTUR
Standard packing **02001**
jar, 900 ml

4-pack **02002**
4 x 900 ml

clinic packing **02003**
1 x 5400 ml

refill bag **02681**
900 ml

silasoft® N
standard
tube-packing **02004**
160 ml

4-pack tubes **02005**
4 x 160 ml

standard
bottle-packing **02320**
160 ml

4-pack bottles **02321**
4 x 160 ml

silasoft® S
standard
tube-packing **02275**
160 ml

4-pack tubes **02276**
4 x 160 ml

sta-seal f
standard-packing **02012**
measuring syringes
3 x 80 ml
catalyst liquid bottle 10 ml

large-packing **02013**
measuring syringes 9 x 80 ml
catalyst liquid bottles 3 x
10 ml

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