Sockelgerät

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Study model base unit

by Prof. Hinz

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Instructions for use

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Mounting of jaw models

The study model base unit by Prof. Dr. Hinz is suitable for locating three-dimensionally oriented jaw models. A uniform shape and height of the jaw models can be determined when using the study model base unit and base trays. The models can be fixed on the base unit with the webs supplied.

The plastic base trays are in conformity with international standards and are different in shape for maxilla and mandibula.

The base trays are available in two sizes for orthodontic treatment and periodontology.

Observe the following when using this unit:

The **median raphe plane** is represented by a sagittal recess that is lasered into the upper part of the base unit. This recess serves as an optical reference to fix the centre of the upper jaw.

After the median raphe plane has been adjusted the **tuber plane** is automatically at right angle.

The **occlusal plane** running horizontally along the base plate of the base unit is obtained from the bicuspid cusps and molars resting on the upper part of base unit.

1. Preparation

Pour the impressions with a flat base and grind the models all round. Trimming of the model base shall be slightly conical. Carve the teeth and vestibular fold. Briefly water the model base and mark the median raphe plane on the maxillary model.

Place the maxillary tray into the profiles of the vertical back wall of the base unit and lower it down to the base plate. Make sure the model is movable in all directions and check the correct model height using the base unit slide.

2. Aligning the maxillary model

Pour plaster into the maxillary tray and allow the plaster to set for a short time. Then place the maxillary model onto the plaster and lower the slide as far as it will go.

Make sure the slide is pushed downwards in the middle and not at the side. The pushing hand should run along the back wall.


Align the maxillary model to the maxillary centre using the base unit slide.


Now, determine the occlusal plane by placing the bicuspid and molar cusps on the upper part of base unit aligning the distal surfaces of the canine teeth to the front edge of the base unit.

Remove the slide and allow the plaster to continue setting.

3. Locating the mandibular model

After plaster curing clean the tray under running water. Then, put together the maxillary model and mandibular model with the situation bite and secure this status using a rubber band.

Then, insert the mandibular tray into the base unit as described on step 2 and fill it with plaster.

Completion

Attach the base unit slide with the fixed models at the back wall of the base unit and slowly lower it before the models are lowered as far as it will go.

Allow to set the plaster for a short time, then remove the rubber band and take off the slide with the maxillary tray.

After that take the located mandibula from the slide and remove any excessive plaster.

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