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Light-curing resin based on (meth)acrylate, for the generative fabrication of dental models for DLP printers with UV-LED 385 nm

## Suitable for the following DLP-Printer

| Asiga<br>385 nm               | MiiCraft<br>385 nm | Rapidshape<br>385 nm | W2P<br>385 nm |
|-------------------------------|--------------------|----------------------|---------------|
| MAX UV<br>Pro2 UV<br>Pico2 UV | MiiCraft y-Serie   | DII-Serie            | Solflex Serie |

Printers may only be operated using material parameters authorised by DETAX.

## Processing

- The properties of the final product depend, among other things, on post-processing. Therefore it must be ensured that the light unit is in an orderly condition and that the models are completely cured (observe process description on page 2).
- After storage, the material in the bottle should be shaken intensively and homogenized with a bottle roller before use.
- Processing temperature 23 °C ± 2 °C.

## Safety Information

- Only to be used by trained specialists for the specified area of application.
- Product contains (meth)acrylates and phosphine oxides. The ingredients of **FREEPRINT® model basic** may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further.
- Avoid direct contact with the liquid material and the components before post-curing. Irritating to eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- In case of contact with skin, wash immediately with plenty of soap and water.
- The processing instructions and precautionary measures must be observed.
- Refer to the relevant safety data sheet for hazard and safety information.

## Notes

- DETAX shall not be held liable for any damage caused by misuse.
- Always keep container tightly closed, carefully close immediately after each use.
- Read and understand the safety data sheet!!

## Storage

FREEPRINT® model basic is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

## Adverse effects

Product may cause allergic reactions.

## Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations..

## Area of application:

3D Print working and presentation models for the entire spectrum of models in dentistry

## Processing:

at 23 °C ± 2 °C

## Storage

15 °C  
59 °F

28 °C  
82 °F

## Ordering information:

### FREEPRINT® model basic 385

1.000 g bottle  
light grey **02111**

### FREEPRINT® model 385

1.000 g bottle  
sand **03778**  
ivory **03780**  
grey **03782**  
caramel **02850**

### FREEPRINT® model 405

1.000 g bottle  
sand **03065**  
ivory **03779**  
grey **03781**

### FREEPRINT® gingiva 385

500 g bottle, gingiva **02820**  
1.000 g bottle, gingiva **02843**

### FREEPRINT® model T 385

1.000 g bottle blue **02332**

### FREEPRINT® model 2.0 385

1.000 g bottle  
sand **02128**  
light-grey **02099**  
grey **02177**  
caramel **02850**  
white **02148**

500 g bottle  
caramel **02844**

## Manufacturing process

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

## Construction process

Generation of a Print Job complying with machine and material parameters

## Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

## Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity  $\geq 98\%$ ) for 3 min. in an ultrasonic bath.

## Cleaning

Then thoroughly clean the openings, cavities and gap areas, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

## Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity  $\geq 98\%$ ) for 3 min. in an ultrasonic bath. Prior to drying, check the openings and additional bore holes for residues.

## Drying

Heat the construction components for 30 min. in an oven to approx. 40 °C to remove the solvent residues from the cleaning process.

## Post-exposure

Post-exposure is performed with a xenon photoflash unit (e.g. Otofash G171) with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.