FREEPRINT®

3D PRINTING MATERIAL

INDIVIDUAL FUNCTIONAL TRY-INS

DETAX

DE	Gebrauchsanweisung.....	3
EN	Instructions for use	6
FR	Mode d'emploi	9
ES	Modo de empleo.....	12
IT	Istruzioni per l'uso.....	15
DA	Brugsanvisning.....	18
EL	Οδηγίες χρήσης.....	21
HR	Upute za uporabu	24
LV	Lietošanas instrukcija.....	27
LT	Naudojimo instrukcija	30
NL	Gebruiksaanwijzing.....	33
NO	Brugsanvisning.....	36
PL	Instrukcja stosowania	39

PT	Manual de instruções.....	42
RO	Instrucțiuni de utilizare	45
SV	Brugsanvisning	48
SK	Návod na použitie.....	51
SL	Navodila za uporabo	54
CS	Návod k použití.....	57
HU	Használati útmutató.....	60
RU	Инструкция по применению	63
TR	Kullanım kılavuzu	66
FI	Käyttöohje.....	69
ET	Kasutusjuhend	72
BG	Инструкция за употреба.....	75

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Resin for dental 3D printing

INDICATION

Individual functional try-ins of complete dentures

PATIENT TARGET GROUP

Persons being treated in the context of a dental procedure.

INTENDED USERS

Dentist, dental technician

SUITABLE FOR THE FOLLOWING DLP PRINTERS

"Annex 1" (enclosed separately)

PROCESSING

- ▶ The properties of the final product depend, among other things, on the finishing process. Correct post-exposure is important for biocompatibility. Therefore, it must be ensured that the exposure unit is in proper condition and the molded parts are completely cured (refer to the process description).
- ▶ Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- ▶ A layer thickness of 100 µm is suitable for the printing process.
- ▶ Polish the surface mechanically if necessary. Pre-polishing is done with rotating brushes and pre-polishing paste, high-gloss polishing with buffing and high-gloss polishing agents for resin.
- ▶ Processing temperature 23 °C ± 2 °C.

SAFETY INFORMATION

- ▶ Only for the specified intended use by trained specialists.
- ▶ Avoid direct contact with the liquid material and the components before post-curing, in particular in pregnant / breastfeeding women. Irritating to airways, eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- ▶ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ In case of contact with skin, wash immediately with plenty of soap and water.
- ▶ Biocompatibility is guaranteed only after complete polymerization.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

NOTES

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly closed, carefully close immediately after each use.
- ▶ Read and understand the safety data sheet

For users and/or patients:

Any serious incident that has occurred in connection with this product must be reported immediately to incident@detax.com and to the competent public authority of the Member State in which the user and/or patient resides.

STORAGE

- Store **FREEPRINT® tryin** in a dry place (at 15 °C – 28 °C) and away from light. Even slight exposure to light can trigger polymerization.
- To protect against contamination, cover the material in the tray with the lid or a glass plate.

CONTRAINDICATION

Contains (meth)acrylates and phosphine oxides.

Ingredients of **FREEPRINT® tryin** may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further. Apply **FREEPRINT® tryin** intraorally only when fully polymerized.

SIDE EFFECTS

Product may cause allergic reactions.

DISPOSAL

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Storage:

Processing:

at 23 °C ± 2 °C

Manufacturing process

Data processing and creation of the support structure according to the specifications of the CAD software manufacturer

Building process

Creation of a print job in compliance with the machine and material parameters.

Finishing process

After raising the platform, a dripping time of approx. 10 minutes is recommended. If possible, the finishing should take place immediately after the building process.

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DETAX GmbH