MUPPY

Perle starr

D
GB
F
I
E

Patienteninformation

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D
GB
F
I
E

Patient information

MUPPY if there is a muscle-weak tongue

MUPPY with a pearl is used in tongue disorders, atypical swallowing, speaking difficulties and salivary flow caused by an open mouth position.

This special **MUPPY** is mostly used by children with Down syndrome and cerebral palsy so the muscle-weak tongue can be guided to the palate by the pearl attached to the wire.

MUPPY-rules

Studies have shown that more than 90 % of children give up damaging sucking habits with the help of **MUPPY**. Your support and a few rules help your child:

- 1 **MUPPY** is a substitute for thumb and dummy. The substitution must be voluntary.
- 2 **MUPPY** must always be worn at night.
- 3 To get used to it easier children should also wear **MUPPY** a few hours during the day (when playing or watching television).
- 4 As a rule, weaning is successful within a week. However, **MUPPY** should also be worn for 6–8 weeks after successful weaning of the dummy ("recidivism").
- 5 **MUPPY** can be stored in the hygiene box.
- 6 Do not forget to regularly brush the teeth and **MUPPY**!

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The movable pearl attached to a wire guides the tongue to the correct position to the palate, in order to bring about the mouth closure afterwards. At the same time, the tongue is guided upwards and to the front to prevent malpositioned teeth.

If your child is having speech therapy treatment, you should take note of the following:

- Follow your speech therapist's individual wearing and exercise instructions exactly.
- Support your child during this stage, as he/she has to feel the importance of mouth closure in order to take the exercises seriously.
- Make your child aware of the "rules of play" and supervise the necessary exercises.

Instructions for use

Before using **MUPPY** for the first time, clean it thoroughly under warm running water. Check **MUPPY** regularly for places that have been bitten through! Replace **MUPPY** when necessary.

Cleaning instructions

Clean thoroughly with a toothbrush under warm running water. Alternatively, you can also use the OROFAN® dental cleaning tablets. **MUPPY** is neither dishwasher-safe nor microwave-safe. Only clean the hygienic **MUPPY**-Box with a damp cloth.

MUPPY designs

MUPPY pearl is available in two different sizes:

- with violet ring for deciduous dentition (size I)
- with blue ring for mixed dentition (size II)

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