core-X[®] flow

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core-X® flow

Dual Cure Core Build-Up Material and Cement for Endodontic Posts

CAUTION: For dental use only.
In the USA: Rx only

Content	Page
1 Product description	2
2 General safety notes	3
3 Step-by-step instructions	4
4 Hygiene	7
5 Lot number and expiration date	7

1 Product description

core-X® flow consists of two-components, base and catalyst, which when mixed form a dual cured, highly filled, composite resin core build-up and post-cementation material.

core-X® flow uses a biocompatible urethane resin and is supplied in a tooth-colored shade which is ideal in situations where esthetics and show-through of the core are of primary concern.

1.1 Delivery forms

core-X® flow is available in a convenient dual-barreled syringe. This delivery system simplifies direct, precise intraoral delivery and minimizes product waste.

1.2 Composition

- Urethane Dimethacrylate
- Di- & Tri-functional Methacrylates
- Barium Boron Fluoroaluminosilicate Glass
- Camphorquinone (CQ) Photoinitiator
- Photoaccelerators
- Silicon Dioxide
- Benzoyl Peroxide

1.3 Indications

- Vital or non-vital tooth core build-up (replacement of existing restorations and/or lost tooth structure) as a base prior to fabricating an indirect restoration.
- Cementation of endodontic fiber-posts.

1.4 Contraindication

- **core-X® flow** is contraindicated in patients with a history of allergic reaction to methacrylate resins or any other of the components.

1.5 Compatible adhesives

core-X® flow is to be used after the application of a suitable dentin and enamel adhesive.

core-X® flow is compatible with numerous dentin/enamel adhesive systems including all Dentsply Sirona adhesives designed for use with dual/self cured composite restoratives (see complete Directions for Use of selected adhesive). Other dentin/enamel adhesive systems may be used at the discretion and sole responsibility of the dental practitioner.

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2.1 Warnings

- **core-X® flow** contains polymerizable methacrylates which may irritate to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, redish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitisation or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitisation of mucosa persists, seek medical attention.

2.2 Precautions

1. This product is intended to be used only as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.
2. Contact with saliva and blood during composite placement may cause failure of the restoration. Use of rubber dam or adequate isolation is recommended.
3. Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
4. Syringe should be tightly closed by replacing the original cap immediately after use.
5. **core-X® flow** is intended to be used as a dual cure product. To obtain optimal curing follow the step-by-step instructions. Self cured, dual cured and light cure methods are integrated into the procedural steps.
6. **core-X® flow** should be mixed in equal volumes directly through the supplied mix tips for optimum performance. Variations from this 1:1 ratio may adversely affect working time and compromise final physical properties. Dispensing and mixing other than through the supplied mixing tips is not recommended.
7. **core-X® flow** should extrude easily. **DO NOT USE EXCESSIVE FORCE.** Excessive pressure may result in unanticipated extrusion of the material or cause syringe rupture.
8. **core-X® flow** when used for post-cementation has a working time of only 40 seconds! Prepare and organize the necessary materials and instruments to allow a fast insertion of the post.

9. Interactions

- Eugenol-containing dental materials should not be used with **core-X® flow** because they may interfere with hardening and cause softening of the polymeric components of the materials.
- Variable in-vitro data exist regarding use of light cured-only adhesives without Self Cure Activator in conjunction with dual cured or self cured composite restorative materials such as **core-X® flow** in limited or no light curing applications. Chemical/product incompatibility may adversely affect product efficacy, leading to premature restoration failure.
- As with any dual cure resin system, the use of a dual cure adhesive system can shorten working time. This effect should be investigated in the laboratory prior to clinical use.
- The **core-X® flow** is a light cured material. Therefore, it should be protected from ambient light. Proceed immediately once material has been placed.

2.3 Storage conditions

core-X® flow syringes should be stored with the original cap only, at temperatures between 10 °C/50 °F and 24 °C/75 °F. Keep out of direct sunlight and protect from moisture. Refrigerated storage is acceptable when not in use. Allow material to reach room temperature prior to use. Do not freeze. Do not use after expiration date.

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2.4 Adverse reactions

Product may irritate the eyes and skin. **Eye contact:** Irritation and possible corneal damage. **Skin contact:** Irritation or possible allergic response. Reddish rashes may be seen on the skin. **Mucous Membranes:** Inflammation (see 2.1 Warnings).

3 Step-by-step instructions

Working and Setting Times for Core Build-Up and Endodontic Post Cementation

Working Time at body temperature (37 °C) – for manipulation/post seating after intraoral placement	40 seconds
Working Time at room temperature (22 °C) – for application after start of mix	1 minute 30 seconds
Setting Time Light Cure	20 seconds¹ (3 mm depth)
Setting Time Self Cure (catalyst initiated)	2-3 minutes

3.1. Core Build-Up

3.1.1 Cleaning and Preparation of Tooth

- Finish preparation by removing existing restorations and/or caries.
- Clean freshly cut enamel and dentin with water spray and air dry. Do not desiccate dentin.
- Place pin, post or matrix as needed.

3.1.2 Pulp protection, tooth conditioning/dentin pretreatment, adhesive application

Refer to applicable manufacturer's directions for pulp protection, tooth conditioning and/or adhesive application. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of **core-X® flow**.

3.1.3 Dispensing and mixing core-X® flow dual-barreled syringe

Danger of injury due to excessive force.

1. Apply slow and steady pressure on the syringe plunger.
2. Do not use excessive force – syringe rupture may result.

1. If refrigerated prior to use, remove from refrigerator and bring to room temperature.
2. Remove syringe cap and dispense and discard a small amount of material from the dual-barreled syringe. Be sure material is flowing freely from both ports. Holding syringe vertically, carefully wipe away excess so base and catalyst do not cross-contaminate and cause obstruction of the ports. Save syringe cap for replacement following use.
3. Install a mixing tip on the cartridge by lining up the v-shaped notch on the outside of the mix tip with the v-shape notch on the syringe flange. Turn colored mix tip cap 90 degrees in a clockwise direction to lock in place on syringe.
4. Gently depress syringe plungers to begin the flow of material. **DO NOT USE EXCESSIVE FORCE.** If force is encountered, remove syringe from operating field, remove and discard mix tip. Check for obstruction and confirm material is flowing from both syringe barrels. Wipe barrels and install new mix tip as outlined above. Dispense a small amount through the mix tip onto a mixing pad and discard.
5. Proceed immediately to placement, **Step 3.1.4. Placement of core-X® flow.**

3.1.4 Placement of core-X® flow

core-X® flow produces a syringable mixed consistency for direct placement using the dual-barreled syringe and a banded matrix retainer or clear crown form. Alternatively, mixed material may be expressed onto a clean paper pad and manually placed with the operator's placement instrument(s) of choice. Working time of mixed material at room temperature (22 °C/72 °F) is approximately 1'30" (one minute thirty seconds) when protected from ambient light. Elevated oral temperatures and/or contact with dual cure adhesives may shorten work time (**see 2.2 Precautions**).

Any of the following placement techniques are acceptable. The preferred method is to use both visible light and chemical self cure (**see Step 3.1.4.3 Dual Cure**).

1. Visible Light Cure

Place mixed material directly in 2-3 mm increments and light cure increment for 20 seconds¹.

2. Chemical Self Cure

Place large increment(s) or load mixed material into a clear crown form and seat onto the preparation. Allow to self cure approximately 2-3 minutes. Visible light curing¹ of external surface is optional, but advisable.

3. Dual Cure

After placing the first 2-3 mm increment, and visible light curing for 20 seconds¹, large increments may then be placed or material loaded into a clear crown form, seated onto the preparation and allowed to self cure. Visible light curing¹ of external surface is optional, but advisable.

3.1.5 Finishing the core-X® flow

After initial self cure set and removal of matrix, a surface light cure of the material for an additional 20 seconds¹ is recommended. Initial crown preparation (gross reduction) may begin immediately using rotary instruments. If used per **step 3.1.4.2 (Chemical Self Cure)** or **3.1.4.3 (Dual Cure)**, utilizing chemical self cure, greatest bulk hardness is obtained 10 minutes after start of mixing at which time the final preparation may be completed.

3.2 Endodontic Post Cementation

Tooth Conditioner/Dentin Pretreatment/Adhesive application.

Refer to applicable manufacturer's directions for tooth conditioning and/or adhesive application. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of **core-X® flow**, as outlined below.

Note: If core build-up will immediately follow post cementation, treat all available tooth structure beyond post space with adhesive for core build-up bonding.

1. Preparation of post:

Treat surface of post according to adhesive manufacturer's and/or dental laboratory's instructions, i.e., adhesive application.

2. Dispensing and mixing **core-X® flow** dual-barreled syringe.

1. Prepare syringe, bleed material and install mix tip as outlined above **3.1.3. (Dispensing and mixing core-X® flow dual-barreled syringe 1-4).**

3. Proceed immediately with placement. Dispense the same material directly from the syringe onto the post surface and to the orifice of the post hole preparation.

4. Immediately place post into final position and stabilize.

Important note: Work time of applied cement within a post space pretreated with dual cure adhesive is 40 seconds maximum from application. Post must be fully seated and stabilized within 40" (**see 2.2 Precautions, 9. Interactions**).

Inadequate insertion of post due to short working time.

1. At body temperature remaining working time is 40 seconds from start of application.
2. Prepare and organize necessary material and instruments to allow a fast insertion of the post.

Once stabilized, light cure, using a visible light curing unit designed to cure CQ initiated methacrylates (spectral output including 470 nm), with a minimum output of 550 mW/cm². Cure all accessible areas of the post for 20 seconds.

Inadequate polymerization due to insufficient curing.

1. Check compatibility of curing light.
2. Check curing cycle.
3. Check curing output before each procedure.

5. Proceed with the core build-up as described in **chapter 3.1 (Core Build-Up)**.

4 Hygiene

4.1 Cleaning the syringe

For dual-barreled syringe, remove used mix tip and discard appropriately. Replace original syringe cap prior to storing. The dual-barreled syringe may be cleaned by scrubbing with a disposable towel soaked with hot water and soap or detergent.

4.2 Disinfecting the syringe

Disinfect the syringe with a hospital-level, tuberculocidal disinfectant solution according to national/local regulations. Iodophors, sodium hypochlorite (5.25%), chlorine dioxide and dual or synergized quaternary ammoniums are approved disinfectants. Some phenolic-based agents and iodophor-based products may cause surface staining. The disinfectant manufacturer's directions should be followed properly for optimum results. Water-based disinfectant solutions are preferred.

5 Lot number () and expiration date ()

1. Do not use after expiration date.
ISO standard is used: "YYYY-MM" or "YYYY-MM-DD".
2. The following numbers should be quoted in all correspondence:
 - Reorder number
 - Lot number
 - Expiration date

© Dentsply Sirona 2007-11-28

[These Directions for Use are based on Master Version 01]

¹ Use suitable visible light curing unit designed to cure materials containing camphorquinone (CQ) initiator, i.e. spectral output containing 470 nm. Minimum light output must be at least 550 mW/cm² exposure. Refer to curing light manufacturer's recommendations for compatibility and curing recommendations.

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Manufactured by
Dentsply Caulk
38 West Clarke Avenue
Milford, DE 19963
USA
www.dentsplysirona.com

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CAUTION: For dental use only.

USA: Rx only.

Content	Page
1 Product description	2
2 Safety notes	3
3 Step-by-step instructions	5
4 Hygiene	12
5 Lot number, expiration date and disposal	13

1 Product description

Prime&Bond® XP adhesive is a universal etch-and-rinse one-bottle dental adhesive designed to bond resin-based materials to enamel and dentin as well as to metals and ceramic.

Prime&Bond® XP adhesive stands for eXtra Performance due to high bond strength to enamel and dentin, easy and comfortable application and a high degree of technique robustness.

Prime&Bond® XP adhesive, when used on its own, is designed to be used with light cured resin-based restorative materials and cements.

Prime&Bond® XP adhesive, when mixed with Self Cure Activator (available separately), is designed to be used with dual cure/self cure resin materials.

1.1 Indications

Self-priming adhesive for:

- Direct, light cured composite and compomer restoratives.
- Indirect restorations; light cured, resin cemented veneers.
- Composite/compomer, ceramic and amalgam repairs.
- Cavity varnish for use with fresh amalgam.
- Direct, dual cure or self cure composite restorations and core build-ups.
- Indirect restorations; dual cured and self cured resin cemented inlays, onlays, crown and bridge retainers.
- Dual cured and self cured resin cemented endodontic post cementation.
- Adhesive bonding of direct amalgam restorations.

1.2 Contraindications

- Use with patients who have a history of severe allergic reactions to acrylate or methacrylate resins or any of the other components.
- Direct application to dental pulp (direct pulp capping).

1.3 Delivery forms¹

- Tri-layer bottle with flip-top-cap
- Unit Dose container

1.4 Composition

Prime&Bond® XP:

- Carboxylic acid modified dimethacrylate (TCB resin)
- Phosphoric acid modified acrylate resin (PENTA)
- Urethane dimethacrylate (UDMA)
- Triethyleneglycol dimethacrylate (TEGDMA)
- 2-Hydroxyethylmethacrylate (HEMA)
- Butylated benzenediol (stabilizer)
- Ethyl-4(dimethylamino)benzoate
- Camphorquinone
- Functionalised amorphous silica
- Tertiary butanol

1.5 Compatible tooth conditioners

Prime&Bond® XP adhesive is used following phosphoric acid conditioning with conventional dentin/enamel conditioners including Dentsply Sirona 34% and 36% conditioners (see Directions for Use of selected conditioner).

1.6 Compatible restorative materials

Prime&Bond® XP adhesive is used in conjunction with suitable restorative or cement material and is chemically compatible with conventional methacrylate-based composite restorative materials including Dentsply Sirona visible light cured composite restorative or cement materials (see Directions for Use of selected restorative or cement material).

2 Safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Directions for Use.

Safety alert symbol.

- This is the safety alert symbol. It is used to alert you to potential personal injury hazards.
- Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

The material contains polymerizable acrylate and methacrylate monomers which may be irritating to skin, eyes and oral mucosa. May cause sensitization by skin contact and allergic contact dermatitis in susceptible persons.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes rinse with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, red-dish rashes may be seen on the skin. If contact with skin occurs, remove material with cotton and alcohol and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, remove material from the tissues. Flush mucosa with plenty of water and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.

¹ Some delivery forms may not be available in all countries.

2.2 Precautions

This product is intended to be used only as specifically outlined in these Directions for Use. Any use of this product inconsistent with these Directions for Use is at the discretion and sole responsibility of the dental practitioner.

- Use protective measures for the dental team and patients such as glasses and rubber dam in accordance with local best practice.
- Contact with saliva, blood and sulcus fluid during application may cause failure of the restoration. Use adequate isolation such as rubber dam.
- Tightly close bottle immediately after use.
- Use in a well-ventilated area. Avoid inhaling vapor.
- Flammable: Keep away from sources of ignition. Take precautionary measures against static discharges.
- Unit Dose containers and Applicator Tips are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- The bottles cannot be reprocessed. To prevent bottles from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the bottle is handled offside the dental unit with clean/disinfected gloves. During treatment clinicians with patient contact should not handle bottle.
- Interactions:
 - Do not use eugenol- and hydrogen peroxide-containing materials in conjunction with this product since they may interfere with hardening of the product.
 - Avoid the product saturating the gingival retraction cord. If it soaks into the cord, it may set hard and bond the cord to the underlying tooth surface making removal difficult.
 - If mineral-impregnated (e.g., ferric compounds) retraction cords and/or hemostatic solutions are used in conjunction with adhesive procedures, marginal seal may be adversely affected, allowing microleakage, subsurface staining and/or restoration failure. If gingival retraction is necessary, use of plain, non-impregnated cord is recommended.
 - Prime&Bond® XP adhesive is a light cured material. Protect from ambient light.
 - Insufficient data exist to support the use of desensitizing agents and/or cavity cleansing agents with Prime&Bond® XP adhesive. Recommended pre-treatment is pumice or prophylaxis paste with a rubber cup (see Step-by-step instructions).
 - If H₂O₂ has been used to clean the cavity, proper rinsing is essential. Higher concentration H₂O₂ may interfere with the setting of polymerizable material and should not be used prior to the application of Prime&Bond® XP adhesive.
 - Various in-vitro data exist regarding use of light cured only adhesives such as Prime&Bond® XP adhesive without Self Cure Activator in conjunction with dual cured or self cured resin restorative or cement materials in limited or no light curing applications. Chemical/product incompatibility may adversely affect product efficacy, leading to premature restoration failure.

2.3 Adverse reactions

- Eye contact: Irritation and possible corneal damage.
- Skin contact: Sensitization. Irritation or possible allergic response. Reddish rashes may be seen on the skin.
- Contact with mucous membranes: Inflammation (see Warnings).

2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Keep the product out of direct sunlight and store in a well-ventilated place.
- Store the product at temperatures between 10 °C and 28 °C (50 °F and 82 °F). Use the product at room temperature.
- Protect from moisture.
- Do not use after expiration date.

3 Step-by-step instructions

3.1 Direct restoration (light cured composite resin and compomers)

3.1.1 Cleaning

1. Clean uninstrumented enamel and dentin with a rubber cup and pumice or a non-fluoride cleaning paste such as Nupro® prophylaxis paste.
2. Wash thoroughly with water spray and air dry.
3. Clean freshly instrumented enamel and dentin with water spray, then air dry.

3.1.2 Pulp protection

1. In deep cavities cover the dentin close to the pulp (less than 1 mm) with a hard-setting calcium hydroxide liner (Dycal® calcium hydroxide liner, see Directions for Use) leaving the rest of the cavity surface free for bonding.

3.1.3 Tooth conditioning/dentin pre-treatment

When used as a bonding agent for direct restorative materials, it is recommended to follow the total etch technique described in this chapter. Prior to use, please see the Directions for Use of the respective conditioner.

Application of 34%-36% phosphoric acid conditioner

1. After application of rubber dam or other suitable isolation technique, apply the respective conditioner.
2. Apply the respective conditioner to the cavity surfaces starting at the enamel margins.
3. For best results, condition enamel for at least 15 seconds and dentin for 15 seconds.

Rinsing and blot drying

1. Remove conditioner with aspirator tube and/or vigorous water spray and rinse conditioned areas thoroughly for at least 15 seconds.
2. Remove rinsing water completely by blowing gently with an air syringe or by blot drying with a cotton pellet. Do not desiccate dentin.
3. Proceed immediately to application of Prime&Bond® XP adhesive.

Contamination.

Restoration failure.

1. Once the surfaces have been properly treated, they must be kept uncontaminated.
2. If salivary contamination occurs, thoroughly clean with vigorous water spray, dry, and repeat conditioning procedure for 5 seconds only. Rinse and dry as described above.

3.1.4 Dosage

Flip-top Bottle:

1. Pull off originality seal.
2. To open, take the bottle in the hand, place a thumb in the thumb recess and push the cap up until it remains in the open position.
3. Hold the bottle upside-down in a vertical position and dispense 1-2 drops of adhesive into a Clixdish™ Mixing Well or standard dappen dish or mixing well.
4. If necessary clean dropper of the bottle using a soft paper tissue.
5. Carefully close the bottle again by applying pressure to the lever from above. The cap snaps into place when pressed closed with a perceptible “click”.

**Technique Tip: CliXdish™.**

Material in a closed CliXdish™ Mixing Well will remain usable for **up to 30 minutes**. Material dispensed to a standard mixing well must be used immediately.

Unit Dose container for single patient use

1. Grasp container at each end, or insert into holder, placing thumb along center score.
2. Firmly apply pressure until container separates.
3. Holder may be placed on the tabletop or held between fingers for convenience.
4. Insert disposable applicator into opening to saturate Applicator Tip.

3.1.5 Application and curing

1. Using the Applicator Tip, apply the product wetting all cavity surfaces uniformly. Avoid pooling. This surface should remain undisturbed for 20 seconds.
2. Evaporate solvent by thoroughly drying with clean, dry air from a dental syringe for at least 5 seconds. Use a vacuum aspirator to prevent spatter of the adhesive onto the mucosa. Surface should have a uniform glossy appearance. If not, repeat application and air dry.
3. Cure Prime&Bond® XP adhesive for 10 seconds² using a curing light.

**Insufficient curing.**

Inadequate polymerization.

1. Check compatibility of curing light.
2. Check curing cycle.
3. Check curing irradiance.

4. Place restorative material over the cured Prime&Bond® XP adhesive as per restorative material manufacturer's Directions for Use.

3.2 Light cured composite, ceramic and amalgam repairs**3.2.1 Preparation**

1. Roughen and clean fractured restoration surface with a fine diamond bur and create mechanical retention, if possible.

3.2.2 Tooth conditioning/dentin pre-treatment

1. Etch fractured tooth surface with the respective conditioner as outlined in 3.1.3 Tooth conditioning/dentin pre-treatment.

3.2.3 Treatment of restoration

1. Etch fractured restoration surface with the respective conditioner for at least 15 seconds as outlined in 3.1.3 Tooth conditioning/dentin pre-treatment.

3.2.4 Application and curing

1. Apply adhesive to the fractured restoration surfaces and light cure² the adhesive [3.1.5 Application and curing].

3.2.5 Completion

1. Complete repair with placement and curing of desired shade(s) of light cured composite restorative following manufacturer's Directions for Use.

² Check curing light for minimum light output of at least 800 mW/cm². Cure for at least 20 seconds if light output is between 550 and 800 mW/cm².

3.3 Cementation of light transmitting indirect restoration with a light cured resin cement

3.3.1 Treatment of restoration

1. Treat internal bonding surface of restoration according to manufacturer's or dental laboratory's Directions for Use, i.e. etching, mechanical roughening and/or Silane Coupling Agent application.
2. Apply a single coat of Prime&Bond® XP adhesive to the internal bonding surface of the restoration. Immediately air dry for 5 seconds.

Light curing of applied Prime&Bond® XP adhesive is not necessary.

3.3.2 Treatment of tooth preparation

1. Clean and condition tooth preparation, apply and light cure the material as described in 3.1 Direct restoration (light cured composite resin and compomers).

3.3.3 Cementation

1. Prepare and apply light cured resin cement according to manufacturer's Directions for Use.

3.4 Cavity varnish for use with fresh amalgam

When used as a cavity varnish, Prime&Bond® XP adhesive is not an amalgam adhesive.

3.4.1 Preparation

1. Finish preparation by removing existing restorations and/or caries.

3.4.2 Cleaning

Clean preparation and place a hard-setting calcium hydroxide liner (Dycal® calcium hydroxide liner, see Directions for Use) if needed, according to:

1. Cleaning [3.1.1 Cleaning].
2. Pulp protection [3.1.2 Pulp protection].

3.4.3 Tooth conditioning/dentin pre-treatment

1. Rinse and carefully air dry cavity preparation, but do not desiccate exposed dentin.
2. Condition enamel and dentin following section 3.1.3 Tooth conditioning/dentin pre-treatment.

3.4.4 Application

1. Apply and light cure Prime&Bond® XP adhesive as described for direct restorations, see section 3.1.3 Tooth conditioning/dentin pretreatment.

3.4.5 Completion

1. Place and condense amalgam (e.g. Dispersalloy® Dispersed Phase Alloy, see Directions for Use) as per manufacturer's directions.

3.5 Direct restorations (dual cure or self cure composite/core build-up)

3.5.1 Preparation

1. Finish preparation by removing existing restorations and/or caries.

3.5.2 Cleaning

Clean preparation and place a hard-setting calcium hydroxide liner (Dycal® calcium hydroxide liner, see Directions for Use) if needed, according to:

1. Cleaning [3.1.1 Cleaning].
2. Pulp protection [3.1.2 Pulp protection].
3. Place pin, post or matrix as needed.

3.5.3 Tooth conditioning/dentin pre-treatment

1. When used as a bonding agent for composite materials, it is recommended to follow the total etch technique described in section 3.1.3 Tooth conditioning/dentin pre-treatment.

3.5.4 Application

1. Place 1-2 drops of adhesive into a clean CliXdish™ Mixing Well or standard plastic mixing well.
2. Close bottle promptly.
3. Place an equal number of drops of Self Cure Activator into the same mixing well.
4. Replace cap promptly.
5. Mix contents for 1-2 seconds with a clean, unused brush tip.
6. Using the disposable brush supplied, apply mixed adhesive/activator to wet all the tooth surfaces uniformly. Avoid pooling. These surfaces should remain undisturbed for 20 seconds.
7. Remove solvent by gently drying with clean, dry air from a dental syringe for at least 5 seconds. Surface should have a uniform glossy appearance. If not, repeat application and air dry. Surface should not show areas of excessive adhesive thickness or pooling. Repeat air drying/evaporation outlined above if necessary.

Light curing note.

Light curing of applied Prime&Bond® XP adhesive and activator mixture is not necessary when used in conjunction with Dentsply Sirona manufactured dual or self cure resin restoratives. If curing is desired, cure mixed adhesive/activator for 10 seconds³ using a curing light.

8. Dispense and mix dual cure or self cure composite per manufacturer's instructions.

3.5.5 Curing

Placement of self cure composite material: Follow section 3.5.5.2 Chemical self cure, below.

Placement of dual cure composite material: Any of the following placement techniques are acceptable. The preferred method is to use both visible light and chemical self cure (see step 3.5.5.3 Dual cure).

3.5.5.1 Visible light cure

1. Place mixed material directly in increments and light cure each increment per manufacturer's instructions.

3.5.5.2 Chemical self cure

1. Place large increment(s) or load mixed material into a clear crown form and seat onto the preparation.
2. Allow to self cure per manufacturer's instructions.
3. Visible light curing of external surface is optional, but advisable.

3.5.5.3 Dual cure

1. After placing the first increment, and light curing, large increments may then be placed or material loaded into a clear crown form and seated onto the preparation.
2. Allow to self cure per manufacturer's instructions.
3. Visible light curing of external surface is optional, but advisable.

³ Check curing light for minimum light output of at least 800 mW/cm². Cure for at least 20 seconds if light output is between 550 and 800 mW/cm².

3.5.6 Finishing

1. After allowing manufacturer's recommended set time, remove matrix, and apply an optional but advisable surface light curing.
2. Gross reduction may begin immediately.

3.6 Indirect restorations (dual cure inlays, onlays, crowns and bridge retainers)

3.6.1 Cleaning

1. See section 3.1.1 Cleaning.

3.6.2 Tooth conditioning/dentin pre-treatment

1. When used as a bonding agent for inlays or onlays, crowns and bridge retainers, it is recommended to follow the total etch technique described in section 3.1.3 Tooth conditioning/dentin pre-treatment.
2. Once the surfaces have been properly treated, they must be kept uncontaminated. If salivary contamination occurs, repeat procedure beginning at section 3.6.1.

3.6.3 Application

1. When used as a bonding agent for inlays, onlays, crowns or bridge retainers, it is recommended to follow the dual cure mixed adhesive/activator application technique described in section 3.5.4.

Light curing note.

Light curing of applied Prime&Bond® XP adhesive and activator mixture is not necessary when used in conjunction with Dentsply Sirona manufactured dual cure resin cement. If curing is desired, cure mixed adhesive/activator for 10 seconds⁴ using a curing light.

3.6.4 Preparation of restoration

1. Treat surface of restoration according to manufacturer's or dental laboratory's instructions, i.e. etching, mechanical roughening and/or Calibra® Silane Coupling Agent (available separately, see complete Directions for Use) application.
2. Apply a single coat of mixed adhesive/activator to the internal bonding surface of the restoration. Immediately air dry for 5 seconds.

Light curing note.

Light curing of applied Prime&Bond® XP adhesive and activator mixture is not necessary when used in conjunction with Dentsply Sirona manufactured dual cure resin cement. If curing is desired, cure mixed adhesive/activator for 10 seconds⁴ using a curing light.

3.6.5 Cementation

1. Prepare and apply dual cured resin cement, according to manufacturer's instructions.

3.7 Endodontic post cementation

3.7.1 Cleaning

1. See section 3.1.1 Cleaning.

⁴ Check curing light for minimum light output of at least 800 mW/cm². Cure for at least 20 seconds if light output is between 550 and 800 mW/cm².

3.7.2 Tooth conditioning/dentin pre-treatment

1. Rinse and thoroughly dry the prepared post preparation space (using air and paper points).
2. Attach disposable needle to end of syringe. Needle tip may be bent for easy access.
3. Gently extrude respective 34%-36% conditioner to the post space and maintain contact for 15 seconds followed by a 15-second rinse.
4. The preparation post space should then be dried with a gentle air blast and paper points to remove residual moisture, but do not desiccate the conditioned dentin surface.
5. Once the surfaces have been properly treated, they must be kept uncontaminated.
6. If salivary contamination occurs, repeat procedure beginning at step 3.7.1.

3.7.3 Application

1. Place 1-2 drops of Prime&Bond® XP adhesive into a clean CliXdish™ Mixing Well or standard plastic mixing well. Close bottle promptly.
2. Place an equal number of drops of Self Cure Activator into the same mixing well. Replace cap promptly. Mix contents for 1-2 seconds with a clean, unused brush tip.
3. Apply mixed adhesive/activator to post preparation with the brush provided, being sure to apply generous amounts to the preparation orifice. A paper point prewetted with the adhesive mixture may aid in bringing the adhesive mixture down to the deepest portion of the preparation. Avoid pooling. Maintain contact of adhesive/activator with tooth structure for at least 20 seconds.
4. Remove solvent by gently drying with clean, dry air from a dental syringe for at least 5 seconds. Surface should have a uniform glossy appearance. If not, repeat application and air dry. Surface should not show areas of excessive adhesive thickness or pooling. Repeat air drying/evaporation outlined above if necessary. Use of clean, dry paper points may aid in thorough removal of solvent/excess adhesive in post space.

Light curing note.

Light curing of applied Prime&Bond® XP adhesive and activator mixture is not necessary when used in conjunction with Dentsply Sirona manufactured dual cure resin cements or restoratives. If curing is desired, cure mixed adhesive/activator for 10 seconds⁵ using a curing light.

3.7.4 Preparation of post

1. Treat surface of post according to manufacturer's or dental laboratory's instructions, i.e. etching, mechanical roughening and/or Calibra® Silane Coupling Agent (available separately, see complete Directions for Use) application.
2. Apply a single coat of mixed adhesive/activator to the post. Immediately air dry for 5 seconds.

Light curing note.

Light curing of applied Prime&Bond® XP adhesive and activator mixture is not necessary when used in conjunction with Dentsply Sirona manufactured dual cure resin cements or restoratives. If curing is desired, cure mixed adhesive/activator for 10 seconds⁵ using a curing light.

⁵ Check curing light for minimum light output of at least 800 mW/cm². Cure for at least 20 seconds if light output is between 550 and 800 mW/cm².

3.7.5 Post Cementation

1. Mix dual cured resin cement components according to manufacturer's directions and spread on surface of post and/or into the post preparation with a syringe tip, Lentulo Spiral, or metal file.
2. Seat post immediately.
3. Clean up excess with appropriate instruments.
4. Light cure the coronal portion of cemented post for 20 seconds to aid in post stabilization once fully seated.
5. Proceed with core build-up and/or preparation as directed by resin cement manufacturer's instructions.

3.8. Adhesive bonding of direct amalgam restorations

3.8.1 Preparation

1. Finish preparation by removing existing restorations and/or caries.

3.8.2 Cleaning

Clean preparation and place a hard-setting calcium hydroxide liner (Dycal® calcium hydroxide liner, see Directions for Use) if needed, according to:

1. Cleaning [3.1.1 Cleaning].
2. Pulp protection [3.1.2 Pulp protection].

3.8.3 Tooth conditioning/dentin pre-treatment

1. When used as a bonding agent for amalgam materials, it is recommended to follow the total etch technique described in section 3.1.3 Tooth conditioning/dentin pre-treatment.
2. Once the surfaces have been properly treated, they must be kept uncontaminated. If salivary contamination occurs, repeat procedure beginning at section 3.8.2.

3.8.4 Application

1. When used as a bonding agent for amalgam materials, it is recommended to follow the dual cure mixed adhesive/activator application technique described in section 3.5.4.

3.8.5 Application of Amalgam Bonding Base and Amalgam Bonding Catalyst

1. Dispense two drops of Amalgam Bonding Base into a separate, clean plastic mixing well. Replace cap securely.
2. Dispense two drops of Amalgam Bonding Catalyst into the same mixing well. Replace cap securely. Mix contents of plastic well for 1 to 2 seconds with new, clean disposable brush tip.
3. Using the disposable brush, apply the mixed Amalgam Bonding Agents by coating the entire preparation.

3.8.6 Completion

1. Immediately begin placing and condensing amalgam (e.g., Dispersalloy® Dispersed Phase Alloy) as per manufacturer's directions.

4 Hygiene

Cross-contamination.

Infection.

1. Do not reuse single use products. Dispose in accordance with local regulations.
2. Reprocess reusable products as described below.

4.1 Prime&Bond® XP bottle

Cross-contamination.

Infection.

1. The bottle cannot be reprocessed.
2. To prevent bottles from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the bottle is handled offside the dental unit with clean/disinfected gloves.
3. During treatment clinicians with patient contact should not handle bottle.

4.2 CliXdish™ Mixing Well

For reprocessing instructions please refer to the Directions for Use of the CliXdish™, which is available on our webpage at www.dentsply.eu/IFU. If requested, we will send you a free printed copy of the Directions for Use in the language you require within 7 days. You can use the order form provided on our webpage for this purpose.

4.3 Unit Dose Holder

For reprocessing instructions please refer to the Directions for Use of the Unit Dose Holder, which is available on our webpage at www.dentsply.eu/IFU. If requested, we will send you a free printed copy of the Directions for Use in the language you require within 7 days. You can use the order form provided on our webpage for this purpose.

5 Lot number () , expiration date () and disposal

1. Do not use after expiration date.
ISO standard is used: "YYYY-MM" or "YYYY-MM-DD".
2. The following numbers should be quoted in all correspondence:
 - Reorder number
 - Lot number
 - Expiration date
3. Dispose in accordance with local regulations.

© Dentsply Sirona 2018-11-26

[These Directions for Use are based on Master Version 06]

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Manufactured by
Dentsply DeTrey GmbH
De-Trey-Str. 1
78467 Konstanz
GERMANY
www.dentsplysirona.com

Distributed in U.S.A. by
Dentsply Caulk
38 West Clarke Avenue
Milford, DE 19963
USA
www.dentsplysirona.com

X-Post[®]

Endodontic Fiber Post

Endodontischer Glasfaserstift
Tenon endodontique fibré
Perno endodontico in fibra
Postes de fibra de endodoncia

Instructions for Use	English	2
Gebrauchsanweisung	Deutsch	4
Mode d'emploi	Français	6
Istruzioni per l'uso	Italiano	8
Instrucciones de uso	Español	10

X-Post®

Endodontic Fiber Post

Caution: For dental use only. **In the USA:** RX only.

Composition

- Quartz fibers 60% volume
- Epoxy resin 40% volume

Indications for use

These instruments are to be used only in a clinical or hospital environment, by qualified users.

Application field

Endodontic fiber post for direct, adhesive core-build up. In case of insufficient residual tooth substance (< 4 mm), the **X-Post®** may be used to support the coronal restoration.

Contraindications

Insufficient residual dentin: At least 2 mm of tooth structure is required around the preparation.

Warnings

None known.

Precautions

The post must be cleaned with alcohol before insertion in the canal.
Avoid touching the posts with your fingers after cleaning.
Shortening the post should be done outside of the mouth.
Dentsply Sirona recommends the use of a rubber dam.

Adverse Reactions

In the present technical state, no adverse reaction has been reported so far.

Step-by-step Instructions

1. Root canal preparation. Ensure that sufficient space exists to place the post. Remove residual gutta percha out of the root canal.
2. Select correct **X-Post®** size, according to anatomical situation of tooth using the radiograph and the information provided below.
Select the Largo® Peeso Reamer and Easypost™ Precision Drill corresponding to the selected **X-Post®** size.

X-Post®	No. 1	No. 2	No. 3	No. 4
Color Code	Yellow	Red	Blue	Green
Largo® Peeso Reamer	No. 1	No. 1	No. 2	No. 2
Easypost™ Precision Drill	No. 1	No. 2	No. 3	No. 4
Tip diameter (mm)	0.80	0.80	0.80	0.80
Head diameter (mm)	1.35	1.47	1.67	1.83

3. Determine the length of the coronal build-up which will account for 1/3 of the final post length. Remove the root canal obturation with the selected Largo® Peeso Reamer (rotation speed 800-1200 rpm) as deep as necessary for the post to be inserted at 2/3 of its length. At least 4 mm of root canal filling should remain in the apical region.
In case of curved canals, the length has to be reduced.
4. Drill the canal with the selected Easypost™ Precision Drill (rotation speed 1000-1200 rpm).
5. Check that the post is properly seated in the canal.
6. Shorten the post to its final length with a diamond disc outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
7. Clean the post with alcohol.
8. Apply etchant in the canal and to the exposed dentin for 15 seconds (e.g. DeTrey® Conditioner 36). Rinse for 10 seconds. Dry gently with paper points but leave surface moist. Do not put the etchant in contact with the gingiva.
Note: phosphoric etching is optional for Prime&Bond active™, Prime&Bond universal™ and Prime&Bond elect®.¹
9. Mix Prime&Bond® adhesives² and Self Cure Activator. Apply 1 coat of the mixture with a brush (e.g. Applicator Tips) in the root canal and leave for 20 seconds. Remove excess using paper points.
10. Gently air dry all surfaces. Enamel and dentin surfaces should have a uniform, glossy appearance. If not, repeat the application.
11. Light cure the adhesive mixture in the root canal for 10 seconds³.
Note: Light cure step is not required when using core-X® flow + Self Cure Activator with either Prime&Bond active™, Prime&Bond universal™, Prime&Bond elect® or Prime&Bond® XP.

12. Apply a single coat of Prime&Bond® adhesive² and Self Cure Activator mixture to the post. Gently air dry for 5 seconds.

For further treatment procedure, please consult Illustrated Technique Guide and/or the Instructions for Use of core-X® flow or Calibra®.

Disinfection, Cleaning and Sterilization

For hygiene and cross contamination caution, all instruments not marked "sterile" must be cleaned, disinfected and sterilized before each usage to prevent any contamination according to ISO 17664. This concerns the first use as well as the subsequent ones.

General recommendation

- Single use marked instruments are not approved for re-use.
- The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments.
- For your own safety, please wear personal protective equipments (gloves, glasses).
- Use only disinfecting solution which is approved for their efficacy (DGHM-listing, CE marking, FDA approval).
- Tungsten Carbide burs, Plastic Support, Hand instruments and NiTi instruments are degraded by Peroxide hydrogen (H₂O₂) solution.
- NiTi Instruments are degraded if immersed more than 5 minutes in a solution of NaOCl at more than 5%.
- Instruments made of aluminium are degraded in presence of caustic soda solutions with mercury salt. Do not use acid (pH < 6) or alkaline (pH > 8) solutions.
- After 5 cycles, degradation of marking may appear.

¹ Some adhesives are not available in all countries.

² Prime&Bond active™, Prime&Bond universal™, Prime&Bond elect®, Prime&Bond® XP (in some countries sold as XP Bond®), prime&bond® one Etch&Rinse, Prime&Bond® NT.

³ Check curing light for minimum light output of at least 800 mW/cm². Cure for at least 20 seconds if light output is between 500 and 800 mW/cm².

Step-by-step procedure

Implantable Devices

Instruments: Following usage

Instruments: First usage

A. Devices

	Operation	Operating mode	Warning			
1.	Pre-Disinfection or Decontamination	- Soak immediately just after usage all instruments in a disinfectant solution combined with proteolytic enzyme if possible.	- Follow instructions and observe concentrations and immersion time given by the manufacturer (an excessive concentration may cause corrosion or other defects on instruments). - The disinfectant solution should be aldehyde free (to avoid blood impurities fixation). - Do not use disinfectant solution containing Phenol or any products which are not compatible with the instruments (see General recommendation). - For visible impurities are observed on instruments a pre-cleaning is recommended by brushing them manually with soft material.		X	
2a.	Automated Cleaning / Disinfection	- Disassemble the devices (silicone stops to be removed). - Place them in a kit, support or container. - Put them in the thermoisinfectant (for at least 10 min at 93 °C or Ao value > 3000).	- Discard any instruments with large obvious defects (broken, bent). - Avoid any contact between instruments or posts when placing in the thermoisinfectant. - Follow instructions and observe concentrations given by the manufacturer (see also General recommendation). - The thermoisinfectant is not recommended for instruments made of aluminium, tungsten carbide or carbon steel. In case of chemical disinfection, there is a danger of remnants of the disinfectant on the instruments. - Sufficient rinsing step should be available in the program. Use purified water for rinsing (max 10 germs/ml and max 0.25 endotoxin units/ml) like Aqua Purificata. - Use filtered air for drying the instruments or posts. - Use only approved thermoisinfectant according to EN ISO 15883, maintain and calibrate it regularly. - If possible, an automated procedure is preferred.	X	X	X
2b.	Manual Cleaning / Disinfection	- Disassemble the devices (silicone stops to be removed). - Immerse them in the disinfectant solution assisted with ultrasonic device if suitable. - Rinse the devices thoroughly with clean, demineralised or distilled water and then dry them with filtered compressed air.	- No visible impurities should be observed on the instruments. - Discard any instruments with large obvious defects (broken, bent, and twisted). - Avoid any contact between instruments or posts when placing in the solution, use kits, supports or container. - Follow instructions and observe concentrations and time given by the manufacturer (see also General recommendation). - Use purified water for rinsing (max 10 germs/ml and max 0.25 endotoxin units/ml) like Aqua Purificata. - If a disinfectant solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before autoclaving.	X	X	X
3.	Inspection	- Inspect devices and sort out those with defects. - Assemble the devices (stops).	- Dirty instruments must be cleaned and disinfected again. - Discard instruments which show any deformations (bent, twisted), damages (broken, corroded) or defects (loss of colour coding or marking) affecting the resistance, the safety or the performance of the instrument or posts. - Protect carbon steel bur with corrosion inhibitor before packaging.		X	X
4.	Packaging	- Pack the devices in "Sterilisation pouches".	- Check the validity period of the pouch given by the manufacturer to determine the shelf life. - Use packaging which are resistant up to a temperature of 141 °C (286 °F) and in accordance with EN ISO 11607.	X	X	X
5.	Sterilization	- Steam sterilisation at: 134 °C (237 °F) during 18 min. for endodontic instruments and root canal filling instruments. - 134 °C (237 °F) during at least 5 min. for all other devices.	- The instruments, posts and the plastic supports must be sterilized according to the packaging labelling. - Use fractionated vacuum or gravity (less preferred) autoclaves (according to EN 13060, EN 285). - Use validated sterilisation procedure according to ISO 17665-1. - Respect maintenance procedure of the autoclave device given by the manufacturer. - Use only the listed sterilization procedures.	X	X	X
6.	Storage	- Keep devices in sterilization packaging in a dry and clean environment.	- Sterility cannot be guaranteed if packaging is open, damaged or wet (check the packaging before using the instruments).	X	X	X

B. Filling material

	Operation	Operating mode	Warning
1.	Disinfection	- Immerse the obturation devices in alcohol during 30 sec at ambient temperature.	- Do not use disinfectant solution containing Phenol or any products which are not compatible with the instruments (see General recommendation).

© Dentsply Sirona 2017-11

[These Instructions for Use are based on Master Version "X POST RADIX FIBER POST_DFU_1117_MASTER_DSE_EN Rev 00"]

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