Equitamp[®]

OXIDISED REGENERATED CELLULOSE

Registered by:

Equimedical BV
Zwanenburgerdijk 349
1161 NN Zwanenburg
The Netherlands
www.equimedical.com

© Equimedical BV, 2018

EQUIMEDICAL[®]

Equitamp®

Intentionally left blank

Equitamp®

Instruction for use.....	4	(UK)
Mode d'emploi	5	(FR)
Gebrauchsinformation	6	(DE)
Istruzioni per l'uso.....	7	(IT)
Instrucciones de uso	8	(SP)
Instruções para utilização	9	(PT)
Gebruiksaanwijzing	10	(NL)
Brugsvejledning.....	11	(DK)
Οδηγίες χρήσης	12	(GR)
Kullanım kılavuzu	13	(TR)
)Instrucțiuni de utilizare.....	14	(RO)
Instrukcja stosowania	15	(PL)
Инструкция по применению.....	16	(RU)
Инструкции за използване	17	(BG)
Használati utasítás	18	(HU)
Příbalová informace	19	(CZ)
Қолдану ережесі	20	(KZ)
Uputstvo za upotrebu	21	(RS)
Käyttöohje.....	22	(FI)
Bruksanvisning	23	(SV)
Návod na použitie.....	24	(SK)
Navodila za uporabo.....	25	(SL)
Upute za uporabu	26	(HR)
تعليمات للاستخدام	27	(AR)

Instructions for use

Equitamp®

Oxidised regenerated cellulose with haemostatic effect.

Equitamp is prepared by oxidising a suitable form of cellulose, pure viscose. This is followed by additional processes in order to obtain a pure and high-quality form of oxidised and regenerated cellulose. It is strong and although a slight discoloration may occur with age, this does not affect performance. Equitamp is immediately available for use in the operating theatre and does not require any sterilisation by dry heat or autoclaving. Equitamp can be sutured or cut without fraying. Equitamp is double sterile packed.

Equitamp can be used in many areas of surgery, e.g. cardiovascular surgery, haemorrhoidectomy, implantation of vascular prostheses, biopsies, lung operations, surgery to the face and jaw, gastric resection, operations to the throat or nose, liver and gall bladder operations, gynaecological operations, thoracic and abdominal sympathectomies, neurosurgery, especially cerebral operations, thyroid operations, skin transplants, treatment of superficial injuries.

Indication:

Equitamp is used adjunctively in surgical procedures to assist in the control of capillary, venous, and small arterial haemorrhages when ligation or other conventional methods of control are impractical or ineffective. Equitamp can be cut to size for use in endoscopic procedures.

Equitamp is indicated also for adjunctive use in dental applications to assist in the control of bleeding in exodontias and oral surgery. It may also be used to help achieve haemostasis after single or multiple tooth extractions, alveoloplasty, gingival haemorrhage, impactions, biopsies, and other procedures in the oral cavity.

Bactericidal:

In addition to its local haemostatic properties, Equitamp has been proven by an independent study to be bactericidal in vitro against a wide range of gram positive and gram negative organisms including the following emerging multi-drug resistant micro-organisms. A log 7 reduction of above mentioned micro-organisms was observed with usage of Equitamp after 24 hours.

Staphylococcus aureus meticylinooporny (MRSA)
Staphylococcus epidermidis meticylinooporny (MRSE)
Enterococcus faecalis oporny na wankomycynę (VRE)
Streptococcus pneumoniae penicylinooporny (PRSP)

However, Equitamp is not an alternative to specially designed systemically applied antimicrobial agents.

Haemostatic mechanism:

When Equitamp comes into contact with blood, it will absorb the blood and gradually swell, eventually dissolving into a gelatinous material. By briefly applying pressure at this point, the material will adhere to the wound, effectively sealing the ends of the venous capillaries and resulting in a mechanical haemostatic effect. The coagulation cascade is activated, transforming soluble fibrinogen into a net of insoluble fibrin which stops the bleeding. When implanted into tissue, Equitamp should be absorbed within a maximum of 30 days. Under normal circumstances, Equitamp should be absorbed in 8 days.

Mode of application:

Equitamp should be used dry... Equitamp must be pressed against the haemorrhage site. The gauze may remain in situ (after strict aseptic measures have been taken) and the wound may be closed since the gauze is completely resorbed... Equitamp does not provoke any harmful tissue reaction or any reaction of hypersensitivity. Equitamp can be discharged by natural means, for instance by the uterus, once haemostasis has been achieved due to the fact that haemostasis has been achieved and Equitamp can not adhere to any tissue as a result.

Warning:

When placed into cavities or closed tissue spaces, minimal preliminary compression is advised and care should be exercised to avoid overpacking (the gauze expands upon absorption of liquid). Equitamp may swell to its original size and absorbed fluids may increase the risk of nerve damage... Due to the fact that Equitamp may form a nidus for infection, it must not be left in infected areas; it must be removed once bleeding has been controlled.

Keep dry

Sterile

Use only once

Store under 30°C

Do not use when package has been opened or has been damaged.

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