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Intended use
Addition-vulcanising silicone for the production of orthodontic elasto-positioners

Indication for use
Manufacturing orthodontic appliances to prevent a patient's teeth from shifting or to move teeth to a final position.

Product description
A clear-transparent, especially tear resistant, addition-vulcanising silicone in double cartridges for the production of orthodontic appliances. Due to the platinum hardener system it shows practically no shrinkage, the appliance is exact and stable and shows an excellent, checked biocompatibility (contains no peroxide). As a result of the treatment with Odontosil positioner the teeth are placed into the correct position. The Odontosil Finisher is at the same time also a positioner and keeps the teeth in the requested position. The conventional indirect method of wax modellation in a flask is also possible as well as the modellation directly at the Set-up. Due to the vulcanisation of Odontosil at room temperature and the thixotropic consistency the faster, direct modellation is possible. Odontosil is a clear transparent material; due to the coating with Odontosil 2K Lack the surface can be sealed. Final hardnesses of 50 and 60 Shore A are available and can be combined among each other within one workstep.

The field of indication comprises elastic appliances for the fine adjustment of the teeth after multiband-treatment. Due to its high elasticity, it is possible to move the teeth also without a previous multiband-technique. For this reason the Odontosil positioner becomes an active appliance with a comprehensive field of application. For grinding we recommend the special Odontosil grinding caps and for the final glasse the Odontosil 2K Lack.

Odontosil vulcanises by mixing the two components. Temperature is not necessary but can accelerate the polymerisation.

Technical Data	
Odontosil 50 / Odontosil 60	
Total processing time	≥ 5min.
Deflasking	≥ 32 min.
Setting end	ca. 120 min.
Consistency	type 2 – medium consistency
Deformation under pressure	< 4,0 %
Retention after deformation	> 99,5 %
Colour	transparent

Preparation of the model
Manufacturing of a Set-up-model. For this purpose teeth resp. jaw segments from the situation model are detached, changed according to the statement of the dentist and afterwards connected with the model in the requested position.

The Set-up-model is fixed together with the counter bite in the articulator or occludator. The production of a positioner requires a disconnection of the jaws among each other for approx. 2.5–3.0 mm, but the habitual occlusion of the patient should not be exceeded.

Processing instruction
Direct Method
Open the articulator with the Set-up-model and apply Odontosil in the requested layer thickness onto both dental archs as well as onto the gingiva parts up to the reflection (oral, buccal, occlusal) by means of the Injector and mixing canulas. Afterwards close the articulator and smooth the Odontosil surface carefully with the finger moistened with soapsud. The polymerisation is effected at room temperature. To avoid bubbles in the silicone the polymerisation can also be effected under pressure (at least 3 bar). After the polymerisation (approx. 30 minutes) the positioner is released form the model, the edges are cut with a scissors and afterwards put into the final form and strength with special Odontosil grinding caps. For the homogenisation of he positioner surface a coating with Odontosil 2K Lack is recommended (see instruction). The coating also takes place at room temperature in approx. 30 minutes.

Indirect method
A plaster duplicate of the modelled Set-up-model is produced and articulated. Onto this duplication models the positioner is modelled into its final form and strength with wax. The wax positioner is embedded with plaster into a brass flask and after hardening boiled out in the Laborat. The hollow form should be coated with an alginate-insulation (Isolat film). Apply Odontosil by means of the Injector directly into the flask. In order to achieve an exact air bubble free result the closed flask should be stored under pressure of 3 bar. This pressure should be hold up until the end of the vulcanisation time (approx. 30 minutes) as e. g. with a flask clamp. After the deflasking the finished positioner and the removal of possible plaster surplus Odontosil 2K Lacquer can be applied for the sealing of the surface. At room temperature the vulcanisation time of the lacquer is about 30 min.

Important advice
• Apply the Odontosil 2K Lack within the setting end (120 min.) of the Odontosil base silicone to ensure a bond between both materials.
• Due to filling reasons an uneven transportation of the components can occur when using new cartridges. For this reason the first transported quantity in the size of a pea should be thrown away.
• Never transport against a hardened mixing canula. This could lead to leakiness in the cartridge system and

therefore to an emission of unmixed, non-hardened single components.
• Leave mixing canula on the cartridge for sealing effect.
• Wear protective clothing! The cleaning of silicone dirtying on the clothes is not possible (also not chemical).
• Sensitive persons can show allergic reactions to one of the product components.
• A disinfection with all commercial materials is possible.
• The vulcanisation can be disturbed by using latex gloves, hand creme, cleaning agents etc..
• The product is exclusively designed for the a. m. field of application and is only to be used by qualified and instructed persons.
• The general information of this product is based on detailed application practice and intensive research work. Due to further development we reserve the right of modification.

Storage conditions
At room temperature, i.e. 18–28°C.
Shelf life: See expiry date.
Odontosil 2K Lack is not to be used after its expiry date has passed.

Delivery forms	
REF D5038150	Odontosil 50 2 x 50 ml double cartridges, incl. 12 mixing canulas, green
REF D5038229	Odontosil 50 8 x 50 ml double cartridges
REF D5038160	Odontosil 60 2 x 50 ml double cartridges, incl. 12 mixing canulas, green
REF D5038227	Odontosil 60 8 x 50 ml double cartridges
REF D38234	Odontosil 2K Lack 2 x 50 ml double cartridges, incl.12 mixing canulas, lightblue

This product has been developed for dental use and must be applied according to the instructions. The manufacturer shall assume no liability in the event of damage resulting from misuse or improper application. Furthermore, the user is responsible for checking that the product is suitable for the intended purpose, especially if this purpose is not mentioned in the instructions.

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