



W&H adaptor set for cleaning units and disinfectors

Tips for preparation in cleaning units and disinfectors

Cleaning adaptor for ultrasonic scaler tips (prophylaxis, periodontology, endodontics, restorative dentistry & prosthetics)

- › All scaler tips with a M3 x 0.5 mm thread (internal thread – viewed from the tip) fit on the cleaning adaptor with the REF 06895000
- › EMS tips with a M3 x 0.5 mm thread (internal thread) e.g., EMS C (EMS tips do not fit the quick-clamping system.)
- › NSK tips with a M3 x 0.5 mm thread (internal thread), compatible with EMS handpiece
All other NSK tips do not fit
- › Mectron tips with a M3 x 0.5 mm thread (internal thread), e.g., Mectron S7
- › Satelec tips are not compatible
(M3 x 0.6 mm thread)

Cleaning adaptors for ultrasonic scaler tips (surgery)

- › Only W&H Piezomed instruments fit the adaptor set for cleaning units and disinfectors
- › Tips like Sirona SIROSON, NSK SG3, Satelec CS3 and Mectron OT 7 are not compatible due to the different thread sizes

Cleaning adaptor for sonic scaler tips (prophylaxis, periodontology, endodontics, restorative dentistry & prosthetics, surgery)

- › All scaler tips with a M3 x 0.5 mm thread (external thread – viewed from the tip) fit on the cleaning adaptor with the REF 06894700
- › Sirona tips with a M3 x 0.5 mm thread (external thread) like SIROAIR
- › KaVo tips with a M3 x 0.5 mm thread (external thread) like KaVo tip 36/10
(KaVo tips do not fit the quick-clamping system.)
- › NSK tips with a M3 x 0.5 mm thread (external thread) like NSK S1
- › Brasseler/Komet tips with a M3 x 0.5 mm thread (external thread) like SF8878KD

**External coolant tubes, spray clips,
Y manifolds, etc., can also be cleaned**

Errors and omissions excepted

W&H Dentalwerk Bürmoos GmbH

Ignaz-Glaser-Straße 53, Postfach 1

5111 Bürmoos, **Austria**

t +43 6274 6236-0

f +43 6274 6236-55

office@wh.com

wh.com