

Alcohol Torch

Economical alcohol torch that produces a needle-point flame



The classic Buffalo Dental Alcohol Torch produces a needle-point flame with only a gentle squeeze of the plastic bottle, allowing for precision control.

The Alcohol Torch now features an optional 3-point base, which reduces the chance of knocking the bottle over.

Features and Benefits:

- 5 1/4" high torch uses denatured alcohol fuel for cost savings and convenience
- Base of bottle is weighted with marbles to further improve safety
- Ideal for use in applying rolled edges to vacuum formed adaptations (flame edge and rub with thumb)
- Versatile torch can also be used as an alcohol lamp
- 3-point base for added stability
- The base has an easy slide fit and can remain with the torch as you use it or on the bench as a receptacle
- Includes plastic cap to help keep fuel from evaporating when torch is not in use (for use after torch is extinguished)

Ordering Information:

Item#	Description
82440	Plastic Alcohol Torch w/ 7" Wick
82410	Plastic Alcohol Torch w/ 7" Wick and Base
57560-1	Replacement Wick, 7" long, Pkg. of 1
57560-1D	Replacement Wicks, 7" long, Pkg. of 12
82400	Fluorescent Base Only
Torchcover Replacement Clear Cap for Torch	



Buffalo Dental Manufacturing Co., Inc.

159 Lafayette Dr., PO Box 678, Syosset, NY 11791-0678 | Tel: 516.496.7200 or 800.828.0203

Fax: 516.496.7751 | Email: sales@mailbdm.com | www.buffalodental.com