

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproCleaner**
Issue/Revision: 20.02.2024

Version: 3.1
Replaces version: 3.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Trade name: **AlproCleaner**
UFI: YDX5-N7K4-C00F-9P5P

1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses: Cleaning agent
Intended purpose: Ready-to-use liquid cleaning agent for the spittoon bowls of dental units.
Uses advised against: None at intended use.
Note: The product is intended for professional users.

1.3. Details of the supplier of the safety data sheet

Manufacturer/Supplier: ALPRO MEDICAL GMBH
Mooswiesenstraße 9
D-78112 St. Georgen (Germany)
Telephone: +49 7725 9392-0
Telefax: +49 7725 9392-91
E-mail: info@alpro-medical.de
Internet: www.alpro-medical.de

E-mail address for the competent person responsible for the safety data sheet: doku@alpro-medical.de

1.4. Emergency telephone number

In-house emergency telephone number: +49 7725 9392-0
Monday – Friday from 08:00 am to 04:30 pm (UTC+1);
for chemical information and legal information on
hazardous substances only
Poison centre: +49 761 19240
Poisoning information centre, Freiburg, Germany
(24 h / 7 d), English is spoken
National Poisons Information Service (UK): +44 344 892 0111
National Poisons Information Service (NPIS)
(24 h / 7 d), Medical Professionals Only

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification in accordance with Regulation (EC) No 1272/2008 [CLP]

Classification	Classification procedure
Eye Irrit. 2; H319	Calculation method

Full text of hazard classes as well as H-phrases: see under SECTION 16.1.

2.2. Label elements

Label elements in accordance with Regulation (EC) No 1272/2008 [CLP]

Hazard pictograms:



Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproCleaner**
Issue/Revision: 20.02.2024

Version: 3.1
Replaces version: 3.0

Signal word:		Warning
Hazard components for labelling:		-
H- phrases:	H319	Causes serious eye irritation.
P- phrases:	P280 P305+P351+P338 P337+P313	Wear eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

2.3. Other hazards

The substances in the mixture do not meet the PBT/vPvB criteria according to Regulation (EC) No 1907/2006, Annex XIII (see SECTION 12.5.)

The substances in the mixture do not have any endocrine disrupting properties according to Regulation (EC) No. 1907/2006, Annex XIV (see Section 11 and Section 12.6.). They are not on the list of substances of very high concern for authorization according to Regulation (EC) No. 1907/2006, Article 59, Paragraph 10).

No further hazards known.

SECTION 3: Composition/information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

Chemical characterisation: Mixture of substances listed below with non-hazardous additions in aqueous solution.

Hazardous ingredients

Chemical name	Identification numbers	Classification in accordance with Regulation (EC) No 1272/2008	Weight %
Propan-2-ol	CAS No: 67-63-0 EC No: 200-661-7 Index No: 603-117-00-0 REACH Registration No: 01-2119457558-25-XXXX	Flam. Liq. 2; H225 Eye Irrit. 2; H319 STOT SE 3; H336	≥ 1 - < 5
Trisodium nitilotriacetate	CAS No: 5064-31-3 EC No: 225-768-6 Index No: 607-620-00-6 REACH Registration No: 01-2119519239-36-XXXX	Carc. 2; H351 Acute Tox. 4; H302 Eye Irrit. 2; H319 <i>Specific concentration limits:</i> Carc. 2; H351: C ≥ 5 %	≥ 1 - < 5
1-methoxypropan-2-ol	CAS No: 107-98-2 EC No: 203-539-1 Index No: 603-064-00-3 REACH Registration No: 01-2119457435-35-XXXX	Flam. Liq. 3; H226 STOT SE 3; H336	≥ 0.1 - < 1

Full text of hazard classes and H-phrases: see SECTION 16.1.

Occupational exposure limits: see SECTION 8.1

SECTION 4: First aid measures

4.1. Description of first aid measures

General information: First aider: Pay attention to self-protection!

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproCleaner**
Issue/Revision: 20.02.2024

Version: 3.1
Replaces version: 3.0

- | | |
|-------------------------|---|
| Following inhalation: | Move affected person into fresh air and keep still and warm. In case of continued complaints seek medical advice. |
| Following skin contact: | Wash skin immediately with plenty of water and soap. In case of skin reactions, consult a physician. |
| Following eye contact: | Flush eyes immediately with flowing water for 10 to 15 minutes holding eyelids apart. Remove contact lenses, if present and easy to do. Consult an ophthalmologist. |
| Following ingestion: | Rinse mouth with water. Let drink plenty of water. Do not induce vomiting. Consult a physician immediately. |

4.2. Most important symptoms and effects, both acute and delayed

Causes serious eye irritation.

4.3. Indication of any immediate medical attention and special treatment needed

No information available.

SECTION 5: Firefighting measures

5.1. Extinguishing media

- | | |
|---------------------------------|--|
| Suitable extinguishing media: | Water spray jet, alcohol resistant foam, extinguishing powder, carbon dioxide (CO ₂) |
| Unsuitable extinguishing media: | Full water jet |

5.2. Special hazards arising from the substance or mixture

Hazardous combustion products: Carbon monoxide (CO), carbon dioxide (CO₂), nitrogen oxides (NO_x)

5.3. Advice for firefighters

- | | |
|-------------------------------|--|
| Special protective equipment: | Wear self-contained breathing apparatus. |
| Further information: | Cool endangered containers with water spray jet. |

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

For non-emergency personnel

- Use personal protective equipment. See SECTION 8.2.
- Avoid skin and eye contact. Special danger of slipping by leaked/spilled product.

For emergency responders

- Use personal protective equipment. See SECTION 8.2.

6.2. Environmental precautions

Do not discharge into drains or surface and ground water.

6.3. Methods and material for containment and cleaning up

Containment

For large spills, dyke spilled material or otherwise contain material to ensure runoff does not reach a waterway. Cover or seal drains.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproCleaner**
Issue/Revision: 20.02.2024

Version: 3.1
Replaces version: 3.0

Cleaning up

Wipe up small amounts with absorbent material (e.g. cloth, fleece). Absorb large amounts with liquid-binding material (sand, diatomaceous earth, universal binder, sawdust). Collect in suitable, closed containers for disposal. Clean contaminated surfaces thoroughly.

Other information

Inappropriate containment and cleaning methods are not known.

6.4. Reference to other sections

Information on safe handling see SECTION 7.1.

Information on personal protective equipment see SECTION 8.2.

Information on disposal see SECTION 13.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions

Avoid contact with skin and eyes. Fill refill packages only in labelled original bottles.

Do not mix with other cleaners or disinfectants.

Advice on general occupational hygiene

When using do not eat, drink or smoke. Wash hands before breaks and at end of work. Keep away from medicines, food, feed, cosmetics and stimulants.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels: Keep only in the original container. Keep container tightly closed and kept upright to prevent any leakage.

Advice on common storage: Not necessary

Further information on storage conditions: Not necessary

7.3. Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific end uses are stipulated.

Industry and sector specific guidance

[DE] TRGS 525 – Hazardous substances in medical care facilities (Section 7 Activities with disinfectants); Issue: September 2014; source: GMBI 2014 page 1294-1307 of 13.10.2014 [No 63]; 10.07.2015 [No 27]; www.baua.de.

[DE] DGUV rules 207-206 – Prevention of chemical risks when handling disinfectants in health service, Issue: 2016.12; source: <https://publikationen.dguv.de>

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproCleaner**
Issue/Revision: 20.02.2024

Version: 3.1
Replaces version: 3.0

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Occupational exposure limits

Country	Limit values				Legal basis	Remarks
	Long term (8 hours)		Short term (15 minutes)			
	ppm	mg/m ³	ppm	mg/m ³		
Propan-2-ol (CAS No: 67-63-0)						
EU						no limit value specified
UK	400	999	500	1250	EH40	
Trisodium nitrilotriacetate (CAS No: 5064-31-3)						
EU						no limit value specified
UK						no limit value specified
1-methoxypropan-2-ol (CAS No: 107-98-2)						
EU	100	375	150	568	2000/39/EG	Skin
UK	100	375	150	560	EH40	Sk

Used abbreviations, symbols, numerals and explanations

Skin A significant uptake of the substance through the skin is possible.
Sk Can be absorbed through the skin. The assigned substances are those for which there are concerns that dermal absorption will lead to systemic toxicity.

Biological limit values

Country	Parameter	Limit value	Test material	Sampling time	Legal basis
Propan-2-ol (CAS No: 67-63-0)					
Germany	Acetone	25 mg/l	Whole blood	End of exposition, resp. end of shift	TRGS 903
	Acetone	25 mg/l	Urine	End of exposition, resp. end of shift	TRGS 903
1-methoxypropan-2-ol (CAS No: 107-98-2)					
Germany	1-methoxypropan-2-ol	15 mg/l	Whole blood	End of exposition, resp. end of shift	TRGS 903
	1-methoxypropan-2-ol	15 mg/l	Urine	End of exposition, resp. end of shift	TRGS 903

Information on monitoring procedures

BS EN 482:2021-04; Title: Workplace exposure - Procedures for the determination of the concentration of chemical agents;
British version of EN 482:2021

BS EN 689:2018; Title: Workplace atmospheres - Guidance for the assessment of exposure by inhalation to chemical agents for comparison with limit values and measurement strategy;
British version of EN 689:2018

BS EN 14042:2003-04-24; Title: Workplace atmospheres - Guide for the application and use of procedures for the assessment of exposure to chemical and biological agents;
British version of EN 14042:2003)

CEN/TR 17055:2017; Title: Workplace exposure. Measurement of chemical agents complying with the requirements given in EN 482 and either one of EN 838, EN 1076, EN 13205, EN 13890 and EN 13936.
Choice of procedures

ISO TR 14294:2011; Title: Workplace atmospheres. Measurement of dermal Exposure. Principles and methods

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproCleaner**
Issue/Revision: 20.02.2024

Version: 3.1
Replaces version: 3.0

8.2. Exposure controls

Appropriate engineering controls

Technical and organisational protective measures

The eyewash station (or eyewash bottle) must be located near the workplace.

Personal protective equipment

Eye/face protection: Safety glasses with side protection according to BS EN 166

Skin protection:

Hand protection: Protective gloves according to BS EN ISO 374-1 and BS EN 21420

Splash guard:

Protective gloves: type C; permeation-resistant at least 10 minutes

Permanent contact (> 480 min):

Protective gloves: type A or B; code letters: G, K, O, H;
permeation-resistant at least 30 minutes

Other skin protection: Not necessary when used as intended.

Respiratory protection: Not necessary when used as intended.

Thermal hazards: No special protective measures necessary.

Environmental exposure controls

Not necessary

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Appearance/physical state:	slightly cloudy, colourless liquid	
Odour:	alcoholic, of lemon	
Odour threshold:	no data available	
Melting point/freezing point:	no data available	
Initial boiling point and boiling range:	no data available	
Flammability (solid, gas):	not applicable	
Lower explosive limit:	not applicable	
Upper explosive limit:	not applicable	
Flash point:	> 60 °C	
Auto-ignition temperature:	no data available	
pH (undiluted):	12.0 – 13.0	(20 °C)
Kinematic viscosity:	no data available	
Solubility in water:	completely soluble	
Partition coefficient: n-octanol/water	not applicable	
Vapour pressure:	no data available	(... °C)
Density:	1.000 – 1.010 g/cm ³	(20 °C)
Relative vapour density:	no data available	
Particle characteristics:	not applicable	

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproCleaner**
Issue/Revision: 20.02.2024

Version: 3.1
Replaces version: 3.0

9.2. Other information

9.2.1. Information with regard to physical hazard classes

Explosive substances/mixtures and products containing explosives:	not applicable
Flammable gases:	not applicable
Aerosols:	not applicable
Oxidising gases:	not applicable
Gases under pressure:	not applicable
Flammable liquids:	not applicable
Flammable solids:	not applicable
Self-reactive substances and mixtures:	not applicable
Pyrophoric liquids:	not applicable
Pyrophoric solids:	not applicable
Self-heating substances and mixtures:	not applicable
Substances and mixtures, which emit flammable gases in contact with water:	not applicable
Oxidising liquids:	not applicable
Oxidising solids:	not applicable
Organic peroxides:	not applicable
Substances and mixtures corrosive to metals:	not applicable
Desensitised explosives:	not applicable

9.2.2. Other safety characteristics

Electrical conductivity (undiluted): 10,000 – 13,000 $\mu\text{S}/\text{cm}$ (20 °C)

SECTION 10: Stability and reactivity

10.1. Reactivity

No hazardous reactions when handled and stored as intended.

10.2. Chemical stability

The product is stable when handled and stored as intended.

10.3. Possibility of hazardous reactions

None known

10.4. Conditions to avoid

None known

10.5. Incompatible materials

No incompatible materials known when used as intended.

10.6. Hazardous decomposition products

Does not decompose when used as intended.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproCleaner**
Issue/Revision: 20.02.2024

Version: 3.1
Replaces version: 3.0

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity

Product

Acute toxicity - oral:	Acute Toxicity Estimate ATE _{mix} > 2000 mg/kg => no classification
Acute toxicity - dermal:	Acute Toxicity Estimate ATE _{mix} > 2000 mg/kg => no classification
Acute toxicity - inhalation:	Acute Toxicity Estimate ATE _{mix} > 20 mg/l => no classification

Ingredients

Trisodium nitrilotriacetate (CAS No: 5064-31-3):

Acute toxicity - oral: 1470 mg/kg bw; species: rat; strain: Wistar method: EU method B.1

Skin corrosion/irritation

Product

No classification. [calculation method]

Serious eye damage/irritation

Product

Causes serious eye irritation. [calculation method]

Respiratory or skin sensitisation

Product

No data available.

Germ cell mutagenicity

Product

No data available.

Carcinogenicity

Product

No classification. [calculation method]

Ingredients

Trisodium nitrilotriacetate (CAS No: 5064-31-3):

Carcinogenicity: NOAEL: 9 mg/kg bw/d; species: rat; strain: Fischer 344; oral uptake; target organ: urogenital, kidney; method: study

Reproductive toxicity

Product

No classification. [calculation method]

STOT-single exposure

Product

No classification. [calculation method]

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproCleaner**
Issue/Revision: 20.02.2024

Version: 3.1
Replaces version: 3.0

Ingredients

Propan-2-ol (CAS No: 67-63-0):

May cause drowsiness or dizziness.

1-methoxypropan-2-ol (CAS No: 107-98-2):

May cause drowsiness or dizziness.

STOT-repeated exposure

Product

No data available.

Aspiration hazard

Product

No data available.

Information on other hazards

Endocrine disrupting properties:

No substances are contained that have endocrine disrupting properties for humans.

SECTION 12: Ecological information

12.1. Toxicity

No classification. [calculation method]

12.2. Persistence and degradability

No data available.

12.3. Bioaccumulative potential

No data available.

12.4. Mobility in soil

No data available.

12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to Regulation (EC) No 1907/2006, Annex XIII.

12.6. Endocrine disrupting properties

No substances are contained that have endocrine disrupting properties for non-target organisms.

12.7. Other adverse effects

No data available.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal of the product

Product residues must be disposed of as hazardous waste in compliance with the Directive 2008/98/EC on waste as well as national and regional regulations. Do not dispose of via the waste water. Leave product in the original container as possible. Do not mix with other waste materials.

Waste codes / waste designations according to EWC

Product residues: 16 10 01* aqueous liquid wastes containing hazardous substances

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproCleaner**
Issue/Revision: 20.02.2024

Version: 3.1
Replaces version: 3.0

Disposal of the packaging

Packaging contaminated with product is considered as hazardous waste and must be disposed of accordingly.

Waste codes / waste designations according to EWC

Contaminated packaging: 15 01 10* packaging containing residues of or contaminated by hazardous substances

Recommendation

Contaminated packaging must be emptied optimally and can be recycled after appropriate cleaning (rinse with water).

From a dilution to 1 % the concentrate is no longer classified as a hazardous substance.

SECTION 14: Transport information

14.0. Transport classification

No dangerous good in sense of the transport regulations in road traffic (ADR), railway traffic (RID), inland waterway traffic (ADN), maritime traffic (IMDG-Code) and air traffic (ICAO-TI/IATA-DGR).

14.1. UN number

-

14.2. UN proper shipping name

-

14.3. Transport hazard class(es)

-

14.4. Packing group

-

14.5. Environmental hazards

-

14.6. Special precautions for user

-

14.7. Maritime transport in bulk according to IMO instruments

-

14.8. Further information

-

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

EU-Regulations

REGULATION (EC) No 1005/2009 on substances that deplete the ozone layer
not applicable

REGULATION (EC) No 2019/1021 on persistent organic pollutants
not applicable

REGULATION (EU) No 649/2012 concerning the export and import of hazardous chemicals

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproCleaner**
Issue/Revision: 20.02.2024

Version: 3.1
Replaces version: 3.0

not applicable

REGULATION (EU) No 648/2004 on detergents

Anionic surfactants: < 5 %
NTA (nitrilotriacetic acid): < 5 %
Disinfectants
Perfumes

DIRECTIVE 2012/18/EU (Seveso III Directive) on the control of major-accident hazards involving dangerous substances, amending and subsequently repealing Council Directive 96/82/EC

not applicable

DIRECTIVE 2010/75/EU on industrial emissions (integrated pollution prevention and control)

VOC content: < 4 %

REACH – List of substances subject to authorisation (Annex XIV)

not applicable

REACH – Restrictions on the manufacture, placing on the market and use of certain dangerous substances, mixtures and articles (Annex XVII)

not applicable

COUNCIL DIRECTIVE 94/33/EC on the protection of young people at work

not applicable

COUNCIL DIRECTIVE 92/85/EEC on the introduction of measures to encourage improvements in the safety and health at work of pregnant workers and workers who have recently given birth or are breastfeeding

not applicable

15.2. Chemical safety assessment

For this mixture no chemical safety assessment has been carried out.

SECTION 16: Other information

16.1. Full text of hazard classes and H-phrases

Hazard classes

Acute Tox.	Acute toxicity
Carc.	Carcinogenicity
Eye Irrit.	Eye irritation
Flam. Liq.	Flammable liquid
STOT SE	Specific target organ toxicity (single exposure)

H-phrases (Hazard statements)

H225	Highly flammable liquid and vapour.
H226	Flammable liquid and vapour.
H302	Harmful if swallowed.
H319	Causes serious eye irritation.
H336	May cause drowsiness or dizziness.
H351	Suspected of causing cancer <state route of exposure if it is conclusively proven that no other routes of exposure cause the hazard>.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproCleaner**
Issue/Revision: 20.02.2024

Version: 3.1
Replaces version: 3.0

16.2. Abbreviations and acronyms

ADN	<u>A</u> ccord européen relatif au transport international des marchandises <u>d</u> angereuses par voie de <u>n</u> avigation intérieure (European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways)
ADR	<u>A</u> ccord européen relatif au transport international des marchandises <u>d</u> angereuses par <u>r</u> oute (European Agreement concerning the International Carriage of Dangerous Goods by Road)
BS	<u>B</u> ritish <u>S</u> tandards
CAS	<u>C</u> hemical <u>A</u> bstracts <u>S</u> ervice
CLP	Regulation on <u>C</u> lassification, <u>L</u> abelling and <u>P</u> ackaging of Substances and Mixtures
[DE]	National German regulations
DGUV	<u>D</u> eutsche <u>G</u> esetzliche <u>U</u> nfall <u>v</u> ersicherung (English: German statutory accident insurance)
EC	<u>E</u> uropean <u>C</u> ommunity
EEC	<u>E</u> uropean <u>E</u> conomic <u>C</u> ommunity
EN	European Standard
EU	<u>E</u> uropean <u>U</u> nion
EWC	<u>E</u> uropean <u>W</u> aste <u>C</u> atalogue
GMBI	<u>G</u> emeinsames <u>M</u> inisterial <u>b</u> latt (English: Joint Ministerial Gazette)
IATA-DGR	<u>I</u> nternational <u>A</u> ir <u>T</u> ransport <u>A</u> ssociation - <u>D</u> angerous <u>G</u> oods <u>R</u> egulations
ICAO-TI	Technical Instructions For The Safe Transport of Dangerous Goods by Air
IMDG-Code	<u>I</u> nternational <u>M</u> aritime Code for <u>D</u> angerous <u>G</u> oods
LC ₅₀	Median lethal concentration
LD ₅₀	Median lethal dose
NOAEL	No Observed Adverse Effect Level
OECD	<u>O</u> rganization for <u>E</u> conomic <u>C</u> o-operation and <u>D</u> evelopment
PBT	<u>P</u> ersistent, <u>b</u> ioaccumulative and <u>t</u> oxic
ppm	<u>P</u> arts <u>p</u> er <u>m</u> illion
REACH	<u>R</u> egistration, <u>E</u> valuation, <u>A</u> uthorisation and Restriction of <u>C</u> hemicals
RID	<u>R</u> èglement concernant le transport <u>I</u> nternational ferroviaire de marchandises <u>D</u> angereuses (Regulations Concerning the International Carriage of Dangerous Goods by Rail)
TRGS	<u>T</u> echnische <u>R</u> egeln für <u>G</u> efahr <u>s</u> toffe (English: Technical Rules for Hazardous Substances)
UN	<u>U</u> nited <u>N</u> ations
UTC	Coordinated Universal Time (French: Temps Universel Coordonné)
VOC	Volatile Organic Compounds
vPvB	<u>V</u> ery persistent and <u>v</u> ery <u>b</u> ioaccumulative

16.3. Key literature references and sources for data

- Regulation (EC) No 1907/2006 (REACH), Annex II
- European Chemicals Agency (ECHA) – Guidance on the compilation of safety data sheets; Version 4.0 (December 2022); <https://echa.europa.eu/documents>
- Regulation (EC) No 1272/2008 (CLP regulation)
- European Chemicals Agency (ECHA) – Guidance on Labelling and Packaging in accordance with Regulation (EC) No 1272/2008; Version 4.2 (03/2021); <https://echa.europa.eu/documents>
- European Chemicals Agency (ECHA), Registered substances; <https://echa.europa.eu/information-on-chemicals/registered-substances>
- European Chemicals Agency (ECHA), C&L Classification and Labelling Inventory; <https://echa.europa.eu/information-on-chemicals/cl-inventory-database>
- Institute for Occupational Safety and Health of the German Social Accident Insurance (IFA): GESTIS database on hazardous substances and GESTIS - International limit values for chemical agents; <https://www.dguv.de/ifa/index.jsp>
- German Environmental Agency (Umweltbundesamt), Section IV 2.4: Office of Documentation and Information on Substances Hazardous to Waters RIGOLETTO (catalogue of Substances Hazardous to Waters); <https://webrigoletto.uba.de/rigoletto>

Safety data sheet
according to Regulation (EC) No 1907/2006



Trade name: **AlproCleaner**
Issue/Revision: 20.02.2024

Version: 3.1
Replaces version: 3.0

16.4. Methods according to Article 9 of Regulation (EC) No 1272/2008 for the evaluation of information for classification purposes

Calculation method according to the criteria in Annex I 1272/2008.
Flash point according to EN ISO 2719:2002.
pH value measurement.
Material compatibility and corrosiveness in practical tests.

16.5. Training advice

Provide adequate information, instructions and training for users.

16.6. Indication of changes

A dash in the left hand margin indicates an amendment from the previous version.

The information given in the safety data sheet only applies to the described product in connection with its intended use. This information is based on the latest state of our knowledge at the time of revision. In particular, it describes our product under the aspect of its hazards and safety measures to be taken. It does not constitute any guarantee of product properties and quality features.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-DD**

Issue/Revision: 09.03.2023

Print date: 09.03.2023

Version: 2.1

Replaces version: 2.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Trade name: **AlproJet-DD**
UFI: V3YH-04WM-800W-8VG0

1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses: Cleaning and disinfecting agent
Intended purpose: Liquid concentrate with enhanced disinfecting effect for the daily cleaning and maintenance of dental aspiration systems with and without amalgam separator.
Uses advised against: None at intended use.
Note: The product is intended for professional users.

1.3. Details of the supplier of the safety data sheet

Manufacturer/Supplier: ALPRO MEDICAL GMBH
Mooswiesenstraße 9
D-78112 St. Georgen (Germany)
Telephone: +49 7725 9392-0
Telefax: +49 7725 9392-91
E-mail: info@alpro-medical.de
Internet: www.alpro-medical.com

E-mail address for the competent person responsible for the safety data sheet: doku@alpro-medical.de

1.4. Emergency telephone number

In-house emergency telephone number: +49 7725 9392-0
Monday – Friday from 08:00 am to 04:30 pm (UTC+1);
for chemical information and legal information on hazardous substances only
Poison centre: +49 761 19240
Poisoning information centre, Freiburg, Germany
(24 h / 7 d), English is spoken

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification in accordance with Regulation (EC) No 1272/2008 [CLP]

Classification	Classification procedure
Skin Corr. 1B; H314	Calculation method
STOT SE 3; H335	Calculation method

Full text of hazard classes as well as H-phrases: see under SECTION 16.1.

2.2. Label elements

Label elements in accordance with Regulation (EC) No 1272/2008 [CLP]

Hazard pictograms:



Safety data sheet
according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-DD**

Issue/Revision: 09.03.2023

Print date: 09.03.2023

Version: 2.1

Replaces version: 2.0

Signal word:	Danger
Hazard components for labelling:	2-Aminoethanol (141-43-5); Benzalkonium chloride (85409-22-9)
H-phrases:	H314 Causes severe skin burns and eye damage. H335 May cause respiratory irritation.
P-phrases:	P280 Wear protective gloves/protective clothing/eye protection/face protection. P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower]. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P310 Immediately call a POISON CENTER/ doctor.

2.3. Other hazards

The substances in the mixture do not meet the PBT/vPvB criteria according to Regulation (EC) No 1907/2006, Annex XIII.

No further hazards known.

SECTION 3: Composition/information on ingredients

3.2. Mixtures

Chemical characterisation: Mixture of substances listed below with non-hazardous additions in aqueous solution.

Hazardous ingredients

<i>Chemical name</i>	<i>Identification numbers</i>	<i>Classification in accordance with Regulation (EC) No 1272/2008</i>	<i>Weight %</i>
2-Aminoethanol	CAS No: 141-43-5 EC No: 205-483-3 Index No: 603-030-00-8 REACH Registration No: 01-2119486455-28-XXXX	Acute Tox. 4; H332 Acute Tox. 4; H312 Acute Tox. 4; H302 Skin Corr. 1B; H314 STOT SE 3; H335 <i>Specific concentration limits:</i> STOT SE 3; H335: C ≥ 5 %	≥ 5 - < 15
Trisodium nitrilotriacetate	CAS No: 5064-31-3 EC No: 225-768-6 Index No: 607-620-00-6 REACH Registration No: 01-2119519239-36-XXXX	<i>Carc. 2; H351</i> Acute Tox. 4; H302 Eye Irrit. 2; H319 <i>Specific concentration limits:</i> <i>Carc. 2; H351: C ≥ 5 %</i>	≥ 1 - < 5
Propan-2-ol	CAS No: 67-63-0 EC No: 200-661-7 Index No: 603-117-00-0 REACH Registration No: 01-2119457558-25-XXXX	Flam. Liq. 2; H225 Eye Irrit. 2; H319 STOT SE 3; H336	≥ 1 - < 5
Benzalkonium chloride	CAS No: 85409-22-9 EC No: 287-089-1	Acute Tox. 4; H302 Skin Corr. 1B; H314 Aquatic Acute 1; H400	≥ 1 - < 5

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-DD**

Issue/Revision: 09.03.2023

Print date: 09.03.2023

Version: 2.1

Replaces version: 2.0

Quaternary ammonium compounds, benzyl-C8-18-alkyldimethyl, chlorides	CAS No: 63449-41-2 EC No: 264-151-6 Index No: 612-140-00-5	Acute Tox. 4; H312 Acute Tox. 4; H302 Skin Corr. 1B; H314 Aquatic Acute 1; H400	< 0,2
--	--	--	-------

Full text of hazard classes and H-phrases: see SECTION 16.1.

SECTION 4: First aid measures

4.1. Description of first aid measures

General information:	First aider: Pay attention to self-protection! Remove contaminated, saturated clothing immediately.
Following inhalation:	Move affected person into fresh air and keep still and warm. Seek medical advice.
Following skin contact:	Wash skin immediately with plenty of water and soap. In case of skin reactions, consult a physician.
Following eye contact:	Flush eyes immediately with flowing water for 10 to 15 minutes holding eyelids apart. Remove contact lenses, if present and easy to do. Consult an ophthalmologist.
Following ingestion:	Rinse mouth with water. Let drink plenty of water. Do not induce vomiting (risk of perforation). Consult a physician immediately.

4.2. Most important symptoms and effects, both acute and delayed

Causes severe skin burns and eye damage. May cause respiratory irritation.

4.3. Indication of any immediate medical attention and special treatment needed

No information available.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media: Water spray jet, alcohol resistant foam, extinguishing powder, carbon dioxide (CO₂)

Unsuitable extinguishing media: Full water jet

5.2. Special hazards arising from the substance or mixture

Hazardous combustion products: Carbon monoxide (CO), carbon dioxide (CO₂), nitrogen oxides (NO_x), hydrogen chloride (HCl)

5.3. Advice for firefighters

Special protective equipment: Wear self-contained breathing apparatus.

Further information: Cool endangered containers with water spray jet.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

For non-emergency personnel

Use personal protective equipment. See SECTION 8.2.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-DD**

Issue/Revision: 09.03.2023

Print date: 09.03.2023

Version: 2.1

Replaces version: 2.0

Avoid skin and eye contact. Do not breathe vapours. Provide adequate ventilation. Special danger of slipping by leaked/spilled product. Evacuate danger area. Observe emergency plans. Consult experts.

For emergency responders

Use personal protective equipment. See SECTION 8.2.

6.2. Environmental precautions

Do not discharge into drains or rivers.

6.3. Methods and material for containment and cleaning up

Containment

For large spills, dyke spilled material or otherwise contain material to ensure runoff does not reach a waterway. Cover or seal drains.

Cleaning up

Wipe up small amounts with absorbent material (e.g. cloth, fleece). Absorb large amounts with liquid-binding material (sand, diatomaceous earth, universal binder, sawdust). Collect in suitable, closed containers for disposal. Clean contaminated surfaces thoroughly.

Other information

Inappropriate containment and cleaning methods are not known.

6.4. Reference to other sections

Information on safe handling see SECTION 7.1.

Information on personal protective equipment see SECTION 8.2.

Information on disposal see SECTION 13.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions

Avoid contact with skin and eyes. Avoid breathing aerosols and vapours. Keep container tightly closed. Fill refill packages only in labelled original bottles.

Advice on general occupational hygiene

When using do not eat, drink or smoke. Wash hands before breaks and at end of work. Keep away from food and drink.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels: Keep only in the original container. Keep container tightly closed and kept upright to prevent any leakage.

Advice on common storage: Keep away from food, drink and animal feedingstuffs.

Further information on storage conditions: Not necessary

Storage class ([DE] TRGS 510): LGK 8B Non-combustible corrosive hazardous substances

7.3. Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific end uses are stipulated.

Industry and sector specific guidance

[DE] TRGS 525 – Hazardous substances in medical care facilities (Section 7 Activities with disinfectants); Issue: September 2014;
Source: GMBI 2014 page 1294-1307 of 13.10.2014 [No 63]; www.baua.de

Safety data sheet
according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-DD**

Issue/Revision: 09.03.2023

Print date: 09.03.2023

Version: 2.1

Replaces version: 2.0

[DE] DGUV rules 107-002 (former BGR 206) - Disinfection works in health service
 Issue: July 1999; Source: www.dguv.de/publikationen

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Occupational exposure limits

Country	Limit values				Legal basis	Remarks
	Long term (8 hours)		Short term (15 minutes)			
	ppm	mg/m ³	ppm	mg/m ³		
2-Aminoethanol (CAS No: 141-43-5)						
EU	1	2.5	3	7.6	2006/15/EC	Skin
UK	1	2.5	3	7.6	EH40	Sk
Propan-2-ol (CAS No: 67-63-0)						
EU						no limit value specified
UK	400	999	500	1250	EH40	

Used abbreviations, symbols, numerals and explanations in column „Remarks“

Skin A significant uptake of the substance through the skin is possible.
 Sk Can be absorbed through the skin. The assigned substances are those for which there are concerns that dermal absorption will lead to systemic toxicity.

Biological limit values

Country	Parameter	Limit value	Test material	Sampling time	Legal basis
Propan-2-ol (CAS No: 67-63-0)					
Germany	Acetone	25 mg/l	Whole blood	End of exposition, resp. end of shift	TRGS 903
	Acetone	25 mg/l	Urine	End of exposition, resp. end of shift	TRGS 903

Information on monitoring procedures

BS EN 482:2012-04-30; Title: Workplace exposure - General requirements for the performance of procedures for the measurement of chemical agents;
 British version of EN 482:2012

BS EN 689:1996-04-15; Title: Workplace atmospheres - Guidance for the assessment of exposure by inhalation to chemical agents for comparison with limit values and measurement strategy;
 British version of EN 689:1995

BS EN 14042:2003-04-24; Title: Workplace atmospheres - Guide for the application and use of procedures for the assessment of exposure to chemical and biological agents;
 British version of EN 14042:2003

8.2. Exposure controls

Appropriate engineering controls

Technical and organisational protective measures

The eyewash station (or eyewash bottle) and emergency shower must be located near the workplace.

Personal protective equipment

Eye/face protection: Safety glasses with side protection according to EN 166

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-DD**

Issue/Revision: 09.03.2023

Print date: 09.03.2023

Version: 2.1
Replaces version: 2.0

Skin protection:

Hand protection: Protective gloves according to EN 374
Splash guard:
Disposable gloves made of nitrile rubber (thickness 0.11 mm)
Permanent contact (> 480 min):
Protective gloves made of nitrile rubber (thickness 0.40 mm)

Other skin protection: Long-sleeved protective clothing (lab coat)

Respiratory protection: Not necessary when used as intended.

Thermal hazards: No special protective measures necessary.

Environmental exposure controls

Do not discharge into drains.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Appearance:	clear, blue-green liquid	
Odour:	of lemon	
Odour threshold:	no data available	
pH (20 g/l H ₂ O):	11.0 – 12.0	(20 °C)
Melting point/freezing point:	no data available	
Initial boiling point and boiling range:	no data available	
Flash point:	> 60 °C	
Evaporation rate:	no data available	
Flammability (solid, gas):	not applicable	
Lower explosive limit:	not applicable	
Upper explosive limit:	not applicable	
Vapour pressure:	no data available	(... °C)
Vapour density:	no data available	
Relative density:	1.030 – 1.035	(20 °C)
Solubility in water:	completely soluble	
Partition coefficient: n-octanol/water	not applicable	
Auto-ignition temperature:	not applicable	
Decomposition temperature:	no data available	
Viscosity:	no data available	
Explosive properties:	none	
Oxidising properties:	none	

9.2. Other information

Refractive index n _D :	1.3638-1.3723	(20 °C)
Electrical conductivity (20 g/l H ₂ O):	1500-2000 µS/cm	(20 °C)

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-DD**

Issue/Revision: 09.03.2023

Print date: 09.03.2023

Version: 2.1

Replaces version: 2.0

SECTION 10: Stability and reactivity

10.1. Reactivity

No hazardous reactions when handled and stored as intended.

10.2. Chemical stability

The product is stable when handled and stored as intended.

10.3. Possibility of hazardous reactions

None known

10.4. Conditions to avoid

None known

10.5. Incompatible materials

None known

10.6. Hazardous decomposition products

Does not decompose when used as intended.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity

Product

Acute toxicity - oral:	Acute Toxicity Estimate ATE _{mix} > 2000 mg/kg => no classification
Acute toxicity - dermal:	Acute Toxicity Estimate ATE _{mix} > 2000 mg/kg => no classification
Acute toxicity - inhalation:	Acute Toxicity Estimate ATE _{mix} > 20 mg/l => no classification

Ingredients

2-Aminoethanol (CAS No: 141-43-5):

Acute toxicity - oral: LD₅₀: 1515 mg/kg; species: rat; method: OECD 401

Acute toxicity - inhalation: LC₅₀: > 1.3 mg/l; species: rat; 6 h; vapour

Trisodium nitrilotriacetate (CAS No: 5064-31-3):

Acute toxicity - oral: LD₅₀: 1000 - 2000 mg/kg; species: rat; method: (BASF-test)

Benzalkonium chloride (CAS No: 85409-22-9):

Acute toxicity - oral: LD₅₀: approx. 344 mg/kg; species: rat

Quaternary ammonium compounds, benzyl-C8-18-alkyldimethyl, chlorides (CAS No: 63449-41-2):

Acute toxicity - oral: LD₅₀: approx. 398 mg/kg; species: rat

Acute toxicity - dermal: LD₅₀: 1420 mg/kg; species: rat

Skin corrosion/irritation

Product

Causes severe skin burns. [calculation method]

Safety data sheet
according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-DD**

Issue/Revision: 09.03.2023

Print date: 09.03.2023

Version: 2.1

Replaces version: 2.0

Serious eye damage/irritation

Product

Causes serious eye damage. [calculation method]

Respiratory or skin sensitisation

Product

No data available.

Germ cell mutagenicity

Product

No data available.

Carcinogenicity

Product

No classification. [calculation method]

Ingredients

Trisodium nitrilotriacetate (CAS No: 5064-31-3):

The substance was shown to have a carcinogenic effect in animal studies with long-term administration of large amounts via the drinking water or via the food. With single or short-term intake of the substance a carcinogenic effect is however practically ruled out.

Reproductive toxicity

Product

No data available.

STOT-single exposure

Product

May cause respiratory irritation. [calculation method]

Ingredients

2-Aminoethanol (CAS No: 141-43-5):

May cause respiratory irritation.

Propan-2-ol (CAS No: 67-63-0):

May cause drowsiness or dizziness.

STOT-repeated exposure

Product

No data available.

Aspiration hazard

Product

No data available.

SECTION 12: Ecological information

12.1. Toxicity

No classification. [calculation method]

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-DD**

Issue/Revision: 09.03.2023

Print date: 09.03.2023

Version: 2.1

Replaces version: 2.0

12.2. Persistence and degradability

Biodegradability:

The product is biodegradable according to OECD criteria. The statement has been derived from the properties of the ingredients.

12.3. Bioaccumulative potential

No data available.

12.4. Mobility in soil

No data available.

12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to Regulation (EC) No 1907/2006, Annex XIII.

12.6. Other adverse effects

No data available.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal of the product

Product residues must be disposed of as hazardous waste in compliance with the Directive 2008/98/EC on waste as well as national and regional regulations. Do not dispose of via the waste water. Leave product in the original container as possible. Do not mix with other waste materials.

Waste codes / waste designations according to EWC

Product residues: 16 10 03* aqueous concentrates containing hazardous substances

Disposal of the packaging

Packaging contaminated with product is considered as hazardous waste and must be disposed of accordingly.

Waste codes / waste designations according to EWC

Contaminated packaging: 15 01 10* packaging containing residues of or contaminated by hazardous substances

Recommendation

Contaminated packaging must be emptied optimally and can be recycled after appropriate cleaning (rinse with water).

SECTION 14: Transport information

14.0. Transport classification

Dangerous good in sense of the transport regulations in road traffic (ADR), railway traffic (RID), inland waterway traffic (ADN), maritime traffic (IMDG-Code) and air traffic (ICAO-TI/IATA-DGR).

14.1. UN number

UN 1903

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-DD**

Issue/Revision: 09.03.2023

Print date: 09.03.2023

Version: 2.1

Replaces version: 2.0

14.2. UN proper shipping name

ADR/RID/ADN

DISINFECTANT, LIQUID, CORROSIVE, N.O.S. (Ethanolamine, Benzalkonium chloride)

IMDG-Code/ICAO-TI/IATA-DGR

DISINFECTANT, LIQUID, CORROSIVE, N.O.S. (Ethanolamine, Benzalkonium chloride)

14.3. Transport hazard class(es)

Class: 8

Subsidiary risk(s): -

14.4. Packing group

III

14.5. Environmental hazards

ADR/RID/ADN

Environmentally Hazardous: No

IMDG-Code

Marine Pollutant: No

14.6. Special precautions for user

Not necessary.

14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

14.8. Further information

Transport category according to ADR section 1.1.3.6: 3

Maximum total quantity per transport unit according to ADR section 1.1.3.6: 1000 L

Limited quantity (Maximum quantity per inner packaging) according to ADR/RID/ADN/IMDG-Code: 5 L

Classification code according to ADR/RID/ADN: C9

Hazard identification number according to ADR/RID: 80

Tunnel restriction code according to ADR/RID: E

Segregation group according to IMDG-Code section 5.4.1.5.11.1: IMDG-Code- Segregation group 18 – alkalis

EmS codes: F-A, S-B

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

EU-Regulations

REGULATION (EC) No 1005/2009 on substances that deplete the ozone layer

not applicable

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-DD**

Issue/Revision: 09.03.2023

Print date: 09.03.2023

Version: 2.1

Replaces version: 2.0

REGULATION (EC) No 850/2004 on persistent organic pollutants and amending Directive 79/117/EEC
not applicable

REGULATION (EU) No 649/2012 concerning the export and import of hazardous chemicals
not applicable

DIRECTIVE 2012/18/EU (Seveso III Directive) on the control of major-accident hazards involving dangerous substances, amending and subsequently repealing Council Directive 96/82/EC
not applicable

DIRECTIVE 2010/75/EU on industrial emissions (integrated pollution prevention and control)
not applicable

REACH – List of substances subject to authorisation (Annex XIV)
not applicable

REACH – Restrictions on the manufacture, placing on the market and use of certain dangerous substances, mixtures and articles (Annex XVII)
not applicable

COUNCIL DIRECTIVE 94/33/EC on the protection of young people at work
Observe employment restrictions for juveniles.

COUNCIL DIRECTIVE 92/85/EEC on the introduction of measures to encourage improvements in the safety and health at work of pregnant workers and workers who have recently given birth or are breastfeeding

Observe employment restrictions for pregnant and nursing mothers.

15.2. Chemical safety assessment

For this mixture no chemical safety assessment has been carried out.

SECTION 16: Other information

16.1. Full text of hazard classes and H-phrases

Hazard classes

Acute Tox.	Acute toxicity
Aquatic Acute	Acute aquatic hazard
Carc.	Carcinogenicity
Eye Irrit.	Eye irritation
Flam. Liq.	Flammable liquid
Skin Corr.	Skin corrosion
STOT SE	Specific target organ toxicity (single exposure)

H-phrases (Hazard statements)

H225	Highly flammable liquid and vapour.
H302	Harmful if swallowed.
H312	Harmful in contact with skin.
H314	Causes severe skin burns and eye damage.
H319	Causes serious eye irritation.
H332	Harmful if inhaled.
H335	May cause respiratory irritation.
H336	May cause drowsiness or dizziness.
H351	Suspected of causing cancer <state route of exposure if it is conclusively proven that no other routes of exposure cause the hazard>.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-DD**

Issue/Revision: 09.03.2023

Print date: 09.03.2023

Version: 2.1

Replaces version: 2.0

H400 Very toxic to aquatic life.

16.2. Abbreviations and acronyms

ADN	<u>A</u> ccord européen relatif au transport international des marchandises <u>d</u> angereuses par voie de <u>n</u> avigation intérieure (European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways)
ADR	<u>A</u> ccord européen relatif au transport international des marchandises <u>d</u> angereuses par <u>r</u> oute (European Agreement concerning the International Carriage of Dangerous Goods by Road)
BGR	<u>B</u> erufsgenossenschaftliche <u>R</u> egeln (English: Employers' liability insurance association rules)
BS	<u>B</u> ritish <u>S</u> tandards
CAS	<u>C</u> hemical <u>A</u> bstracts <u>S</u> ervice
CLP	Regulation on <u>C</u> lassification, <u>L</u> abelling and <u>P</u> ackaging of Substances and Mixtures
[DE]	National German regulations
DGUV	<u>D</u> eutsche <u>G</u> esetzliche <u>U</u> nfall <u>v</u> ersicherung (English: German statutory accident insurance)
EC	<u>E</u> uropean <u>C</u> ommunity
EEC	<u>E</u> uropean <u>E</u> conomic <u>C</u> ommunity
EmS	<u>E</u> mergency <u>S</u> chedules (Emergency response procedures for ships carrying dangerous goods)
EN	European Standard
EU	<u>E</u> uropean <u>U</u> nion
EWC	<u>E</u> uropean <u>W</u> aste <u>C</u> atalogue
GHS	<u>G</u> lobally <u>H</u> armonized <u>S</u> ystem of Classification, Labelling and Packaging of Chemicals
GMBI	<u>G</u> emeinsames <u>M</u> inisterial <u>b</u> latt (English: Joint Ministerial Gazette)
IATA-DGR	<u>I</u> nternational <u>A</u> ir <u>T</u> ransport <u>A</u> ssociation - <u>D</u> angerous <u>G</u> oods <u>R</u> egulations
IBC-Code	International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk
ICAO-TI	Technical Instructions For The Safe Transport of Dangerous Goods by Air
IMDG-Code	International <u>M</u> aritime Code for <u>D</u> angerous <u>G</u> oods
LC ₅₀	Median lethal concentration
LD ₅₀	Median lethal dose
LGK	<u>L</u> ager <u>k</u> lasse (English: Storage class)
MARPOL	International Convention for the Prevention of Marine Pollution from Ships
N.O.S.	Not otherwise specified
NOAEL	<u>N</u> o <u>O</u> bserved <u>A</u> dverse <u>E</u> ffect <u>L</u> evel (dose at which no adverse effect is found)
OECD	<u>O</u> rganization for <u>E</u> conomic <u>C</u> o-operation and <u>D</u> evelopment
PBT	<u>P</u> ersistent, <u>b</u> ioaccumulative and <u>t</u> oxic
ppm	<u>P</u> arts <u>p</u> er <u>m</u> illion
REACH	<u>R</u> egistration, <u>E</u> valuation, <u>A</u> uthorisation and Restriction of <u>C</u> hemicals
RID	<u>R</u> èglement concernant le transport <u>I</u> nternational ferroviaire de marchandises <u>D</u> angereuses (Regulations Concerning the International Carriage of Dangerous Goods by Rail)
TRGS	<u>T</u> echnische <u>R</u> egeln für <u>G</u> efahr <u>s</u> toffe (English: Technical Rules for Hazardous Substances)
UN	<u>U</u> nited <u>N</u> ations
UTC	Coordinated Universal Time (French: Temps Universel Coordonné)
vPvB	<u>V</u> ery <u>p</u> ersistent and <u>v</u> ery <u>b</u> ioaccumulative

16.3. Key literature references and sources for data

- Regulation (EC) No 1907/2006 (REACH), Annex II
- European Chemicals Agency (ECHA) – Guidance on the compilation of safety data sheets; Version 2.1 (February 2014); http://echa.europa.eu/documents/10162/13643/sds_en.pdf
- GISBAU (Hazardous substances information system of the BG BAU) – course „safety data sheet“; <http://www.bgbau.de/gisbau/SDB/lehrgang/lehrgang.htm>
- Regulation (EC) No 1272/2008 (CLP regulation)

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-DD**

Issue/Revision: 09.03.2023

Print date: 09.03.2023

Version: 2.1

Replaces version: 2.0

- European Chemicals Agency (ECHA) – Guidance on Labelling and Packaging in accordance with Regulation (EC) No 1272/2008 (04/2011);
http://echa.europa.eu/documents/10162/13562/clp_labelling_en.pdf
- European Chemicals Agency (ECHA), Registered substances;
<http://echa.europa.eu/information-on-chemicals/registered-substances>
- European Chemicals Agency (ECHA), C&L Classification and Labelling Inventory;
<http://echa.europa.eu/information-on-chemicals/cl-inventory-database>
- Institute for Occupational Safety and Health of the German Social Accident Insurance (IFA):
GESTIS database on hazardous substances and GESTIS - International limit values for chemical agents;
<http://www.dguv.de/dguv/ifa/index.jsp>
- German Environmental Agency (Umweltbundesamt), Section IV 2.4: Office of Documentation and Information on Substances Hazardous to Waters RIGOLETTO (catalogue of Substances Hazardous to Waters); <http://webrigoletto.uba.de/rigoletto>

16.4. Training advice

Provide adequate information, instructions and training for users.

16.5. Indication of changes

A dash in the left hand margin indicates an amendment from the previous version.

The information given in the safety data sheet only applies to the described product in connection with its intended use. This information is based on the latest state of our knowledge at the time of revision. In particular, it describes our product under the aspect of its hazards and safety measures to be taken. It does not constitute any guarantee of product properties and quality features.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-W**
Issue/Revision: 27.04.2023

Version: 3.0
Replaces version: 2.1

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Trade name: **AlproJet-W**
UFI: V9JD-X52K-M00T-H3NC

1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses: Cleaning and preservative agent
Intended purpose: Liquid concentrate for the weekly acid contra-cleaning and maintenance of dental aspiration systems with and without amalgam separator.
Uses advised against: None at intended use.
Note: The product is intended for professional users.

1.3. Details of the supplier of the safety data sheet

Manufacturer/Supplier: ALPRO MEDICAL GMBH
Mooswiesenstraße 9
D-78112 St. Georgen (Germany)
Telephone: +49 7725 9392-0
Telefax: +49 7725 9392-91
E-mail: info@alpro-medical.de
Internet: www.alpro-medical.de

E-mail address for the competent person responsible for the safety data sheet: doku@alpro-medical.de

1.4. Emergency telephone number

In-house emergency telephone number: +49 7725 9392-0
Monday – Friday from 08:00 am to 04:30 pm (UTC+1);
for chemical information and legal information on hazardous substances only
Poison centre: +49 761 19240
Poisoning information centre, Freiburg, Germany
(24 h / 7 d), English is spoken
National Poisons Information Centre (UK): +44 344 892 0111
National Poisons Information Service (NPIS)
(24 h / 7 d), Medical Professionals Only

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification in accordance with Regulation (EC) No 1272/2008 [CLP]

<i>Classification</i>	<i>Classification procedure</i>
Eye Irrit. 2; H319	Calculation method
Skin Irrit. 2; H315	Calculation method

Full text of hazard classes as well as H-phrases: see under SECTION 16.1.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-W**
Issue/Revision: 27.04.2023

Version: 3.0
Replaces version: 2.1

2.2. Label elements

Label elements in accordance with Regulation (EC) No 1272/2008 [CLP]

Hazard pictograms:



Signal word: Warning

Hazard components for labelling: -

H-phrases: H315 Causes skin irritation.
H319 Causes serious eye irritation.

P-phrases: P280 Wear protective gloves/eye protection/face protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313 If eye irritation persists: Get medical advice/attention.

2.3. Other hazards

The substances in the mixture do not meet the PBT/vPvB criteria according to Regulation (EC) No 1907/2006, Annex XIII (see SECTION 12.5.)

The substances in the mixture do not have any endocrine disrupting properties according to Regulation (EC) No. 1907/2006, Annex XIV (see Section 11 and Section 12.6.). They are not on the list of substances of very high concern for authorization according to Regulation (EC) No. 1907/2006, Article 59, Paragraph 10).

No further hazards known.

SECTION 3: Composition/information on ingredients

3.1. Substances

Not applicable.

3.2. Mixtures

Chemical characterisation: Mixture of substances listed below with non-hazardous additions in aqueous solution.

Hazardous ingredients

Chemical name	Identification numbers	Classification in accordance with Regulation (EC) No 1272/2008	Weight %
Sulphamidic acid	CAS No: 5329-14-6 EC No: 226-218-8 Index No: 016-026-00-0 REACH Registration No: 01-2119488633-28-XXXX	Eye Irrit. 2; H319 Skin Irrit. 2; H315 Aquatic Chronic 3; H412	≥ 1 - < 5

Full text of hazard classes and H-phrases: see SECTION 16.1.

Occupational exposure limits: see SECTION 8.1.

SECTION 4: First aid measures

4.1. Description of first aid measures

General information: First aider: Pay attention to self-protection!

Following inhalation: Move affected person into fresh air and keep still and warm. In case of continued complaints seek medical advice.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-W**
Issue/Revision: 27.04.2023

Version: 3.0
Replaces version: 2.1

- Following skin contact: Wash skin immediately with plenty of water and soap. In case of skin reactions, consult a physician.
- Following eye contact: Flush eyes immediately with flowing water for 10 to 15 minutes holding eyelids apart. Remove contact lenses, if present and easy to do. Consult an ophthalmologist.
- Following ingestion: Rinse mouth with water. Let drink plenty of water. Do not induce vomiting. Consult a physician immediately.

4.2. Most important symptoms and effects, both acute and delayed

Causes skin irritation. Causes serious eye irritation.

4.3. Indication of any immediate medical attention and special treatment needed

No information available.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media: Water spray jet, alcohol resistant foam, extinguishing powder, carbon dioxide (CO₂)

Unsuitable extinguishing media: Full water jet

5.2. Special hazards arising from the substance or mixture

Hazardous combustion products: Carbon monoxide (CO), carbon dioxide (CO₂), nitrogen oxides (NO_x), sulfur oxides (S_xO_y)

5.3. Advice for firefighters

Special protective equipment: Wear self-contained breathing apparatus.

Further information: Cool endangered containers with water spray jet.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

For non-emergency personnel

Use personal protective equipment. See SECTION 8.2.

Avoid skin and eye contact. Special danger of slipping by leaked/spilled product.

For emergency responders

Use personal protective equipment. See SECTION 8.2.

6.2. Environmental precautions

Do not discharge into drains or surface and ground water.

6.3. Methods and material for containment and cleaning up

Containment

For large spills, dyke spilled material or otherwise contain material to ensure runoff does not reach a waterway. Cover or seal drains.

Cleaning up

Wipe up small amounts with absorbent material (e.g. cloth, fleece). Absorb large amounts with liquid-binding material (sand, diatomaceous earth, universal binder, sawdust). Collect in suitable, closed containers for disposal. Clean contaminated surfaces thoroughly.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-W**
Issue/Revision: 27.04.2023

Version: 3.0
Replaces version: 2.1

Other information

Inappropriate containment and cleaning methods are not known.

6.4. Reference to other sections

Information on safe handling see SECTION 7.1.
Information on personal protective equipment see SECTION 8.2.
Information on disposal see SECTION 13.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions

Avoid contact with skin and eyes. Fill refill packages only in labelled original bottles.
Do not mix with alkaline solutions or other cleaners or disinfectants.

Advice on general occupational hygiene

When using do not eat, drink or smoke. Wash hands before breaks and at end of work. Keep away from medicines, food, feed, cosmetics and stimulants.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels: Keep only in the original container. Keep container tightly closed and kept upright to prevent any leakage.

Advice on common storage: Not necessary

Further information on storage conditions: store between -5 °C and +30 °C.

7.3. Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific end uses are stipulated.

Industry and sector specific guidance

[DE] TRGS 525 – Hazardous substances in medical care facilities (Section 7 Activities with disinfectants); Issue: September 2014; Source: GMBI 2014 page 1294-1307 of 13.10.2014 [No 63]; 10.07.2015 [No 27]; www.baua.de

[DE] DGUV rules 207-206 – Prevention of chemical risks when handling disinfectants in health service, Issue: 2016.12; Source: www.dguv.de/publikationen

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Occupational exposure limits

Does not contain substances above concentration limits fixing an occupational exposure limit.

Biological limit values

Does not contain substances above concentration limits fixing a biological limit value.

Information on monitoring procedures

BS EN 482:2021-04; Title: Workplace exposure - Procedures for the determination of the concentration of chemical agents;
British version of EN 482:2021

BS EN 689:2018; Title: Workplace atmospheres - Guidance for the assessment of exposure by inhalation to chemical agents for comparison with limit values and measurement strategy;

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-W**
Issue/Revision: 27.04.2023

Version: 3.0
Replaces version: 2.1

British version of EN 689:2018

BS EN 14042:2003-04-24; Title: Workplace atmospheres - Guide for the application and use of procedures for the assessment of exposure to chemical and biological agents; British version of EN 14042:2003)

CEN/TR 17055:2017; Title: Workplace exposure. Measurement of chemical agents complying with the requirements given in EN 482 and either one of EN 838, EN 1076, EN 13205, EN 13890 and EN 13936. Choice of procedures

ISO TR 14294:2011; Title: Workplace atmospheres. Measurement of dermal Exposure. Principles and methods

8.2. Exposure controls

Appropriate engineering controls

Technical and organisational protective measures

The eyewash station (or eyewash bottle) must be located near the workplace.

Personal protective equipment

Eye/face protection: Safety glasses with side protection according to BS EN 166

Skin protection:

Hand protection: Protective gloves according to BS EN ISO 374-1 and BS EN 21420

Splash guard:

Protective gloves: type C; permeation-resistant at least 10 minutes

Permanent contact (> 480 min):

Protective gloves: type A or B; code letters: S, G; permeation-resistant at least 30 minutes

Other skin protection: Not necessary when used as intended.

Respiratory protection: Not necessary when used as intended.

Thermal hazards: No special protective measures necessary.

Environmental exposure controls

Do not discharge into drains.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Appearance/physical state: clear, orange liquid

Odour: slightly aromatic

Odour threshold: no data available

Melting point/freezing point: no data available

Initial boiling point and boiling range: no data available

Flammability: not applicable

Lower explosive limit: not applicable

Upper explosive limit: not applicable

Flash point: > 60 °C

Auto-ignition temperature: no data available

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-W**
Issue/Revision: 27.04.2023

Version: 3.0
Replaces version: 2.1

pH (50 g/l H ₂ O):	1.5 – 2.5	(20 °C)
Kinematic viscosity:	no data available	
Solubility in water:	completely soluble	
Partition coefficient: n-octanol/water	not applicable	
Vapour pressure:	no data available	(... °C)
Density:	1.050 – 1.060 g/cm ³	(20 °C)
Relative vapour density:	no data available	
Particle characteristics:	not applicable	

9.2. Other information

9.2.1. Information with regard to physical hazard classes

Explosive substances/mixtures and products containing explosives:	not applicable
Flammable gases:	not applicable
Aerosols:	not applicable
Oxidising gases:	not applicable
Gases under pressure:	not applicable
Flammable liquids:	not applicable
Flammable solids:	not applicable
Self-reactive substances and mixtures:	not applicable
Pyrophoric liquids:	not applicable
Pyrophoric solids:	not applicable
Self-heating substances and mixtures:	not applicable
Substances and mixtures, which emit flammable gases in contact with water:	not applicable
Oxidising liquids:	not applicable
Oxidising solids:	not applicable
Organic peroxides:	not applicable
Substances and mixtures corrosive to metals:	May be corrosive to metals. [Literature research]
Desensitised explosives:	not applicable

9.2.2. Other safety characteristics

Electrical conductivity (50 g/l H ₂ O):	6000-7000 µS/cm	(20 °C)
--	-----------------	---------

SECTION 10: Stability and reactivity

10.1. Reactivity

Reacts with bases.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-W**
Issue/Revision: 27.04.2023

Version: 3.0
Replaces version: 2.1

10.2. Chemical stability

The product is stable when handled and stored as intended.

10.3. Possibility of hazardous reactions

Reacts with bases.

10.4. Conditions to avoid

None known

10.5. Incompatible materials

Stainless steel can oxidize.

10.6. Hazardous decomposition products

Does not decompose when used as intended.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity

Product

Acute toxicity - oral:	Acute Toxicity Estimate ATE _{mix} > 2000 mg/kg => no classification
Acute toxicity - dermal:	Acute Toxicity Estimate ATE _{mix} > 2000 mg/kg => no classification
Acute toxicity - inhalation:	Acute Toxicity Estimate ATE _{mix} > 20 mg/l => no classification

Skin corrosion/irritation

Product

Causes skin irritation. [calculation method]

Serious eye damage/irritation

Product

Causes serious eye irritation. [calculation method]

Respiratory or skin sensitisation

Product

No data available.

Germ cell mutagenicity

Product

No data available.

Carcinogenicity

Product

No data available.

Reproductive toxicity

Product

No data available.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-W**
Issue/Revision: 27.04.2023

Version: 3.0
Replaces version: 2.1

STOT-single exposure

Product

No classification. [calculation method]

STOT-repeated exposure

Product

No data available.

Aspiration hazard

Product

No data available.

Information on other hazards

Endocrine disrupting properties:

No substances are contained that have endocrine disrupting properties for humans.

SECTION 12: Ecological information

12.1. Toxicity

No classification. [calculation method]

12.2. Persistence and degradability

Biodegradability:

The product is biodegradable according to OECD criteria. The statement has been derived from the properties of the ingredients.

12.3. Bioaccumulative potential

No data available.

12.4. Mobility in soil

No data available.

12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to Regulation (EC) No 1907/2006, Annex XIII.

12.6. Endocrine disrupting properties

No substances are contained that have endocrine disrupting properties for non-target organisms.

12.7. Other adverse effects

No data available.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal of the product

Product residues must be disposed of as hazardous waste in compliance with the Directive 2008/98/EC on waste as well as national and regional regulations. Do not dispose of via the waste water. Leave product in the original container as possible. Do not mix with other waste materials.

Waste codes / waste designations according to EWC

Product residues: 16 10 03* aqueous concentrates containing hazardous substances

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-W**
Issue/Revision: 27.04.2023

Version: 3.0
Replaces version: 2.1

Disposal of the packaging

Packaging contaminated with product is considered as hazardous waste and must be disposed of accordingly.

Waste codes / waste designations according to EWC

Contaminated packaging: 15 01 10* packaging containing residues of or contaminated by hazardous substances

Recommendation

Contaminated packaging must be emptied optimally and can be recycled after appropriate cleaning (rinse with water).

From a dilution to 1 % the concentrate is no longer classified as a hazardous substance.

SECTION 14: Transport information

14.0. Transport classification

No dangerous good in sense of the transport regulations in road traffic (ADR), railway traffic (RID), inland waterway traffic (ADN), maritime traffic (IMDG-Code) and air traffic (ICAO-TI/IATA-DGR).

14.1. UN number

-

14.2. UN proper shipping name

-

14.3. Transport hazard class(es)

-

14.4. Packing group

-

14.5. Environmental hazards

-

14.6. Special precautions for user

-

14.7. Maritime transport in bulk according to IMO instruments

-

14.8. Further information

-

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

EU-Regulations

REGULATION (EC) No 1005/2009 on substances that deplete the ozone layer

not applicable

REGULATION (EC) No 2019/1021 on persistent organic pollutants

not applicable

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-W**
Issue/Revision: 27.04.2023

Version: 3.0
Replaces version: 2.1

REGULATION (EU) No 649/2012 concerning the export and import of hazardous chemicals
not applicable

REGULATION (EU) No 648/2004 on detergents
nonionic surfactants: < 5 %
phosphates: < 5 %

DIRECTIVE 2012/18/EU (Seveso III Directive) on the control of major-accident hazards involving dangerous substances, amending and subsequently repealing Council Directive 96/82/EC
not applicable

DIRECTIVE 2010/75/EU on industrial emissions (integrated pollution prevention and control)
VOC content: < 1 %

REACH – List of substances subject to authorisation (Annex XIV)
not applicable

REACH – Restrictions on the manufacture, placing on the market and use of certain dangerous substances, mixtures and articles (Annex XVII)
not applicable

COUNCIL DIRECTIVE 94/33/EC on the protection of young people at work
not applicable

COUNCIL DIRECTIVE 92/85/EEC on the introduction of measures to encourage improvements in the safety and health at work of pregnant workers and workers who have recently given birth or are breastfeeding
not applicable

15.2. Chemical safety assessment

For this mixture no chemical safety assessment has been carried out.

SECTION 16: Other information

16.1. Full text of hazard classes and H-phrases

Hazard classes

Aquatic Chronic	Long-term aquatic hazard
Eye Irrit.	Eye irritation
Skin Irrit.	Skin irritation

H-phrases (Hazard statements)

H315	Causes skin irritation.
H319	Causes serious eye irritation.
H412	Harmful to aquatic life with long lasting effects.

16.2. Abbreviations and acronyms

ADN	<u>A</u> ccord européen relatif au transport international des marchandises <u>d</u> angereuses par voie de <u>n</u> avigation intérieure (European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways)
ADR	<u>A</u> ccord européen relatif au transport international des marchandises <u>d</u> angereuses par <u>r</u> oute (European Agreement concerning the International Carriage of Dangerous Goods by Road)
BS	<u>B</u> ritish <u>S</u> tandards
CAS	<u>C</u> hemical <u>A</u> bstracts <u>S</u> ervice
CLP	Regulation on <u>C</u> lassification, <u>L</u> abelling and <u>P</u> ackaging of Substances and Mixtures

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-W**
Issue/Revision: 27.04.2023

Version: 3.0
Replaces version: 2.1

[DE]	National German regulations
DGUV	<u>D</u> eutsche <u>G</u> esetzliche <u>U</u> nfall <u>v</u> ersicherung (English: German statutory accident insurance)
EC	<u>E</u> uropean <u>C</u> ommunity
EEC	<u>E</u> uropean <u>E</u> conomic <u>C</u> ommunity
EN	European Standard
EU	<u>E</u> uropean <u>U</u> nion
EWC	<u>E</u> uropean <u>W</u> aste <u>C</u> atalogue
GMBI	<u>G</u> emeinsames <u>M</u> inisterial <u>b</u> latt (English: Joint Ministerial Gazette)
IATA-DGR	<u>I</u> nternational <u>A</u> ir <u>T</u> ransport <u>A</u> ssociation - <u>D</u> angerous <u>G</u> oods <u>R</u> egulations
ICAO-TI	Technical Instructions For The Safe Transport of Dangerous Goods by Air
IMDG-Code	<u>I</u> nternational <u>M</u> aritime Code for <u>D</u> angerous <u>G</u> oods
LC ₅₀	Median lethal concentration
LD ₅₀	Median lethal dose
OECD	<u>O</u> rganization for <u>E</u> conomic <u>C</u> o-operation and <u>D</u> evelopment
PBT	<u>P</u> ersistent, <u>b</u> ioaccumulative and <u>t</u> oxic
ppm	<u>P</u> arts <u>p</u> er <u>m</u> illion
REACH	<u>R</u> egistration, <u>E</u> valuation, <u>A</u> uthorisation and Restriction of <u>C</u> hemicals
RID	<u>R</u> èglement concernant le transport <u>I</u> nternational ferroviaire de marchandises <u>D</u> angereuses (Regulations Concerning the International Carriage of Dangerous Goods by Rail)
TRGS	<u>T</u> echnische <u>R</u> egeln für <u>G</u> efahr <u>s</u> toffe (English: Technical Rules for Hazardous Substances)
UN	<u>U</u> nited <u>N</u> ations
UTC	Coordinated Universal Time (French: Temps Universel Coordonné)
VOC	Volatile Organic Compounds
vPvB	<u>V</u> ery persistent and <u>v</u> ery <u>b</u> ioaccumulative

16.3. Key literature references and sources for data

- Regulation (EC) No 1907/2006 (REACH), Annex II
- European Chemicals Agency (ECHA) – Guidance on the compilation of safety data sheets; Version 4.0 (December 2022); <https://echa.europa.eu/documents>
- Regulation (EC) No 1272/2008 (CLP regulation)
- European Chemicals Agency (ECHA) – Guidance on Labelling and Packaging in accordance with Regulation (EC) No 1272/2008; Version 4.2 (03/2021); <https://echa.europa.eu/documents>
- European Chemicals Agency (ECHA), Registered substances; <https://echa.europa.eu/information-on-chemicals/registered-substances>
- European Chemicals Agency (ECHA), C&L Classification and Labelling Inventory; <https://echa.europa.eu/information-on-chemicals/cl-inventory-database>
- Institute for Occupational Safety and Health of the German Social Accident Insurance (IFA): GESTIS database on hazardous substances and GESTIS - International limit values for chemical agents; <https://www.dguv.de/ifa/index.jsp>
- German Environmental Agency (Umweltbundesamt), Section IV 2.4: Office of Documentation and Information on Substances Hazardous to Waters RIGOLETTO (catalogue of Substances Hazardous to Waters); <https://webrigoletto.uba.de/rigoletto>

16.4. Methods according to Article 9 of Regulation (EC) No 1272/2008 for the evaluation of information for classification purposes

Calculation method according to the criteria in Annex I 1272/2008.
Flash point according to EN ISO 2719:2002.
pH value measurement.
Material compatibility and corrosiveness in practical tests

16.5. Training advice

Safety data sheet
according to Regulation (EC) No 1907/2006



Trade name: **AlproJet-W**
Issue/Revision: 27.04.2023

Version: 3.0
Replaces version: 2.1

Provide adequate information, instructions and training for users.

16.6. Indication of changes

A dash in the left hand margin indicates an amendment from the previous version.

The information given in the safety data sheet only applies to the described product in connection with its intended use. This information is based on the latest state of our knowledge at the time of revision. In particular, it describes our product under the aspect of its hazards and safety measures to be taken. It does not constitute any guarantee of product properties and quality features.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **SteriCleaner**

Issue/Revision: 09.04.2020

Print date: 23.10.2020

Version: 2.0

Replaces version: 1.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Trade name: **SteriCleaner**

1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses: Cleaning agent

Intended purpose: Solution ready for use for cleaning of sterilizing chambers, instrument trays and steel surfaces.

Uses advised against: None at intended use.

Note: The product is intended for professional users.

1.3. Details of the supplier of the safety data sheet

Manufacturer/Supplier: ALPRO MEDICAL GMBH
Mooswiesenstraße 9
D-78112 St. Georgen (Germany)
Telephone: +49 7725 9392-0
Telefax: +49 7725 9392-91
E-mail: alpro@alpro-medical.de
Internet: www.alpro-medical.com

E-mail address for the competent person responsible for the safety data sheet: doku@alpro-medical.de

1.4. Emergency telephone number

In-house emergency telephone number: +49 7725 9392-0
Monday – Friday from 08:00 am to 04:30 pm (UTC+1);
for chemical information and legal information on
hazardous substances only

Poison centre: +49 761 19240
Poisoning information centre, Freiburg, Germany
(24 h / 7 d), English is spoken

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification in accordance with Regulation (EC) No 1272/2008 [CLP]

<i>Classification</i>	<i>Classification procedure</i>
Flam. Liq. 3; H226	On basis of test data (Flash point)
Eye Irrit. 2; H319	Calculation method

Full text of hazard classes as well as H-phrases: see under SECTION 16.1.

2.2. Label elements

Label elements in accordance with Regulation (EC) No 1272/2008 [CLP]

Hazard pictograms:



Signal word: Warning

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **SteriCleaner**

Issue/Revision: 09.04.2020

Print date: 23.10.2020

Version: 2.0

Replaces version: 1.0

Hazard components
for labelling:

-

H-phrases:	H226	Flammable liquid and vapour.
	H319	Causes serious eye irritation.
P-phrases:	P210	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
	P233	Keep container tightly closed.
	P280	Wear eye protection/face protection.
	P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
	P337+P313	If eye irritation persists: Get medical advice/attention.

2.3. Other hazards

The substances in the mixture do not meet the PBT/vPvB criteria according to Regulation (EC) No 1907/2006, Annex XIII.

No further hazards known.

SECTION 3: Composition/information on ingredients

3.2. Mixtures

Chemical characterisation: Mixture of substances listed below with non-hazardous additions in aqueous solution.

Hazardous ingredients

Chemical name	Identification numbers	Classification in accordance with Regulation (EC) No 1272/2008	Weight %
Propan-2-ol	CAS No: 67-63-0 EC No: 200-661-7 Index No: 603-117-00-0 REACH Registration No: 01-2119457558-25-XXXX	Flam. Liq. 2; H225 Eye Irrit. 2; H319 STOT SE 3; H336	≥ 15 - < 20

Full text of hazard classes and H-phrases: see SECTION 16.1.

SECTION 4: First aid measures

4.1. Description of first aid measures

General information:	First aider: Pay attention to self-protection!
Following inhalation:	Move affected person into fresh air and keep still and warm. In case of continued complaints seek medical advice.
Following skin contact:	Wash skin immediately with plenty of water and soap. In case of skin reactions, consult a physician.
Following eye contact:	Flush eyes immediately with flowing water for 10 to 15 minutes holding eyelids apart. Remove contact lenses, if present and easy to do. Consult an ophthalmologist.
Following ingestion:	Rinse mouth with water. Let drink plenty of water. Do not induce vomiting. Consult a physician immediately.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **SteriCleaner**

Issue/Revision: 09.04.2020

Print date: 23.10.2020

Version: 2.0

Replaces version: 1.0

4.2. Most important symptoms and effects, both acute and delayed

Causes serious eye irritation.

4.3. Indication of any immediate medical attention and special treatment needed

No information available.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media: Water spray jet, alcohol resistant foam, extinguishing powder, carbon dioxide (CO₂)

Unsuitable extinguishing media: Full water jet

5.2. Special hazards arising from the substance or mixture

Hazardous combustion products: Carbon monoxide (CO), carbon dioxide (CO₂)

5.3. Advice for firefighters

Special protective equipment: Wear self-contained breathing apparatus.

Further information: Cool endangered containers with water spray jet.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

For non-emergency personnel

Use personal protective equipment. See SECTION 8.2.

Avoid skin and eye contact. Do not breathe vapours. Remove all sources of ignition. Provide adequate ventilation. Special danger of slipping by leaked/spilled product. Evacuate danger area. Observe emergency plans. Consult experts.

For emergency responders

Use personal protective equipment. See SECTION 8.2.

6.2. Environmental precautions

Do not discharge into drains or rivers.

6.3. Methods and material for containment and cleaning up

Containment

For large spills, dyke spilled material or otherwise contain material to ensure runoff does not reach a waterway. Cover or seal drains.

Cleaning up

Wipe up small amounts with absorbent material (e.g. cloth, fleece). Absorb large amounts with liquid-binding material (sand, diatomaceous earth, universal binder, sawdust). Collect in suitable, closed containers for disposal. Clean contaminated surfaces thoroughly.

Other information

Inappropriate containment and cleaning methods are not known.

6.4. Reference to other sections

Information on safe handling see SECTION 7.1.

Information on personal protective equipment see SECTION 8.2.

Safety data sheet
according to Regulation (EC) No 1907/2006



Trade name: **SteriCleaner**

Issue/Revision: 09.04.2020

Print date: 23.10.2020

Version: 2.0

Replaces version: 1.0

Information on disposal see SECTION 13.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions

Avoid contact with skin and eyes. Avoid breathing aerosols and vapours. Keep away from sources of ignition. Provide adequate ventilation. Keep container tightly closed. Fill refill packages only in labelled original bottles.

Advice on general occupational hygiene

When using do not eat, drink or smoke. Wash hands before breaks and at end of work. Keep away from food and drink.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels: Keep only in the original container. Keep container tightly closed and kept upright to prevent any leakage.

Advice on common storage: Not necessary

Further information on storage conditions: Not necessary

Storage class ([DE] TRGS 510): LGK 3 Flammable liquids

7.3. Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific end uses are stipulated.

Industry and sector specific guidance

[DE] TRGS 525 – Hazardous substances in medical care facilities (Section 7 Activities with disinfectants); Issue: September 2014;

Source: GMBI 2014 page 1294-1307 of 13.10.2014 [No 63]; www.baua.de

[DE] DGUV rules 107-002 (former BGR 206) - Disinfection works in health service

Issue: July 1999; Source: www.dguv.de/publikationen

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Occupational exposure limits

Country	Limit values				Legal basis	Remarks
	Long term (8 hours)		Short term (15 minutes)			
	ppm	mg/m ³	ppm	mg/m ³		
Propan-2-ol (CAS No: 67-63-0)						
EU						no limit value specified
UK	400	999	500	1250	EH40	

Used abbreviations, symbols, numerals and explanations in column „Remarks“

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **SteriCleaner**

Issue/Revision: 09.04.2020

Print date: 23.10.2020

Version: 2.0

Replaces version: 1.0

Biological limit values

Country	Parameter	Limit value	Test material	Sampling time	Legal basis
Propan-2-ol (CAS No: 67-63-0)					
Germany	Acetone	25 mg/l	Whole blood	End of exposition, resp. end of shift	TRGS 903
	Acetone	25 mg/l	Urine	End of exposition, resp. end of shift	TRGS 903

Information on monitoring procedures

BS EN 482:2012-04-30; Title: Workplace exposure - General requirements for the performance of procedures for the measurement of chemical agents;
British version of EN 482:2012

BS EN 689:1996-04-15; Title: Workplace atmospheres - Guidance for the assessment of exposure by inhalation to chemical agents for comparison with limit values and measurement strategy;
British version of EN 689:1995

BS EN 14042:2003-04-24; Title: Workplace atmospheres - Guide for the application and use of procedures for the assessment of exposure to chemical and biological agents;
British version of EN 14042:2003

8.2. Exposure controls

Appropriate engineering controls

Technical and organisational protective measures

The eyewash station (or eyewash bottle) must be located near the workplace.

Personal protective equipment

Eye/face protection:	Safety glasses with side protection according to EN 166
Skin protection:	
Hand protection:	Protective gloves according to EN 374 are recommended <u>Splash guard:</u> Disposable gloves made of nitrile rubber (thickness 0.11 mm) <u>Permanent contact (> 480 min):</u> Protective gloves made of nitrile rubber (thickness 0.40 mm)
Other skin protection:	Not necessary when used as intended.
Respiratory protection:	Not necessary when used as intended.
Thermal hazards:	No special protective measures necessary.

Environmental exposure controls

Do not discharge into drains.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Appearance:	clear, colourless liquid
Odour:	alcoholic
Odour threshold:	no data available
pH (undiluted):	9.5 – 10.5 (20 °C)
Melting point/freezing point:	no data available

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **SteriCleaner**

Issue/Revision: 09.04.2020

Print date: 23.10.2020

Version: 2.0
Replaces version: 1.0

Initial boiling point and boiling range:	no data available	
Flash point:	28 °C	(EN ISO 2719:2002)
Evaporation rate:	no data available	
Flammability (solid, gas):	not applicable	
Lower explosive limit:	Propan-2-ol: 2 vol%	
Upper explosive limit:	Propan-2-ol: 13 vol%	
Vapour pressure:	no data available	(... °C)
Vapour density:	no data available	
Relative density:	0.965 – 0.975	(20 °C)
Solubility in water:	completely soluble	
Partition coefficient: n-octanol/water	not applicable	
Auto-ignition temperature:	not applicable	
Decomposition temperature:	no data available	
Viscosity:	no data available	
Explosive properties:	The product is not explosive. However, formation of explosive vapour/air mixtures is possible.	
Oxidising properties:	none	

9.2. Other information

Refractive index nD:	1.3509-1.3540	(20 °C)
Electrical conductivity (undiluted):	1100-1300 µS/cm	(20 °C)

SECTION 10: Stability and reactivity

10.1. Reactivity

No hazardous reactions when handled and stored as intended.

10.2. Chemical stability

The product is stable when handled and stored as intended.

10.3. Possibility of hazardous reactions

Vapours may form explosive mixtures with air.

10.4. Conditions to avoid

Keep away from heat and sources of ignition.

10.5. Incompatible materials

None known

10.6. Hazardous decomposition products

Does not decompose when used as intended.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **SteriCleaner**

Issue/Revision: 09.04.2020

Print date: 23.10.2020

Version: 2.0

Replaces version: 1.0

Acute toxicity

Product

Acute toxicity - oral: Acute Toxicity Estimate $ATE_{mix} > 2000$ mg/kg
=> no classification

Acute toxicity - dermal: Acute Toxicity Estimate $ATE_{mix} > 2000$ mg/kg
=> no classification

Acute toxicity - inhalation: Acute Toxicity Estimate $ATE_{mix} > 20$ mg/l
=> no classification

Skin corrosion/irritation

Product

No classification. [calculation method]

Serious eye damage/irritation

Product

Causes serious eye irritation. [calculation method]

Respiratory or skin sensitisation

Product

No data available.

Germ cell mutagenicity

Product

No data available.

Carcinogenicity

Product

No classification. [calculation method]

Reproductive toxicity

Product

No data available.

STOT-single exposure

Product

No classification. [calculation method]

Ingredients

Propan-2-ol (CAS No: 67-63-0):
May cause drowsiness or dizziness.

STOT-repeated exposure

Product

No data available.

Aspiration hazard

Product

No data available.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **SteriCleaner**

Issue/Revision: 09.04.2020

Print date: 23.10.2020

Version: 2.0

Replaces version: 1.0

SECTION 12: Ecological information

12.1. Toxicity

No classification. [calculation method]

12.2. Persistence and degradability

Biodegradability:

The product is biodegradable according to OECD criteria. The statement has been derived from the properties of the ingredients.

12.3. Bioaccumulative potential

No data available.

12.4. Mobility in soil

No data available.

12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to Regulation (EC) No 1907/2006, Annex XIII.

12.6. Other adverse effects

No data available.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal of the product

Product residues must be disposed of as hazardous waste in compliance with the Directive 2008/98/EC on waste as well as national and regional regulations. Do not dispose of via the waste water. Leave product in the original container as possible. Do not mix with other waste materials.

Waste codes / waste designations according to EWC

Product residues: 07 06 04* other organic solvents, washing liquids and mother liquors

Disposal of the packaging

Packaging contaminated with product is considered as hazardous waste and must be disposed of accordingly.

Waste codes / waste designations according to EWC

Contaminated packaging: 15 01 10* packaging containing residues of or contaminated by hazardous substances

Recommendation

Contaminated packaging must be emptied optimally and can be recycled after appropriate cleaning (rinse with water).

SECTION 14: Transport information

14.0. Transport classification

Dangerous good in sense of the transport regulations in road traffic (ADR), railway traffic (RID), inland waterway traffic (ADN), maritime traffic (IMDG-Code) and air traffic (ICAO-TI/IATA-DGR).

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **SteriCleaner**

Issue/Revision: 09.04.2020

Print date: 23.10.2020

Version: 2.0
Replaces version: 1.0

14.1. UN number

UN 1987

14.2. UN proper shipping name

ADR/RID/ADN

ALCOHOLS, N.O.S. (Isopropanol)

IMDG-Code/ICAO-TI/IATA-DGR

ALCOHOLS, N.O.S. (Isopropanol)

14.3. Transport hazard class(es)

Class: 3

Subsidiary risk(s): -

14.4. Packing group

III

14.5. Environmental hazards

ADR/RID/ADN

Environmentally Hazardous: No

IMDG-Code

Marine Pollutant: No

14.6. Special precautions for user

Not necessary.

14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

14.8. Further information

Transport category according to ADR section 1.1.3.6: 3

Maximum total quantity per transport unit
according to ADR section 1.1.3.6: 1000 L

Limited quantity (Maximum quantity per inner
packaging) according to ADR/RID/ADN/IMDG-Code: 5 L

Classification code according to ADR/RID/ADN: F1

Hazard identification number according to
ADR/RID: 30

Tunnel restriction code according to ADR/RID: D/E

Segregation group according to IMDG-Code
section 5.4.1.5.11.1: -

EmS codes: F-E, S-D

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

EU-Regulations

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **SteriCleaner**

Issue/Revision: 09.04.2020

Print date: 23.10.2020

Version: 2.0

Replaces version: 1.0

REGULATION (EC) No 1005/2009 on substances that deplete the ozone layer

not applicable

REGULATION (EC) No 850/2004 on persistent organic pollutants and amending Directive 79/117/EEC

not applicable

REGULATION (EU) No 649/2012 concerning the export and import of hazardous chemicals

not applicable

DIRECTIVE 2012/18/EU (Seveso III Directive) on the control of major-accident hazards involving dangerous substances, amending and subsequently repealing Council Directive 96/82/EC

<i>Hazard category</i>	<i>Qualifying quantity (tonnes) (lower-tier establishment)</i>	<i>Qualifying quantity (tonnes) (upper-tier establishment)</i>
P5c FLAMMABLE LIQUIDS	5.000	50.000

DIRECTIVE 2010/75/EU on industrial emissions (integrated pollution prevention and control)

not applicable

REACH – List of substances subject to authorisation (Annex XIV)

not applicable

REACH – Restrictions on the manufacture, placing on the market and use of certain dangerous substances, mixtures and articles (Annex XVII)

not applicable

COUNCIL DIRECTIVE 94/33/EC on the protection of young people at work

not applicable

COUNCIL DIRECTIVE 92/85/EEC on the introduction of measures to encourage improvements in the safety and health at work of pregnant workers and workers who have recently given birth or are breastfeeding

not applicable

15.2. Chemical safety assessment

For this mixture no chemical safety assessment has been carried out.

SECTION 16: Other information

16.1. Full text of hazard classes and H-phrases

Hazard classes

Eye Irrit.	Eye irritation
Flam. Liq.	Flammable liquid
STOT SE	Specific target organ toxicity (single exposure)

H-phrases (Hazard statements)

H225	Highly flammable liquid and vapour.
H226	Flammable liquid and vapour.
H319	Causes serious eye irritation.
H336	May cause drowsiness or dizziness.

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **SteriCleaner**

Issue/Revision: 09.04.2020

Print date: 23.10.2020

Version: 2.0

Replaces version: 1.0

16.2. Abbreviations and acronyms

ADN	<u>A</u> ccord européen relatif au transport international des marchandises <u>d</u> angereuses par voie de <u>n</u> avigation intérieure (European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways)
ADR	<u>A</u> ccord européen relatif au transport international des marchandises <u>d</u> angereuses par <u>r</u> oute (European Agreement concerning the International Carriage of Dangerous Goods by Road)
BGR	<u>B</u> erufsgenossenschaftliche <u>R</u> egeln (English: Employers' liability insurance association rules)
BS	<u>B</u> ritish <u>S</u> tandards
CAS	<u>C</u> hemical <u>A</u> bstracts <u>S</u> ervice
CLP	Regulation on <u>C</u> lassification, <u>L</u> abelling and <u>P</u> ackaging of Substances and Mixtures
[DE]	National German regulations
DGUV	<u>D</u> eutsche <u>G</u> esetzliche <u>U</u> nfall <u>v</u> ersicherung (English: German statutory accident insurance)
EC	<u>E</u> uropean <u>C</u> ommunity
EEC	<u>E</u> uropean <u>E</u> conomic <u>C</u> ommunity
EmS	<u>E</u> mergency <u>S</u> chedules (Emergency response procedures for ships carrying dangerous goods)
EN	European Standard
EU	<u>E</u> uropean <u>U</u> nion
EWC	<u>E</u> uropean <u>W</u> aste <u>C</u> atalogue
GHS	<u>G</u> lobally <u>H</u> armonized <u>S</u> ystem of Classification, Labelling and Packaging of Chemicals
GMBI	<u>G</u> emeinsames <u>M</u> inisterial <u>bl</u> att (English: Joint Ministerial Gazette)
IATA-DGR	<u>I</u> nternational <u>A</u> ir <u>T</u> ransport <u>A</u> ssociation - <u>D</u> angerous <u>G</u> oods <u>R</u> egulations
IBC-Code	International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk
ICAO-TI	Technical Instructions For The Safe Transport of Dangerous Goods by Air
IMDG-Code	International <u>M</u> aritime Code for <u>D</u> angerous <u>G</u> oods
LGK	<u>L</u> ager <u>k</u> lasse (English: Storage class)
MARPOL	International Convention for the Prevention of Marine Pollution from Ships
N.O.S.	Not otherwise specified
OECD	<u>O</u> rganization for <u>E</u> conomic <u>C</u> o-operation and <u>D</u> evelopment
PBT	<u>P</u> ersistent, <u>b</u> ioaccumulative and <u>t</u> oxic
ppm	<u>P</u> arts <u>p</u> er <u>m</u> illion
REACH	<u>R</u> egistration, <u>E</u> valuation, <u>A</u> uthorisation and Restriction of <u>C</u> hemicals
RID	<u>R</u> èglement concernant le transport <u>I</u> nternational ferroviaire de marchandises <u>D</u> angereuses (Regulations Concerning the International Carriage of Dangerous Goods by Rail)
TRGS	<u>T</u> echnische <u>R</u> egeln für <u>G</u> efahr <u>s</u> toffe (English: Technical Rules for Hazardous Substances)
UN	<u>U</u> nited <u>N</u> ations
UTC	Coordinated Universal Time (French: Temps Universel Coordonné)
vPvB	<u>V</u> ery <u>p</u> ersistent and <u>v</u> ery <u>b</u> ioaccumulative
WGK	<u>W</u> assergefährdungs <u>k</u> lasse (English: water hazard class)

16.3. Key literature references and sources for data

- Regulation (EC) No 1907/2006 (REACH), Annex II
- European Chemicals Agency (ECHA) – Guidance on the compilation of safety data sheets; Version 2.1 (February 2014); http://echa.europa.eu/documents/10162/13643/sds_en.pdf
- GISBAU (Hazardous substances information system of the BG BAU) – course „safety data sheet“; <http://www.bgbau.de/gisbau/SDB/lehrgang/lehrgang.htm>
- Regulation (EC) No 1272/2008 (CLP regulation)
- European Chemicals Agency (ECHA) – Guidance on Labelling and Packaging in accordance with Regulation (EC) No 1272/2008 (04/2011); http://echa.europa.eu/documents/10162/13562/clp_labelling_en.pdf
- European Chemicals Agency (ECHA), Registered substances; <http://echa.europa.eu/information-on-chemicals/registered-substances>

Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **SteriCleaner**

Issue/Revision: 09.04.2020

Print date: 23.10.2020

Version: 2.0

Replaces version: 1.0

- European Chemicals Agency (ECHA), C&L Classification and Labelling Inventory; <http://echa.europa.eu/information-on-chemicals/cl-inventory-database>
- Institute for Occupational Safety and Health of the German Social Accident Insurance (IFA): GESTIS database on hazardous substances and GESTIS - International limit values for chemical agents; <http://www.dguv.de/dguv/ifa/index.jsp>
- German Environmental Agency (Umweltbundesamt), Section IV 2.4: Office of Documentation and Information on Substances Hazardous to Waters RIGOLETTO (catalogue of Substances Hazardous to Waters); <http://webrigoletto.uba.de/rigoletto>

16.4. Training advice

Provide adequate information, instructions and training for users.

16.5. Indication of changes

A dash in the left hand margin indicates an amendment from the previous version.

The information given in the safety data sheet only applies to the described product in connection with its intended use. This information is based on the latest state of our knowledge at the time of revision. In particular, it describes our product under the aspect of its hazards and safety measures to be taken. It does not constitute any guarantee of product properties and quality features.
