

Innovative Construction

First pressure-moulding material with two different degrees of hardness.

Superior wearing comfort

This custom sports mouth guard has the perfect fit. During sports, the slight impression worked into the mouth guard stabilises the teeth in the lower jaw.



Can be combined

Your braces treatment isn't affected.

Highly effective

Wearing a sports mouth guard significantly reduces the risk of injury to lips, teeth, jaws, and temporomandibular joint.

Impact absorbing

The impact is absorbed elastically and distributed with less force over a time delay.

Hygienic

CETRON® cleansing powder, toothbrush and toothpaste are the best choice for gentle care.

Precision-Fit Sports Mouth Guard

in just a few steps.

1 Consult with your dentist

Depending on the types of sports and recreational activities practiced, the dentist will advise you or your child on selecting the optimal sports mouth guard.

2 Impression

An impression of the upper jaw is the basis for sports mouth guards.

3 Producing the sports mouth guard

It is custom-made with precision in the laboratory. It covers the upper row of teeth and has a slight impression for the lower teeth.

4 Fitting and examination

The fit will be tested regularly by your dentist and adjusted as needed on an individual basis.

Provided by:



BIOPLAST® – the safe, ideal sports mouth guard for any application.



BIOPLAST®

What is it?

BIOPLAST® is a high-quality ethylene vinyl acetate (EVA) distinguished by its first-class impact absorption and outstanding elasticity despite its extremely light weight.

BIOPLAST® sports mouth guards for active sports are custom-made just for you by dental technicians, who work together with your dentist. This guarantees that your personal mouth guard will suit your particular needs for all types of sports.

For any challenge

Three sports mouth guard types



BIOPLAST® MASTER

- // proven, two-layer sports mouth guard
- // for all recreational sports
- // double-layer laminate construction for the best protection and comfort
- // 10 colour variations are available



BIOPLAST® XTREME

- // with integrated reinforcement section in the front teeth area
- // suited for contact sports, as it protects against hard and concentrated hits
- // double-layer laminate construction for the best protection and comfort



BIOPLAST® XTREME PRO

- // with integrated reinforcement area and supplemental hardelastic insert around the front teeth area
- // three-layer material structure for maximum protection and comfort
- // for professional contact sports with expected concentrated, hard hits (impact is distributed across a greater area)



A large number of sports accidents include injuries to the teeth and the periodontium.

