



UNIGLOVES® BLACK LATEX


SIZE	ITEM-NO
------	---------

S	2902
---	------

M	2903
---	------

L	2904
---	------

XL	2905
----	------



EIGENSCHAFTEN

MATERIAL	<u>Latex</u>
COLOUR	<u>black</u>
LÄNGE	<u>240 mm</u>
SIZE	<u>XS (5.5-6.0), S (6.5-7.0), M (7.5-8.0), L (8.5-9), XL (9.5-10)</u>
PRODUCT FEATURES	<u>beaded cuff, micro-roughened palm, non-sterile, powder-free</u>
AQL	<u>1.5</u>
NORM GLOVES	<u>CE, EN 455, T. 1-4</u>
VERPACKUNGSMENGE	<u>box of 100 pieces</u>

VERSANDEINHEIT

10 boxes