

**Trade name:** BEGO PMMA Splint E (REF 511069)

**Product no.:** REF511069

**Current version :** 2.0.0, issued: 20.10.2021

**Replaced version:** 1.0.1, issued: 26.04.2021

**Region:** GB

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1 Product identifier

**Trade name**

**BEGO PMMA Splint E (REF 511069)**

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses of the substance or mixture**

Dental Splint with thermo effect

**Uses advised against**

No data available.

### 1.3 Details of the supplier of the safety data sheet

**Address**

BEGO Medical GmbH  
Wilhelm-Herbst-Str. 1  
28359 Bremen

Telephone no. +49 421 2028 - 0  
Fax no. +49 421 2028 - 115  
e-mail msds@bego.com

**Advice on Safety Data Sheet**

msds@bego.com

### 1.4 Emergency telephone number

For medical advice (in German and English):  
+49 (0)551 192 40 (Giftinformationszentrum Nord)

## SECTION 2: Hazards identification

### 2.1 Classification of the substance or mixture

**Classification information**

This product does not meet the classification and labelling criteria given in the Regulation (EC) No 1272/2008 (CLP).

### 2.2 Label elements

Not relevant

### 2.3 Other hazards

PBT assessment  
No data available.

vPvB assessment  
No data available.

## SECTION 3: Composition/information on ingredients

### 3.1 Substances

Not applicable. The product is not a substance.

### 3.2 Mixtures

**Hazardous ingredients**

This product does not contain substances to be mentioned according to EU Regulation No 1907/2006 (REACH), annex II.

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

**General information**

with 1907/2006/EC

**Trade name:** BEGO PMMA Splint E (REF 511069)

**Product no.:** REF511069

**Current version :** 2.0.0, issued: 20.10.2021

**Replaced version:** 1.0.1, issued: 26.04.2021

**Region:** GB

In case of persisting adverse effects, consult a physician. Remove contaminated clothing and shoes and launder thoroughly before reusing.

**After inhalation**

Ensure supply of fresh air. Remove affected person from the immediate area.

**After skin contact**

When in contact with the skin, clean with soap and water.

**After eye contact**

Remove contact lenses. Rinse eye thoroughly under running water keeping eyelids wide open and protecting the unaffected eye (at least 10 to 15 minutes). In case of irritation consult an ophthalmologist.

**After ingestion**

Rinse the mouth thoroughly with water. Never give anything by mouth to an unconscious person. Do not induce vomiting.

**4.2 Most important symptoms and effects, both acute and delayed**

No data available.

**4.3 Indication of any immediate medical attention and special treatment needed**

No data available.

## SECTION 5: Firefighting measures

**5.1 Extinguishing media**

**Suitable extinguishing media**

Carbon dioxide; Extinguishing powder; Water spray jet; Foam

**Unsuitable extinguishing media**

High power water jet

**5.2 Special hazards arising from the substance or mixture**

No data available.

**5.3 Advice for firefighters**

Wear protective clothing. Use self-contained breathing apparatus.

## SECTION 6: Accidental release measures

**6.1 Personal precautions, protective equipment and emergency procedures**

**For non-emergency personnel**

Refer to protective measures listed in sections 7 and 8. Avoid dust formation.

**For emergency responders**

No data available. Personal protective equipment (PPE) - see Section 8.

**6.2 Environmental precautions**

Do not discharge into the drains/surface waters/groundwater.

**6.3 Methods and material for containment and cleaning up**

Collect mechanically. Avoid raising dust. When collected, handle material as described under the section heading "Disposal considerations".

**6.4 Reference to other sections**

No data available.

## SECTION 7: Handling and storage

**7.1 Precautions for safe handling**

**Advice on safe handling**

Avoid the formation and deposition of dust. Provide good ventilation at the work area (local exhaust ventilation, if necessary).

**General protective and hygiene measures**

**Trade name:** BEGO PMMA Splint E (REF 511069)

**Product no.:** REF511069

**Current version :** 2.0.0, issued: 20.10.2021

**Replaced version:** 1.0.1, issued: 26.04.2021

**Region:** GB

Wash hands before breaks and after work. Do not eat, drink or smoke during work time. Keep away from foodstuffs and beverages. Do not inhale dust. Remove soiled or soaked clothing immediately. Provide eye wash fountain in work area. Avoid contact with eyes and skin.

**Advice on protection against fire and explosion**

Dust can form an explosive mixture with air.

**7.2 Conditions for safe storage, including any incompatibilities**

**Technical measures and storage conditions**

Keep container tightly closed and dry in a cool, well-ventilated place. Protect from heat and direct sunlight.

**Requirements for storage rooms and vessels**

Containers which are opened must be carefully closed and kept upright to prevent leakage. Keep only in the original container.

**Incompatible products**

None known

**7.3 Specific end use(s)**

No data available.

**SECTION 8: Exposure controls/personal protection**

**8.1 Control parameters**

**Occupational exposure limit values**

No	Substance name	CAS no.	EC no.
1	Dust		
	<b>List of approved workplace exposure limits (WELs) / EH40</b>		
	Dust respirable		
	WEL long-term (8-hr TWA reference period)	4	mg/m <sup>3</sup>
	Comments	see Definition 44 "Dust"	
	<b>List of approved workplace exposure limits (WELs) / EH40</b>		
	Dust inhalable		
	WEL long-term (8-hr TWA reference period)	10	mg/m <sup>3</sup>
	Comments	see Definition 44 "Dust"	

**8.2 Exposure controls**

**Appropriate engineering controls**

Ensure adequate ventilation, local exhaust at the work station if necessary.

**Personal protective equipment**

**Respiratory protection**

If workplace exposure limits are exceeded, a respiration protection approved for this particular job must be worn. In case of dust formation, take appropriate measures for breathing protection in the event that workplace threshold values are not specified.

**Eye / face protection**

Safety glasses with side protection shield (EN 166)

**Hand protection**

In case of intensive contact, wear protective gloves (EN 374). Before use, the protective gloves should be tested in any case for its specific work-station suitability (i.e. mechanical resistance, product compatibility and antistatic properties). Adhere to the manufacturer's instructions and information relating to the use, storage, care and replacement of protective gloves. Protective gloves shall be replaced immediately when physically damaged or worn. Design operations thus to avoid permanent use of protective gloves.

**Other**

Chemical-resistant work clothes.

**Environmental exposure controls**

No data available.

with 1907/2006/EC

Trade name: BEGO PMMA Splint E (REF 511069)

Product no.: REF511069

Current version : 2.0.0, issued: 20.10.2021

Replaced version: 1.0.1, issued: 26.04.2021

Region: GB

## SECTION 9: Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

<b>State of aggregation</b>	
solid	
<b>Form/Colour</b>	
solid	
Various, depending on coloration	
<b>Odour</b>	
No data available	
<b>pH value</b>	
No data available	
<b>Boiling point / boiling range</b>	
No data available	
<b>Melting point/freezing point</b>	
No data available	
<b>Decomposition temperature</b>	
No data available	
<b>Flash point</b>	
Value	> 250 °C
<b>Ignition temperature</b>	
No data available	
<b>Flammability</b>	
No data available	
<b>Lower explosion limit</b>	
No data available	
<b>Upper explosion limit</b>	
No data available	
<b>Vapour pressure</b>	
No data available	
<b>Relative vapour density</b>	
No data available	
<b>Relative density</b>	
No data available	
<b>Density</b>	
No data available	
<b>Solubility</b>	
No data available	
<b>Partition coefficient n-octanol/water (log value)</b>	
No data available	
<b>Viscosity</b>	
No data available	
<b>Particle characteristics</b>	
No data available	

### 9.2 Other information

<b>Other information</b>
--------------------------

with 1907/2006/EC

**Trade name:** BEGO PMMA Splint E (REF 511069)**Product no.:** REF511069**Current version :** 2.0.0, issued: 20.10.2021**Replaced version:** 1.0.1, issued: 26.04.2021**Region:** GB

No data available.

**SECTION 10: Stability and reactivity****10.1 Reactivity**

No data available.

**10.2 Chemical stability**

Stable under recommended storage and handling conditions (See section 7).

**10.3 Possibility of hazardous reactions**

None, when used as directed.

**10.4 Conditions to avoid**

Heat, naked flames and other ignition sources. Temperatures &gt; 250°C

**10.5 Incompatible materials**

None known.

**10.6 Hazardous decomposition products**

In case of thermal decomposition, combustible, eyes and respiratory system irritant vapours will be formed, consisting mainly of methyl methacrylate

**SECTION 11: Toxicological information****11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008****Acute oral toxicity**

No data available

**Acute dermal toxicity**

No data available

**Acute inhalational toxicity**

No data available

**Skin corrosion/irritation**

No data available

**Serious eye damage/irritation**

No data available

**Respiratory or skin sensitisation**

No data available

**Germ cell mutagenicity**

No data available

**Reproduction toxicity**

No data available

**Carcinogenicity**

No data available

**STOT - single exposure**

No data available

**STOT - repeated exposure**

No data available

**Aspiration hazard**

No data available

**Delayed and immediate effects as well as chronic effects from short and long-term exposure**

Inhalation of dusts may irritate the respiratory tract. Contact with the skin and eyes may cause mechanical irritation.

with 1907/2006/EC

**Trade name:** BEGO PMMA Splint E (REF 511069)

**Product no.:** REF511069

**Current version :** 2.0.0, issued: 20.10.2021

**Replaced version:** 1.0.1, issued: 26.04.2021

**Region:** GB

**11.2 Information on other hazards**

**Endocrine disrupting properties**

No data available.

**Other information**

No data available.

**SECTION 12: Ecological information**

**12.1 Toxicity**

**Toxicity to fish (acute)**

No data available

**Toxicity to fish (chronic)**

No data available

**Toxicity to Daphnia (acute)**

No data available

**Toxicity to Daphnia (chronic)**

No data available

**Toxicity to algae (acute)**

No data available

**Toxicity to algae (chronic)**

No data available

**Bacteria toxicity**

No data available

**12.2 Persistence and degradability**

No data available.

**12.3 Bioaccumulative potential**

No data available.

**12.4 Mobility in soil**

No data available.

**12.5 Results of PBT and vPvB assessment**

**Results of PBT and vPvB assessment**

PBT assessment

No data available.

vPvB assessment

No data available.

**12.6 Endocrine disrupting properties**

No data available.

**12.7 Other adverse effects**

No data available.

**12.8 Other information**

**Other information**

Do not discharge product unmonitored into the environment.

**SECTION 13: Disposal considerations**

**13.1 Waste treatment methods**

**Product**

Allocation of a waste code number, according to the European Waste Catalogue, should be carried out in agreement with the regional waste disposal company.

**Packaging**

with 1907/2006/EC

**Trade name:** BEGO PMMA Splint E (REF 511069)

**Product no.:** REF511069

**Current version :** 2.0.0, issued: 20.10.2021

**Replaced version:** 1.0.1, issued: 26.04.2021

**Region:** GB

Residues must be removed from packaging and when emptied completely disposed of in accordance with the regulations for waste removal. Incompletely emptied packaging must be disposed of in the form of disposal specified by the regional disposer.

## SECTION 14: Transport information

### 14.1 Transport ADR/RID/ADN

The product is not subject to ADR/RID/ADN regulations.

### 14.2 Transport IMDG

The product is not subject to IMDG regulations.

### 14.3 Transport ICAO-TI / IATA

The product is not subject to ICAO-TI / IATA regulations.

### 14.4 Other information

No data available.

### 14.5 Environmental hazards

Information on environmental hazards, if relevant, please see 14.1 - 14.3.

### 14.6 Special precautions for user

No data available.

### 14.7 Maritime transport in bulk according to IMO instruments

Not relevant

## SECTION 15: Regulatory information

### 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

#### EU regulations

##### **Regulation (EC) No 1907/2006 (REACH) Annex XIV (List of substances subject to authorisation)**

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substances considered as substances requiring authorisation as listed on Annex XIV of the REACH regulation (EC) 1907/2006.

##### **REACH candidate list of substances of very high concern (SVHC) for authorisation**

According to available data and the information provided by preliminary suppliers, the product does not contain substances that are considered substances meeting the criteria for inclusion in annex XIV (List of Substances Subject to Authorisation) as laid down in Article 57 and article 59 of REACH (EC) 1907/2006.

##### **Regulation (EC) No 1907/2006 (REACH) Annex XVII: RESTRICTIONS ON THE MANUFACTURE, PLACING ON THE MARKET AND USE OF CERTAIN DANGEROUS SUBSTANCES, MIXTURES AND ARTICLES**

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substances subject to restriction as listed in Annex XVII of the REACH regulation (EC) 1907/2006.

##### **Directive 2012/18/EU on the control of major-accident hazards involving dangerous substances**

This product is not subject to Part 1 or 2 of Annex I.

##### **Other regulations**

Adhere to the national sanitary and occupational safety regulations when using this product.

### 15.2 Chemical safety assessment

A chemical safety assessment has not been carried out for this mixture.

## SECTION 16: Other information

#### **Sources of key data used to compile the data sheet:**

Regulation (EC) No 1907/2006 (REACH), 1272/2008 (CLP) as amended in each case.

Directives 2000/39/EC, 2006/15/EC, 2009/161/EU, (EU) 2017/164.

National Threshold Limit Values of the corresponding countries as amended in each case.

**Trade name:** BEGO PMMA Splint E (REF 511069)

**Product no.:** REF511069

**Current version :** 2.0.0, issued: 20.10.2021

**Replaced version:** 1.0.1, issued: 26.04.2021

**Region:** GB

---

Transport regulations according to ADR, RID, IMDG, IATA as amended in each case.

The data sources used to determine physical, toxic and ecotoxic data, are indicated directly in the corresponding section.

**Creation of the safety data sheet**

UMCO GmbH - D-21107 Hamburg, Georg-Wilhelm-Strasse 187, Tel.: +49(40)555 546 300, Fax: +49(40)555 546 357, e-mail: [umco@umco.de](mailto:umco@umco.de)

This information is based on our present knowledge and experience.

The safety data sheet describes products with a view to safety requirements.

It does not however, constitute a guarantee for any specific product properties and shall not establish a legally valid contractual relationship.

Alterations/supplements:

Alterations to the previous edition are marked in the left-hand margin.

Document protected by copyright. Alterations or reproductions require the express written permission of UMCO GmbH.

Prod-ID 701680