

Safety Data Sheet



according Regulation (EC) No. 2015/830
Date of issue: 02.12.2022
Revision date: - / Version: 1.0

Seite 1 von 6
Druckdatum: 02.09.2025

Pumice paste

SECTION 1: Identification of the substance/mixture and of the company/undertaking:

1.1 Product identifier

Product name: Pumice paste

1.2 Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses: Industrial use, Professional use

1.3 Details of the supplier of the safety data sheet

Manufacturer/Supplier: ERNST HINRICHS Dental GmbH
Street / mailbox: Borsigstr. 1
Country code. / postal code / city: D - 38644 Goslar
Phone: 0 53 21 / 5 06 24
Fax: 0 53 21 / 5 08 81
E-mail / Website: info@hinrichs-dental.de / www.hinrichs-dental.de
Further information obtainable from: ERNST HINRICHS Dental GmbH

1.4 Emergency telephone number

ERNST HINRICHS Dental GmbH: +49 (0) 53 21 / 5 06 24 (Mon-Fri. 8 a.m. – 4 p.m.)

SECTION 2: Hazards identification:

2.1 Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]: Not classified.

2.2 Label elements:

Labelling according to Regulation (EC) No. 1272/2008 [CLP]:
EUH210 Safety data sheet available on request

2.3 Other hazards

Results of PBT and vPvB assessment

PBT: Not applicable.
vPvB: Not applicable.

SECTION 3: Composition/information on ingredients:

3.1 Substance Not applicable

3.2 Mixture

Dangerous components:

(CAS no) 1332-09-8	Pumice	70 - 90%
(EC no) 603-719-3	Substance with an occupational exposure limit	

Full text of H-statements: see section 16

SECTION 4: First aid measures:

4.1 Description of first aid measures:

First-aid measures general:

Never give anything by mouth to an unconscious person. If unconscious, place in the recovery position and seek medical advice.

First-aid measures after inhalation:

Allow breathing of fresh air. If breathing stops, give artificial respiration. Get medical advice/attention if you feel unwell.

Safety Data Sheet



according Regulation (EC) No. 2015/830
Date of issue: 02.12.2022
Revision date: - / Version: 1.0

Seite 2 von 6
Druckdatum: 02.09.2025

Pumice paste

First-aid measures after skin contact:	Take off contaminated clothing and wash it before reuse. Wash with plenty of soap and water. If skin irritation occurs: Get medical advice/attention.
First-aid measures after eye contact:	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
First-aid measures after ingestion:	Rinse mouth. Get medical advice/attention if you feel unwell.
4.2 Most important symptoms and effects, both acute and delayed:	Not expected to present a significant hazard under anticipated conditions of normal use.
4.3 Indication of any immediate medical attention and special treatment needed:	Treat symptomatically.

SECTION 5: Firefighting measures:

5.1 Extinguishing media: Suitable extinguishing media:	Making extinguishing agents environment-friendly. Foam. Fire-extinguishing powder. Carbon dioxide. Water spray. Sand.
Unsuitable extinguishing media:	Do not use a heavy water stream.
5.2 Special hazards arising from the substance or mixture:	The product itself does not burn.
5.3 Advice for firefighters: Firefighting instructions	Use water spray or fog for cooling exposed containers. Prevent fire-fighting water from entering environment.
Protection during firefighting:	Do not enter fire area without proper protective equipment, including respiratory protection.

SECTION 6: Accidental release measures:

6.1 Personal precautions, protective equipment and emergency procedures:	Provide adequate ventilation. Evacuate unnecessary personnel. Use personal protective equipment as required. Avoid formation of aerosols.
6.2 Environmental precautions:	Prevent entry to sewers and public waters. Notify authorities if product enters sewers or public waters.
6.3 Methods and material for containment and cleaning up:	Send for recovery or disposal in suitable containers. Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders). Collect in suitable container for disposal.
6.4 Reference to other sections:	Concerning personal protective equipment to use, see section 8. Concerning disposal elimination after cleaning, see section 13.

SECTION 7: Handling and storage:

7.1 Precautions for safe handling: Additional hazards when processed:	Ensure good ventilation of the work station. Avoid formation of aerosols.
--	---

Safety Data Sheet



according Regulation (EC) No. 2015/830
Date of issue: 02.12.2022
Revision date: - / Version: 1.0

Seite 3 von 6
Druckdatum: 02.09.2025

Pumice paste

Hygiene measures:

Handle in accordance with good industrial hygiene and safety procedures. When using do not eat, drink or smoke. Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Contaminated work clothing should not be allowed out of the workplace.

7.2 Conditions for safe storage, including any incompatibilities:

Storage conditions:

Store in original container. Store tightly closed in a dry and cool place.

Storage class:

12 - Non-flammable liquid

7.3 Specific end use(s):

No additional information available.

SECTION 8: Exposure controls/personal protection

8.1 Control parameters:

No additional information available.

DNEL:

No further relevant information available.

PNEC:

No further relevant information available.

8.2 Exposure controls

Appropriate engineering controls:

Use adequate ventilation to keep vapour concentrations below applicable standard.

Personal protective equipment:

Hand protection:

Wear suitable gloves. Nitrile rubber. Latex. 0.4 mm. The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

Eye protection:

Chemical goggles or safety glasses.

Respiratory protection:

Where exposure through inhalation may occur from use, respiratory protection equipment is recommended. Breathing apparatus with filter P.

Skin and body protection:

Wear suitable protective clothing.

SECTION 9: Physical and chemical properties:

9.1 Information on basic physical and chemical properties:

Appearance:

Colorless to whitish liquid

Physical state:

Liquid.

Colour:

Colorless to whitish.

Odour:

Almost odourless.

Odour threshold:

No data available.

pH:

< 8

Melting point/Freezing point:

No data available.

Boiling point

No data available.

Flash point:

No data available.

Flammability (solid, gas):

Not applicable.

Auto-ignition temperature:

No data available.

Decomposition temperature:

No data available.

Safety Data Sheet



according Regulation (EC) No. 2015/830
Date of issue: 02.12.2022
Revision date: - / Version: 1.0

Seite 4 von 6
Druckdatum: 02.09.2025

Pumice paste

Oxidising properties:	No.
Explosive properties:	No data available.
Explosive limits:	No data available.
Vapour pressure:	No data available.
Density:	No data available.
Relative density:	No data available.
Relative vapour density at 20 °C:	No data available.
Relative evaporation rate (butylacetate=1):	No data available.
Solubility:	No data available.
Log Pow:	No data available.
Viscosity, dynamic:	No data available.
Viscosity, kinematic:	No data available.

9.2 Other information: No additional information available.

SECTION 10: Stability and reactivity:

10.1 Reactivity:	No dangerous reactions known under normal conditions of use.
10.2 Chemical stability:	Stable under use and storage conditions as recommended in section 7. No decomposition if stored and applied as directed.
10.3 Possibility of hazardous reactions:	None under normal use.
10.4 Conditions to avoid:	No information available.
10.5 Incompatible materials:	None under normal use.
10.6 Hazardous decomposition products:	No information available.

SECTION 11: Toxicological information:

11.1 Information on toxicological effects	
Acute toxicity:	Not classified
Skin corrosion/irritation:	Not classified Based on available data, the classification criteria are not met
Serious eye damage/irritation:	Not classified Based on available data, the classification criteria are not met
Respiratory or skin sensitisation	Not classified Based on available data, the classification criteria are not met
Germ cell mutagenicity:	Not classified Based on available data, the classification criteria are not met

Safety Data Sheet



according Regulation (EC) No. 2015/830
Date of issue: 02.12.2022
Revision date: - / Version: 1.0

Seite 5 von 6
Druckdatum: 02.09.2025

Pumice paste

Carcinogenicity:	Not classified Based on available data, the classification criteria are not met
Reproductive toxicity:	Not classified Based on available data, the classification criteria are not met
Specific target organ toxicity (single exposure):	Not classified Based on available data, the classification criteria are not met
Specific target organ toxicity (repeated exposure):	Not classified Based on available data, the classification criteria are not met
Aspiration hazard:	Not classified Based on available data, the classification criteria are not met
Potential adverse human health effects and symptoms:	Based on available data, the classification criteria are not met

SECTION 12: Ecological information:

12.1 Toxicity:	
Aquatic toxicity:	No further relevant information available.
12.2 Persistence and degradability:	No information available.
12.3 Bioaccumulative potential:	No information available.
12.4 Mobility in soil:	No information available.
12.5 Results of PBT and vPvB assessment:	
PBT:	Not applicable.
vPvB:	Not applicable.
12.6 Other adverse effects:	No other effects known.

SECTION 13: Disposal considerations:

13.1 Waste treatment methods	
Regional legislation (waste):	Dispose in a safe manner in accordance with local/national regulations.
Waste disposal recommendations:	When totally empty, containers are recyclable like any other packing.
Waste code:	The waste code number according to the Ordinance on the European Waste Catalogue (AVV) depends on the waste producer and can therefore vary for any given product. The waste code number is therefore to be gleaned separately from each waste producer.

Safety Data Sheet



according Regulation (EC) No. 2015/830
Date of issue: 02.12.2022
Revision date: - / Version: 1.0

Seite 6 von 6
Druckdatum: 02.09.2025

Pumice paste

SECTION 14: Transport information:

14.1 UN number:	
ADR, IMDG, IATA:	Not applicable
14.2 UN proper shipping name	
ADR, IMDG, IATA:	Not applicable
14.3 Transport hazard class(es)	
ADR, IMDG, IATA:	Not applicable
14.4 Packing group:	
ADR, IMDG, IATA:	Not applicable
14.5 Environmental hazards:	No
Marine pollutant:	No
Other information:	No supplementary information available
14.6 Special precautions for user:	Not applicable.
14.7 Transport in bulk according to Annex II of MARPOL and the IBC Code:	Not applicable.

SECTION 15: Regulatory information:

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture	
EU-Regulations	This mixture is classified as not hazardous according to regulation (EC) 1272/2008.
National regulations	No additional information available.
15.2 Chemical safety assessment:	A chemical safety assessment has not been carried out.

SECTION 16: Other information:

These data are based on our present knowledge. However, they shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.
The data given here only applies when product used for proper application(s). The product is not sold as suitable for other applications - usage in such may cause risks not mentioned in this sheet. Do not use for other application(s) without seeking advice from manufacturer.

Abbreviations and acronyms:

ADR:	Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
IMDG:	International Maritime Code for Dangerous Goods
IATA:	International Air Transport Association
GHS:	Globally Harmonised System of Classification and Labelling of Chemicals
EINECS:	European Inventory of Existing Commercial Chemical Substances
ELINCS:	European List of Notified Chemical Substances
CAS:	Chemical Abstracts Service (division of the American Chemical Society)
DNEL:	Derived No-Effect Level (REACH)
PNEC:	Predicted No-Effect Concentration (REACH)
LC50:	Lethal concentration, 50 percent
LD50:	Lethal dose, 50 percent