

according to Regulation (EC) No 1907/2006

UC31 BioSonic General Purpose Cleaner and Super Rust Inhibitor

Revision date: 13.05.2016

Print date: 21.03.2017

Page 1 of 5

SECTION 1: Identification of the substance/mixture and of the company/undertaking**1.1. Product identifier**

UC31 BioSonic General Purpose Cleaner and Super Rust Inhibitor

Further trade namesCat. No UC31
48062**1.2. Relevant identified uses of the substance or mixture and uses advised against****1.3. Details of the supplier of the safety data sheet****Manufacturer**

Company name: Coltène/Whaledent Inc.
Street: 235 Ascot Parkway
Place: Cuyahoga Falls, OH 44223/USA
Telephone: +1 (330) 916 8800
Telefax: +1 (330) 916 7077
e-mail: info.us@coltene.com
Internet: www.coltene.com

Authorized Representative

Company name: Coltène/Whaledent GmbH & Co. KG
Street: Raiffeisenstraße 30
Place: D-89129 Langenau
Telephone: +49 (7345) 805 0
Telefax: +49 (7345) 805 201
e-mail: info.de@coltene.com
Internet: www.coltene.com

SECTION 2: Hazards identification**2.1. Classification of the substance or mixture****Regulation (EC) No. 1272/2008**

Hazard categories:
Serious eye damage/eye irritation: Eye Dam. 1
Reproductive toxicity: Repr. 1B
Hazard Statements:
Causes serious eye damage.
May damage fertility. May damage the unborn child.

2.2. Label elements**Regulation (EC) No. 1272/2008****Hazard components for labelling**

disodium tetraborate pentahydrate; borax pentahydrate
Alkyl Polyglycosides

Signal word: Danger**Pictograms:****Hazard statements**

H318 Causes serious eye damage.
H360FD May damage fertility. May damage the unborn child.

according to Regulation (EC) No 1907/2006

UC31 BioSonic General Purpose Cleaner and Super Rust Inhibitor

Revision date: 13.05.2016

Print date: 21.03.2017

Page 2 of 5

Precautionary statements

P201	Obtain special instructions before use.
P202	Do not handle until all safety precautions have been read and understood.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P264	Wash hands thoroughly after handling.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a doctor.
P308+P313	IF exposed or concerned: Get medical advice/attention.
P405	Store locked up.

SECTION 3: Composition/information on ingredients**3.2. Mixtures****Hazardous components**

CAS No	Chemical name			Quantity
	EC No	Index No	REACH No	
	Classification according to Regulation (EC) No. 1272/2008 [CLP]			
12179-04-3	disodium tetraborate pentahydrate; borax pentahydrate			5 - < 10 %
	215-540-4	005-011-02-9		
	Repr. 1B; H360FD			
68515-73-1	Alkyl Polyglycosides			5 - < 10 %
	500-220-1			
	Eye Dam. 1, Aquatic Acute 3; H318 H402			

Full text of H and EUH statements: see section 16.

SECTION 4: First aid measures**4.1. Description of first aid measures****After inhalation**

Provide fresh air.
In all cases of doubt, or when symptoms persist, seek medical advice.

After contact with skin

After contact with skin, wash immediately with: Water and soap.
Change contaminated clothing.
In case of skin irritation, seek medical treatment.

After contact with eyes

After contact with the eyes, rinse with water with the eyelids open for a sufficient length of time, then consult an ophthalmologist immediately.

After ingestion

Do NOT induce vomiting. Let water be drunken in little sips (dilution effect). Seek medical advice.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**

Water spray. Foam. Extinguishing powder. Carbon dioxide (CO₂). Sand.

Additional information

Co-ordinate fire-fighting measures to the fire surroundings.

SECTION 6: Accidental release measures

according to Regulation (EC) No 1907/2006

UC31 BioSonic General Purpose Cleaner and Super Rust Inhibitor

Revision date: 13.05.2016

Print date: 21.03.2017

Page 3 of 5

6.1. Personal precautions, protective equipment and emergency procedures

See protective measures under point 7 and 8.

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up

Wipe up with absorbent material (eg. cloth, fleece). Clear contaminated areas thoroughly.

SECTION 7: Handling and storage**7.1. Precautions for safe handling****Advice on safe handling**

Avoid contact with skin and eyes.

for professional use only Keep out of the reach of children.

Advice on protection against fire and explosion

Keep away from sources of ignition - No smoking.

Further information on handling

Always close containers tightly after the removal of product.

7.2. Conditions for safe storage, including any incompatibilities**Requirements for storage rooms and vessels**

Keep/Store only in original container.

Further information on storage conditions

Recommended storage temperature: >0 °C

SECTION 8: Exposure controls/personal protection**8.1. Control parameters****8.2. Exposure controls****Protective and hygiene measures**

When using do not eat, drink or smoke.

Eye/face protection

Framed glasses.

Hand protection

Wear suitable gloves.

Skin protection

Lab apron.

Respiratory protection

Provide adequate ventilation.

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

Physical state:	liquid
Colour:	green
Odour:	characteristic

pH-Value:	7 - 9
-----------	-------

Changes in the physical state

Initial boiling point and boiling range:	100 °C
--	--------

Flash point:	> 100 °C
--------------	----------

Test method

according to Regulation (EC) No 1907/2006

UC31 BioSonic General Purpose Cleaner and Super Rust Inhibitor

Revision date: 13.05.2016

Print date: 21.03.2017

Page 4 of 5

Density: 1.08 g/cm³
Water solubility: easily soluble.
Viscosity / dynamic: 60 mPa·s Brookfield

SECTION 10: Stability and reactivity**10.2. Chemical stability**

The substance is chemically stable under recommended conditions of storage, use and temperature.

10.4. Conditions to avoid

heat.

10.5. Incompatible materials

Oxidizing agents.

SECTION 11: Toxicological information**11.1. Information on toxicological effects****Acute toxicity**

Based on available data, the classification criteria are not met.

Irritation and corrosivity

Causes serious eye damage.

Skin corrosion/irritation: Based on available data, the classification criteria are not met.

Sensitising effects

Based on available data, the classification criteria are not met.

Carcinogenic/mutagenic/toxic effects for reproduction

May damage fertility. May damage the unborn child.

Germ cell mutagenicity: Based on available data, the classification criteria are not met.

Carcinogenicity: Based on available data, the classification criteria are not met.

STOT-single exposure

Based on available data, the classification criteria are not met.

STOT-repeated exposure

Based on available data, the classification criteria are not met.

Aspiration hazard

Based on available data, the classification criteria are not met.

Additional information on tests

SECTION 12: Ecological information**12.2. Persistence and degradability**

Readily biodegradable (according to OECD criteria).

12.3. Bioaccumulative potential

No indication of bioaccumulation potential.

Further information

According to the criteria of the EC-classification and labelling "dangerous for the environment" (93/21/EEC) the material/product is not to be classified as dangerous to the environment.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Advice on disposal**

This material and its container must be disposed of in a safe way.

according to Regulation (EC) No 1907/2006

UC31 BioSonic General Purpose Cleaner and Super Rust Inhibitor

Revision date: 13.05.2016

Print date: 21.03.2017

Page 5 of 5

Dispose of waste according to applicable legislation.

SECTION 14: Transport information**Land transport (ADR/RID)****14.1. UN number:** ---**Marine transport (IMDG)****14.1. UN number:** ---**Air transport (ICAO-TI/IATA-DGR)****14.1. UN number:** ---**Other applicable information**

No dangerous good in sense of these transport regulations.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****EU regulatory information****Additional information**

Safety Data Sheet according to Regulation (EC) No. 1907/2006 (REACH)

National regulatory information**Additional information**

Canada WHMIS: This product has been classified according to the hazard criteria of the CPR and the MSDS contains all the information required by the CPR.

SECTION 16: Other information**Relevant H and EUH statements (number and full text)**

H318	Causes serious eye damage.
H360FD	May damage fertility. May damage the unborn child.
H402	Harmful to aquatic life.

Further Information

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)