

Name : CAVITY LINER LIGHT CURE - ART.11035



## SAFETY DATA SHEET

(REACH regulation (EC) n° 1907/2006 - n° 453/2010)

### SECTION 1 : IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

#### 1.1. Product identifier

Product name : CAVITY LINER LIGHT CURE

Product code : ART.11035.

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

For dental professional use only.

#### 1.3. Details of the supplier of the safety data sheet

Registered company name : PRODUITS DENTAIRES SA.

Address : Rue des Bosquets 18 .1800 .Vevey .Switzerland.

Telephone : +41 21 921 26 31. Fax : +41 21 921 39 79.

info@pdsa.ch

http://www.pdsa.ch

#### 1.4. Emergency telephone number : +41 21 921 26 31.

Association/Organisation : .

### SECTION 2 : HAZARDS IDENTIFICATION

#### 2.1. Classification of the substance or mixture

##### In compliance with EC regulation No. 1272/2008 and its amendments.

Skin sensitisation, Category 1 (Skin Sens. 1, H317).

This mixture does not present a physical hazard. Refer to the recommendations regarding the other products present on the site.

This mixture does not present an environmental hazard. No known or foreseeable environmental damage under standard conditions of use.

##### In compliance with directives 67/548/EEC, 1999/45/EC and their amendments.

Skin sensitisation (Xi, R 43).

Aquatic environmental hazard, chronic toxicity: harmful (R 52/53).

This mixture does not present a physical hazard. Refer to the recommendations regarding the other products present on the site.

#### 2.2. Label elements

##### In compliance with EC regulation No. 1272/2008 and its amendments.

Hazard pictograms :



GHS07

Signal Word :

WARNING

Product identifiers :

607-132-00-3

DIMETHYLAMINOETHYL METHACRYLATE 2-

Hazard statements :

H317

May cause an allergic skin reaction.

Precautionary statements - Prevention :

P280

Wear protective gloves/protective clothing/eye protection/face protection.

**Name : CAVITY LINER LIGHT CURE - ART.11035**

**Precautionary statements - Response :**

P302 + P352 IF ON SKIN: Wash with plenty of soap and water.  
P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

**Precautionary statements - Disposal :**

P501 Dispose of contents/container to ...

**2.3. Other hazards**

The mixture does not contain any substances classified as 'Substances of Very High Concern' (SVHC) by the European Chemical Agency (ECHA) under article 57 of REACH: <http://echa.europa.eu/fr/candidate-list-table>

The mixture satisfies neither the PBT nor the vPvB criteria for mixtures in accordance with annexe XIII of the REACH regulations EC 1907/2006.

**SECTION 3 : COMPOSITION/INFORMATION ON INGREDIENTS**

**3.2. Mixtures**

**Composition :**

| Identification                                                                                      | (EC) 1272/2008                                                                                                              | 67/548/EEC                       | Note | %               |
|-----------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|----------------------------------|------|-----------------|
| INDEX: 607-132-00-3<br>CAS: 2867-47-2<br>EC: 220-688-8<br><br>DIMETHYLAMINOETHYL<br>METHACRYLATE 2- | GHS07<br>Wng<br>Acute Tox. 4, H312<br>Acute Tox. 4, H302<br>Eye Irrit. 2, H319<br>Skin Irrit. 2, H315<br>Skin Sens. 1, H317 | Xn<br>Xn;R21/22<br>Xi;R36/38-R43 |      | 2.5 <= x % < 10 |

**SECTION 4 : FIRST AID MEASURES**

As a general rule, in case of doubt or if symptoms persist, always call a doctor.

NEVER induce swallowing by an unconscious person.

**4.1. Description of first aid measures**

**In the event of exposure by inhalation :**

In the event of massive inhalation, remove the person exposed to fresh air. Keep warm and at rest.

**In the event of splashes or contact with eyes :**

Wash thoroughly with soft, clean water for 15 minutes holding the eyelids open.

If necessary refer the patient to an ophtalmologist and show the composition.

**In the event of splashes or contact with skin :**

In the event of an allergic reaction, seek medical attention.

If the contaminated area is widespread and/or there is damage to the skin, a doctor must be consulted or the patient transferred to hospital.

Wash the skin thoroughly with soap and fresh water.

Be careful with product which can remain between the skin and the clothes, the watch, the jewels.

**In the event of swallowing :**

Do not give the patient anything orally.

In the event of swallowing, if the quantity is small (no more than one mouthful), rinse the mouth with water and consult a doctor.

**4.2. Most important symptoms and effects, both acute and delayed**

No data available.

**4.3. Indication of any immediate medical attention and special treatment needed**

No data available.

Name : CAVITY LINER LIGHT CURE - ART.11035

---

---

## SECTION 5 : FIREFIGHTING MEASURES

Non-flammable.

### 5.1. Extinguishing media

#### Suitable methods of extinction

In the event of a fire, use :

- carbon dioxide (CO<sub>2</sub>)
- powder
- sprayed water or water mist
- stable foams in alcohols

#### Unsuitable methods of extinction

In the event of a fire, do not use :

- water jet

### 5.2. Special hazards arising from the substance or mixture

A fire will often produce a thick black smoke. Exposure to decomposition products may be hazardous to health.

Do not breathe in smoke.

In the event of a fire, the following may be formed :

- carbon monoxide (CO)
- carbon dioxide (CO<sub>2</sub>)

### 5.3. Advice for firefighters

Wear self-contained breathing apparatus.

---

## SECTION 6 : ACCIDENTAL RELEASE MEASURES

### 6.1. Personal precautions, protective equipment and emergency procedures

Consult the safety measures listed under headings 7 and 8.

#### For non first aid worker

Avoid any contact with the skin and eyes.

#### For first aid worker

First aid workers will be equipped with suitable personal protective equipment (See section 8).

Avoid any contact with the skin and eyes.

### 6.2. Environmental precautions

Prevent any material from entering drains or waterways.

### 6.3. Methods and material for containment and cleaning up

No data available.

### 6.4. Reference to other sections

No data available.

---

## SECTION 7 : HANDLING AND STORAGE

Individuals with a history of skin sensitisation should not, under any circumstance, handle this mixture.

### 7.1. Precautions for safe handling

Always wash hands after handling.

Avoid any contact with skin and eyes.

#### Fire prevention :

Prevent access by unauthorised personnel.

**Recommended equipment and procedures :**

For personal protection, see section 8.

Packages which have been opened must be reclosed carefully.

**7.2. Conditions for safe storage, including any incompatibilities**

Keep in the original container, respect the storage temperatures mentioned on the container, the leaflet and the box.

**Storage**

Respect storage precautions mentioned in the instructions for use.

**Packaging**

Always conserve in the original packaging.

**7.3. Specific end use(s)**

Consult instructions for use enclosed.

---

**SECTION 8 : EXPOSURE CONTROLS/PERSONAL PROTECTION**

**8.1. Control parameters**

No data available.

**8.2. Exposure controls**

**Personal protection measures, such as personal protective equipment**

Pictogram(s) indicating the obligation of wearing personal protective equipment (PPE) :



**- Eye / face protection**

Avoid contact with eyes.

**- Hand protection**

Wear gloves according to the hygienic rules in force in dental surgeries.

**- Body protection**

Avoid skin contact.

After contact with the product, all parts of the body that have been soiled must be washed.

---

**SECTION 9 : PHYSICAL AND CHEMICAL PROPERTIES**

**9.1. Information on basic physical and chemical properties**

**General information :**

|                  |                |
|------------------|----------------|
| Physical state : | Paste.         |
| Colour:          | white          |
| Odour:           | characteristic |

**Important health, safety and environmental information**

|                                           |                                |
|-------------------------------------------|--------------------------------|
| pH :                                      | Not stated.<br>Slightly basic. |
| Boiling point/boiling range :             | Not relevant.                  |
| Flash point interval :                    | Not relevant.                  |
| Vapour pressure (50°C) :                  | Not relevant.                  |
| Density :                                 | < 1                            |
| Water solubility :                        | Insoluble.                     |
| Melting point/melting range :             | Not relevant.                  |
| Self-ignition temperature :               | Not relevant.                  |
| Decomposition point/decomposition range : | Not relevant.                  |

Name : CAVITY LINER LIGHT CURE - ART.11035

---

## 9.2. Other information

No data available.

---

## SECTION 10 : STABILITY AND REACTIVITY

### 10.1. Reactivity

No data available.

### 10.2. Chemical stability

This mixture is stable under the recommended handling and storage conditions in section 7.

Do not use if the expiry date is exceeded.

### 10.3. Possibility of hazardous reactions

No data available.

### 10.4. Conditions to avoid

Avoid :

- exposure to light

### 10.5. Incompatible materials

### 10.6. Hazardous decomposition products

The thermal decomposition may release/form :

- carbon monoxide (CO)
- carbon dioxide (CO<sub>2</sub>)

---

## SECTION 11 : TOXICOLOGICAL INFORMATION

### 11.1. Information on toxicological effects

May cause an allergic reaction by skin contact.

#### 11.1.1. Substances

No toxicological data available for the substances.

#### 11.1.2. Mixture

No toxicological data available for the mixture.

---

## SECTION 12 : ECOLOGICAL INFORMATION

### 12.1. Toxicity

#### 12.1.2. Mixtures

No aquatic toxicity data available for the mixture.

### 12.2. Persistence and degradability

No data available.

### 12.3. Bioaccumulative potential

No data available.

### 12.4. Mobility in soil

No data available.

### 12.5. Results of PBT and vPvB assessment

No data available.

### 12.6. Other adverse effects

No data available.

Name : CAVITY LINER LIGHT CURE - ART.11035

---

### SECTION 13 : DISPOSAL CONSIDERATIONS

Proper waste management of the mixture and/or its container must be determined in accordance with Directive 2008/98/EC.

#### 13.1. Waste treatment methods

Do not pour into drains or waterways.

#### Waste :

Waste management is carried out without endangering human health, without harming the environment and, in particular without risk to water, air, soil, plants or animals.

Recycle or dispose of waste in compliance with current legislation, preferably via a certified collector or company.

Do not contaminate the ground or water with waste, do not dispose of waste into the environment.

#### Soiled packaging :

Empty container completely. Keep label(s) on container.

Give to a certified disposal contractor.

---

### SECTION 14 : TRANSPORT INFORMATION

Exempt from transport classification and labelling.

Transport product in compliance with provisions of the ADR for road, RID for rail, IMDG for sea and ICAO/IATA for air transport (ADR 2013 - IMDG 2012 - ICAO/IATA 2014).

---

### SECTION 15 : REGULATORY INFORMATION

#### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

##### - Classification and labelling information included in section 2:

The following regulations have been used:

- Directive 67/548/EEC and its adaptations
- Directive 1999/45/EC and its adaptations
- Regulation EC 1272/2008 modified by regulation EC 618/2012
- EU Regulation No. 1272/2008 amended by EU Regulation No. 758/2013.

##### - Container information:


The mixture is contained in packaging that does not exceed 125 ml.

##### - Particular provisions :

No data available.

##### - Standardised American system for the identification of hazards presented by the product in view of emergency procedures (NFPA 704) :

NFPA 704, Labelling: Health=2 Inflammability=1 Instability/Reactivity=1 Specific Risk=none



#### 15.2. Chemical safety assessment

No data available.

## SECTION 16 : OTHER INFORMATION

Since the user's working conditions are not known by us, the information supplied on this safety data sheet is based on our current level of knowledge and on national and community regulations.

The mixture must not be used for other uses than those specified in section 1 without having first obtained written handling instructions.

It is at all times the responsibility of the user to take all necessary measures to comply with legal requirements and local regulations.

The information in this safety data sheet must be regarded as a description of the safety requirements relating to the mixture and not as a guarantee of the properties thereof.

### In compliance with directives 67/548/EEC, 1999/45/EC and their amendments.

Hazard symbols :



Irritant

Contains :

607-132-00-3

DIMETHYLAMINOETHYL METHACRYLATE 2-

Risk phrase :

R 52/53

Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

R 43

May cause sensitisation by skin contact.

Safety phrase :

S 24

Avoid contact with skin.

S 37

Wear suitable gloves.

S 60

This material and its container must be disposed of as hazardous waste.

S 61

Avoid release to the environment. Refer to special instructions/Safety data sheets.

### Title for H, EUH and R indications mentioned in section 3 :

H302

Harmful if swallowed.

H312

Harmful in contact with skin.

H315

Causes skin irritation.

H317

May cause an allergic skin reaction.

H319

Causes serious eye irritation.

R 21/22

Harmful in contact with skin and if swallowed.

R 36/38

Irritating to eyes and skin.

R 43

May cause sensitisation by skin contact.

### Abbreviations :

ADR : European agreement concerning the international carriage of dangerous goods by Road.

IMDG : International Maritime Dangerous Goods.

IATA : International Air Transport Association.

ICAO : International Civil Aviation Organisation

RID : Regulations concerning the International carriage of Dangerous goods by rail.

GHS07 : Exclamation mark