



UNIGLOVES COMFORT®

SIZE	ITEM-NO
------	---------



XS	1201
----	------

S	1202
---	------

M	1203
---	------

L	1204
---	------

XL	1205
----	------



EIGENSCHAFTEN

MATERIAL	<u>Latex</u>
COLOUR	<u>natural latex</u>
LÄNGE	<u>240 mm</u>
SIZE	<u>XS (5.5-6.0), S (6.5-7.0), M (7.5-8.0), L (8.5-9), XL (9.5-10)</u>
AQL	<u>1.5</u>
PRODUCT FEATURES	<u>beaded cuff, chlorinated on the inside & outside, non-sterile, powder-free, smooth palm</u>
NORM GLOVES	<u>CE, EN 374, EN 455, T. 1-4</u>

VERPACKUNGSGEHALT

box of 100 pieces

VERSANDEINHEIT

10 boxes