

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Mixture
 Product name : Cementspacer blue, grey and red

1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses

Main use category : Professional use
 Use of the substance/mixture : Preparation intended for dental technical applications

Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Kerr Italia S.r.l.
 Via Passanti, 174
 Scafati (SA) 84018 (Italy)
 T +39 081 850 3511 / Information Phone Number: 1-800-841-1428 (Customer Service)
safety@envistaco.com

1.4. Emergency telephone number

Emergency number : CHEMTREC® Emergency Call Center. Emergency Telephone Number (for USA only) 001-800-424-9300
 International and Maritime Telephone Number +1 (703) 527-3887

Country	Official advisory body	Address	Emergency number	Comment
Gibraltar	GHA Call Centre Zone 2, Level3, St Bernard's Hospital	Harbour Views Road	+350 200 79700 +350 200 72266	
Ireland	National Poisons Information Centre Beaumont Hospital	PO Box 1297 Beaumont Road 9 Dublin	+353 1 809 2566 (Healthcare professionals-24/7) +353 1 809 2166 (public, 8am - 10pm, 7/7)	
Malta	Medicines & Poisons Info Office	Mater Dei Hospital MSD 2090 Msida	+356 2545 6508	
United Kingdom	National Poisons Information Service (Newcastle Unit)	Claremont Place Newcastle-upon-Tyne, Newcastle	+44 191 2606182 +44 191 2606180	Hours of operation: 24hrs

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Flam. Liq. 2 H225
 Eye Irrit. 2 H319
 STOT SE 3 H335
 STOT SE 3 H336

Full text of hazard classes and H-statements : see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :



Safety Data Sheet

Cementspacer blue, grey and red

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

	GHS02	GHS07
Signal word (CLP)	: Danger	
Hazard statements (CLP)	: H225 - Highly flammable liquid and vapour. H319 - Causes serious eye irritation. H335 - May cause respiratory irritation. H336 - May cause drowsiness or dizziness.	
Precautionary statements (CLP)	: P210 - Keep away from heat/sparks/open flames/hot surfaces. – No smoking. P243 - Take precautionary measures against static discharge. P261 - Avoid breathing vapours. P262 - Do not get in eyes, on skin, or on clothing. P280 - Wear eye protection, protective gloves. P403+P233 - Store in a well-ventilated place. Keep container tightly closed. P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.	
EUH-statements	: EUH066 - Repeated exposure may cause skin dryness or cracking.	

2.3. Other hazards

Other hazards which do not result in classification : None under normal conditions.

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII

This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

The mixture does not contain substance(s) included in the list established in accordance with Article 59(1) of REACH for having endocrine disrupting properties, or is not identified as having endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at a concentration equal to or greater than 0,1 %

SECTION 3: Composition/information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

Name	Product identifier	Conc. (% w/w)	Classification according to Regulation (EC) No. 1272/2008 [CLP]
pentan-3-one; diethyl ketone	(CAS-No.) 96-22-0 (EC-No.) 202-490-3 (EC Index-No.) 606-006-00-5 (REACH-no) 01-2119531111-60	25 – 50	Flam. Liq. 2, H225 STOT SE 3, H335 STOT SE 3, H336
butanone; ethyl methyl ketone	(CAS-No.) 78-93-3 (EC-No.) 201-159-0 (EC Index-No.) 606-002-00-3 (REACH-no) 01-2119457290-43	25 – 50	Flam. Liq. 2, H225 Eye Irrit. 2, H319 STOT SE 3, H336
n-butyl acetate	(CAS-No.) 123-86-4 (EC-No.) 204-658-1 (EC Index-No.) 607-025-00-1 (REACH-no) 01-2119485493-29	10 – 25	Flam. Liq. 3, H226 STOT SE 3, H336

Full text of H- and EUH-statements: see section 16

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general	: Never give anything by mouth to an unconscious person. If you feel unwell, seek medical advice (show the label where possible).
First-aid measures after inhalation	: Remove person to fresh air and keep comfortable for breathing. Get medical advice/attention if you feel unwell.
First-aid measures after skin contact	: Generally the product does not irritate the skin.
First-aid measures after eye contact	: Rinse cautiously with water for several minutes. Obtain medical attention if pain, blinking or redness persists.
First-aid measures after ingestion	: Rinse mouth. If you feel unwell, seek medical advice.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects	: In all cases of doubt, or when symptoms persist, seek medical attention.
Symptoms/effects after inhalation	: May cause respiratory irritation. May cause drowsiness or dizziness.

Safety Data Sheet

Cementspacer blue, grey and red

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

Symptoms/effects after eye contact : Causes serious eye irritation.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically. In all cases of doubt, or when symptoms persist, seek medical attention.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Dry powder. Carbon dioxide. Water spray. Fight larger fires with water spray or alcohol resistant foam.
Unsuitable extinguishing media : Do not use a heavy water stream.

5.2. Special hazards arising from the substance or mixture

Fire hazard : Highly flammable liquid and vapour.
Hazardous decomposition products in case of fire : Carbon dioxide. Carbon monoxide.

5.3. Advice for firefighters

Firefighting instructions : Use water spray or fog for cooling exposed containers.
Protection during firefighting : Do not enter fire area without proper personal protective equipment, including respiratory protection (EN137).

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures : Remove ignition sources. No open flames. No smoking. Use special care to avoid static electric charges.
Do not eat, drink or smoke in areas where product is used. Avoid contact with skin and eyes.

6.1.1. For non-emergency personnel

Emergency procedures : Evacuate unnecessary personnel.

6.1.2. For emergency responders

Protective equipment : Equip cleanup and emergency crew with proper protection.
Emergency procedures : Ventilate area.

6.2. Environmental precautions

Avoid release to the environment. Prevent entry to sewers and public waters.

6.3. Methods and material for containment and cleaning up

For containment : Collect all waste in suitable and labelled containers and dispose according to local legislation.
Methods for cleaning up : Small quantities of liquid spill: take up in non-combustible absorbent material and shovel into container for disposal. Take up large spills with pump or vacuum and finish with dry chemical absorbent.

6.4. Reference to other sections

For further information refer to section 13. See Section 8. Exposure controls and personal protection.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Do not eat, drink or smoke when using this product. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Take precautionary measures against static discharge.
Hygiene measures : Do not eat, drink or smoke when using this product. Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep only in original container. Keep container tightly closed. Store in a dry place. Store in a well-ventilated place. Keep cool.
Incompatible materials : Refer to Section 10 on Incompatible Materials.
Storage area : Store in a well-ventilated place.

7.3. Specific end use(s)

No additional data.

Safety Data Sheet

Cementspacer blue, grey and red

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

8.1.1 National occupational exposure and biological limit values

n-butyl acetate (123-86-4)	
United Kingdom - Occupational Exposure Limits	
Local name	Butyl acetate
WEL TWA (OEL TWA) [1]	724 mg/m ³
WEL TWA (OEL TWA) [2]	150 ppm
WEL STEL (OEL STEL)	966 mg/m ³
WEL STEL (OEL STEL) [ppm]	200 ppm
Regulatory reference	EH40/2005 (Fourth edition, 2020). HSE

pentan-3-one; diethyl ketone (96-22-0)	
United Kingdom - Occupational Exposure Limits	
Local name	Pentan-3-one
WEL TWA (OEL TWA) [1]	716 mg/m ³
WEL TWA (OEL TWA) [2]	200 ppm
WEL STEL (OEL STEL)	895 mg/m ³
WEL STEL (OEL STEL) [ppm]	250 ppm
Regulatory reference	EH40/2005 (Fourth edition, 2020). HSE

butanone; ethyl methyl ketone (78-93-3)	
Gibraltar - Occupational Exposure Limits	
Local name	Butanone
OEL TWA	600 mg/m ³
OEL TWA [ppm]	200 ppm
OEL STEL	900 mg/m ³
OEL STEL [ppm]	300 ppm
Regulatory reference	Factories (Control of Chemical Agents at Work) Regulations 2003 (LN. 2018/181)
Ireland - Occupational Exposure Limits	
Local name	Methyl ethyl ketone (MEK)
OEL TWA [1]	600 mg/m ³
OEL TWA [2]	200 ppm
OEL STEL	900 mg/m ³
OEL STEL [ppm]	300 ppm
Remark	Sk (Substances which have the capacity to penetrate intact skin when they come in contact with it, and be absorbed into the body), IOELV (Indicative Occupational Exposure Limit Values)
Regulatory reference	Chemical Agents Code of Practice 2021
Malta - Occupational Exposure Limits	
Local name	Butanone
OEL TWA	600 mg/m ³
OEL TWA [ppm]	200 ppm
OEL STEL	900 mg/m ³
OEL STEL [ppm]	300 ppm

Safety Data Sheet

Cementspacer blue, grey and red

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

Regulatory reference	S.L.424.24 - Chemical Agents at Work Regulations (L.N.57 of 2018)
United Kingdom - Occupational Exposure Limits	
Local name	Butan-2-one (methyl ethyl ketone)
WEL TWA (OEL TWA) [1]	600 mg/m ³
WEL TWA (OEL TWA) [2]	200 ppm
WEL STEL (OEL STEL)	899 mg/m ³
WEL STEL (OEL STEL) [ppm]	300 ppm
Remark	Sk (Can be absorbed through the skin. The assigned substances are those for which there are concerns that dermal absorption will lead to systemic toxicity)
Regulatory reference	EH40/2005 (Fourth edition, 2020). HSE
United Kingdom - Biological limit values	
Local name	Butan-2-one (methyl ethyl ketone)
BMGV	70 µmol/l Parameter: butan-2-one - Medium: urine - Sampling time: Post shift
Regulatory reference	EH40/2005 (Fourth edition, 2020). HSE

8.1.2. Recommended monitoring procedures

No additional information available

8.1.3. Air contaminants formed

No additional information available

8.1.4. DNEL and PNEC

No additional information available

8.1.5. Control banding

No additional information available

8.2. Exposure controls

8.2.1. Appropriate engineering controls

Appropriate engineering controls:

Ensure good ventilation of the work station.

8.2.2. Personal protection equipment

Personal protective equipment:

Gloves. Safety glasses.



8.2.2.1. Eye and face protection

Eye protection:

Wear security glasses which protect from splashes. STANDARD EN 166.

8.2.2.2. Skin protection

Skin and body protection:

In case of possible repeated skin contact wear protective clothing

Hand protection:

Wear suitable gloves. Nitrile rubber gloves. Butyl rubber. Viton. Layer thickness : 0,09mm. Breakthrough time : >480 min. STANDARD EN 374.

8.2.2.3. Respiratory protection

Respiratory protection:

No special respiratory protection equipment is recommended under normal conditions of use with adequate ventilation. In case of inadequate ventilation wear respiratory protection. Extra personal protection: A/P2 filter respirator for organic vapour and harmful dust. Filter AX (brown). STANDARD EN 149:2001 + A1:2009

8.2.2.4. Thermal hazards

No additional information available

Safety Data Sheet

Cementspacer blue, grey and red

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

8.2.3. Environmental exposure controls

Other information:

Do not eat, drink or smoke during use. Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Personal protective equipment should be chosen according to the CEN standards and in discussion with the supplier of the protective equipment.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	: Liquid
Colour	: Various colours.
Appearance	: Liquid.
Odour	: Characteristic.
Odour threshold	: Not available
Melting point	: Not determined.
Freezing point	: Not determined.
Boiling point	: 79 °C
Flammability	: Not applicable
Explosive properties	: Not explosive.
Explosive limits	: Not available
Lower explosive limit (LEL)	: 1.2 vol %
Upper explosive limit (UEL)	: 11.5 vol %
Flash point	: -4 °C
Auto-ignition temperature	: Not self-igniting.
Decomposition temperature	: Not determined.
pH	: Not determined.
Viscosity, kinematic	: Not determined.
Viscosity, dynamic	: Not determined.
Solubility	: Insoluble in water.
Partition coefficient n-octanol/water (Log Kow)	: Not available
Vapour pressure	: 105 hPa
Vapour pressure at 50 °C	: Not available
Density	: 0.9 g/cm ³
Relative density	: Not determined.
Relative vapour density at 20 °C	: Not determined.
Particle size	: Not applicable
Particle size distribution	: Not applicable
Particle shape	: Not applicable
Particle aspect ratio	: Not applicable
Particle aggregation state	: Not applicable
Particle agglomeration state	: Not applicable
Particle specific surface area	: Not applicable
Particle dustiness	: Not applicable

9.2. Other information

9.2.1. Information with regard to physical hazard classes

No additional information available

9.2.2. Other safety characteristics

Relative evaporation rate (butylacetate=1)	: Not determined.
VOC content	: 67 %

SECTION 10: Stability and reactivity

10.1. Reactivity

None to our knowledge.

10.2. Chemical stability

Stable under normal conditions of use.

10.3. Possibility of hazardous reactions

No polymerization.

10.4. Conditions to avoid

No flames, no sparks. Eliminate all sources of ignition. Keep away from heat and direct sunlight.

Safety Data Sheet

Cementspacer blue, grey and red

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

10.5. Incompatible materials

None to our knowledge.

10.6. Hazardous decomposition products

No decomposition if stored and used normally.

SECTION 11: Toxicological information

11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

Acute toxicity (oral) : Not classified
Acute toxicity (dermal) : Not classified
Acute toxicity (inhalation) : Not classified
Additional information : Based on available data, the classification criteria are not met

n-butyl acetate (123-86-4)

LD50 oral rat	10768 mg/kg
LD50 dermal rabbit	> 17600 mg/kg
LC50 Inhalation - Rat (Vapours)	40 mg/l/4h

pentan-3-one; diethyl ketone (96-22-0)

LD50 oral rat	2140 mg/kg
LD50 dermal rabbit	16000 mg/kg
LC50 Inhalation - Rat (Vapours)	> 29 mg/l/4h

butanone; ethyl methyl ketone (78-93-3)

LD50 oral rat	> 2000 mg/kg
LD50 dermal rabbit	> 5000 mg/kg
LC50 Inhalation - Rat	34.5 mg/l/4h

Skin corrosion/irritation : Not classified
pH: Not determined.
Additional information : Based on available data, the classification criteria are not met
Serious eye damage/irritation : Causes serious eye irritation.
pH: Not determined.
Respiratory or skin sensitisation : Not classified
Additional information : Based on available data, the classification criteria are not met
Germ cell mutagenicity : Not classified
Additional information : Based on available data, the classification criteria are not met
Carcinogenicity : Not classified
Additional information : Based on available data, the classification criteria are not met
Reproductive toxicity : Not classified
Additional information : Based on available data, the classification criteria are not met
STOT-single exposure : May cause respiratory irritation. May cause drowsiness or dizziness.

n-butyl acetate (123-86-4)

STOT-single exposure	May cause drowsiness or dizziness.
----------------------	------------------------------------

pentan-3-one; diethyl ketone (96-22-0)

STOT-single exposure	May cause respiratory irritation. May cause drowsiness or dizziness.
----------------------	--

butanone; ethyl methyl ketone (78-93-3)

STOT-single exposure	May cause drowsiness or dizziness.
----------------------	------------------------------------

STOT-repeated exposure : Not classified

Safety Data Sheet

Cementspacer blue, grey and red

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

Additional information : Based on available data, the classification criteria are not met

Aspiration hazard : Not classified

Additional information : Based on available data, the classification criteria are not met

Cementspacer blue, grey and red

Viscosity, kinematic	Not determined.
----------------------	-----------------

11.2. Information on other hazards

11.2.1. Endocrine disrupting properties

Adverse health effects caused by endocrine disrupting properties : The mixture does not contain substance(s) included in the list established in accordance with Article 59(1) of REACH for having endocrine disrupting properties, or is not identified as having endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at a concentration equal to or greater than 0,1 %

11.2.2 Other information

Potential adverse human health effects and symptoms : For further information see section 4

SECTION 12: Ecological information

12.1. Toxicity

Hazardous to the aquatic environment, short-term (acute) : Not classified

Hazardous to the aquatic environment, long-term (chronic) : Not classified

n-butyl acetate (123-86-4)

LC50 - Fish [1]	18 mg/l (96 hours - Pimephales promelas)
-----------------	--

EC50 - Crustacea [1]	10 – 100 mg/l (48 hours - Daphnia magna)
----------------------	--

EC50 - Other aquatic organisms [1]	24 mg/l EC50 24h - Daphnia magna [mg/l]
------------------------------------	---

ErC50 algae	674.7 mg/l (IC50, 72 hours -Selenastrum)
-------------	--

pentan-3-one; diethyl ketone (96-22-0)

LC50 - Fish [1]	< 1540 mg/l (96 hours - Pimephales promelas)
-----------------	--

LC50 - Other aquatic organisms [1]	> 5000 mg/l (EC50, 72 hours - Selenastrum capricornutum)
------------------------------------	--

EC50 - Crustacea [1]	> 500 mg/l (48 hours - Daphnia magna)
----------------------	---------------------------------------

ErC50 algae	> 500 mg/l (IC50, 72 hours -Selenastrum)
-------------	--

butanone; ethyl methyl ketone (78-93-3)

LC50 - Fish [1]	> 400 mg/l (96 hours - Lepomis macrochirus)
-----------------	---

EC50 - Crustacea [1]	> 300 mg/l (48 hours - Daphnia magna)
----------------------	---------------------------------------

Threshold limit - Algae [1]	110 mg/l
-----------------------------	----------

12.2. Persistence and degradability

Cementspacer blue, grey and red

Persistence and degradability	No data available.
-------------------------------	--------------------

n-butyl acetate (123-86-4)

BOD (% of ThOD)	> 0.5 % ThOD BOD5/COD
-----------------	-----------------------

Biodegradation	98 % (OECD 301B method)
----------------	-------------------------

pentan-3-one; diethyl ketone (96-22-0)

BOD (% of ThOD)	0.58 % ThOD BOD5/COD
-----------------	----------------------

Safety Data Sheet

Cementspacer blue, grey and red

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

12.3. Bioaccumulative potential

Cementspacer blue, grey and red	
Bioaccumulative potential	No data.

n-butyl acetate (123-86-4)	
Bioconcentration factor (BCF REACH)	14
Partition coefficient n-octanol/water (Log Pow)	1.81

pentan-3-one; diethyl ketone (96-22-0)	
Partition coefficient n-octanol/water (Log Pow)	0.99

butanone; ethyl methyl ketone (78-93-3)	
Partition coefficient n-octanol/water (Log Pow)	0.61

12.4. Mobility in soil

Cementspacer blue, grey and red	
Ecology - soil	Insoluble in water.

12.5. Results of PBT and vPvB assessment

Cementspacer blue, grey and red	
This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII	
This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII	

12.6. Endocrine disrupting properties

No additional information available

12.7. Other adverse effects

Other adverse effects : None to our knowledge.
Additional information : No other effects known

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods : Must not be disposed together with household garbage. Do not allow product to reach the sewage system.
Product/Packaging disposal recommendations : Dispose in a safe manner in accordance with local/national regulations. Empty container completely. After cleaning, recycle or dispose of at an authorised site.
European List of Waste (LoW) code : 18 01 06* - chemicals consisting of or containing dangerous substances

SECTION 14: Transport information

In accordance with ADR / IMDG / IATA / ADN / RID /

ADR	IMDG	IATA	ADN	RID
14.1. UN number or ID number				
UN 1263				
14.2. UN proper shipping name				
PAINT (butanone; ethyl methyl ketone ; pentan-3-one; diethyl ketone ; n-butyl acetate)	PAINT (butanone; ethyl methyl ketone ; pentan-3-one; diethyl ketone ; n-butyl acetate)	Paint (butanone; ethyl methyl ketone ; pentan-3-one; diethyl ketone ; n-butyl acetate)	PAINT (butanone; ethyl methyl ketone ; pentan-3-one; diethyl ketone ; n-butyl acetate)	PAINT (butanone; ethyl methyl ketone ; pentan-3-one; diethyl ketone ; n-butyl acetate)

Safety Data Sheet

Cementspacer blue, grey and red

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

Transport document description				
UN 1263 PAINT (butanone; ethyl methyl ketone ; pentan-3-one; diethyl ketone ; n-butyl acetate), 3, II, (D/E)	UN 1263 PAINT (butanone; ethyl methyl ketone ; pentan-3-one; diethyl ketone ; n-butyl acetate), 3, II	UN 1263 Paint (butanone; ethyl methyl ketone ; pentan-3-one; diethyl ketone ; n-butyl acetate), 3, II	UN 1263 PAINT (butanone; ethyl methyl ketone ; pentan-3-one; diethyl ketone ; n-butyl acetate), 3, II	UN 1263 PAINT (butanone; ethyl methyl ketone ; pentan-3-one; diethyl ketone ; n-butyl acetate), 3, II
14.3. Transport hazard class(es)				
3	3	3	3	3
				
14.4. Packing group				
II	II	II	II	II
14.5. Environmental hazards				
Dangerous for the environment : No	Dangerous for the environment : No Marine pollutant : No	Dangerous for the environment : No	Dangerous for the environment : No	Dangerous for the environment : No
No supplementary information available				

14.6. Special precautions for user

Overland transport

Classification code (ADR) : F1
 Special provisions (ADR) : 163, 367, 640C, 650
 Limited quantities (ADR) : 5I
 Excepted quantities (ADR) : E2
 Orange plates :



EAC code : •3YE

Transport by sea

Special provisions (IMDG) : 163, 367
 Limited quantities (IMDG) : 5 L
 Excepted quantities (IMDG) : E2
 EmS-No. (Fire) : F-E
 EmS-No. (Spillage) : S-E

Air transport

PCA Excepted quantities (IATA) : E2
 PCA Limited quantities (IATA) : Y341
 Special provisions (IATA) : A3, A72, A192

Inland waterway transport

Classification code (ADN) : F1
 Special provisions (ADN) : 163, 367, 640C, 650
 Limited quantities (ADN) : 5 L
 Excepted quantities (ADN) : E2

Rail transport

Special provisions (RID) : 163, 367, 640C, 650
 Limited quantities (RID) : 5L
 Excepted quantities (RID) : E2
 Hazard identification number (RID) : 33

14.7. Maritime transport in bulk according to IMO instruments

Not applicable

Safety Data Sheet

Cementspacer blue, grey and red

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions

Contains no substance on the REACH candidate list

Contains no REACH Annex XIV substances

Contains no substance subject to Regulation (EU) No 649/2012 of the European Parliament and of the Council of 4 July 2012 concerning the export and import of hazardous chemicals.

Contains no substance subject to Regulation (EU) No 2019/1021 of the European Parliament and of the Council of 20 June 2019 on persistent organic pollutants

VOC content : 67 %

15.1.2. National regulations

REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006

15.2. Chemical safety assessment

No chemical safety assessment has been carried out

SECTION 16: Other information

Data sources : REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

Other information : None.

Full text of H- and EUH-statements:

EUH066	Repeated exposure may cause skin dryness or cracking.
Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
Flam. Liq. 2	Flammable liquids, Category 2
Flam. Liq. 3	Flammable liquids, Category 3
H225	Highly flammable liquid and vapour.
H226	Flammable liquid and vapour.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
H336	May cause drowsiness or dizziness.
STOT SE 3	Specific target organ toxicity – Single exposure, Category 3, Respiratory tract irritation

The information in this safety data sheet is based on information from the manufacturer/supplier, present European and national legislation, and presupposes that the product is used within the specified area of application.