






BASIC SHADE IMAGES



- INCISAL
- TRANSLUCENT
- BODY
- CERVICAL
- OPAQUE
- METAL



LIGHT-CURING TIME

		Solidilite V	Solidilite EX
SHOFU Universal Opaque	Pre-Opaque	60	60
	Opaque	180	180
Flowable Composite (Pre-light-cure)		60	60
Final		180	300

(Seconds)

COLOUR	W0	W1	W2	W3	A1	A2	A3	A3.5	A4	rootA
UNIVERSAL OPAQUE	W0	W1	W2	W3	A1	A2	A3	A3.5	A4	rootA
CERVICAL	-				AC1			AC2		-
OPAQUE DENTIN	-				ODA1	ODA2	ODA3	ODA3.5	ODA4	ODrootA
BODY	W0B	W1B	W2B	W3B	A1B	A2B	A3B	A3.5B	A4B	rootAB
INCISAL	56		57		58		59		60	

COLOUR	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
UNIVERSAL OPAQUE	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
CERVICAL	-	BC1		BC2	-	CC1	CC2	DC1		DC2	
OPAQUE DENTIN	ODB1	ODB2	ODB3	ODB4	ODC1	ODC2	ODC3	ODC4	ODD2	ODD3	ODD4
BODY	B1B	B2B	B3B	B4B	C1B	C2B	C3B	C4B	D2B	D3B	D4B
INCISAL	57	58	59	60	58	59	60	59	60	59	

UNIVERSAL OPAQUE	MO	BGO	GO	WO	GUM-O	MO	VPO	VPO	PO	FreeO	
TRANSLUCENT	T	HVT	LVT	T-Glass	BG	GT	CT-A	CT-B			
INCISAL	56	57	58	59	60						
CONCENTRATE	MY	WE	OC	AM-Y	W						
GUM COLOR	GUM-L	GUM-D	GUM-Or	GUM-Br	GUM-V	GUM-R	GUM-LP	GUM-DP	GUM-SP	GUM-P	GUM-T

The printed colours might be different from the actual colours in this system.