




DIAGNOcam 2170 This is the basic equipment with DIAGNOcam, accessories and software for single place option. The software works with an integrated database with patient data.

Mehrplatzoption The software and data is stored on the server of the practice. Like this you have access to the data from different computers.

VDDS DIAGNOcam pictures can be saved in the patients' files of the practice management software. It can only be used with patient management software with VDDS interface. This interface is only relevant in Germany, Austria and Switzerland.

TWAIN Enables the transfer of the pictures to your x-ray software or any image editing software (If this function is supported by the connected software). The transfer of the data can be inconvenient, depending on your software.



Additional Questions:

Which practice management software do you use?

Do you have your own IT specialist?

☐ Yes - Phone number: _____

☐ No

Where would you like to illustrate the pictures?

☐ Monitor on the treatment unit – Which treatment unit do you use?

☐ Computer – Which computer / operating system do you use?

☐ Apple: _____

☐ Microsoft: _____

☐ Others: _____



KaVo. Dental Excellence.