

Hardwareanforderungen Intraoralscanner

CPU	Desktop: Intel Core i7 9700 / AMD Ryzen7 2700; Recommended: Intel Core i7 12700K / AMD Ryzen 7 5800X; Laptop: Intel Core i7 10750H / AMD Ryzen 7-Pro 4750U; Recommended: Intel Core i7 12700H / AMD Ryzen 7 6800H	System	Windows 10/11 64 Bit, Home/Pro/Corporate
Graphics-Card	NVIDIA GeForce 1660GTX or above, more than 6GB; Recommended: NVIDIA GeForce 3060 or above graphics with 6GB or more memory. Note: AMD Graphics cards are not supported.	Monitor resolution	1920 x 1080 or above
Memory	More than 16 GB	Port	Type USB-C, or USB 3.0
Harddrive	More than 512 GB SSD	Power consumption	25 V/A