



Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 16.10.2012

Version number 1

Revision: 16.10.2012

SECTION 1: Identification of the substance/mixture and of the company/undertaking

- **1.1 Product identifier**
- **Trade name:** DETARTRINE 150 Z
- **Article number / formula:** 274S
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**
No further relevant information available.
- **Application of the substance / the preparation :** dental product
- **1.3 Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
SEPTODONT.
58 rue du Pont de Creteil
94100 Saint Maur des Fosses
France
Téléphone: +33(0)149767000
Fax: +33(0)148855401
msds@septodont.com
<http://www.septodont.com>
- **1.4 Emergency telephone number:**
France (ORFILA 24h/24) - Tel : +33 (0)1 45 42 59 59
Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24)
EU Tel : 112

SECTION 2: Hazards identification

- **2.1 Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008 :**
The product is not classified according to the CLP regulation.
- **Classification according to Directive 67/548/EEC or Directive 1999/45/EC :** Not applicable.
- **Information concerning particular hazards for human and environment:**
The product does not have to be labelled due to the calculation procedure of the "General Classification guideline for preparations of the EU" in the latest valid version.
- **Classification system:**
The classification is according to the latest editions of the EU-lists, and extended by company and literature data.
- **2.2 Label elements**
- **Labelling according to EU guidelines:**
Observe the general safety regulations when handling chemicals.
The product has not been classified and marked in accordance with EU Directives / respective national laws. 2008/58/CE (30eme ATP) ; 2009/2/CE (31eme ATP) ; 2006/8/CE
- **Special labelling of certain preparations:**
Safety data sheet available for professional user on request.
- **2.3 Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

SECTION 3: Composition/information on ingredients

- **3.2 Chemical characterization: Mixtures**
- **Description:** Mixture: consisting of the following components.

(Contd. on page 2)

GB



Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 16.10.2012

Version number 1

Revision: 16.10.2012

Trade name: DETARTRINE 150 Z

(Contd. of page 1)

· Dangerous components:

CAS: 10101-52-7	Zirconium(IV) silicate ☒ Xn R20; ☒ Xi R36/37/38 ⚠ Acute Tox. 4, H332; Skin Irrit. 2, H315; Eye Irrit. 2, H319; STOT SE 3, H335	10-25%
CAS: 64-17-5 EINECS: 200-578-6 Reg.nr.: 01-2119457610-43-xxxx	ethanol ☒ F R11 ⚠ Flam. Liq. 2, H225; ⚠ Eye Irrit. 2, H319	≤2.5%
CAS: 89-83-8 EINECS: 201-944-8	thymol ☒ C R34; ☒ Xn R22; ☒ N R51/53 ⚠ Skin Corr. 1B, H314; ⚠ Aquatic Chronic 2, H411; ⚠ Acute Tox. 4, H302	≤1%

· Additional information: For the wording of the listed risk phrases, refer to section 16.

SECTION 4: First aid measures

· 4.1 Description of first aid measures

- General information:** No special measures required.
- After excessive inhalation:** Supply fresh air; consult a doctor in case of complaints.
- After skin contact:** Generally the product does not irritate the skin.
- After eye contact:** Rinse opened eye for several minutes under running water.
- After swallowing:** If symptoms persist consult a doctor.

· 4.2 Most important symptoms and effects, both acute and delayed :

No further relevant information available.

· 4.3 Indication of any immediate medical attention and special treatment needed :

No further relevant information available.

SECTION 5: Firefighting measures

· 5.1 Suitable extinguishing agents:CO₂, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.**· 5.2 Special hazards arising from the substance or mixture :** No further relevant information available.**· 5.3 Advice for firefighters****· Protective equipment:** No special measures required.

SECTION 6: Accidental release measures

· 6.1 Personal precautions, protective equipment and emergency procedures : Not required.**· 6.2 Environmental precautions:** Do not allow to enter sewers/ surface or ground water.**· 6.3 Methods and material for containment and cleaning up:** Pick up mechanically.**· 6.4 Reference to other sections :**

No dangerous substances are released.

See Section 7 for information on safe handling.

See Section 8 for information on personal protective equipment.

See Section 13 for disposal information.

SECTION 7: Handling and storage

· 7.1 Precautions for safe handling : No special measures required.

(Contd. on page 3)



Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 16.10.2012

Version number 1

Revision: 16.10.2012

Trade name: DETARTRINE 150 Z

(Contd. of page 2)

- **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **7.3 Specific end use(s) :** No further relevant information available.

SECTION 8: Exposure controls/personal protection

- **Additional information about design of technical facilities:** No further data; see section 7.

· 8.1 Control parameters

- **Ingredients with limit values that require monitoring at the workplace:**

56-81-5 glycerol (10-25%)WEL (Great Britain) Long-term value: 10 mg/m³**64-17-5 ethanol (≤2.5%)**WEL (Great Britain) Long-term value: 1920 mg/m³, 1000 ppm

- **Additional information:**

As a basis for the production of this document, the most current valid lists were used.

· 8.2 Exposure controls

- **Personal protective equipment:**

- **General protective and hygienic measures:**

The usual precautionary measures are to be adhered to when handling chemicals.

- **Respiratory protection:** Not required.

- **Protection of hands:**

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

- **Gloves material :**

The selection of suitable gloves does not only depend on the material, but also on further marks of quality which may vary from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked before use.

- **Eye protection:** Not required.

SECTION 9: Physical and chemical properties

- **9.1 Information on basic physical and chemical properties**

- **General Information**

- **Appearance:**

Form: Pasty**Colour:** Pink

- **Odour:** Characteristic

- **Odour threshold:** Not determined.

- **pH-value:** Not applicable.

- **Change in condition**

Melting point/Melting range: Undetermined.**Boiling point/Boiling range:** Undetermined.

(Contd. on page 4)

GB

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 16.10.2012

Version number 1

Revision: 16.10.2012

Trade name: DETARTRINE 150 Z

(Contd. of page 3)

· Flash point:	Not applicable.
· Flammability (solid, gaseous):	Not determined.
· Ignition temperature:	400°C
· Decomposition temperature:	Not determined.
· Self-igniting:	Product is not self igniting.
· Danger of explosion:	Product does not present an explosion hazard.
· Explosion limits:	
Lower:	0.9 Vol %
Upper:	Not determined.
· Vapour pressure at 20°C:	0.1 hPa
· Density:	Not determined.
· Relative density :	Not determined.
· Vapour density :	Not applicable.
· Evaporation rate :	Not applicable.
· Solubility in / Miscibility with water:	Insoluble.
· Partition coefficient (n-octanol/water):	Not determined.
· Viscosity:	
Dynamic:	Not applicable.
Kinematic:	Not applicable.
· Solvent content:	
VOC (EC)	0.00 %
· 9.2 Other information :	No further relevant information available.

SECTION 10: Stability and reactivity

- **10.1 Reactivity**
- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions :** No dangerous reactions known.
- **10.4 Conditions to avoid :** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

SECTION 11: Toxicological information

- **11.1 Information on toxicological effects**
- **Acute toxicity:**
- **Primary irritant effect:**
- **On the skin:** No irritant effect.
- **On the eyes:** No irritant effect.
- **Sensitization:** No sensitizing effects known.

(Contd. on page 5)



Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 16.10.2012

Version number 1

Revision: 16.10.2012

Trade name: DETARTRINE 150 Z

(Contd. of page 4)

Additional toxicological information:

The product is not subject to classification according to the calculation method of the General EU Classification Guidelines for Preparations as issued in the latest version.
When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us.

SECTION 12: Ecological information

12.1 Toxicity

• **Aquatic toxicity:** No further relevant information available.

• **12.2 Persistence and degradability :** No further relevant information available.

• **12.3 Bioaccumulative potential :** No further relevant information available.

• **12.4 Mobility in soil :** No further relevant information available.

Additional ecological information:**General notes:**

Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water.

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

12.5 Results of PBT and vPvB assessment

• **PBT:** Not applicable.

• **vPvB:** Not applicable.

• **12.6 Other adverse effects :** No further relevant information available.

SECTION 13: Disposal considerations

13.1 Waste treatment methods

• **Recommendation :** Smaller quantities can be disposed of with household waste.

Uncleaned packaging:

• **Recommendation:** Disposal must be made according to official regulations.

SECTION 14: Transport information

14.1 UN-Number

• **ADR, ADN, IMDG, IATA** Void

14.2 UN proper shipping name

• **ADR, ADN, IMDG, IATA** Void

14.3 Transport hazard class(es)

• **ADR, ADN, IMDG, IATA**

• **Class** Void

14.4 Packing group

• **ADR, IMDG, IATA** Void

14.5 Environmental hazards:

• **Marine pollutant:** No

14.6 Special precautions for user

Not applicable.

(Contd. on page 6)

GB



Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 16.10.2012

Version number 1

Revision: 16.10.2012

Trade name: DETARTRINE 150 Z

(Contd. of page 5)

- | | |
|---|-----------------|
| · 14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code | Not applicable. |
| · UN "Model Regulation": | - |

SECTION 15: Regulatory information

- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Relevant phrases**

- H225 Highly flammable liquid and vapour.
H302 Harmful if swallowed.
H314 Causes severe skin burns and eye damage.
H315 Causes skin irritation.
H319 Causes serious eye irritation.
H332 Harmful if inhaled.
H335 May cause respiratory irritation.
H411 Toxic to aquatic life with long lasting effects.

R11 Highly flammable.

R20 Harmful by inhalation.

R22 Harmful if swallowed.

R34 Causes burns.

R36/37/38 Irritating to eyes, respiratory system and skin.

R51/53 Toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

VOC: Volatile Organic Compounds (USA, EU)