


Directa Crowns Molars Mould Guide


Translucent/Translucide/Translucida/Translucenti

DIRECTA Directa Temporary Crown


Upper Right/Höger ÖK/Haut droit/OK-Rechts/Superior Derecha/Superiore Destro




Upper Left/Vänster ÖK/Haut gauche/OK-Links/Superior Izquierda/Superiore Sinistro



Lower Right/Höger UK/Bas droite/UK-Rechts/Inferior Derecha/Inferiore Destro



Lower Left/Vänster UK/Bas gauche/UK-Links/Inferior Izquierda/Inferio Sinistro



Opaque/Opak/Opaca/Opache
Translucent/Translucide/Translucida/Translucenti

Centrals/Centraler/Centrales/OK-Front/Centrales/Centrali

mm	↑	14.0	13.0	12.4	12.0	11.5	11.0	10.6	10.6	11.0	11.5	12.0	12.4	13.0	14.0
mm	↔	10.0	9.2	8.7	8.4	8.2	7.8	7.5	7.5	7.8	8.2	8.4	8.7	9.2	10.0
		100	101	10	11	12	13	14	15	16	17	18	19	102	103
		centrals													

Laterals/Lateraler/Latérales/OK-Front/Laterales/Laterali

Teeth	2	20	21	22	23	24	25	26	27	28	29	200
Height (mm)	11.0	10.4	10.1	9.5	8.8	8.2	8.2	8.8	9.5	10.1	10.4	11.0
Width (mm)	7.3	6.8	6.5	6.2	5.9	5.5	5.5	5.9	6.2	6.5	6.8	7.3

Cuspids/Hörntänder/Canines/Eckzähne/Caninos/Cuspidi

mm ↑	13.0	12.4	11.6	11.4	11.1	10.8	10.5	10.5	10.8	11.1	11.4	11.6	12.4	13.0
mm ↔	9.0	8.5	8.1	7.9	7.7	7.5	7.3	7.3	7.5	7.7	7.9	8.1	8.5	9.0
	300	301	30	31	32	33	34	35	36	37	38	39	302	303
	cuspid													

**Upper bicuspids/ÖK-premolarer/Prémolaires supérieures/
OK Prämolaren/Premolares superiores/Bicuspidi superiori**

mm	↕	9.8	9.4	8.9	8.6	8.4		7.5	8.1	8.5	9.0	9.5
mm	↔	7.4	6.9	6.7	6.3	6.0		6.0	6.3	6.7	7.2	7.4
												
		40	41	42	43	44		45	46	47	48	49
<i>bicusplids upper</i>												

**Lower bicuspid/UK-premolarer/Prémolaires inférieures/
UK Prämolaren/Premolares inferiores/Bicuspidi inferiori**

mm	9.8	9.5	9.2	8.8	8.2	7.1	7.3	8.0	8.1	8.4
mm	7.4	7.2	6.8	6.6	6.3	6.4	6.8	7.1	7.5	8.0

50 51 52 53 54 55 56 57 58 59

bicuspids lower

**Lower incisors/UK-incisiver/Incisives inférieures/
UK Front/Incisivos inferiores/Incisivi inferiori**

mm	12.0	11.5	11.0	10.5	10.0	8.0	8.5	9.1	9.6	10.3
mm	6.2	5.6	5.3	5.1	4.9	4.5	4.9	5.2	5.5	5.8

Incisors lower