



Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 2 / GB

Date revised: 09.10.2025

Replaces Version: 1 / GB

Print date: 03.11.2025

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Drufosoft / Drufosoft pro

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/preparation

Dental splinting material long-term

1.3. Details of the supplier of the safety data sheet

Address/Manufacturer

Dreve Dentamid GmbH

Max-Planck-Straße 31

DE-59423 Unna

Telephone no. +49 2303 8807-0

Fax no. +49 2303 8807-29

Information provided by / telephone Department Research & Development: Fax: +49 2303 8807-562

E-mail address of person responsible sicherheitsdatenblatt@dreve.com

for this SDS

1.4. Emergency telephone number

Henkel Fire Department / 24h-Emergency-Contact-No.: +49 211 797-3350

SECTION 2: Hazards identification ***

2.1. Classification of the substance or mixture

This product is not classified hazardous in accordance with Regulation (EC) No 1272/2008.

2.2. Label elements

Labelling according to regulation (EC) No 1272/2008

Supplemental information

EUH210

Safety data sheet available on request.

2.3. Other hazards

No special hazards have to be mentioned.

The Substance does not meet PBT-criteria. This substance does not meet the vPvB-criteria. This substance does not have endocrine disrupting properties with respect to humans. This substance does not have endocrine disrupting properties with respect to non-target organisms.

SECTION 3: Composition/information on ingredients ***

3.2. Mixtures



Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 2 / GB

Date revised: 09.10.2025

Replaces Version: 1 / GB

Print date: 03.11.2025

Chemical characterization

Ethylene vinyl acetate copolymer (EVA)

Hazardous ingredients *****vinyl acetate**

| | | | | | |
|--|--------------|-----|---|------|---|
| CAS No. | 108-05-4 | | | | |
| EINECS no. | 203-545-4 | | | | |
| Concentration | >= | 0,1 | < | 1,0 | % |
| Classification (Regulation (EC) No. 1272/2008) | | | | | |
| | Flam. Liq. 2 | | | H225 | |
| | Acute Tox. 4 | | | H332 | |
| | Carc. 2 | | | H351 | |
| | STOT SE 3 | | | H335 | |

Additional remarks:

CLP

Regulation (EC) No 1272/2008, Annex VI, Note D

SECTION 4: First aid measures**4.1. Description of first aid measures****General information**

No special measures required

After inhalation

Ensure supply of fresh air.

After skin contact

After contact with molten product cool quickly with cold water. Consult a doctor if skin irritation persists.

After eye contact

Separate eyelids, wash the eyes thoroughly with water (15 min.). In case of irritation consult an oculist.

After ingestion

Rinse out mouth and give plenty of water to drink.

Adhere to personal protective measures when giving first aid

First aider: Pay attention to self-protection!

4.2. Most important symptoms and effects, both acute and delayed

Until now no symptoms known so far.

4.3. Indication of any immediate medical attention and special treatment needed**Hints for the physician / hazards**

In the case of swallowing with subsequent vomiting, aspiration of the lungs can occur which can lead to chemical pneumonia or asphyxiation.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**Recommended: alcohol resistant foam, CO₂, powders, water spray/mist, Extinguishing measures to suit surroundings

Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 2 / GB

Date revised: 09.10.2025

Replaces Version: 1 / GB

Print date: 03.11.2025

5.3. Advice for firefighters

Special protective equipment for fire-fighting

In case of combustion use a suitable breathing apparatus.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Avoid contact with skin, eyes and clothing.

6.2. Environmental precautions

Do not allow to enter drains or waterways.

6.3. Methods and material for containment and cleaning up

Pick up mechanically. Clean contaminated floors and objects thoroughly, observing environmental regulations.

6.4. Reference to other sections

Refer to protective measures listed in Sections 7 and 8.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on safe handling

Observe the usual precautions for handling chemicals.

Advice on protection against fire and explosion

No special measures required.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

Store product in closed containers.

Hints on storage assembly

Do not store together with foodstuffs.

Further information on storage conditions

Keep container tightly closed and dry.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Other information

Contains no substances with occupational exposure limit values.

8.2. Exposure controls

General protective and hygiene measures

Observe the usual precautions for handling chemicals.

Respiratory protection

Do not inhale vapours



Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 2 / GB

Date revised: 09.10.2025

Replaces Version: 1 / GB

Print date: 03.11.2025

Hand protection

Heat-resistant gloves when using hot material

Eye protection

Safety glasses

Body protection

Clothing as usual in the chemical industry.

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

| | | | |
|---|----------------------------------|----|-------------------|
| Physical state | Thermoforming blank | | |
| Colour | Various, depending on coloration | | |
| Odour | characteristic | | |
| Melting point | | | |
| Value | 60 | to | 90 °C |
| Freezing point | | | |
| Remarks | not determined | | |
| Boiling point or initial boiling point and boiling range | | | |
| Remarks | not determined | | |
| Flammability | | | |
| evaluation | not determined | | |
| Upper and lower explosive limits | | | |
| Remarks | not determined | | |
| Flash point | | | |
| Remarks | Not applicable | | |
| Auto-ignition temperature | | | |
| Remarks | not determined | | |
| Decomposition temperature | | | |
| Remarks | not determined | | |
| pH value | | | |
| Remarks | not determined | | |
| Viscosity | | | |
| Remarks | not determined | | |
| Solubility(ies) | | | |
| Remarks | not determined | | |
| Partition coefficient n-octanol/water (log value) | | | |
| Remarks | not determined | | |
| Vapour pressure | | | |
| Remarks | not determined | | |
| Density and/or relative density | | | |
| Value | 0,95 | | g/cm ³ |
| Temperature | 20 | °C | |
| Relative vapour density | | | |
| Remarks | not determined | | |



Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 2 / GB

Date revised: 09.10.2025

Replaces Version: 1 / GB

Print date: 03.11.2025

9.2. Other information**Odour threshold**

Remarks not determined

Evaporation rate (ether = 1) :

Remarks not determined

Solubility in water

Remarks virtually insoluble

Explosive properties

evaluation not determined

Oxidising properties

Remarks not determined

Other information

None known

SECTION 10: Stability and reactivity**10.1. Reactivity**

No hazardous reactions when stored and handled according to prescribed instructions.

10.2. Chemical stability

No hazardous reactions known.

10.3. Possibility of hazardous reactions

No hazardous reactions known.

10.4. Conditions to avoid

Protect from heat and direct sunlight

10.5. Incompatible materials

None known

SECTION 11: Toxicological information**11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008****Acute oral toxicity**

Remarks Based on available data, the classification criteria are not met.

Acute oral toxicity (Components)**vinyl acetate**Species rat (male)
LD50 3470 mg/kg**Acute dermal toxicity**

Remarks Based on available data, the classification criteria are not met.

Acute dermal toxicity (Components)**vinyl acetate**Species rabbit
LD50 7440 mg/kg



Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 2 / GB

Date revised: 09.10.2025

Replaces Version: 1 / GB

Print date: 03.11.2025

Acute inhalational toxicity

Remarks Based on available data, the classification criteria are not met.

Acute inhalative toxicity (Components)**vinyl acetate**

| | | | |
|----------------------|--------|---|--------|
| Species | rat | | |
| LC50 | 4000 | | ppm(V) |
| Duration of exposure | 4 | h | |
| Administration/Form | Vapors | | |

Skin corrosion/irritation

Remarks Based on available data, the classification criteria are not met.

Serious eye damage/irritation

Remarks Based on available data, the classification criteria are not met.

Sensitization

Remarks Based on available data, the classification criteria are not met.

Subacute, subchronic, chronic toxicity

Remarks not determined

Mutagenicity

Remarks Based on available data, the classification criteria are not met.

Reproductive toxicity

Remarks Based on available data, the classification criteria are not met.

Carcinogenicity

Remarks Based on available data, the classification criteria are not met.

Carcinogenicity (Components)

| | |
|----------------------|------------------------------|
| vinyl acetate | |
| evaluation | Suspected of causing cancer. |

Specific Target Organ Toxicity (STOT)**Single exposure**

Remarks Based on available data, the classification criteria are not met.

Repeated exposure

Remarks Based on available data, the classification criteria are not met.

Specific Target Organ Toxicity (STOT) (Components)**vinyl acetate**

| | |
|------------------------|-----------------------------------|
| Single exposure | |
| evaluation | May cause respiratory irritation. |

Aspiration hazard

Based on available data, the classification criteria are not met.

11.2. Information on other hazards**Endocrine disrupting properties with respect to humans**

This substance does not have endocrine disrupting properties with respect to humans.

Other information

No toxicological data are available.

SECTION 12: Ecological information

Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 2 / GB

Date revised: 09.10.2025

Replaces Version: 1 / GB

Print date: 03.11.2025

12.1. Toxicity

General information

not determined

Fish toxicity (Components)

vinyl acetate

| | | | |
|----------------------|---|---|------|
| Species | Fathead minnow (<i>Pimephales promelas</i>) | | |
| NOEC | 0,551 | | mg/l |
| Duration of exposure | 34 | d | |
| Method | OECD 210 | | |

Daphnia toxicity (Components)

vinyl acetate

| | | | |
|----------------------|---------------|---|------|
| Species | Daphnia magna | | |
| EC50 | 12,6 | | mg/l |
| Duration of exposure | 48 | h | |
| Method | OECD 202 | | |

Algae toxicity (Components)

vinyl acetate

| | | | |
|----------------------|---------------------------------|---|------|
| Species | Pseudokirchneriella subcapitata | | |
| ErC50 | 12,7 | | mg/l |
| Duration of exposure | 72 | h | |
| Method | OECD 201 | | |

12.2. Persistence and degradability

General information

not determined

Biodegradability (Components)

vinyl acetate

| | | | | |
|--|----|----|----|---|
| Value | 82 | to | 98 | % |
| Duration of test evaluation | 14 | d | | |
| Readily biodegradable (according to OECD criteria) | | | | |

12.3. Bioaccumulative potential

General information

not determined

Partition coefficient n-octanol/water (log value)

Remarks not determined

Octanol/water partition coefficient (log Pow) (Components)

vinyl acetate

| | | | |
|-------------|------|----|--|
| log Pow | 0,73 | | |
| Temperature | 25 | °C | |

12.4. Mobility in soil

General information

not determined

12.5. Results of PBT and vPvB assessment

General information

not determined

Results of PBT and vPvB assessment



Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 2 / GB

Date revised: 09.10.2025

Replaces Version: 1 / GB

Print date: 03.11.2025

The Substance does not meet PBT-criteria.
This substance does not meet the vPvB-criteria.

12.6 Endocrine disrupting properties**Endocrine disrupting properties with respect to the environment**

This substance does not have endocrine disrupting properties with respect to non-target organisms.

12.7. Other adverse effects**General information**

not determined

General information / ecology

Do not discharge product unmonitored into the environment.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Disposal recommendations for the product**

Must not be disposed together with household garbage.
Dispose of waste according to applicable legislation.

Disposal recommendations for packaging

Packaging that cannot be cleaned should be disposed off in agreement with the regional waste disposal company.

SECTION 14: Transport information

| | Land transport ADR/RID | Marine transport IMDG/GGVSee | Air transport ICAO/IATA |
|----------------------------------|--|---|---|
| 14.1. UN number or ID number | The product does not constitute a hazardous substance in land transport. | The product does not constitute a hazardous substance in sea transport. | The product does not constitute a hazardous substance in air transport. |
| 14.2. UN proper shipping name | - | - | - |
| 14.3. Transport hazard class(es) | | - | - |
| Label | | | |
| 14.4. Packing group | | - | - |

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****Other information**

Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 2 / GB

Date revised: 09.10.2025

Replaces Version: 1 / GB

Print date: 03.11.2025

All components are contained in the TSCA inventory or exempted.

15.2. Chemical safety assessment

For this preparation a chemical safety assessment has not been carried out.

SECTION 16: Other information

Hazard statements listed in Chapter 2/3

| | |
|------|-------------------------------------|
| H225 | Highly flammable liquid and vapour. |
| H332 | Harmful if inhaled. |
| H335 | May cause respiratory irritation. |
| H351 | Suspected of causing cancer. |

CLP categories listed in Chapter 2/3

| | |
|--------------|--|
| Acute Tox. 4 | Acute toxicity, Category 4 |
| Carc. 2 | Carcinogenicity, Category 2 |
| Flam. Liq. 2 | Flammable liquid, Category 2 |
| STOT SE 3 | Specific target organ toxicity - single exposure, Category 3 |

Supplemental information

Relevant changes compared with the previous version of the safety data sheet are marked with: ***
This information is based on our present state of knowledge. However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.