

**Renfert**

making work easy




EASY blank wax

**NEW** CAD/CAM wax blank

# EASY blank wax

## CAD/CAM wax blank

CAD/CAM wax blank for the production of digitally designed crowns and bridges for press and casting techniques.



Milled result: Chip formation  
Renfert EASYblank wax

Milled result: Chip formation  
competitor product

### Technical Data

|                    | No. 7700014                              | No. 7700020                              |
|--------------------|--|--|
| Dimensions (Ø x H) | Ø 98.5 mm x 14 mm<br>(with 10 mm collar) | Ø 98.5 mm x 20 mm<br>(with 10 mm collar) |
| Dropping point     | 100–120°C                                | 100–120°C                                |
| Weight             | approx. 99 g                             | approx. 138 g                            |
| Color              | blue                                     | blue                                     |

### Ordering information

|   |             |
|---|-------------|
| 1 EASY blank wax, Ø 98.5 x 14 mm, 1 piece | No. 7700014 |
| 2 EASY blank wax, Ø 98.5 x 20 mm, 1 piece | No. 7700020 |

### Advantages

- Precisely milled results with extremely smooth surfaces and strong edge stability due to the high degree of material homogeneity.
- Especially suitable for full ceramic restorations as the wax will burn-out without residue.
- Reliable milling process due to small milled chips which do not stick to the cutter or the object surface.

### Making work easy

The Renfert EASY blank wax is the optimal CAD/CAM wax blank for milled restorations in the customarily Renfert quality "made in Germany".

„Making work easy“ stands for outstanding product characteristics, e.g. high material homogeneity. This allows, amongst other things, a perfect fit and surface quality to be achieved in the final restoration.

