

Safety Data Sheet EXTRUDE Base & Catalyst (Cartridges & tubes) Viscosities: Medium; Extra; Wash; MPV

SECTION 1: IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

1.1. Product identifier

Product name : EXTRUDE Base & Catalyst (Cartridges & tubes) Viscosities: Medium; Extra; Wash; MPV

1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses

Main use category : Professional use Function or use category : Dental materials.

Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier
Kerr Italia S.r.I.
Via Passanti, 332
84018 Scafati (SA) - Italy
T +39-081-850-8311
E-mail: safety@kerrhawe.com

Manufacturer
Kerr Corporation
1717 West Collins Avenue
92867 Orange – CALIFORNIA (U.S.A.)
T 00-800-41-050-505
safety@kerrhawe.com

Contact person: safety@kerrhawe.com - tel. 00-800-41-050-505 (08.00-17.00)

1.4. Emergency telephone number

Emergency number : CHEMTREC® Emergency Call Center. Emergency Telephone Number (for USA only) 001-

800-424-9300 International and Maritime Telephone Number +1 (703) 527-3887

Country	Organisation/Company	Address	Emergency number
United Kingdom	National Poisons Information Service (Newcastle Unit)	Claremont Place Newcastle-upon-Tyne, Newcastle	+44 191 2606182/+44 191 2606180 24H

SECTION 2: HAZARDS IDENTIFICATION

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Extra phrases : The product is seen as a medical device and therefore not subject to labelling (EU-

regulation 1907/2006, article 2, paragraph 6c).

2.3. Other hazards

Other hazards not contributing to the : None under normal conditions.

classification

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

3.1. Substance

Not applicable

3.2. Mixture



Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Wollastonite	(CAS No) 13983-17-0 (EC no) 237-772-5 (REACH-no) N/A		Not classified
silica, vitreous	(CAS No) 60676-86-0 (EC no) 262-373-8 (REACH-no) N/A		Not classified

Full text of H-statements: see section 16

SECTION 4: FIRST AID MEASURES

4.1. Description of first aid measures

First-aid measures general : If medical advice is needed, have product container or label at hand.

First-aid measures after inhalation : Remove victim to fresh air and keep at rest in a position comfortable for breathing. Get

medical advice/attention if you feel unwell.

First-aid measures after skin contact : Gently wash with plenty of soap and water. If skin irritation occurs: Get medical

advice/attention.

First-aid measures after eye contact : IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if

present and easy to do. Continue rinsing. Get medical advice/attention if you feel unwell.

: If swallowed, rinse mouth with water (only if the person is conscious). Call a POISON

CENTER or doctor/physician if you feel unwell.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/injuries : Not expected to present a significant hazard under anticipated conditions of normal use.

4.3. Indication of any immediate medical attention and special treatment needed

No specific measures identified.

First-aid measures after ingestion

SECTION 5: FIREFIGHTING MEASURES

5.1. Extinguishing media

Suitable extinguishing media : Use extinguishing media appropriate for surrounding fire. Foam, carbon dioxide (CO2) and

powder.

Unsuitable extinguishing media : Do not use a heavy water stream.

5.2. Special hazards arising from the substance or mixture

Fire hazard : Non flammable. Hazardous polymerization may occur if exposed to high temperature.

Explosion hazard : Product is not explosive.

Hazardous decomposition products in case of

fire

Carbon dioxide. Carbon monoxide.

5.3. Advice for firefighters

Firefighting instructions : Use water spray or fog for cooling exposed containers.

Protection during firefighting : Do not enter fire area without proper protective equipment, including respiratory protection.

SECTION 6: ACCIDENTAL RELEASE MEASURES

6.1. Personal precautions, protective equipment and emergency procedures

General measures : Avoid contact with skin and eyes.

For non-emergency personnel

Protective equipment : See Heading 8.

Emergency procedures : Evacuate unnecessary personnel.

For emergency responders

No additional information available

6.2. Environmental precautions

Avoid release to the environment.

6.3. Methods and material for containment and cleaning up

For containment : Collect all waste in suitable and labelled containers and dispose according to local

legislation.

Methods for cleaning up : Small quantities of liquid spill: take up in non-combustible absorbent material and shovel

into container for disposal. Large spills: scoop solid spill into closing containers.

6.4. Reference to other sections

For further information refer to section 13.



SECTION 7: HANDLING AND STORAGE

7.1. Precautions for safe handling

Precautions for safe handling : Do not eat, drink or smoke when using this product. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Avoid contact with skin and

eyes.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No

smoking. Keep container tightly closed. Store in original container. Store in dry protected location to prevent any moisture contact. Store at room temperature (> 11°C) to prevent

solidification

Incompatible materials : Oxidizing substances. Organic peroxides. Amines.

Storage area : Store in a well-ventilated place.

7.3. Specific end use(s)

Consult the supplier for further information.

SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1. Control parameters

silica, vitreous (60676-86-0)		
United Kingdom	Local name	Silica, fused
United Kingdom	WEL TWA (mg/m³)	0.08 mg/m³ respirable dust

8.2. Exposure controls

Appropriate engineering controls : Ensure good ventilation of the work station.

Personal protective equipment : Gloves

Hand protection : Wear suitable gloves. Nitrile rubber gloves. Layer thickness: 0,09mm. Breakthrough time:

>480 min. STANDARD EN 374.

Eye protection : Use splash goggles when eye contact due to splashing is possible. STANDARD EN 166.

Skin and body protection : Wear suitable protective clothing

Respiratory protection : In case of insufficient ventilation, wear suitable respiratory equipment. Extra personal

protection: P2 filter respirator for harmful particles. Standard EN 143.



Other information : Do not eat, drink or smoke during use. Wash hands and other exposed areas with mild

soap and water before eating, drinking or smoking and when leaving work.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

9.1. Information on basic physical and chemical properties

Physical state : Liquid
Appearance : Paste.
Colour : Various.
Odour : Fruity.

Odour threshold : No data available pH : No data available Relative evaporation rate (butylacetate=1) : No data available Melting point : No data available Freezing point : No data available Boiling point : No data available

Flash point : 142 °C

Auto-ignition temperature : No data available Decomposition temperature : No data available Flammability (solid, gas) : No data available Vapour pressure : No data available Relative vapour density at 20 °C : No data available Relative density : No data available Density : ≈ 1.25 g/cm³

Solubility : Material insoluble in water.



300916 EXTRUDE Base & Catalyst (Cartridges & tubes) Viscosities: Medium; Extra; Wash; MPV

Log Pow : No data available
Viscosity, kinematic : No data available
Viscosity, dynamic : No data available
Explosive properties : Product is not explosive.

Oxidising properties : Non flammable.

Explosive limits : No data available

9.2. Other information

No additional information available

SECTION 10: STABILITY AND REACTIVITY

10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No polymerization.

10.4. Conditions to avoid

No flames, no sparks. Eliminate all sources of ignition.

10.5. Incompatible materials

Oxidizing agent.

10.6. Hazardous decomposition products

No decomposition if stored normally.

SECTION 11: TOXICOLOGICAL INFORMATION

11.1. Information on toxicological effects

Acute toxicity : Not classified Skin corrosion/irritation : Not classified

Contact during a long period may cause slight irritation

Serious eye damage/irritation : Not classified

Liquid splashes in the eye may cause irritation

Respiratory or skin sensitisation : Not classified
Germ cell mutagenicity : Not classified
Carcinogenicity : Not classified
Reproductive toxicity : Not classified
Specific target organ toxicity (single exposure) : Not classified

By prolonged exposure: May cause minor irritation to the respiratory tract and to other

mucous membranes

May cause minor irritation to the respiratory tract and to other mucous membranes

Specific target organ toxicity (repeated

exposure)

: Not classified

Aspiration hazard : Not classified

SECTION 12: ECOLOGICAL INFORMATION

12.1. Toxicity

Ecology - general : The product is not considered harmful to aquatic organisms or to cause long-term adverse

effects in the environment.

12.2. Persistence and degradability

EXTRUDE Base & Catalyst (Cartridges & tubes) Viscosities: Medium; Extra; Wash; MPV		
Persistence and degradability	No data available.	

12.3. Bioaccumulative potential

EXTRUDE Base & Catalyst (Cartridges & tubes) Viscosities: Medium; Extra; Wash; MPV		
Bioaccumulative potential	No data.	

12.4. Mobility in soil

EXTRUDE Base & Catalyst (Cartridges & tubes) Viscosities: Medium; Extra; Wash; MPV		
Ecology - soil	Material insoluble in water.	

06/11/2015

12.5. Results of PBT and vPvB assessment

EXTRUDE Base & Catalyst (Cartridges & tubes) Viscosities: Medium; Extra; Wash; MPV

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII

This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

12.6. Other adverse effects

Other adverse effects : None to our knowledge.

Additional information : No other effects known

SECTION 13: DISPOSAL CONSIDERATIONS

13.1. Waste treatment methods

Regional legislation (waste) : Product is not hazardous waste.

Waste treatment methods : Recover the product with absorbent material. Dispose of contents/container in accordance

with licensed collector's sorting instructions.

Waste disposal recommendations : Dispose in a safe manner in accordance with local/national regulations.

Ecology - waste materials : Avoid release to the environment.

European List of Waste (LoW) code : 18 01 07 - chemicals other than those mentioned in 18 01 06

SECTION 14: TRANSPORT INFORMATION

In accordance with ADR / RID / IMDG / IATA / ADN

14.1. UN number

Not regulated for transport

14.2. UN proper shipping name

14.3. Transport hazard class(es)

14.4. Packing group

14.5. Environmental hazards

Dangerous for the environment : No

Other information : No supplementary information available

14.6. Special precautions for user

14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

SECTION 15: REGULATORY INFORMATION

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

EU-Regulations

Contains no REACH substances with Annex XVII restrictions

Contains no substance on the REACH candidate list

Contains no REACH Annex XIV substances

National regulations

EC-regulation 453/2010/EC, 1907/2006/EC (REACH), 1272/2008/EC (CLP), 790/2009/EC. Transport of dangerous goods (ADR/RID, IMDG, IATA/ICAO). Workplace exposure limits.

15.2. Chemical safety assessment

No chemical safety assessment has been carried out for the substance or the mixture by the supplier

SECTION 16: OTHER INFORMATION

Indication of changes:

Regulatory information.

		_	
2.1	Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD]	Removed	
3.2	Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD]	Removed	

Data sources : EC-regulation 453/2010/EC, 1907/2006/EC (REACH), 1272/2008/EC (CLP), 790/2009/EC.

Transport of dangerous goods (ADR/RID, IMDG, IATA/ICAO). Workplace exposure limits.



300916 EXTRUDE Base & Catalyst (Cartridges & tubes) Viscosities: Medium; Extra; Wash; MPV 06/11/2015

 Date of issue
 : 15/04/2013

 Revision date
 : 06/11/2015

 Supersedes
 : 20/05/2014

Version : 3.0

Signature : A. Åsebø Murel

The information in this safety data sheet is based on information from the manufacturer/supplier, present European and national legislation, and presupposes that the product is used within the specified area of application.