

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 19 Jun 2026

Print date: 19 Jun 2026

Version: 1.1



Page 1/11

FOKALDRY

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Trade name/designation:

FOKALDRY

Other means of identification:

Article No. : 0032305 (80 ml), 0032306 (500 ml)

Additional information:

Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH). As a medical device this product is not subject to Regulation (EC) no. 1272/2008 [CLP] and do not require labeling according this regulation.

1.2 Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/mixture:

Medical device for dental application. Observe instructions for use.
for professional use only.

1.3 Details of the supplier of the safety data sheet

Manufacturer/Supplier:

lege artis Pharma GmbH + Co. KG

Breitwasenring 1

72135 Dettenhausen

GERMANY

Telephone: +49 (0) 71 57 / 56 45 - 0

Telefax: +49 (0) 71 57 / 56 45 - 50

E-mail: info@legeartis.de

Website: www.legeartis.de

E-mail (competent person): sicherheitsdatenblaetter@legeartis.de

1.4 Emergency phone number

24h: +49 (0) 551 / 192 40, +49 (0) 71 57 / 56 45 - 0 (Only available during office hours.)

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Hazard classes and hazard categories	Hazard statements	Classification procedure
Serious eye damage/eye irritation (<i>Eye Irrit. 2</i>)	H319: Causes serious eye irritation.	Calculation
Flammable liquids (<i>Flam. Liq. 2</i>)	H225: Highly flammable liquid and vapour.	Test data
Aspiration hazard (<i>Asp. Tox. 1</i>)	H304: May be fatal if swallowed and enters airways.	Test data
Skin corrosion/irritation (<i>Skin Irrit. 2</i>)	H315: Causes skin irritation.	Calculation
STOT-single exposure (<i>STOT SE 3</i>)	H336: May cause drowsiness or dizziness.	Calculation
Hazardous to the aquatic environment (<i>Aquatic Chronic 3</i>)	H412: Harmful to aquatic life with long lasting effects.	Calculation

Additional information:

Contains: Isopropyl ether, Solvent (C6-C7), Acetone, isopropanol

May form explosive peroxides. Repeated exposure may cause skin dryness or cracking.

Precautionary statements:

Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Avoid breathing dust/fume/gas/mist/vapours/spray. Wear protective gloves/protective clothing/eye protection/face protection. IF SWALLOWED: Immediately call a POISON CENTER. IF ON SKIN: Gently wash with plenty

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 19 Jun 2026

Print date: 19 Jun 2026

Version: 1.1



Page 2/11

FOKALDRY

of soap and water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

2.2 Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

According to EC directives or the corresponding national regulations the product does not have to be labelled.

Additional information:

As a medical device this product is not subject to Regulation (EC) no. 1272/2008 [CLP] and do not require labeling according this regulation.

2.3 Other hazards

Adverse human health effects and symptoms:

Special danger of slipping by leaking/spilling product.

SECTION 3: Composition/information on ingredients

3.2 Mixtures

Description:

A mixture of the substances listed below, containing non-hazardous additives.

Ingredients:

Product identifiers	Substance name Classification according to Regulation (EC) No. 1272/2008 [CLP]	Content
CAS No.: 108-20-3 EC No.: 203-560-6 Index No.: 603-045-00-X	diisopropyl ether Flam. Liq. 2 (H225), STOT SE 3 (H336) Danger EUH019-EUH066	40 – 60 weight-%
CAS No.: 67-63-0 EC No.: 200-661-7 Index No.: 603-117-00-0 REACH No.: 01-2119457558-25	propan-2-ol Eye Irrit. 2 (H319), Flam. Liq. 2 (H225), STOT SE 3 (H336) Danger	10 – 20 weight-%
CAS No.: 64742-49-0 EC No.: 265-151-9 Index No.: 649-328-00-1 REACH No.: 01-2119475133-43	Naphtha (petroleum), hydrotreated light Aquatic Chronic 2 (H411), Asp. Tox. 1 (H304), Flam. Liq. 2 (H225), STOT SE 3 (H336), Skin Irrit. 2 (H315) Danger	10 – 20 weight-%
CAS No.: 67-64-1 EC No.: 200-662-2 Index No.: 606-001-00-8 REACH No.: 01-2119471330-49	acetone Eye Irrit. 2 (H319), Flam. Liq. 2 (H225), STOT SE 3 (H336) Danger EUH066	10 – 20 weight-%
CAS No.: 109-66-0 EC No.: 203-692-4 Index No.: 601-006-00-1 REACH No.: 01-2119459286-30	pentane Aquatic Chronic 2 (H411), Asp. Tox. 1 (H304), Flam. Liq. 2 (H225), STOT SE 3 (H336) Danger EUH066	5 – 15 weight-%

Full text of H- and EUH-phrases: see section 16.

SECTION 4: First aid measures

4.1 Description of first aid measures

General information:

In case of accident or unwellness, seek medical advice immediately (show directions for use or safety data sheet if possible). Take off immediately all contaminated clothing. If unconscious but breathing normally, place in recovery position and seek medical advice. Do not leave affected person unattended.

Following inhalation:

Provide fresh air. Get medical advice/attention if you feel unwell.

In case of skin contact:

Gently wash with plenty of soap and water. In case of skin irritation, consult a physician.

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 19 Jun 2026

Print date: 19 Jun 2026

Version: 1.1



Page 3/11

FOKALDRY

After eye contact:

Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

Following ingestion:

Rinse mouth. Do NOT induce vomiting. Observe risk of aspiration if vomiting occurs. Get immediate medical advice/attention.

4.2 Most important symptoms and effects, both acute and delayed

Causes skin irritation. Causes serious eye irritation. May cause drowsiness or dizziness. May be fatal if swallowed and enters airways.

4.3 Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1 Extinguishing media

Suitable extinguishing media:

Water spray jet, alcohol resistant foam, Extinguishing powder, Carbon dioxide (CO₂).

5.2 Special hazards arising from the substance or mixture

Highly flammable liquid and vapour. The vapour of the product is heavier than air and may accumulate below ground level, in pits, channels and basements in higher concentration. Beware of reignition.

Hazardous combustion products:

In case of fire may be liberated: Gases/vapours, toxic

5.3 Advice for firefighters

Wear a self-contained breathing apparatus and chemical protective clothing. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

6.1.1 For non-emergency personnel

Personal precautions:

Do not breathe gas/vapour/aerosol. Avoid contact with eyes and skin. Special danger of slipping by leaking/spilling product. Remove persons to safety.

Protective equipment:

Wear protective gloves/protective clothing/eye protection/face protection.

6.1.2 For emergency responders

Personal protection equipment:

Personal protection equipment: see section 8

6.2 Environmental precautions

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

6.3 Methods and material for containment and cleaning up

For containment:

Wipe up with absorbent material (eg. cloth, fleece). Handling larger quantities: Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Collect in closed and suitable containers for disposal.

For cleaning up:

Water (with cleaning agent)

6.4 Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 19 Jun 2026

Print date: 19 Jun 2026

Version: 1.1



Page 4/11

FOKALDRY

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Protective measures

Advices on safe handling:

Do not breathe gas/vapour/aerosol. Avoid contact with eyes and skin. Take off immediately all contaminated clothing. Wear personal protection equipment (refer to section 8).

Fire prevent measures:

Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep cool. Protect from sunlight.

Environmental precautions:

Discharge into the environment must be avoided.

Advices on general occupational hygiene

When using do not eat, drink or smoke. Wash hands and face before breaks and after work and take a shower if necessary. Apply skin care products after work.

7.2 Conditions for safe storage, including any incompatibilities

Technical measures and storage conditions:

Keep container tightly closed in a cool, well-ventilated place.

Protect from sunlight.

Recommended storage temperature: 15 - 25 °C

Requirements for storage rooms and vessels:

Floors should be impervious, resistant to liquids and easy to clean. Keep/Store only in original container.

Hints on storage assembly:

Keep away from food, drink and animal feedingstuffs.

7.3 Specific end use(s)

Recommendation:

Medical device for dental application.

Observe instructions for use.

SECTION 8: Exposure controls / Personal protection

8.1 Control Parameters

8.1.1 Occupational exposure limit values

Limit value type (country)	Substance name	① Long-term occupational exposure limit value ② Short-term occupational exposure limit value ③ Instantaneous value ④ Monitoring and observation processes ⑤ Remark
IOELV (EU)	acetone CAS No.: 67-64-1 EC No.: 200-662-2	① 500 ppm (1,210 mg/m ³)
IOELV (EU)	pentane CAS No.: 109-66-0 EC No.: 203-692-4	① 1,000 ppm (3,000 mg/m ³)

8.1.2 Biological limit values

No data available.

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 19 Jun 2026

Print date: 19 Jun 2026

Version: 1.1



Page 5/11

FOKALDRY

8.1.3 DNEL-/PNEC-values

Substance name	DNEL value	① DNEL type ② Exposure route
propan-2-ol CAS No.: 67-63-0 EC No.: 200-661-7	500 mg/m ³	① DNEL worker ② Long-term - inhalation, systemic effects
propan-2-ol CAS No.: 67-63-0 EC No.: 200-661-7	89 mg/m ³	① DNEL consumer ② Long-term - inhalation, systemic effects
propan-2-ol CAS No.: 67-63-0 EC No.: 200-661-7	888 mg/kg bw/day	① DNEL worker ② Long-term - dermal, systemic effects
propan-2-ol CAS No.: 67-63-0 EC No.: 200-661-7	319 mg/kg	① DNEL consumer ② Long-term - dermal, systemic effects
propan-2-ol CAS No.: 67-63-0 EC No.: 200-661-7	26 mg/kg	① DNEL consumer ② Long-term - oral, systemic effects
acetone CAS No.: 67-64-1 EC No.: 200-662-2	1,210 mg/m ³	① DNEL worker ② Long-term - inhalation, systemic effects
acetone CAS No.: 67-64-1 EC No.: 200-662-2	200 mg/m ³	① DNEL consumer ② Long-term - inhalation, systemic effects
acetone CAS No.: 67-64-1 EC No.: 200-662-2	2,420 mg/m ³	① DNEL worker ② Acute - inhalation, systemic effects
acetone CAS No.: 67-64-1 EC No.: 200-662-2	186 mg/kg	① DNEL worker ② Long-term - dermal, systemic effects
acetone CAS No.: 67-64-1 EC No.: 200-662-2	62 mg/kg	① DNEL consumer ② Long-term - dermal, systemic effects
acetone CAS No.: 67-64-1 EC No.: 200-662-2	62 mg/kg	① DNEL consumer ② Long-term - oral, systemic effects

Substance name	PNEC value	① PNEC type
propan-2-ol CAS No.: 67-63-0 EC No.: 200-661-7	140.9 mg/L	① PNEC aquatic, freshwater
propan-2-ol CAS No.: 67-63-0 EC No.: 200-661-7	140.9 mg/L	① PNEC aquatic, marine water
propan-2-ol CAS No.: 67-63-0 EC No.: 200-661-7	2,251 mg/L	① PNEC sewage treatment plant
propan-2-ol CAS No.: 67-63-0 EC No.: 200-661-7	552 mg/kg	① PNEC sediment, freshwater
propan-2-ol CAS No.: 67-63-0 EC No.: 200-661-7	160 mg/kg	① PNEC secondary poisoning
propan-2-ol CAS No.: 67-63-0 EC No.: 200-661-7	140.9 mg/L	① PNEC aquatic, intermittent release
propan-2-ol CAS No.: 67-63-0 EC No.: 200-661-7	28 mg/kg	① PNEC soil, freshwater

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 19 Jun 2026

Print date: 19 Jun 2026

Version: 1.1



Page 6/11

FOKALDRY

Substance name	PNEC value	① PNEC type
acetone CAS No.: 67-64-1 EC No.: 200-662-2	10.6 mg/L	① PNEC aquatic, freshwater
acetone CAS No.: 67-64-1 EC No.: 200-662-2	1.06 mg/L	① PNEC aquatic, marine water
acetone CAS No.: 67-64-1 EC No.: 200-662-2	100 mg/L	① PNEC sewage treatment plant
acetone CAS No.: 67-64-1 EC No.: 200-662-2	30.4 mg/kg	① PNEC sediment, freshwater
acetone CAS No.: 67-64-1 EC No.: 200-662-2	3.04 mg/kg	① PNEC sediment, marine water
acetone CAS No.: 67-64-1 EC No.: 200-662-2	21 mg/L	① PNEC aquatic, intermittent release
acetone CAS No.: 67-64-1 EC No.: 200-662-2	33.3 mg/kg	① PNEC soil, freshwater

8.2 Exposure controls

8.2.1 Appropriate engineering controls

Technical measures and the application of suitable work processes have priority over personal protection equipment.

8.2.2 Personal protection equipment



Eye/face protection:

Eye glasses with side protection (EN 166).

Skin protection:

Tested protective gloves must be worn (EN ISO 374).

Suitable material: IIR (Butyl rubber)

Thickness of the glove material: $\geq 0,7$ mm

Breakthrough time: ≥ 480 min

The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. The exact breakthrough time must be requested from the protective glove manufacturer and must be observed. Check leak tightness/impermeability prior to use. In the case of wanting to use the gloves again, clean them before taking off and air them well.

Respiratory protection:

Usually no personal respiratory protection necessary. Provide adequate ventilation.

Wear breathing apparatus if exposed to vapours/dusts/aerosols.

Filtering device (full mask or mouthpiece) with filter: AX

Other protection measures:

Wear suitable protective clothing.

8.2.3 Environmental exposure controls

No data available.

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Physical state: liquid

Form: liquid

Colour: colourless

Odour: like Isopropyl ether

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 19 Jun 2026

Print date: 19 Jun 2026

Version: 1.1



Page 7/11

FOKALDRY

Flammability: Yes

Safety relevant basic data

Parameter	Value	at	① Method ② Remark
pH	<i>Not applicable.</i>		
Melting point	<i>No data available.</i>		
Freezing point	<i>No data available.</i>		
Initial boiling point and boiling range	$\geq 44\text{ °C}$		
Flash point	$\approx -10 - 0\text{ °C}$		② (estimated)
Evaporation rate	<i>No data available.</i>		
Auto-ignition temperature	$\geq 280\text{ °C}$		
Upper/lower Flammability or explosive limits	1 - 21 Vol-%		② (Isopropyl ether)
Vapour pressure	<i>No data available.</i>		
Vapour density	<i>No data available.</i>		
Density	0.7 g/mL	20 °C	
Bulk density	<i>Not applicable.</i>		
Water solubility	practically insoluble		
Dynamic viscosity	<i>No data available.</i>		
Kinematic viscosity	<i>No data available.</i>		

9.2 Other information

No data available.

SECTION 10: Stability and reactivity

10.1 Reactivity

Highly flammable liquid and vapour.

10.2 Chemical stability

The product is chemically stable under recommended conditions of storage, use and temperature.

10.3 Possibility of hazardous reactions

No hazardous reaction when handled and stored according to provisions.

Exothermic reaction with: Amines; Acids; Oxidising agent, strong

10.4 Conditions to avoid

Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep cool. Protect from sunlight.

10.5 Incompatible materials

Amines; Acids; Oxidising agent, strong; plastic

10.6 Hazardous decomposition products

No known hazardous decomposition products. In case of fire may be liberated: Gases/vapours, toxic

SECTION 11: Toxicological information

11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008

propan-2-ol CAS No.: 67-63-0 EC No.: 200-661-7
LD₅₀ oral: 5,840 mg/kg (Rat)
LD₅₀ dermal: 13,900 mg/kg (Rabbit)
LC₅₀ Acute inhalation toxicity (gas): 25 mg/L 4 h (Rat)
pentane CAS No.: 109-66-0 EC No.: 203-692-4
LD₅₀ oral: 400 mg/kg (Rat)
LC₅₀ Acute inhalation toxicity (vapour): 364 mg/L 4 h (Rat)

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 19 Jun 2026

Print date: 19 Jun 2026

Version: 1.1



Page 8/11

FOKALDRY

acetone CAS No.: 67-64-1 EC No.: 200-662-2

LC₅₀ Acute inhalation toxicity (vapour): 76 mg/L 4 h (Rat)

Acute oral toxicity:

Based on available data, the classification criteria are not met.

Acute dermal toxicity:

Based on available data, the classification criteria are not met.

Acute inhalation toxicity:

Based on available data, the classification criteria are not met.

Skin corrosion/irritation:

Causes skin irritation.

Serious eye damage/irritation:

Causes serious eye irritation.

Respiratory or skin sensitisation:

Based on available data, the classification criteria are not met.

Germ cell mutagenicity:

Based on available data, the classification criteria are not met.

Carcinogenicity:

Based on available data, the classification criteria are not met.

Reproductive toxicity:

Based on available data, the classification criteria are not met.

STOT-single exposure:

Vapours may cause drowsiness and dizziness.

STOT-repeated exposure:

Based on available data, the classification criteria are not met.

Aspiration hazard:

May be fatal if swallowed and enters airways.

11.2 Information on other hazards

Endocrine disrupting properties:

This product does not contain a substance that has endocrine disrupting properties with respect to humans as no components meets the criteria.

SECTION 12: Ecological information

12.1 Toxicity

Aquatic toxicity:

Harmful to aquatic life with long lasting effects.

12.2 Persistence and degradability

Biodegradation:

There are no data available on the mixture itself. Some of the components are poorly biodegradable.

12.3 Bioaccumulative potential

Accumulation / Evaluation:

No indication of bioaccumulation potential.

12.4 Mobility in soil

No data available.

12.5 Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII.

12.6 Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

12.7 Other adverse effects

No data available.

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 19 Jun 2026

Print date: 19 Jun 2026

Version: 1.1



Page 9/11

FOKALDRY

SECTION 13: Disposal considerations

13.1 Waste treatment methods

Waste treatment options

Appropriate disposal / Product:

Dispose of waste according to applicable legislation. Consult the appropriate local waste disposal expert about waste disposal.

Appropriate disposal / Package:

Completely emptied packages can be recycled.

Other disposal recommendations:

The allocation of waste code numbers / waste names must be carried out in accordance with the European Waste Catalogue (EWC). Do not allow to enter into surface water or drains. Collect in closed and suitable containers for disposal. Waste for disposal is to be classified and labelled.

SECTION 14: Transport information

Land transport (ADR/RID)	Inland waterway craft (ADN)	Sea transport (IMDG)	Air transport (ICAO-TI-/IATA-DGR)
14.1 UN number or ID number			
UN 1993	UN 1993	UN 1993	UN 1993
14.2 UN proper shipping name			
Flammable liquid, n.o.s. (DIISOPROPYL ETHER, ISOPROPANOL)	Flammable liquid, n.o.s. (DIISOPROPYL ETHER, ISOPROPANOL)	Flammable liquid, n.o.s. (DIISOPROPYL ETHER, ISOPROPANOL)	Flammable liquid, n.o.s. (DIISOPROPYL ETHER, ISOPROPANOL)
14.3 Transport hazard class(es)			
 3	 3	 3	 3
14.4 Packing group			
II	II	II	II
14.5 Environmental hazards			
No	No	No	No
14.6 Special precautions for user			
No data available.	No data available.	No data available.	No data available.

14.7 Maritime transport in bulk according to IMO instruments

Not applicable.

14.8. Additional information

Small amounts: Transport as "Limited Quantity" according to chapter 3.4 ADR/RID

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1 EU legislation

Authorisations:

Contains no substances listed in REACH Annex XIV.

Restrictions on use:

Contains no substances subject to the restrictions of Annex XVII of the REACH Regulation.

Other regulations (EU):

Directive 98/24/EC of 7 April 1998 on the protection of the health and safety of workers from the risks related to chemical agents at work Observe restrictions to employment for juveniles according to the 'juvenile work protection guideline' (94/33/EC). Observe employment restrictions under the Maternity Protection Directive (92/85/EEC) for expectant or nursing mothers. Directive 93/42/EEC concerning medical devices.

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 19 Jun 2026

Print date: 19 Jun 2026

Version: 1.1



Page 10/11

FOKALDRY

15.1.2 National regulations

No data available.

15.2 Chemical Safety Assessment

No data available.

SECTION 16: Other information

16.1 Indication of changes

Changes executed in version 1.1:

General revision of all sections

16.2 Abbreviations and acronyms

ADN = Accord européen relatif au transport international des marchandises dangereuses par voie de navigation intérieure

ADR = Accord européen relatif au transport international des marchandises Dangereuses par Route

ATE = acute toxicity estimate

CAS = Chemical Abstracts Service

CLP = Classification, Labelling and Packaging

DMEL = Derived Minimum Effect Level

DNEL = Derived No Effect Level

EC₅₀ = Median effective concentration

ECB = European Chemicals Bureau

EEC = European Economic Community

EINECS = European Inventory of Existing Commercial Chemical Substances

ELINCS = European List of Notified Chemical Substances

GHS = Globally Harmonized System of Classification and Labelling of Chemicals

IATA = International Air Transport Association

IBC-Code = International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk

IC₅₀ = Inhibition concentration, 50%

IMDG = International Maritime Code for Dangerous Goods

IUCLID = International Uniform Chemical Information Database

LC₅₀ = Lethal concentration, 50%

LD₅₀ = Median lethal dose

LC0 = lethal concentration, 0%

LOAEL = lowest-observed-adverse-effect level

MARPOL = International Convention for the Prevention of Marine Pollution from Ships

NOAEL = No Observed Adverse Effect Level

NOEC = No Observed Effect Content

PBT = Persistent, Bioaccumulative and Toxic substance

PNEC = Predicted No-Effect Content

RID = Règlement concernant le transport international ferroviaire de marchandises dangereuses

REACH = Registration, Evaluation, Authorisation and Restriction of Chemicals

STP = Sewage Treatment Plant

TLV®/TWA = Threshold limit value - time-weighted average

TLV®STEL = Threshold limit value - short-time exposure limit

VOC = Volatile Organic Compounds

vPvB = very Persistent and very Bioaccumulative

16.3 Key literature references and sources for data

European Chemicals Agency (ECHA): <https://www.echa.europa.eu>

ECHA, C&L Inventory: <https://echa.europa.eu/information-on-chemicals/cl-inventory-database>

ECHA, Registered substances: <https://echa.europa.eu/information-on-chemicals/registered-substances>

GESTIS, Information system on hazardous substances: <https://www.gestis.dguv.de/search>

GESTIS, International Limit Values: <https://limitvalue.ifa.dguv.de>

Hörsch Gefährliche Stoffe und Gemische, 8. Auflage, Dr. Angela Schulz

Safety data sheets of the manufacturers

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 19 Jun 2026

Print date: 19 Jun 2026

Version: 1.1



Page 11/11

FOKALDRY

16.4 Classification for mixtures and used evaluation method according to regulation (EC) No. 1272/2008 [CLP]

Hazard classes and hazard categories	Hazard statements	Classification procedure
Serious eye damage/eye irritation (<i>Eye Irrit. 2</i>)	H319: Causes serious eye irritation.	Calculation
Flammable liquids (<i>Flam. Liq. 2</i>)	H225: Highly flammable liquid and vapour.	Test data
Aspiration hazard (<i>Asp. Tox. 1</i>)	H304: May be fatal if swallowed and enters airways.	Test data
Skin corrosion/irritation (<i>Skin Irrit. 2</i>)	H315: Causes skin irritation.	Calculation
STOT-single exposure (<i>STOT SE 3</i>)	H336: May cause drowsiness or dizziness.	Calculation
Hazardous to the aquatic environment (<i>Aquatic Chronic 3</i>)	H412: Harmful to aquatic life with long lasting effects.	Calculation

16.5 List of relevant hazard statements and/or precautionary statements from sections 2 to 15

Hazard statements	
H225	Highly flammable liquid and vapour.
H304	May be fatal if swallowed and enters airways.
H315	Causes skin irritation.
H319	Causes serious eye irritation.
H336	May cause drowsiness or dizziness.
H411	Toxic to aquatic life with long lasting effects.

Supplemental hazard information	
EUH019	May form explosive peroxides.
EUH066	Repeated exposure may cause skin dryness or cracking.

16.6 Training advice

No data available.

16.7 Additional information

The information in this safety data sheet has been established to our best knowledge and was up-to-date at time of revision. The information is intended to give you advice about the safe handling of the product for storage, processing, transport and disposal. It does not represent a guarantee for the properties of the product described in terms of the legal warranty regulations.