


**Safety Data Sheet**Trade Name: **FORENDO™ PASTE**, Calcium Hydroxide with Iodoform

<b>1.0 Commercial Product Name and Supplier</b>																								
1.1	Commercial product name / designation	<b>FORENDO™ PASTE</b>																						
1.2	Application / Use	Dental material for use by dental professionals.																						
1.2.2	SIC	851 Human health activity																						
1.2.3	Use Category	55																						
1.3	Manufacturer	<b>Pulpdent Corporation</b> 80 Oakland Street, P.O. Box 780 Watertown, MA 02472 USA Telephone: 1 617 926-6666 / Fax: 1 617 926-6262 Email: <a href="mailto:Pulpdent@pulpdent.com">Pulpdent@pulpdent.com</a>																						
1.4	Emergency Telephone Number	1-800-535-5053 (24 Hour / USA)																						
1.5	Authorized European Representative	Advena Limited Tower Business Centre, 2nd Floor, Tower Street, Swatar, BKR 4013 Malta																						
	UK Responsible Person	Advena Limited Pure Offices, Plato Close Warwick, CV34 6WE United Kingdom																						
	CH Authorized Representative	MedEnvoy Switzerland Gotthardstrasse 28, 6302 Zug, Switzerland																						
<b>2.0 Hazards Identification</b>																								
2.1	Classification																							
2.1.1	Classification according to Regulation (EC) No. 1272/2008 [CLP]	<table border="1"> <thead> <tr> <th><u>Hazard Class</u></th> <th><u>Hazard Category</u></th> <th><u>Hazard Statement</u></th> </tr> </thead> <tbody> <tr> <td>Skin irritation</td> <td>2</td> <td>H315</td> </tr> <tr> <td>Eye irritation</td> <td>2</td> <td>H319</td> </tr> <tr> <td>STOT SE Respiratory tract irritant</td> <td>3</td> <td>H335</td> </tr> <tr> <td>Acute toxicity (oral)</td> <td>4</td> <td>H302</td> </tr> <tr> <td>Acute toxicity (dermal)</td> <td>4</td> <td>H312</td> </tr> <tr> <td>Harmful if inhaled</td> <td>4</td> <td>H332</td> </tr> </tbody> </table>	<u>Hazard Class</u>	<u>Hazard Category</u>	<u>Hazard Statement</u>	Skin irritation	2	H315	Eye irritation	2	H319	STOT SE Respiratory tract irritant	3	H335	Acute toxicity (oral)	4	H302	Acute toxicity (dermal)	4	H312	Harmful if inhaled	4	H332	
<u>Hazard Class</u>	<u>Hazard Category</u>	<u>Hazard Statement</u>																						
Skin irritation	2	H315																						
Eye irritation	2	H319																						
STOT SE Respiratory tract irritant	3	H335																						
Acute toxicity (oral)	4	H302																						
Acute toxicity (dermal)	4	H312																						
Harmful if inhaled	4	H332																						
2.2	<b>GHS Label Elements</b>																							
	<b>Hazard Pictogram</b>																							
																								
	<b>Signal Word: Warning</b>																							
	<b>Restricted to use by dental professional only</b>																							
	<b>Hazard Statements:</b>																							
	H315: Causes skin irritation. Category 2																							
	H319: Causes serious eye irritation. Category 2.																							
	H335: STOT SE Respiratory tract irritant. Category 3.																							
	H302: Harmful if swallowed. Category 4																							
	H312: Harmful in contact with skin. Category 4																							
	H332: Harmful if inhaled. Category 4																							
	<b>Precautionary Statements</b>																							

**Safety Data Sheet**Trade Name: **FORENDO™ PASTE**, Calcium Hydroxide with Iodoform

P261: Avoid breathing vapors.

P264: Wash hands thoroughly after handling.

P280: Wear protective gloves, lab coat and eye/face protection.

P305+P351+P338: If in eyes, rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do so. Continue rinsing.

P301+P316: If swallowed, get emergency medical help immediately.

P330: Rinse mouth.

P302+P352: If on skin, wash with plenty of soap and water.

P363: Wash contaminated clothing before reuse.

P304+P340: If inhaled, remove victim to fresh air and keep at rest in a position comfortable for breathing.

<b>3.0 Composition</b>													
3.1	Chemical characterization of the preparation      Calcium hydroxide and Iodoform in a silicone-based paste.												
3.2	Hazardous ingredients												
	<table border="1"> <thead> <tr> <th>CAS Number</th> <th>Name of the Ingredient</th> <th>Concentration</th> <th>Classification per Regulation (EC) No.1272/2008 (CLP)</th> </tr> </thead> <tbody> <tr> <td>1305-62-0</td> <td>Calcium Hydroxide</td> <td>20-30 %</td> <td>Skin irritation, 2, H315 Eye irritation, 2, H319</td> </tr> <tr> <td>75-47-8</td> <td>Iodoform</td> <td>30-40 %</td> <td>STOT SE Respiratory tract irritant, 4, H335 Acute toxicity (oral), 4, H302 Acute toxicity (dermal), 4, H312 Harmful if inhaled, 4, H332</td> </tr> </tbody> </table>	CAS Number	Name of the Ingredient	Concentration	Classification per Regulation (EC) No.1272/2008 (CLP)	1305-62-0	Calcium Hydroxide	20-30 %	Skin irritation, 2, H315 Eye irritation, 2, H319	75-47-8	Iodoform	30-40 %	STOT SE Respiratory tract irritant, 4, H335 Acute toxicity (oral), 4, H302 Acute toxicity (dermal), 4, H312 Harmful if inhaled, 4, H332
CAS Number	Name of the Ingredient	Concentration	Classification per Regulation (EC) No.1272/2008 (CLP)										
1305-62-0	Calcium Hydroxide	20-30 %	Skin irritation, 2, H315 Eye irritation, 2, H319										
75-47-8	Iodoform	30-40 %	STOT SE Respiratory tract irritant, 4, H335 Acute toxicity (oral), 4, H302 Acute toxicity (dermal), 4, H312 Harmful if inhaled, 4, H332										
<b>4.0 First Aid Measures</b>													
4.1	Special Instructions      Calcium hydroxide is a corrosive. Iodoform is considered harmful by inhalation, ingestion or contact with eyes or skin. However, when combined in a silicone-based paste, as in this product, they present a moderate hazard of irritation with prolonged exposure. Show this safety data sheet to medical personnel. Get medical attention in case of uncertainty.												
4.2	Eye Contact      Move eyelids apart; flush with running water for 15+ minutes. Get medical attention.												
4.3	Skin Contact      Remove any contaminated clothing. Wipe off excess from skin. Immediately flush skin with running water. Get medical attention if irritation persists.												
4.4	Ingestion      Rinse mouth with water. Do not induce vomiting unless directed by medical personnel. Give large amounts of water to dilute. Get immediate medical attention for ingestion of large amount. Never give anything by mouth to an unconscious person.												
4.5	Inhalation      Move to fresh air. If necessary, administer artificial respiration and/or oxygen and seek immediate medical attention.												
4.6	Precautions for first responders      Wear safety glasses and gloves to prevent contact.												
4.7	Information for physician												
	Symptoms      Irritation, pain or redness in eyes, on skin, on mucous membranes.												
	Hazards      May be irritating to eyes, skin, respiratory system, gastrointestinal tract or mucous membranes. May cause burns, especially with prolonged exposure. There is a risk of serious damage to eyes.												

**Safety Data Sheet**Trade Name: **FORENDO™ PASTE**, Calcium Hydroxide with Iodoform

Treatment	Same as First Aid above.	
<b>5.0</b>	<b>Fire Fighting Measures</b>	
5.1	Suitable extinguishing media	Not a fire hazard. Extinguish fire with agent suitable for surrounding fire.
5.2	Extinguishing media to avoid	None
5.3	Special exposure hazards in a fire	Toxic gases and vapors may be released under fire conditions.
5.4	Special protective equipment for fire-fighters	Self-contained breathing apparatus; full, protective clothing.
<b>6.0</b>	<b>Accidental Release Measures</b>	
6.1	Personal precautions.	Wear chemical splash goggles, impervious gloves, lab coat with long sleeves and shoe covers.
6.2	Environmental precautions	Avoid release into the environment.
6.3	Method for clean up	Scrape / wipe up with damp paper towels. Wash spill area with detergent and water. Place spill and absorbent towels in covered container for disposal.
<b>7.0</b>	<b>Handling and Storage</b>	
7.1	Handling	For use by dental professionals only. Wear gloves, safety glasses; wash hands after use. Avoid unnecessary exposure. Follow good hygiene practices.
7.2	Storage	Remove applicator tip immediately after use. Keep tightly capped in original container. Store at cool room temperature (<25°C) in a well-ventilated area. Avoid extremes of temperature (>27°C/80°F, <5°C/40°F) and acids.
7.3	Specific uses	Dental preparation
<b>8.0</b>	<b>Exposure Controls / Personal Protection</b>	
8.1	Exposure limit values	Calcium hydroxide powder: PEL/TLV 5 mg/m <sup>3</sup> Iodoform: 0.6 ppm (10 mg/m <sup>3</sup> )
8.2	Exposure controls	.
8.2.1	Occupational exposure controls	Wear safety glasses, lab coat and gloves. No other equipment required under normal conditions of use of this product in the quantity provided.
8.2.1.1	Respiratory protection	Good general ventilation is sufficient to control any airborne vapors.
8.2.1.2	Hand protection	No special requirements other than the usual surgical gloves.
8.2.1.3	Eye protection	No special requirements other than the usual safety glasses.
8.2.1.4	Skin protection	Good personal hygiene and safety practices; wearing a lab coat.
8.2.1.5	Other controls	Emergency eye wash fountain. Wash hands after use.
8.2.2	Environmental exposure controls	Avoid release into environment as calcium hydroxide may cause pH variation and iodoform may be harmful.
<b>9.0</b>	<b>Physical and Chemical Properties</b>	
9.1	Characteristics	
9.1.1	Appearance /Color / Physical state	Smooth, yellow paste
9.1.2	Odor	Characteristic, unpleasant, iodine odor

**Safety Data Sheet**Trade Name: **FORENDO™ PASTE**, Calcium Hydroxide with Iodoform

9.2	Important health, safety and environmental information	
9.2.1	pH	pH > 12
9.2.2	Boiling point	Decomposes
9.2.3	Flash point	Not combustible
9.2.4	Flammability (solid, gas)	Not combustible
9.2.5	Explosive properties	Not applicable
9.2.6	Oxidizing properties	Not determined
9.2.7	Vapor pressure	Not determined
9.2.8	Specific gravity	Not determined
9.2.9	Solubility in water	Nil
9.2.10	Partition coefficient	Not determined
9.2.11	Viscosity	Not determined
9.2.12	Vapor density	Not determined
9.2.13	Evaporation rate	Not determined
<b>10.0</b>	<b>Stability and reactivity</b>	
10.1	Conditions to avoid	Open flame, other ignition sources, temperature over 27°C, intense light.
10.2	Materials to avoid	Strong oxidizers, acetone.
10.3	Hazardous decomposition products	Toxic iodine vapors, caustic fumes of calcium oxide, carbon monoxide, carbon dioxide, hydrogen iodide.
10.4	Further information	Stable under normal conditions.
<b>11.0</b>	<b>Toxicological information</b>	
11.1	Acute toxicity	Iodoform, Oral rat LD <sub>50</sub> : 355 mg/Kg Calcium hydroxide, Oral Rat LD <sub>50</sub> : 7340 mg/Kg
11.2	Irritation and corrosiveness	Forendo presents a minimal health hazard under normal conditions of dental practice. Prolonged exposure may cause burns/irritation of mucous membranes, eyes, skin, gastrointestinal tract, respiratory tract, dermatitis.
11.3	Sensitization	Not applicable.
11.4	Sub-acute, sub-chronic and prolonged toxicity	Chronic exposure may cause skin rash, headache and skin eruptions in severe cases. Persons with pre-existing skin conditions may be more susceptible to the irritating effects of this material.
11.5	Carcinogenicity, Mutagenicity, Reproductive Toxicity	Not considered a carcinogen, mutagen, teratogen or a cause of reproductive toxicity.
11.6	Empirical data	Not available
11.7	Clinical Experience	Similar calcium hydroxide/iodoform preparations have been marketed to and used safely and effectively by dental professionals for many years.
<b>12.0</b>	<b>Ecological Information</b>	
12.1	Ecotoxicity	Avoid release into environment as calcium hydroxide may cause pH variation and iodoform may be harmful. Follow government regulations.

**Safety Data Sheet**Trade Name: **FORENDO™ PASTE**, Calcium Hydroxide with Iodoform

<b>13.0 Disposal Considerations</b>		
13.1	Regulations	Follow all local and national government regulations in disposing material or contaminated packaging.
<b>14.0 Transport Information</b>		
14.1	Restrictions	None. Not regulated.
14.2	UN Number	None
14.3	Technical name	Not applicable
14.4	Packing group	Not applicable
	IATA class	Not applicable
<b>15.0 Regulatory Information</b>		
15.1	EU	Class IIa medical device under MDD 93/42/EEC.
15.2	US FDA	Class II medical device
15.3	Health Canada	Class III medical device
<b>16.0 Other information</b>		
16.1	Hazard Statements	H315: Causes skin irritation. Category 2 H319: Causes serious eye irritation. Category 2. H335: STOT SE Respiratory tract irritant. Category 3. H302: Harmful if swallowed. Category 4 H312: Harmful in contact with skin. Category 4 H332: Harmful if inhaled. Category 4
16.2	Precautionary Statements	P261: Avoid breathing vapors. P264: Wash hands thoroughly after handling. P280: Wear protective gloves, lab coat and eye/face protection. P305+P351+P338: If in eyes, rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do so. Continue rinsing. P301+P316: If swallowed, get emergency medical help immediately. P330: Rinse mouth. P302+P352: If on skin, wash with plenty of soap and water. P363: Wash contaminated clothing before reuse. P304+P340: If inhaled, remove victim to fresh air and keep at rest in a position comfortable for breathing.
16.3	Restrictions on use	Calcium hydroxide/iodoform preparations, like Forendo, are to be used only by dental professionals who have been trained in their use.
16.4	Further information	The information presented herein is believed to be factual as it has been derived from the works of persons believed to be qualified experts. However, nothing contained in this information is to be taken as a warranty or representation for which Pulpdent Corporation bears legal responsibility. The user should review any recommendations in the specific context of the intended use to determine whether they are appropriate.
16.5	Sources of key data	National Institute for Occupational Safety (NIOSH) US Occupational Safety and Health Administration (OSHA)

## Safety Data Sheet

Trade Name: **FORENDO™ PASTE**, Calcium Hydroxide with Iodoform

Eur-Lex European Union Law: Regulation (EC) No. 1272/2008 (CLP), Regulation (EU) No 487/2013 and Regulation (EC) No. 1907/2006 (REACH).

European Chemicals Agency Guidance on the compilation of safety data sheets Version 4.0 April 2020

16.6 Information which has been added, deleted or revised.

This Safety Data Sheet has been revised to meet the requirements of the GHS SDS format, Regulations (EC) No. 487/2013 (CLP), , (EC) No. 1907/2006 (REACH) and ECHA Version 4.0 April 2020.

---