



Trade name: GumQuick scan Comp. A

Substance number: 40921

Version: 1 / GB

Date revised: 10.02.2023

Replaces Version: - / GB

Print date: 24.03.2023

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

GumQuick scan Comp. A

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/preparation

Gingiva mask material

### 1.3. Details of the supplier of the safety data sheet

#### Address/Manufacturer

Dreve Dentamid GmbH

Max-Planck-Straße 31

59423 Unna

Telephone no. +49 2303 8807-0

Fax no. +49 2303 8807-29

Information provided by / telephone Department Research &amp; Development: Fax: +49 2303 8807-562

E-mail address of person responsible sicherheitsdatenblatt@dreve.com

for this SDS

### 1.4. Emergency telephone number

Henkel Fire Department / 24h-Emergency-Contact-No.: +49 211 797-3350

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

This product is not classified hazardous in accordance with Regulation (EC) No 1272/2008.

### 2.2. Label elements

#### Labelling according to regulation (EC) No 1272/2008

#### Supplemental information

EUH210

Safety data sheet available on request.

### 2.3. Other hazards

No special hazards have to be mentioned.

The product contains no PBT substances. The product contains no vPvB substances. This product does not contain a substance that has endocrine disrupting properties with respect to human. The product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms.

## SECTION 3: Composition/information on ingredients

### 3.2. Mixtures

#### Chemical characterization

Addition-vulcanising 2-component silicone



Trade name: GumQuick scan Comp. A

Substance number: 40921

Version: 1 / GB

Date revised: 10.02.2023

Replaces Version: - / GB

Print date: 24.03.2023

## Hazardous ingredients

### Cristobalite

CAS No.	14464-46-1				
EINECS no.	238-455-4				
Concentration	>= 10	<	25	%	
Classification (Regulation (EC) No. 1272/2008)	STOT RE 1		H372		

### White mineral oil

CAS No.	8042-47-5				
EINECS no.	232-455-8				
Registration no.	01-2119487078-27				
Concentration	>= 1	<	10	%	
Classification (Regulation (EC) No. 1272/2008)	Asp. Tox. 1		H304		

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

#### General information

No special measures required

#### After inhalation

Ensure supply of fresh air. In the event of symptoms take medical treatment.

#### After skin contact

In case of contact with skin wash off with warm water. Consult a doctor if skin irritation persists.

#### After eye contact

Separate eyelids, wash the eyes thoroughly with water (15 min.). In case of irritation consult an oculist.

#### After ingestion

Do not induce vomiting - aspiration hazard. Summon a doctor immediately.

#### Adhere to personal protective measures when giving first aid

First aider: Pay attention to self-protection!

### 4.2. Most important symptoms and effects, both acute and delayed

Until now no symptoms known so far.

### 4.3. Indication of any immediate medical attention and special treatment needed

#### Hints for the physician / hazards

In the case of swallowing with subsequent vomiting, aspiration of the lungs can occur which can lead to chemical pneumonia or asphyxiation.

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

#### Suitable extinguishing media

Recommended: alcohol resistant foam, CO<sub>2</sub>, powders, water spray/mist, Extinguishing measures to suit surroundings



Trade name: GumQuick scan Comp. A

Substance number: 40921

Version: 1 / GB

Date revised: 10.02.2023

Replaces Version: - / GB

Print date: 24.03.2023

**Non suitable extinguishing media**

Full water jet

**5.2. Special hazards arising from the substance or mixture**

In case of combustion evolution of dangerous gases possible.

**5.3. Advice for firefighters****Special protective equipment for fire-fighting**

In case of combustion use a suitable breathing apparatus.

**Other information**

Fire residues and contaminated fire-fighting water must be disposed of in accordance with the local regulations.

**SECTION 6: Accidental release measures****6.1. Personal precautions, protective equipment and emergency procedures**

Avoid contact with skin, eyes and clothing.

**6.2. Environmental precautions**

Do not allow to enter drains or waterways.

**6.3. Methods and material for containment and cleaning up**

Pick up mechanically. Clean contaminated floors and objects thoroughly, observing environmental regulations.

**6.4. Reference to other sections**

Refer to protective measures listed in Sections 7 and 8.

**SECTION 7: Handling and storage****7.1. Precautions for safe handling****Advice on safe handling**

Observe the usual precautions for handling chemicals. For personal protection see Section 8.

**Advice on protection against fire and explosion**

No special measures required.

**7.2. Conditions for safe storage, including any incompatibilities****Requirements for storage rooms and vessels**

Store product in closed containers.

**Hints on storage assembly**

Do not store together with foodstuffs.

**Further information on storage conditions**

Keep container tightly closed and dry.

**SECTION 8: Exposure controls/personal protection****8.1. Control parameters**

Trade name: GumQuick scan Comp. A

Substance number: 40921

Version: 1 / GB

Date revised: 10.02.2023

Replaces Version: - / GB

Print date: 24.03.2023

## Exposure limit values

### White mineral oil

List	TRGS 900	
Type	AGW	
A		
Value	5	mg/m <sup>3</sup>
Short term exposure limit	20	mg/m <sup>3</sup>
Maximum limit value: 4(II); Pregnancy group: Y; Status: Sept 2015; Remarks: DGF		

### Other information

Contains no substances with occupational exposure limit values.

### Derived No/Minimal Effect Levels (DNEL/DMEL)

#### White mineral oil

Type of value	Derived No Effect Level (DNEL)	
Reference group	Worker	
Duration of exposure	Repeated exposure	
Route of exposure	inhalative	
Mode of action	Systemic effects	
Concentration	164,6	mg/m <sup>3</sup>
Type of value	Derived No Effect Level (DNEL)	
Reference group	Worker	
Duration of exposure	Repeated exposure	
Route of exposure	dermal	
Mode of action	Systemic effects	
Concentration	217,1	mg/kg/d
Type of value	Derived No Effect Level (DNEL)	
Reference group	Consumer	
Duration of exposure	Repeated exposure	
Route of exposure	inhalative	
Mode of action	Systemic effects	
Concentration	34,78	mg/m <sup>3</sup>
Type of value	Derived No Effect Level (DNEL)	
Reference group	Consumer	
Duration of exposure	Repeated exposure	
Route of exposure	dermal	
Mode of action	Systemic effects	
Concentration	93,02	mg/kg/d
Type of value	Derived No Effect Level (DNEL)	
Reference group	Consumer	
Duration of exposure	Repeated exposure	
Route of exposure	oral	
Mode of action	Systemic effects	
Concentration	25	mg/kg/d

## 8.2. Exposure controls

### General protective and hygiene measures

Observe the usual precautions for handling chemicals. Do not eat, drink or smoke during work time.

Wash hands before breaks and after work. Avoid contact with eyes.

### Respiratory protection

Not necessary.

Trade name: GumQuick scan Comp. A

Substance number: 40921

Version: 1 / GB

Date revised: 10.02.2023

Replaces Version: - / GB

Print date: 24.03.2023

**Hand protection**

Not necessary.

**Eye protection**

Not necessary.

**Body protection**

Clothing as usual in the chemical industry.

## SECTION 9: Physical and chemical properties

**9.1. Information on basic physical and chemical properties**

<b>Physical state</b>	liquid, viscous		
<b>Colour</b>	white		
<b>Odour</b>	characteristic		
<b>Melting point</b>			
Remarks	not determined		
<b>Freezing point</b>			
Remarks	not determined		
<b>Boiling point or initial boiling point and boiling range</b>			
Value	> 300		°C
<b>Flammability</b>			
not determined			
<b>Upper and lower explosive limits</b>			
Remarks	not determined		
<b>Flash point</b>			
Value	> 130		°C
Method	closed cup		
<b>Ignition temperature</b>			
Remarks	not determined		
<b>Decomposition temperature</b>			
Remarks	not determined		
<b>pH value</b>			
Remarks	not determined		
<b>Viscosity</b>			
Remarks	not determined		
<b>Solubility(ies)</b>			
Remarks	not determined		
<b>Partition coefficient n-octanol/water (log value)</b>			
Remarks	not determined		
<b>Vapour pressure</b>			
Remarks	not determined		
<b>Density and/or relative density</b>			
Value	1,34		g/cm <sup>3</sup>
Temperature	20	°C	
<b>Relative vapour density</b>			
Remarks	not determined		

Trade name: GumQuick scan Comp. A

Substance number: 40921

Version: 1 / GB

Date revised: 10.02.2023

Replaces Version: - / GB

Print date: 24.03.2023

**9.2. Other information****Odour threshold**

Remarks not determined

**Evaporation rate (ether = 1) :**

Remarks not determined

**Solubility in water**

Remarks virtually insoluble

**Explosive properties**

evaluation not determined

**Oxidising properties**

Remarks not determined

**Other information**

None known

## SECTION 10: Stability and reactivity

**10.1. Reactivity**

No hazardous reactions when stored and handled according to prescribed instructions.

**10.2. Chemical stability**

No hazardous reactions known.

**10.3. Possibility of hazardous reactions**

No hazardous reactions known.

**10.4. Conditions to avoid**

No hazardous reactions known.

**10.5. Incompatible materials**

None known

**10.6. Hazardous decomposition products**

No hazardous decomposition products known.

## SECTION 11: Toxicological information

**11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008****Acute oral toxicity**

Remarks Based on available data, the classification criteria are not met.

**Acute oral toxicity (Components)****White mineral oil**

Species	rat	
LD50	> 5000	mg/kg
Method	OECD 401	

**Acute dermal toxicity**

Remarks Based on available data, the classification criteria are not met.

**Acute dermal toxicity (Components)**

Trade name: GumQuick scan Comp. A

Substance number: 40921

Version: 1 / GB

Date revised: 10.02.2023

Replaces Version: - / GB

Print date: 24.03.2023

**White mineral oil**

Species	rabbit	
LD50	> 2000	mg/kg
Method	OECD 402	

**Acute inhalational toxicity**

ATE	> 20	mg/l
Administration/Form	Dust/Mist	
Method	calculated value (Regulation (EC) No. 1272/2008)	
Remarks	Based on available data, the classification criteria are not met.	

**Acute inhalative toxicity (Components)****White mineral oil**

Species	rat	
LC50	> 5	mg/l
Duration of exposure	4	h
Administration/Form	Dust/Mist	
Method	OECD 403	

**Skin corrosion/irritation**

Remarks	Based on available data, the classification criteria are not met.
---------	-------------------------------------------------------------------

**Serious eye damage/irritation**

Remarks	Based on available data, the classification criteria are not met.
---------	-------------------------------------------------------------------

**Sensitization**

Remarks	Based on available data, the classification criteria are not met.
---------	-------------------------------------------------------------------

**Subacute, subchronic, chronic toxicity**

Remarks	not determined
---------	----------------

**Mutagenicity**

Remarks	Based on available data, the classification criteria are not met.
---------	-------------------------------------------------------------------

**Reproductive toxicity**

Remarks	Based on available data, the classification criteria are not met.
---------	-------------------------------------------------------------------

**Carcinogenicity**

Remarks	Based on available data, the classification criteria are not met.
---------	-------------------------------------------------------------------

**Specific Target Organ Toxicity (STOT)****Single exposure**

Remarks	Based on available data, the classification criteria are not met.
---------	-------------------------------------------------------------------

**Repeated exposure**

Remarks	Based on available data, the classification criteria are not met.
---------	-------------------------------------------------------------------

**Specific Target Organ Toxicity (STOT) (Components)****Cristobalite****Repeated exposure**

evaluation	Causes damage to organs through prolonged or repeated exposure Route of exposure inhalative
------------	------------------------------------------------------------------------------------------------

**Aspiration hazard**

Remarks	Based on available data, the classification criteria are not met.
---------	-------------------------------------------------------------------

**Aspiration hazard (Components)****White mineral oil**

Harmful:	may cause lung damage if swallowed.
----------	-------------------------------------

**11.2 Information on other hazards**

Trade name: GumQuick scan Comp. A

Substance number: 40921

Version: 1 / GB

Date revised: 10.02.2023

Replaces Version: - / GB

Print date: 24.03.2023

**Endocrine disrupting properties with respect to humans**

The product does not contain a substance that has endocrine disrupting properties with respect to humans.

**Other information**

No toxicological data are available.

## SECTION 12: Ecological information

**12.1. Toxicity****General information**

not determined

**Fish toxicity (Components)****White mineral oil**

Species	golden orfe ( <i>Leuciscus idus</i> )		
LC50	>	10000	mg/l
Duration of exposure	96	h	
Method	OECD 203		

**Daphnia toxicity (Components)****White mineral oil**

Species	Daphnia magna		
LL50	>	100	mg/l
Duration of exposure	48	h	
Method	OECD 202		

**White mineral oil**

Species	Daphnia magna		
NOEC		10	mg/l
Duration of exposure	21	d	
Method	OECD 211		

**Algae toxicity (Components)****White mineral oil**

Species	Pseudokirchneriella subcapitata		
LOEC	>=	100	mg/l
Duration of exposure	72	h	
Method	OECD 201		

**12.2. Persistence and degradability****General information**

not determined

**Biodegradability (Components)****White mineral oil**

Value	31	%
Duration of test evaluation	28	d
	biodegradable	

**12.3. Bioaccumulative potential****General information**

not determined

**Partition coefficient n-octanol/water (log value)**





Trade name: GumQuick scan Comp. A

Substance number: 40921

Version: 1 / GB

Date revised: 10.02.2023

Replaces Version: - / GB

Print date: 24.03.2023

Remarks not determined

**Octanol/water partition coefficient (log Pow) (Components)**

**White mineral oil**

log Pow	<=	4,3	to	18.2
Temperature		20	°C	
Source	calculated value			

**12.4. Mobility in soil**

**General information**

not determined

**12.5. Results of PBT and vPvB assessment**

**General information**

not determined

**Results of PBT and vPvB assessment**

The product contains no PBT substances  
The product contains no vPvB substances.

**12.6 Endocrine disrupting properties**

**Endocrine disrupting properties with respect to the environment**

The product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms.

**12.7. Other adverse effects**

**General information**

not determined

**General information / ecology**

Do not allow to enter soil, waterways or waste water canal. Avoid release into the atmosphere.

**SECTION 13: Disposal considerations**

**13.1. Waste treatment methods**

**Disposal recommendations for the product**

Must not be disposed together with household garbage.  
Dispose of waste according to applicable legislation.

**Disposal recommendations for packaging**

Packaging that cannot be cleaned should be disposed off in agreement with the regional waste disposal company.

**SECTION 14: Transport information**

	Land transport ADR/RID	Marine transport IMDG/GGVSee	Air transport ICAO/IATA
<b>14.1. UN number</b>	The product does not constitute a hazardous substance in land transport.	The product does not constitute a hazardous substance in sea transport.	The product does not constitute a hazardous substance in air transport.



Trade name: GumQuick scan Comp. A

Substance number: 40921

Version: 1 / GB

Date revised: 10.02.2023

Replaces Version: - / GB

Print date: 24.03.2023

14.2. UN proper shipping name	-	-	-
14.3. Transport hazard class(es)		-	-
Label			
14.4. Packing group		-	-
14.5. Environmental hazards		no	-

## SECTION 15: Regulatory information

### 15.2. Chemical safety assessment

For this preparation a chemical safety assessment has not been carried out.

## SECTION 16: Other information

### Hazard statements listed in Chapter 3

H304 May be fatal if swallowed and enters airways.  
H372 Causes damage to organs through prolonged or repeated exposure.

### CLP categories listed in Chapter 3

Asp. Tox. 1 Aspiration hazard, Category 1  
STOT RE 1 Specific target organ toxicity - repeated exposure, Category 1

### Supplemental information

Relevant changes compared with the previous version of the safety data sheet are marked with: \*\*\*  
This information is based on our present state of knowledge. However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.



Trade name: GumQuick scan Comp. B

Substance number: 40922

Version: 1 / GB

Date revised: 13.02.2023

Replaces Version: - / GB

Print date: 13.10.2023

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

GumQuick scan Comp. B

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/preparation

Gingiva mask material

### 1.3. Details of the supplier of the safety data sheet

#### Address/Manufacturer

Dreve Dentamid GmbH

Max-Planck-Straße 31

59423 Unna

Telephone no. +49 2303 8807-0

Fax no. +49 2303 8807-29

Information provided by / telephone Department Research &amp; Development: Fax: +49 2303 8807-562

E-mail address of person responsible sicherheitsdatenblatt@dreve.com

for this SDS

### 1.4. Emergency telephone number

Henkel Fire Department / 24h-Emergency-Contact-No.: +49 211 797-3350

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

This product is not classified hazardous in accordance with Regulation (EC) No 1272/2008.

### 2.2. Label elements

#### Labelling according to regulation (EC) No 1272/2008

#### Supplemental information

EUH210 Safety data sheet available on request.

### 2.3. Other hazards

No special hazards have to be mentioned.

The product contains no PBT substances. The product contains no vPvB substances. This product does not contain a substance that has endocrine disrupting properties with respect to human. The product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms.

## SECTION 3: Composition/information on ingredients

### 3.2. Mixtures

#### Chemical characterization

Addition-vulcanising 2-component silicone

Trade name: GumQuick scan Comp. B

Substance number: 40922

Version: 1 / GB

Date revised: 13.02.2023

Replaces Version: - / GB

Print date: 13.10.2023

**Hazardous ingredients****Cristobalite**

CAS No.	14464-46-1				
EINECS no.	238-455-4				
Concentration	>= 10	<	25	%	
Classification (Regulation (EC) No. 1272/2008)	STOT RE 1		H372		

**SECTION 4: First aid measures****4.1. Description of first aid measures****General information**

No special measures required

**After inhalation**

Ensure supply of fresh air. In the event of symptoms take medical treatment.

**After skin contact**

In case of contact with skin wash off with warm water. Consult a doctor if skin irritation persists.

**After eye contact**

Separate eyelids, wash the eyes thoroughly with water (15 min.). In case of irritation consult an oculist.

**After ingestion**

Do not induce vomiting - aspiration hazard. Summon a doctor immediately.

**Adhere to personal protective measures when giving first aid**

First aider: Pay attention to self-protection!

**4.2. Most important symptoms and effects, both acute and delayed**

Until now no symptoms known so far.

**4.3. Indication of any immediate medical attention and special treatment needed****Hints for the physician / hazards**

In the case of swallowing with subsequent vomiting, aspiration of the lungs can occur which can lead to chemical pneumonia or asphyxiation.

**SECTION 5: Firefighting measures****5.1. Extinguishing media****Suitable extinguishing media**Recommended: alcohol resistant foam, CO<sub>2</sub>, powders, water spray/mist, Extinguishing measures to suit surroundings**Non suitable extinguishing media**

Full water jet

**5.2. Special hazards arising from the substance or mixture**

In case of combustion evolution of dangerous gases possible.

**5.3. Advice for firefighters****Special protective equipment for fire-fighting**

Trade name: GumQuick scan Comp. B

Substance number: 40922

Version: 1 / GB

Date revised: 13.02.2023

Replaces Version: - / GB

Print date: 13.10.2023

In case of combustion use a suitable breathing apparatus.

**Other information**

Fire residues and contaminated fire-fighting water must be disposed of in accordance with the local regulations.

**SECTION 6: Accidental release measures****6.1. Personal precautions, protective equipment and emergency procedures**

Avoid contact with skin, eyes and clothing.

**6.2. Environmental precautions**

Do not allow to enter drains or waterways.

**6.3. Methods and material for containment and cleaning up**

Pick up mechanically. Clean contaminated floors and objects thoroughly, observing environmental regulations.

**6.4. Reference to other sections**

Refer to protective measures listed in Sections 7 and 8.

**SECTION 7: Handling and storage****7.1. Precautions for safe handling****Advice on safe handling**

Observe the usual precautions for handling chemicals. For personal protection see Section 8.

**Advice on protection against fire and explosion**

No special measures required.

**7.2. Conditions for safe storage, including any incompatibilities****Requirements for storage rooms and vessels**

Store product in closed containers.

**Hints on storage assembly**

Do not store together with foodstuffs.

**Further information on storage conditions**

Keep container tightly closed and dry.

**SECTION 8: Exposure controls/personal protection****8.1. Control parameters****Other information**

Contains no substances with occupational exposure limit values.

**8.2. Exposure controls****General protective and hygiene measures**

Observe the usual precautions for handling chemicals. Do not eat, drink or smoke during work time. Wash hands before breaks and after work. Avoid contact with eyes.

**Respiratory protection**

Trade name: GumQuick scan Comp. B

Substance number: 40922

Version: 1 / GB

Date revised: 13.02.2023

Replaces Version: - / GB

Print date: 13.10.2023

Not necessary.

**Hand protection**

Not necessary.

**Eye protection**

Not necessary.

**Body protection**

Clothing as usual in the chemical industry.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

<b>Physical state</b>	liquid, viscous		
<b>Colour</b>	pink		
<b>Odour</b>	characteristic		
<b>Melting point</b>			
Remarks	not determined		
<b>Freezing point</b>			
Remarks	not determined		
<b>Boiling point or initial boiling point and boiling range</b>			
Value	> 300		°C
<b>Flammability</b>			
not determined			
<b>Upper and lower explosive limits</b>			
Remarks	not determined		
<b>Flash point</b>			
Value	> 130		°C
Method	closed cup		
<b>Ignition temperature</b>			
Remarks	not determined		
<b>Decomposition temperature</b>			
Remarks	not determined		
<b>pH value</b>			
Remarks	not determined		
<b>Viscosity</b>			
Remarks	not determined		
<b>Solubility(ies)</b>			
Remarks	not determined		
<b>Partition coefficient n-octanol/water (log value)</b>			
Remarks	not determined		
<b>Vapour pressure</b>			
Remarks	not determined		
<b>Density and/or relative density</b>			
Value	1,34		g/cm <sup>3</sup>
Temperature	20	°C	
<b>Relative vapour density</b>			



Trade name: GumQuick scan Comp. B

Substance number: 40922

Version: 1 / GB

Date revised: 13.02.2023

Replaces Version: - / GB

Print date: 13.10.2023

Remarks not determined

**9.2. Other information****Odour threshold**

Remarks not determined

**Evaporation rate (ether = 1) :**

Remarks not determined

**Solubility in water**

Remarks virtually insoluble

**Explosive properties**

evaluation not determined

**Oxidising properties**

Remarks not determined

**Other information**

None known

**SECTION 10: Stability and reactivity****10.1. Reactivity**

No hazardous reactions when stored and handled according to prescribed instructions.

**10.2. Chemical stability**

No hazardous reactions known.

**10.3. Possibility of hazardous reactions**

No hazardous reactions known.

**10.4. Conditions to avoid**

No hazardous reactions known.

**10.5. Incompatible materials**

None known

**10.6. Hazardous decomposition products**

No hazardous decomposition products known.

**SECTION 11: Toxicological information****11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008****Acute oral toxicity**

Remarks Based on available data, the classification criteria are not met.

**Acute dermal toxicity**

Remarks Based on available data, the classification criteria are not met.

**Acute inhalational toxicity**

ATE &gt; 20 mg/l

Administration/Form Dust/Mist

Method calculated value (Regulation (EC) No. 1272/2008)

Remarks Based on available data, the classification criteria are not met.

**Skin corrosion/irritation**

Trade name: GumQuick scan Comp. B

Substance number: 40922

Version: 1 / GB

Date revised: 13.02.2023

Replaces Version: - / GB

Print date: 13.10.2023

Remarks Based on available data, the classification criteria are not met.

**Serious eye damage/irritation**

Remarks Based on available data, the classification criteria are not met.

**Sensitization**

Remarks Based on available data, the classification criteria are not met.

**Subacute, subchronic, chronic toxicity**

Remarks not determined

**Mutagenicity**

Remarks Based on available data, the classification criteria are not met.

**Reproductive toxicity**

Remarks Based on available data, the classification criteria are not met.

**Carcinogenicity**

Remarks Based on available data, the classification criteria are not met.

**Specific Target Organ Toxicity (STOT)****Single exposure**

Remarks Based on available data, the classification criteria are not met.

**Repeated exposure**

Remarks Based on available data, the classification criteria are not met.

**Specific Target Organ Toxicity (STOT) (Components)****Cristobalite****Repeated exposure**

evaluation Causes damage to organs through prolonged or repeated exposure  
Route of exposure inhalative

**Aspiration hazard**

Based on available data, the classification criteria are not met.

**11.2 Information on other hazards****Endocrine disrupting properties with respect to humans**

The product does not contain a substance that has endocrine disrupting properties with respect to humans.

**Other information**

No toxicological data are available.

## SECTION 12: Ecological information

**12.1. Toxicity****General information**

not determined

**12.2. Persistence and degradability****General information**

not determined

**12.3. Bioaccumulative potential****General information**

not determined





Trade name: GumQuick scan Comp. B

Substance number: 40922

Version: 1 / GB

Date revised: 13.02.2023

Replaces Version: - / GB

Print date: 13.10.2023

**Partition coefficient n-octanol/water (log value)**

Remarks not determined

**12.4. Mobility in soil**

**General information**

not determined

**12.5. Results of PBT and vPvB assessment**

**General information**

not determined

**Results of PBT and vPvB assessment**

The product contains no PBT substances  
The product contains no vPvB substances.

**12.6 Endocrine disrupting properties**

**Endocrine disrupting properties with respect to the environment**

The product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms.

**12.7. Other adverse effects**

**General information**

not determined

**General information / ecology**

Do not allow to enter soil, waterways or waste water canal. Avoid release into the atmosphere.

**SECTION 13: Disposal considerations**

**13.1. Waste treatment methods**

**Disposal recommendations for the product**

Must not be disposed together with household garbage.  
Dispose of waste according to applicable legislation.

**Disposal recommendations for packaging**

Packaging that cannot be cleaned should be disposed off in agreement with the regional waste disposal company.

**SECTION 14: Transport information**

	Land transport ADR/RID	Marine transport IMDG/GGVSee	Air transport ICAO/IATA
14.1. UN number	The product does not constitute a hazardous substance in land transport.	The product does not constitute a hazardous substance in sea transport.	The product does not constitute a hazardous substance in air transport.
14.2. UN proper shipping name	-	-	-
14.3. Transport hazard class(es)	-	-	-



Trade name: GumQuick scan Comp. B

Substance number: 40922

Version: 1 / GB

Date revised: 13.02.2023

Replaces Version: - / GB

Print date: 13.10.2023

Label			
14.4. Packing group		-	-
14.5. Environmental hazards		no	-

## SECTION 15: Regulatory information

### 15.2. Chemical safety assessment

For this preparation a chemical safety assessment has not been carried out.

## SECTION 16: Other information

### Hazard statements listed in Chapter 3

H372 Causes damage to organs through prolonged or repeated exposure.

### CLP categories listed in Chapter 3

STOT RE 1 Specific target organ toxicity - repeated exposure, Category 1

### Supplemental information

Relevant changes compared with the previous version of the safety data sheet are marked with: \*\*\*  
This information is based on our present state of knowledge. However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.



Trade name: DuoSep

Substance number: 26600

Version: 2 / GB

Date revised: 16.05.2023

Replaces Version: 1 / GB

Print date: 16.05.2023

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

DuoSep

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/preparation

Isolating agent between silicone and polyether materials

### 1.3. Details of the supplier of the safety data sheet

#### Address/Manufacturer

Dreve Dentamid GmbH

Max-Planck-Straße 31

59423 Unna

Telephone no.

+49 2303 8807-0

Fax no.

+49 2303 8807-29

Information provided

Department Research &amp; Development: Fax: +49 2303 8807-562

by / telephone

E-mail address of

sicherheitsdatenblatt@dreve.com

person responsible

for this SDS

### 1.4. Emergency telephone number

Henkel Fire Department / 24h-Emergency-Contact-No.: +49 211 797-3350

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

Voluntary product information following the Safety Data Sheet format

This product is not classified hazardous in accordance with Regulation (EC) No 1272/2008.

### 2.2. Label elements

#### Labelling according to regulation (EC) No 1272/2008

The product does not require a hazard warning label in accordance with Regulation (EC) No 1272/2008.

### 2.3. Other hazards

No special hazards have to be mentioned.

The product contains no PBT substances. The product contains no vPvB substances. This product does not contain a substance that has endocrine disrupting properties with respect to human. The product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms.

## SECTION 3: Composition/information on ingredients

### 3.2. Mixtures



Trade name: DuoSep

Substance number: 26600

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 16.05.2023

Print date: 16.05.2023

**Chemical characterization**

Water-based insulating agent

**SECTION 4: First aid measures****4.1. Description of first aid measures****General information**

Remove contaminated, soaked clothing immediately and dispose of safely.

**After inhalation**

Ensure supply of fresh air. In the event of symptoms take medical treatment.

**After skin contact**

In case of contact with skin wash off with warm water. Consult a doctor if skin irritation persists.

**After eye contact**

Separate eyelids, wash the eyes thoroughly with water (15 min.). In case of irritation consult an oculist.

**After ingestion**

Rinse mouth thoroughly with water. Do not induce vomiting. Summon a doctor immediately.

**Adhere to personal protective measures when giving first aid**

First aider: Pay attention to self-protection!

**4.2. Most important symptoms and effects, both acute and delayed**

Until now no symptoms known so far.

**4.3. Indication of any immediate medical attention and special treatment needed****Hints for the physician / hazards**

In the case of swallowing with subsequent vomiting, aspiration of the lungs can occur which can lead to chemical pneumonia or asphyxiation.

**SECTION 5: Firefighting measures****5.1. Extinguishing media****Suitable extinguishing media**Recommended: alcohol resistant foam, CO<sub>2</sub>, powders, water spray/mist, Extinguishing measures to suit surroundings**5.3. Advice for firefighters****Special protective equipment for fire-fighting**

In case of combustion use a suitable breathing apparatus.

**SECTION 6: Accidental release measures****6.1. Personal precautions, protective equipment and emergency procedures**

Refer to protective measures listed in Sections 7 and 8.

**6.2. Environmental precautions**

Trade name: DuoSep

Substance number: 26600

Version: 2 / GB

Date revised: 16.05.2023

Replaces Version: 1 / GB

Print date: 16.05.2023

Do not allow to enter drains or waterways.

### 6.3. Methods and material for containment and cleaning up

Pick up with absorbent material. Clean contaminated floors and objects thoroughly, observing environmental regulations. Dispose of as prescribed.

### 6.4. Reference to other sections

Refer to protective measures listed in Sections 7 and 8.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

#### Advice on safe handling

Provide good ventilation of working area (local exhaust ventilation if necessary).

#### Advice on protection against fire and explosion

No special measures required.

### 7.2. Conditions for safe storage, including any incompatibilities

#### Requirements for storage rooms and vessels

Store product in closed containers.

#### Hints on storage assembly

Do not store together with foodstuffs.

#### Further information on storage conditions

Keep container tightly closed and dry.

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

#### Other information

Contains no substances with occupational exposure limit values.

### 8.2. Exposure controls

#### General protective and hygiene measures

Observe the usual precautions for handling chemicals. Do not eat, drink or smoke during work time. Keep away from foodstuffs and beverages.

#### Respiratory protection

Not necessary, but do not inhale vapours.

#### Hand protection

Not necessary.

#### Eye protection

Safety glasses

#### Body protection

Clothing as usual in the chemical industry.

## SECTION 9: Physical and chemical properties



Trade name: DuoSep

Substance number: 26600

Version: 2 / GB

Date revised: 16.05.2023

Replaces Version: 1 / GB

Print date: 16.05.2023

**9.1. Information on basic physical and chemical properties**

<b>Physical state</b>	liquid		
<b>Colour</b>	colourless		
<b>Odour</b>	characteristic		
<b>Melting point</b>			
Remarks	not determined		
<b>Freezing point</b>			
Remarks	not determined		
<b>Boiling point or initial boiling point and boiling range</b>			
Value	100		°C
<b>Flammability</b>			
evaluation	not determined		
<b>Upper and lower explosive limits</b>			
Remarks	not determined		
<b>Flash point</b>			
Value	80		°C
Method	closed cup		
<b>Ignition temperature</b>			
Value	360		°C
<b>Decomposition temperature</b>			
Remarks	not determined		
<b>pH value</b>			
Remarks	not determined		
<b>Viscosity</b>			
Remarks	not determined		
<b>Solubility(ies)</b>			
Remarks	not determined		
<b>Partition coefficient n-octanol/water (log value)</b>			
Remarks	not determined		
<b>Vapour pressure</b>			
Value	23		hPa
<b>Density and/or relative density</b>			
Value	appr. 1,04		g/cm <sup>3</sup>
Temperature	20	°C	
<b>Relative vapour density</b>			
Remarks	not determined		

**9.2. Other information**

<b>Odour threshold</b>			
Remarks	not determined		
<b>Evaporation rate (ether = 1) :</b>			
Remarks	not determined		
<b>Solubility in water</b>			
Remarks	miscible		
<b>Explosive properties</b>			
evaluation	no		

Trade name: DuoSep

Substance number: 26600

Version: 2 / GB

Date revised: 16.05.2023

Replaces Version: 1 / GB

Print date: 16.05.2023

**Oxidising properties**

Remarks not determined

**Other information**

None known

**SECTION 10: Stability and reactivity****10.1. Reactivity**

No hazardous reactions when stored and handled according to prescribed instructions.

**10.2. Chemical stability**

No hazardous reactions known.

**10.3. Possibility of hazardous reactions**

No hazardous reactions known.

**10.4. Conditions to avoid**

No hazardous reactions known.

**10.5. Incompatible materials**

None known

**10.6. Hazardous decomposition products**

No hazardous decomposition products known.

**SECTION 11: Toxicological information****11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008****Acute oral toxicity**

Remarks Based on available data, the classification criteria are not met.

**Acute dermal toxicity**

Remarks Based on available data, the classification criteria are not met.

**Acute inhalational toxicity**

Remarks Based on available data, the classification criteria are not met.

**Skin corrosion/irritation**

Remarks Based on available data, the classification criteria are not met.

**Serious eye damage/irritation**

Remarks Based on available data, the classification criteria are not met.

**Sensitization**

Remarks Based on available data, the classification criteria are not met.

**Subacute, subchronic, chronic toxicity**

Remarks not determined

**Mutagenicity**

Remarks Based on available data, the classification criteria are not met.

**Reproductive toxicity**

Remarks Based on available data, the classification criteria are not met.

**Carcinogenicity**

Trade name: DuoSep

Substance number: 26600

Version: 2 / GB

Date revised: 16.05.2023

Replaces Version: 1 / GB

Print date: 16.05.2023

Remarks Based on available data, the classification criteria are not met.

### Specific Target Organ Toxicity (STOT)

#### Single exposure

Remarks Based on available data, the classification criteria are not met.

#### Repeated exposure

Remarks Based on available data, the classification criteria are not met.

### Aspiration hazard

Based on available data, the classification criteria are not met.

## 11.2 Information on other hazards

### Endocrine disrupting properties with respect to humans

The product does not contain a substance that has endocrine disrupting properties with respect to humans.

### Other information

No toxicological data are available.

## SECTION 12: Ecological information

### 12.1. Toxicity

#### General information

not determined

### 12.2. Persistence and degradability

#### General information

not determined

### 12.3. Bioaccumulative potential

#### General information

not determined

#### Partition coefficient n-octanol/water (log value)

Remarks not determined

### 12.4. Mobility in soil

#### General information

not determined

### 12.5. Results of PBT and vPvB assessment

#### General information

not determined

#### Results of PBT and vPvB assessment

The product contains no PBT substances  
The product contains no vPvB substances.

### 12.6 Endocrine disrupting properties

#### Endocrine disrupting properties with respect to the environment

The product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms.

### 12.7. Other adverse effects



Trade name: DuoSep

Substance number: 26600

Version: 2 / GB

Date revised: 16.05.2023

Replaces Version: 1 / GB

Print date: 16.05.2023

**General information**

not determined

**General information / ecology**

Do not discharge product unmonitored into the environment.

## SECTION 13: Disposal considerations

**13.1. Waste treatment methods****Disposal recommendations for the product**

Must not be disposed together with household garbage.  
Dispose of waste according to applicable legislation.

**Disposal recommendations for packaging**

Packaging that cannot be cleaned should be disposed off in agreement with the regional waste disposal company.

## SECTION 14: Transport information

	Land transport ADR/RID	Marine transport IMDG/GGVSee	Air transport ICAO/IATA
<b>14.1. UN number or ID number</b>	The product does not constitute a hazardous substance in land transport.	The product does not constitute a hazardous substance in sea transport.	The product does not constitute a hazardous substance in air transport.
<b>14.2. UN proper shipping name</b>	-	-	-
<b>14.3. Transport hazard class(es)</b>		-	-
Label			
<b>14.4. Packing group</b>		-	-

## SECTION 15: Regulatory information

**15.2. Chemical safety assessment**

For this preparation a chemical safety assessment has not been carried out.

## SECTION 16: Other information

**Supplemental information**

Relevant changes compared with the previous version of the safety data sheet are marked with: \*\*\*  
This information is based on our present state of knowledge. However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.



Trade name: DuoSep

Substance number: 26600

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 16.05.2023

Print date: 16.05.2023