

**Safety Information Sheet for Medical Devices
following 1907/2006/EG, Article 31**

Date of issue: 13.06.2024
Printing date: 24.06.2024

Version: 1

This product is a medical device as defined in Directive 93/42/EEC on medical devices (MDD) respectively Regulation (EU) 2017/745 (MDR), which is invasive or used in direct physical contact with the human body. For this product a safety data sheet is not legally required. The provision of this document is on a voluntary basis.

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Trade name:

Harvard TEMP Cem

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses:

Dental product.

Restrictions on use:

This product should only be supplied or handled to dentists or dental technicians or upon their instructions.

1.3 Details of the supplier of the safety information sheet

Manufacturer/supplier:

Harvard Dental International GmbH
Margaretenstrasse 2-4
15366 Hoppegarten, Germany
info@harvard-dental-international.de

1.4 Emergency telephone number

+49 30 99 28 978-0

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Medical devices as defined in Directive 93/42/EEC on medical devices (MDD) respectively Regulation (EU) 2017/245, which are invasive or used in direct physical contact with the human body are exempt from labeling requirements of Regulation (EC) No. 1272/2008 (CLP/GHS). Although not required, classification and labeling are indicated as follows:

Classification:

Acute aquatic toxicity, Category 1, H400 Very toxic to aquatic life.

Chronic aquatic toxicity, Category 1, H410 Very toxic to aquatic life with long lasting effects.

Skin irritation, Category 2; H315 Causes skin irritation.

Eye irritation, Category 2; H319 Causes serious eye irritation.

Specific target organ toxicity – single exposure (STOT SE), Category 3; H335 May cause respiratory Irritation.

2.2 Label elements

CLP Regulation (EC) No. 1272/2008

Signal word:

Warning.

**Safety Information Sheet for Medical Devices
following 1907/2006/EG, Article 31**

Date of issue: 13.06.2024
Printing date: 24.06.2024

Version: 1

Hazard pictograms:



GHS07



GHS09

Hazard statements:

H400 Very toxic to aquatic life.
H315 Causes skin irritation.
H319 Causes serious eye irritation.
H335 May cause respiratory irritation.

Precautionary statements:

P280 Wear protective gloves/protective clothing/eye protection/face protection.
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P302 + P352 IF ON SKIN: Wash with plenty of soap and water.
P273 Avoid release to the environment

2.3 Other hazards

Results of PBT and vPvB assessment:

PBT: not applicable.
vPvB: not applicable.

SECTION 3: Composition/information on ingredients

3.1 Substances

Not applicable.

3.2 Mixtures

Chemical characterization (description):

Organic acids, zinc oxide and methacrylates.

Hazardous components:

CAS No. 1314-13-2	Zinc oxide	GHS09, H410, H411
CAS No. 688-84-6	Methacrylate	GHS07, H315, H319, H335
	Mixture of organic acids	GHS07, H315, H319, H335

SECTION 4: First aid measures

4.1 Description of first aid measures

General information:

Seek medical advice, if symptoms occur, which can be caused by the product.

after skin contact:

Wash with plenty of water and soap.

after inhalation:

Fresh air.

after eye contact:

Rinse with plenty of water and consult a doctor.

**Safety Information Sheet for Medical Devices
following 1907/2006/EG, Article 31**

Date of issue: 13.06.2024
Printing date: 24.06.2024

Version: 1

after swallowing:

Consult a doctor immediately.

SECTION 5: Firefighting measures

5.1 Extinguishing media

Extinguishing media:

CO₂, extinguishing powder, water spray. Use extinguishing method and media suitable to surrounding fire.

5.2 Special hazards arising from the substance or mixture

No further relevant information available.

Protective equipment:

Wear self-contained respiratory protective equipment.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Personal precautions:

Not required.

6.2 Environmental precautions

Do not allow to enter drains/surface water/ ground-water.

6.3 Methods and material for containment and cleaning up

Pick up mechanically. Clean up affected area.

6.4 Reference to other sections

See Section 7: Handling and storage

See Section 8: Exposure controls/personal protection

See Section 13: Disposal considerations

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Handling:

This product should only be supplied or handled to dentists or dental technicians or upon their instructions.

Information for safe handling:

Observe normal care for working with chemicals.

Information about fire and explosion protection:

No special measures required.

7.2 Conditions for safe storage, including any incompatibilities

Storage:

Store at a dry place. Store tightly closed. Storage temperature 8 - 25°C.

Requirements at storerooms and containers:

Store only in the original containers.

**Safety Information Sheet for Medical Devices
following 1907/2006/EG, Article 31**

Date of issue: 13.06.2024
Printing date: 24.06.2024

Version: 1

Information about storage in one common storage facility:

Store away from food stuffs.

Further information on storage conditions:

None.

Storage class:

-

7.3 Specific end use(s)

No further relevant information available.

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

Additional information on the design of technical systems:

-

Components with limits of values to be supervised at the workplace:

The product contains no relevant quantities of components with limits of values to be supervised at the workplace.

Additional information:

The lists valid during the making were used as basis.

8.2 Exposure controls

Personal protective equipment:

General protective and hygienic measures:

Normal hygienic measures. The usual precautionary measures are to be adhered to when handling chemicals. Keep away from foodstuffs and beverages. Avoid contact with eyes or skin.

Respiratory protection:

Required, when aerosols are generated.

Hand protection:

Protective gloves. Due to missing tests no recommendation to the glove material can be given for the product.

Eye protection:

Safety glasses.

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Appearance:

Form:	paste
Color:	white yellowish
Odor:	odorless

Change in condition:

Melting point/melting range:	not applicable
Boiling point/boiling range:	not applicable
Flash point:	not applicable

**Safety Information Sheet for Medical Devices
following 1907/2006/EG, Article 31**

Date of issue: 13.06.2024
Printing date: 24.06.2024

Version: 1

Ignition temperature:	product is not selfigniting
Danger of explosion:	none
Density:	1.71 (20°C) g/cm ³
Vapor pressure:	not applicable (nonvolatile)
Solubility in/miscibility with:	slightly soluble in toluene, organic solvents
Water:	insoluble
Content of solvents:	none
Organic solvents:	none

SECTION 10: Stability and reactivity

10.1 Reactivity

Expected to be non-reactive under normal conditions of use.

10.2 Chemical stability

The product is stable by proper use and storage.

10.4 Conditions to avoid

Heat, sunlight, peroxides, metal powder, strong oxidizing agents, radical initiators.

10.6 Hazardous decomposition products

No information available.

SECTION 11: Toxicological information

11.1 Information on toxicological effects

Acute toxicity:

No data available.

Primary irritation:

Skin: slight irritant.

Eye: irritant.

Sensitization: no information available.

Additional toxicological information:

Allergic reactions to the product after long-term exposure are possible. When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us.

SECTION 12: Ecological information

General information:

Avoid transfer into environment.

Classification of water endangerment:

WGK=2 (polluting substance).

**Safety Information Sheet for Medical Devices
following 1907/2006/EG, Article 31**

Date of issue: 13.06.2024
Printing date: 24.06.2024

Version: 1

SECTION 13: Disposal considerations

13.1 Waste treatment methods

Recommendation:

Use proper landfill disposal or incineration in accordance with local, state and federal regulations.

Uncleaned packaging:

Recommendation:

Disposal must be made according to official regulations.

SECTION 14: Transport information

Land transport ADR/RID and GGVSE (Germany):

ADR, RID

UN 3077 UMWELTGEFÄHRDENDER STOFF, FEST N.A.G (Zinkoxid), III, 9.

Sea transport IMDG Code:

IMDG-Code

UN 3077 ENVIRONMENTALLY HAZARDOUS SUBSTANCE, SOLID, N.O.S. (zinc oxide), 9, III

EmS: F-A S-F.

Air transport ICAO-TI and IATA-DGR:

CAO, PAX

UN 3077 ENVIRONMENTALLY HAZARDOUS SUBSTANCE, SOLID, N.O.S. (zinc oxide), 9, III.

The transport regulations are cited according to international regulations and in the form applicable in Germany. Possible national deviations in other countries are not considered.

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

National regulations:

Water endangerment WGK=2 (polluting substance).

15.2 Chemical safety assessment

A chemical safety assessment has not been carried out.

SECTION 16: Other information

For invasive medical devices or medical devices in direct contact with the human body, the Medical Device Regulation does not require a safety data sheet. The safe use of the product is indicated in the instructions for use and/or the labeling.

The above information is based on our present day knowledge and relates solely to the safety requirements of the product. The data do not signify any warranty with regards to products properties. Users of the product should satisfy themselves that the information given is sufficient and correct for their specific circumstances of use.