



**Safety data sheet  
according to 1907/2006/EC, Article 31**

Printing date 02.08.2022

Version number 4 (replaces version 3)

Revision: 02.08.2022

**SECTION 1: Identification of the substance/mixture and of the company/undertaking**

· **1.1 Product identifier**

· Trade name: **HeraCeram Stainfluid MF  
Butane-1,2-diol**

· CAS Number:

584-03-2

· EC number:

209-527-2

· **1.2 Relevant identified uses of the substance or mixture and uses advised against**

No further relevant information available.

· Application of the substance / the mixture Manufacture of dental prothesis

· **1.3 Details of the supplier of the safety data sheet**

· Manufacturer/Supplier:

Kulzer GmbH

Leipziger Straße 2, 63450 Hanau (Germany)

Tel.: +49 (0)800 4372522

· Informing department: E-Mail: [msds@kulzer-dental.com](mailto:msds@kulzer-dental.com)

· **1.4 Emergency telephone number:** Emergency CONTACT (24-Hour-Number): +49 (0)6132-84463

**SECTION 2: Hazards identification**

· **2.1 Classification of the substance or mixture**

· Classification according to Regulation (EC) No 1272/2008

Eye Irrit. 2 H319 Causes serious eye irritation.

· **2.2 Label elements**

· Labelling according to Regulation (EC) No 1272/2008

The substance is classified and labelled according to the GB CLP regulation.

· Hazard pictograms



GHS07

· Signal word Warning

· Hazard statements

H319 Causes serious eye irritation.

· Precautionary statements

P280

Wear protective gloves / eye protection.

P337+P313 If eye irritation persists: Get medical advice/attention.

· **2.3 Other hazards -**

· Results of PBT and vPvB assessment

· PBT: Not applicable.

· vPvB: Not applicable.

**SECTION 3: Composition/information on ingredients**

· **3.1 Substances**

· CAS No. Designation:

584-03-2 butane-1,2-diol

· Identification number(s):

· EC number: 209-527-2

· Description: -

(Contd. on page 2)



**Safety data sheet  
according to 1907/2006/EC, Article 31**

Printing date 02.08.2022

Version number 4 (replaces version 3)

Revision: 02.08.2022

**Trade name: HeraCeram Stainfluid MF  
Butane-1,2-diol**

(Contd. of page 1)

· **Additional information** For the wording of the listed hazard phrases refer to section 16.

#### **SECTION 4: First aid measures**

- **4.1 Description of first aid measures**
  - **After inhalation** Supply fresh air; consult doctor in case of symptoms.
  - **After skin contact**  
Instantly wash with water and soap and rinse thoroughly.  
If skin irritation continues, consult a doctor.
  - **After eye contact**  
Rinse opened eye for several minutes under running water. Then consult doctor.
  - **After swallowing**  
Rinse out mouth and then drink plenty of water.  
In case of persistent symptoms consult doctor.
- **4.2 Most important symptoms and effects, both acute and delayed**  
No further relevant information available.
- **4.3 Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

#### **SECTION 5: Firefighting measures**

- **5.1 Extinguishing media**
  - **Suitable extinguishing agents**  
CO<sub>2</sub>, extinguishing powder or water jet. Fight larger fires with water jet or alcohol-resistant foam.  
Use fire fighting measures that suit the environment.
- **5.2 Special hazards arising from the substance or mixture**  
Carbon dioxide (CO<sub>2</sub>)  
Carbon monoxide (CO)  
Formation of toxic gases is possible during heating or in case of fire.
- **5.3 Advice for firefighters**
  - **Protective equipment:** Wear self-contained breathing apparatus.
  - **Additional information** -

#### **SECTION 6: Accidental release measures**

- **6.1 Personal precautions, protective equipment and emergency procedures**  
Wear protective equipment. Keep unprotected persons away.
- **6.2 Environmental precautions:** Dilute with much water.
- **6.3 Methods and material for containment and cleaning up:**  
Absorb with liquid-binding material (diatomite, universal binders, for small amounts tissues).  
Send for recovery or disposal in suitable containers.
- **6.4 Reference to other sections**  
See Section 8 for information on personal protection equipment.  
See Section 13 for information on disposal.

#### **SECTION 7: Handling and storage**

- **7.1 Precautions for safe handling** No special measures required.
  - **Information about protection against explosions and fires:** No special measures required.

(Contd. on page 3)



**Safety data sheet  
according to 1907/2006/EC, Article 31**

Printing date 02.08.2022

Version number 4 (replaces version 3)

Revision: 02.08.2022

**Trade name: HeraCeram Stainfluid MF  
Butane-1,2-diol**

(Contd. of page 2)

- **7.2 Conditions for safe storage, including any incompatibilities**
  - **Storage**
    - **Requirements to be met by storerooms and containers:** No special requirements.
    - **Information about storage in one common storage facility:** Not required.
    - **Further information about storage conditions:** None.
- **7.3 Specific end use(s)** No further relevant information available.

### **SECTION 8: Exposure controls/personal protection**

- **8.1 Control parameters**
  - **Components with critical values that require monitoring at the workplace:**

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.  
Not required.

    - **Additional information:** The lists that were valid during the compilation were used as basis.
- **8.2 Exposure controls**
  - **Appropriate engineering controls** No further data; see item 7.
  - **Individual protection measures, such as personal protective equipment**
    - **General protective and hygienic measures**

The usual precautionary measures should be adhered to in handling the chemicals.
    - **Breathing equipment:** Use breathing protection in case of insufficient ventilation.
    - **Hand protection** Check protective gloves prior to each use for their proper condition.
      - **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.
      - **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
      - **For the permanent contact of a maximum of 15 minutes gloves made of the following materials are suitable:**

Butyl rubber, BR  
Nitrile rubber, NBR
    - **Eye/face protection** Tightly sealed safety glasses.
    - **Body protection:** Light weight protective clothing

### **SECTION 9: Physical and chemical properties**

- **9.1 Information on basic physical and chemical properties**
  - **General Information**

· <b>Physical state</b>	Fluid
· <b>Colour:</b>	Colourless
· <b>Smell:</b>	Odourless
· <b>Odour threshold:</b>	Not determined.
· <b>Melting point/freezing point:</b>	Not determined
· <b>Boiling point or initial boiling point and boiling range</b>	230 °C
· <b>Flammability</b>	Not applicable.

(Contd. on page 4)



**Safety data sheet**  
**according to 1907/2006/EC, Article 31**

Printing date 02.08.2022

Version number 4 (replaces version 3)

Revision: 02.08.2022

**Trade name: HeraCeram Stainfluid MF**  
**Butane-1,2-diol**

(Contd. of page 3)

· <b>Lower and upper explosion limit</b>	
· Lower:	Not determined.
· Upper:	Not determined.
· <b>Flash point:</b>	134 °C
· <b>Ignition temperature:</b>	385 °C
· <b>Decomposition temperature:</b>	Not determined.
· <b>SADT</b>	
· pH	Not determined.
· <b>Viscosity:</b>	
· Kinematic viscosity	Not determined.
· dynamic:	Not determined.
· <b>Solubility</b>	
· Water:	Fully miscible
· <b>Partition coefficient n-octanol/water (log value)</b>	Not determined.
· <b>Steam pressure:</b>	Not determined.
· <b>Density and/or relative density</b>	
· Density at 20 °C	1.020 g/cm <sup>3</sup>
· Relative density	Not determined.
· Vapour density	Not determined.

· <b>9.2 Other information</b>	No further relevant information available.
· <b>Appearance:</b>	
· Form:	Fluid
· <b>Important information on protection of health and environment, and on safety.</b>	
· Self-inflammability:	Not determined.
· Explosive properties:	Product is not explosive. Not determined.
· <b>Solvent content:</b>	
· VOC EU	g/l
· <b>Change in condition</b>	
· Evaporation rate	Not determined.

· <b>Information with regard to physical hazard classes</b>	
· Explosives	Void
· Flammable gases	Void
· Aerosols	Void
· Oxidising gases	Void
· Gases under pressure	Void
· Flammable liquids	Void
· Flammable solids	Void
· Self-reactive substances and mixtures	Void
· Pyrophoric liquids	Void
· Pyrophoric solids	Void
· Self-heating substances and mixtures	Void
· Substances and mixtures, which emit flammable gases in contact with water	Void
· Oxidising liquids	Void
· Oxidising solids	Void
· Organic peroxides	Void
· Corrosive to metals	Void

(Contd. on page 5)



**Safety data sheet  
according to 1907/2006/EC, Article 31**

Printing date 02.08.2022

Version number 4 (replaces version 3)

Revision: 02.08.2022

**Trade name: HeraCeram Stainfluid MF  
Butane-1,2-diol**

(Contd. of page 4)

· **Desensitised explosives** Void

**SECTION 10: Stability and reactivity**

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
  - **Conditions to be avoided:** No decomposition if used and stored according to specifications.
- **10.3 Possibility of hazardous reactions** No dangerous reactions known
- **10.4 Conditions to avoid** Heat, flames and sparks.
- **10.5 Incompatible materials:** Strong oxidant, acid chloride, acid anhydride, reductant
- **10.6 Hazardous decomposition products:** None

**SECTION 11: Toxicological information**

- **11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**
  - **Acute toxicity** Based on available data, the classification criteria are not met.

· **LD/LC50 values that are relevant for classification:**

**584-03-2 butane-1,2-diol**

Oral	LD50	>2,000 mg/kg (rat) (OECD 420)
Inhalative	LC50/4 h	370 mg/l (rat)

- **Skin corrosion/irritation** Based on available data, the classification criteria are not met.
  - **Serious eye damage/irritation**  
Causes serious eye irritation.
  - **Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.
  - **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
  - **Carcinogenicity** Based on available data, the classification criteria are not met.
  - **Reproductive toxicity** Based on available data, the classification criteria are not met.
  - **STOT-single exposure** Based on available data, the classification criteria are not met.
  - **STOT-repeated exposure** Based on available data, the classification criteria are not met.
  - **Aspiration hazard** Based on available data, the classification criteria are not met.
- **11.2 Information on other hazards**

· **Endocrine disrupting properties**

Substance is not listed.

**SECTION 12: Ecological information**

- **12.1 Toxicity**

· **Aquatic toxicity:**

**584-03-2 butane-1,2-diol**

EC50/48h	500 mg/l (daphnia)
LC50/96h	>10,000 mg/l (fish)

- **12.2 Persistence and degradability** No further relevant information available.
- **12.3 Bioaccumulative potential** No further relevant information available.
- **12.4 Mobility in soil** No further relevant information available.
- **12.5 Results of PBT and vPvB assessment**
  - **PBT:** Not applicable.
  - **vPvB:** Not applicable.

(Contd. on page 6)



**Safety data sheet  
according to 1907/2006/EC, Article 31**

Printing date 02.08.2022

Version number 4 (replaces version 3)

Revision: 02.08.2022

**Trade name: HeraCeram Stainfluid MF  
Butane-1,2-diol**

(Contd. of page 5)

- **12.6 Endocrine disrupting properties**  
The product does not contain substances with endocrine disrupting properties.
- **12.7 Other adverse effects** No further relevant information available.

### SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
  - **Recommendation** Disposal must be made according to official regulations.
- **Uncleaned packagings:**
  - **Recommendation:** Disposal must be made according to official regulations.
  - **Recommended cleaning agent:** Water, if necessary with cleaning agent.

### SECTION 14: Transport information

- |   |                 |
|---|-----------------|
| · <b>14.1 UN number or ID number</b><br>· ADR, ADN, IMDG, IATA                | Void            |
| · <b>14.2 UN proper shipping name</b><br>· ADR, ADN, IMDG, IATA               | Void            |
| · <b>14.3 Transport hazard class(es)</b><br>· ADR, ADN, IMDG, IATA<br>· Class | Void            |
| · <b>14.4 Packing group</b><br>· ADR, IMDG, IATA                              | Void            |
| · <b>14.5 Environmental hazards:</b><br>· <b>Marine pollutant:</b>            | No              |
| · <b>14.6 Special precautions for user</b>                                    | Not applicable. |
| · <b>14.7 Maritime transport in bulk according to IMO instruments</b>         | Not applicable. |
| · <b>UN "Model Regulation":</b>   | Void            |

### SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**  
No further relevant information available.
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

### SECTION 16: Other information

These data are based on our present knowledge. However, they shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Abbreviations and acronyms:**

SADT: Self Accelerating Decomposition Temperature

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

(Contd. on page 7)



**Safety data sheet**  
**according to 1907/2006/EC, Article 31**

Printing date 02.08.2022

Version number 4 (replaces version 3)

Revision: 02.08.2022

**Trade name: HeraCeram Stainfluid MF**  
**Butane-1,2-diol**

(Contd. of page 6)

*IATA: International Air Transport Association*

*GHS: Globally Harmonised System of Classification and Labelling of Chemicals*

*EINECS: European Inventory of Existing Commercial Chemical Substances*

*CAS: Chemical Abstracts Service (division of the American Chemical Society)*

*VOCV: Lenkungsabgabe auf flüchtigen organischen Verbindungen, Schweiz (Swiss Ordinance on volatile organic compounds)*

*LC50: Lethal concentration, 50 percent*

*LD50: Lethal dose, 50 percent*

*PBT: Persistent, Bioaccumulative and Toxic*

*vPvB: very Persistent and very Bioaccumulative*

*Eye Irrit. 2: Serious eye damage/eye irritation – Category 2*

**\* Data compared to the previous version altered.**

GB