



## Safety Data Sheet

Copyright, 2014, 3M Company All rights reserved. Copying and/or downloading of this information for the purpose of properly utilising 3M products is allowed provided that: (1) the information is copied in full with no changes unless prior written agreement is obtained from 3M, and (2) neither the copy nor the original is resold or otherwise distributed with the intention of earning a profit thereon.

<b>Document group:</b>	31-6679-0	<b>Version number:</b>	1.00
<b>Revision date:</b>	21/05/2014	<b>Supersedes date:</b>	Initial issue.
<b>Transportation version number:</b>	1.00 (21/05/2014)		

This Safety Data Sheet has been prepared in accordance with the REACH Regulation (EC) 1907/2006 and its modifications.

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

3M™ ESPE™IMPRINT™ 4 PENTA™ SUPER QUICK HEAVY BASE

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

##### Identified uses

Dental Product

#### 1.3. Details of the supplier of the substance or mixture

**Address:** 3M United Kingdom PLC, 3M Centre, Cain Road, Bracknell, Berkshire, RG12 8HT.

**Telephone:** +44 (0)1344 858 000

**E Mail:** tox.uk@mmm.com

**Website:** www.3M.com/uk

#### 1.4. Emergency telephone number

+44 (0)1344 858 000

### SECTION 2: Hazard identification

#### 2.1. Classification of the substance or mixture

CLP REGULATION (EC) No 1272/2008

##### CLASSIFICATION:

This material is exempt from hazard classification according to Regulation (EC) No. 1272/2008, as amended, on classification, labelling, and packaging of substances and mixtures.

**Dangerous substances(67/548/EEC)/preparations(1999/45/EC) directive**

#### 2.2. Label elements

CLP REGULATION (EC) No 1272/2008

Not applicable

1% of the mixture consists of components of unknown acute oral toxicity.

**Dangerous substances(67/548/EEC)/preparations(1999/45/EC) directive**

**3M™ ESPE™IMPRINT™ 4 PENTA™ SUPER QUICK HEAVY BASE**

Not applicable

**Notes on labelling**

This product is exempt from labelling per Directive 1999/45/EC as it is defined as a medical device according to Directive 93/42/EEC and is invasive or comes into contact with the human body.

**2.3. Other hazards**

None known.

**SECTION 3: Composition/information on ingredients**

Ingredient	CAS Nbr	EU Inventory	% by Wt	Classification
Silane treated quartz	None		50 - 60	
Siloxanes and silicones, Di-Me, vinyl group-terminated	68083-19-2		20 - 30	
Dimethyl methyl hydrogen silicone fluid	68037-59-2		5 - 15	Xn:R20 (Self Classified) Acute Tox. 4, H332 (Self Classified)
Dimethyl siloxane, reaction product with silica	67762-90-7		1 - 10	
Glycols,polyethylene,methyl 3-[1,3,3,3-tetramethyl-1-(trimethylsiloxy)disiloxanyl]propyl ether	27306-78-1		< 2	Xn:R20; Xi:R36 (Self Classified) Acute Tox. 4, H332; Eye Irrit. 2, H319 (Self Classified)
Aluminium Oxide	1344-28-1	EINECS 215-691-6	< 2	
Silane, trimethyl-2-propenyl-	762-72-1	EINECS 212-104-5	< 2	
Titanium dioxide	13463-67-7	EINECS 236-675-5	< 1.0	
Oils, mint, Mentha arvensis piperascensis, var. piperascens, Labiatae.	68917-18-0		< 0.5	Xn:R22; Xi:R36; R43 (Self Classified) Acute Tox. 4, H302; Eye Irrit. 2, H319; Skin Sens. 1B, H317 (Self Classified)

Please see section 16 for the full text of any R phrases and H statements referred to in this section

Please refer to section 15 for the any applicable Notas that have been applied to the above components

For information on ingredient occupational exposure limits or PBT or vPvB status, see sections 8 and 12 of this SDS

**SECTION 4: First aid measures****4.1. Description of first aid measures****Inhalation**

Remove person to fresh air. If you feel unwell, get medical attention.

**Skin contact**

Immediately wash with soap and water. Remove contaminated clothing and wash before reuse. If signs/symptoms develop, get medical attention.

**Eye contact**

Flush with large amounts of water. Remove contact lenses if easy to do. Continue rinsing. If signs/symptoms persist, get medical attention.

**If swallowed**

Rinse mouth. If you feel unwell, get medical attention.

**4.2. Most important symptoms and effects, both acute and delayed**

See Section 11.1 Information on toxicological effects

**4.3. Indication of any immediate medical attention and special treatment required**

Not applicable

**SECTION 5: Fire-fighting measures**

**5.1. Extinguishing media**

In case of fire: Use a fire fighting agent suitable for ordinary combustible material such as water or foam to extinguish.

**5.2. Special hazards arising from the substance or mixture**

None inherent in this product.

**Hazardous Decomposition or By-Products**

**Substance**

Carbon monoxide.  
Carbon dioxide.  
Hydrogen Chloride  
Irritant vapours or gases.  
Oxides of sulphur.

**Condition**

During combustion.  
During combustion.  
During combustion.  
During combustion.  
During combustion.

**5.3. Advice for fire-fighters**

No unusual fire or explosion hazards are anticipated.

**SECTION 6: Accidental release measures**

**6.1. Personal precautions, protective equipment and emergency procedures**

Ventilate the area with fresh air. Refer to other sections of this SDS for information regarding physical and health hazards, respiratory protection, ventilation, and personal protective equipment.

**6.2. Environmental precautions**

Avoid release to the environment.

**6.3. Methods and material for containment and cleaning up**

Collect as much of the spilled material as possible. Place in a closed container approved for transportation by appropriate authorities. Clean up residue. Seal the container. Dispose of collected material as soon as possible.

**6.4. Reference to other sections**

Refer to Section 8 and Section 13 for more information

**SECTION 7: Handling and storage**

**7.1. Precautions for safe handling**

Avoid prolonged or repeated skin contact. Do not get in eyes, on skin, or on clothing. Do not eat, drink or smoke when using this product. Wash thoroughly after handling. Contaminated work clothing should not be allowed out of the workplace. Wash contaminated clothing before reuse.

**7.2. Conditions for safe storage including any incompatibilities**

Store away from heat.

**7.3. Specific end use(s)**

See information in Section 7.1 and 7.2 for handling and storage recommendations. See Section 8 for exposure controls and personal protection recommendations.

**SECTION 8: Exposure controls/personal protection****8.1 Control parameters****Occupational exposure limits**

<b>Ingredient</b>	<b>CAS Nbr</b>	<b>Agency</b>	<b>Limit type</b>	<b>Additional comments</b>
Aluminium Oxide	1344-28-1	Health and Safety Comm. (UK)	TWA(as inhalable dust):10 mg/m <sup>3</sup> ;TWA(as respirable dust):4 mg/m <sup>3</sup>	
Titanium dioxide	13463-67-7	Health and Safety Comm. (UK)	TWA(Inhalable):10 mg/m <sup>3</sup> ;TWA(respirable):4 mg/m <sup>3</sup>	

Health and Safety Comm. (UK) : UK Health and Safety Commission

TWA: Time-Weighted-Average

STEL: Short Term Exposure Limit

CEIL: Ceiling

**Biological limit values**

No biological limit values exist for any of the components listed in Section 3 of this safety data sheet.

**8.2. Exposure controls****8.2.1. Engineering controls**

Use in a well-ventilated area.

**8.2.2. Personal protective equipment (PPE)****Eye/face protection**

Select and use eye/face protection to prevent contact based on the results of an exposure assessment. The following eye/face protection(s) are recommended:

Safety glasses with side shields.

**Skin/hand protection**

No chemical protective gloves are required. See Section 7.1 for additional information on skin protection.

**Respiratory protection**

Respiratory protection is not required.

**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

<b>Physical state</b>	Solid.
<b>Specific Physical Form:</b>	Paste
<b>Appearance/Odour</b>	slight minty odour;yellow colour paste

## 3M™ ESPE™ IMPRINT™ 4 PENTA™ SUPER QUICK HEAVY BASE

Odour threshold	No data available.
pH	No data available.
Boiling point/boiling range	Not applicable.
Melting point	Not applicable.
Flammability (solid, gas)	Not classified
Explosive properties	Not classified
Oxidising properties	Not classified
Flash point	No flash point
Autoignition temperature	No data available.
Flammable Limits(LEL)	Not applicable.
Flammable Limits(UEL)	Not applicable.
Vapour pressure	No data available.
Relative density	1.5 - 1.6 [Ref Std:WATER=1]
Water solubility	Negligible
Solubility- non-water	No data available.
Partition coefficient: n-octanol/water	Not applicable.
Evaporation rate	No data available.
Vapour density	No data available.
Decomposition temperature	No data available.
Viscosity	No data available.
Density	1.5 g/cm <sup>3</sup> - 1.6 g/cm <sup>3</sup>
<b>9.2. Other information</b>	
Volatile organic compounds (VOC)	Not applicable.
Percent volatile	Not applicable.
VOC less H <sub>2</sub> O & exempt solvents	Not applicable.

## SECTION 10: Stability and reactivity

### 10.1 Reactivity

This material may be reactive with certain agents under certain conditions - see the remaining headings in this section

### 10.2 Chemical stability

Stable.

### 10.3 Possibility of hazardous reactions

Hazardous polymerisation will not occur.

### 10.4 Conditions to avoid

Heat.

### 10.5 Incompatible materials

None known.

### 10.6 Hazardous decomposition products

#### Substance

None known.

#### Condition

Refer to section 5.2 for hazardous decomposition products during combustion.

## SECTION 11: Toxicological information

The information below may not agree with the EU material classification in Section 2 and/or the ingredient classifications in Section 3 if specific ingredient classifications are mandated by a competent authority. In addition, statements and data presented in Section 11 are based on UN GHS calculation rules and classifications derived from 3M assessments.

### 11.1 Information on Toxicological effects

#### Signs and Symptoms of Exposure

Based on test data and/or information on the components, this material may produce the following health effects:

#### Inhalation

This product may have a characteristic odour; however, no adverse health effects are anticipated.

#### Skin contact

Contact with the skin during product use is not expected to result in significant irritation. Allergic skin reaction (non-photo induced): Signs/symptoms may include redness, swelling, blistering, and itching.

#### Eye contact

Contact with the eyes during product use is not expected to result in significant irritation.

#### Ingestion

Gastrointestinal irritation: Signs/symptoms may include abdominal pain, stomach upset, nausea, vomiting and diarrhoea.

#### Carcinogenicity:

Exposures needed to cause the following health effect(s) are not expected during normal, intended use:

Contains a chemical or chemicals which can cause cancer.

#### Toxicological Data

If a component is disclosed in section 3 but does not appear in a table below, either no data are available for that endpoint or the data are not sufficient for classification.

#### Acute Toxicity

Name	Route	Species	Value
Overall product	Dermal		No data available; calculated ATE >5,000 mg/kg
Overall product	Ingestion		No data available; calculated ATE >5,000 mg/kg
Siloxanes and silicones, Di-Me, vinyl group-terminated	Dermal	Rabbit	LD50 > 15,440 mg/kg
Siloxanes and silicones, Di-Me, vinyl group-terminated	Ingestion	Rat	LD50 > 15,440 mg/kg
Dimethyl methyl hydrogen silicone fluid	Dermal	Rabbit	LD50 > 2,000 mg/kg
Dimethyl methyl hydrogen silicone fluid	Inhalation-Dust/Mist (4 hours)	Rat	LC50 4.2 mg/l
Dimethyl methyl hydrogen silicone fluid	Ingestion	Rat	LD50 > 2,000 mg/kg
Dimethyl siloxane, reaction product with silica	Dermal	Rabbit	LD50 > 5,000 mg/kg
Dimethyl siloxane, reaction product with silica	Inhalation-Dust/Mist (4 hours)	Rat	LC50 > 0.691 mg/l
Dimethyl siloxane, reaction product with silica	Ingestion	Rat	LD50 > 5,110 mg/kg
Glycols,polyethylene,methyl 3-[1,3,3,3-tetramethyl-1-(trimethylsiloxy)disiloxanyl]propyl ether	Dermal	Rabbit	LD50 > 2,000 mg/kg
Glycols,polyethylene,methyl 3-[1,3,3,3-tetramethyl-1-(trimethylsiloxy)disiloxanyl]propyl ether	Inhalation-Dust/Mist (4 hours)	Rat	LC50 2 mg/l
Glycols,polyethylene,methyl 3-[1,3,3,3-tetramethyl-1-(trimethylsiloxy)disiloxanyl]propyl ether	Ingestion	Rat	LD50 > 2,000 mg/kg
Aluminium Oxide	Dermal		LD50 estimated to be > 5,000 mg/kg
Aluminium Oxide	Inhalation-	Rat	LC50 > 2.3 mg/l

**3M™ ESPE™IMPRINT™ 4 PENTA™ SUPER QUICK HEAVY BASE**

	Dust/Mist (4 hours)		
Aluminium Oxide	Ingestion	Rat	LD50 > 5,000 mg/kg
Titanium dioxide	Dermal	Rabbit	LD50 > 10,000 mg/kg
Titanium dioxide	Inhalation- Dust/Mist (4 hours)	Rat	LC50 > 6.82 mg/l
Titanium dioxide	Ingestion	Rat	LD50 > 10,000 mg/kg
Oils, mint, Mentha arvensis piperascensis, var. piperascens, Labiatae.	Dermal	Rabbit	LD50 > 5,000 mg/kg
Oils, mint, Mentha arvensis piperascensis, var. piperascens, Labiatae.	Ingestion	Rat	LD50 1,240 mg/kg

ATE = acute toxicity estimate

**Skin Corrosion/Irritation**

Name	Species	Value
Siloxanes and silicones, Di-Me, vinyl group-terminated	Rabbit	No significant irritation
Dimethyl siloxane, reaction product with silica	Rabbit	No significant irritation
Glycols,polyethylene,methyl 3-[1,3,3,3-tetramethyl-1-(trimethylsiloxy)disiloxanyl]propyl ether	Rabbit	No significant irritation
Aluminium Oxide	Rabbit	No significant irritation
Titanium dioxide	Rabbit	No significant irritation
Oils, mint, Mentha arvensis piperascensis, var. piperascens, Labiatae.	Rabbit	Mild irritant

**Serious Eye Damage/Irritation**

Name	Species	Value
Siloxanes and silicones, Di-Me, vinyl group-terminated	Rabbit	Mild irritant
Dimethyl siloxane, reaction product with silica	Rabbit	No significant irritation
Glycols,polyethylene,methyl 3-[1,3,3,3-tetramethyl-1-(trimethylsiloxy)disiloxanyl]propyl ether	Rabbit	Severe irritant
Aluminium Oxide	Rabbit	No significant irritation
Titanium dioxide	Rabbit	No significant irritation
Oils, mint, Mentha arvensis piperascensis, var. piperascens, Labiatae.	In vitro data	Severe irritant

**Skin Sensitisation**

Name	Species	Value
Dimethyl siloxane, reaction product with silica	Human and animal	Not sensitizing
Glycols,polyethylene,methyl 3-[1,3,3,3-tetramethyl-1-(trimethylsiloxy)disiloxanyl]propyl ether	Guinea pig	Not sensitizing
Titanium dioxide	Human and animal	Not sensitizing
Oils, mint, Mentha arvensis piperascensis, var. piperascens, Labiatae.	Guinea pig	Sensitising

**Respiratory Sensitisation**

Name	Species	Value

**Germ Cell Mutagenicity**

Name	Route	Value
Dimethyl siloxane, reaction product with silica	In Vitro	Not mutagenic
Glycols,polyethylene,methyl 3-[1,3,3,3-tetramethyl-1-(trimethylsiloxy)disiloxanyl]propyl ether	In Vitro	Not mutagenic
Glycols,polyethylene,methyl 3-[1,3,3,3-tetramethyl-1-(trimethylsiloxy)disiloxanyl]propyl ether	In vivo	Not mutagenic
Aluminium Oxide	In Vitro	Not mutagenic
Titanium dioxide	In Vitro	Not mutagenic
Titanium dioxide	In vivo	Not mutagenic

**Carcinogenicity**

**3M™ ESPE™IMPRINT™ 4 PENTA™ SUPER QUICK HEAVY BASE**

Name	Route	Species	Value
Dimethyl siloxane, reaction product with silica	Not specified.	Mouse	Some positive data exist, but the data are not sufficient for classification
Aluminium Oxide	Inhalation	Rat	Not carcinogenic
Titanium dioxide	Ingestion	Multiple animal species	Not carcinogenic
Titanium dioxide	Inhalation	Rat	Carcinogenic.

**Reproductive Toxicity****Reproductive and/or Developmental Effects**

Name	Route	Value	Species	Test result	Exposure Duration
Dimethyl siloxane, reaction product with silica	Ingestion	Not toxic to female reproduction	Rat	NOAEL 509 mg/kg/day	1 generation
Dimethyl siloxane, reaction product with silica	Ingestion	Not toxic to male reproduction	Rat	NOAEL 497 mg/kg/day	1 generation
Dimethyl siloxane, reaction product with silica	Ingestion	Not toxic to development	Rat	NOAEL 1,350 mg/kg/day	during organogenesis
Glycols,polyethylene,methyl 3-[1,3,3,3-tetramethyl-1-(trimethylsiloxy)disiloxanyl]propyl ether	Ingestion	Some positive reproductive/developmental data exist, but the data are not sufficient for classification	Rat	NOAEL 450 mg/kg/day	pre mating & during gestation

**Target Organ(s)****Specific Target Organ Toxicity - single exposure**

Name	Route	Target Organ(s)	Value	Species	Test result	Exposure Duration

**Specific Target Organ Toxicity - repeated exposure**

Name	Route	Target Organ(s)	Value	Species	Test result	Exposure Duration
Dimethyl siloxane, reaction product with silica	Inhalation	respiratory system   silicosis	All data are negative	Human	NOAEL Not available	occupational exposure
Aluminium Oxide	Inhalation	pneumoconiosis   pulmonary fibrosis	Some positive data exist, but the data are not sufficient for classification	Human	NOAEL Not available	occupational exposure
Titanium dioxide	Inhalation	respiratory system	Some positive data exist, but the data are not sufficient for classification	Rat	LOAEL 0.010 mg/l	2 years
Titanium dioxide	Inhalation	pulmonary fibrosis	All data are negative	Human	NOAEL Not available	occupational exposure

**Aspiration Hazard**

Name	Value

Please contact the address or phone number listed on the first page of the SDS for additional toxicological information on this material and/or its components.

**SECTION 12: Ecological information**

The information below may not agree with the EU material classification in Section 2 and/or the ingredient classifications in Section 3 if specific ingredient classifications are mandated by a competent authority. In addition, statements and data presented in Section 12 are based on UN GHS calculation rules and classifications derived from 3M assessments.

**12.1. Toxicity**



**3M™ ESPE™IMPRINT™ 4 PENTA™ SUPER QUICK HEAVY BASE**

No product test data available.

Material	CAS Nbr	Organism	Type	Exposure	Test endpoint	Test result
Aluminium Oxide	1344-28-1	Fish	Experimental	96 hours	LC50	>100 mg/l
Aluminium Oxide	1344-28-1	Water flea	Experimental	48 hours	EC50	>100 mg/l
Aluminium Oxide	1344-28-1	Green algae	Experimental	72 hours	EC50	>100 mg/l
Titanium dioxide	13463-67-7	Sheepshead Minnow	Experimental	96 hours	LC50	>240 mg/l
Titanium dioxide	13463-67-7	Water flea	Experimental	48 hours	EC50	>100 mg/l
Aluminium Oxide	1344-28-1	Green algae	Experimental	72 hours	NOEC	>100 mg/l
Titanium dioxide	13463-67-7	Fish	Experimental	30 days	NOEC	>100 mg/l
Titanium dioxide	13463-67-7	Water flea	Experimental	30 days	NOEC	3 mg/l
Dimethyl methyl hydrogen silicone fluid	68037-59-2		Data not available or insufficient for classification			
Dimethyl siloxane, reaction product with silica	67762-90-7		Data not available or insufficient for classification			
Glycols,polyethylene,methyl 3-[1,3,3,3-tetramethyl-1-(trimethylsilyloxy)disiloxanyl]propyl ether	27306-78-1		Data not available or insufficient for classification			
Oils, mint, Mentha arvensis piperascensis, var. piperascens, Labiatae.	68917-18-0		Data not available or insufficient for classification			
Silane, trimethyl-2-propenyl-	762-72-1		Data not available or insufficient for classification			
Siloxanes and silicones, Di-Me, vinyl group-terminated	68083-19-2		Data not available or insufficient for classification			

**12.2. Persistence and degradability**

No test data available.

**3M™ ESPE™IMPRINT™ 4 PENTA™ SUPER QUICK HEAVY BASE**

Material	CAS Nbr	Test type	Duration	Study Type	Test result	Protocol
Dimethyl methyl hydrogen silicone fluid	68037-59-2	Data not available or insufficient for classification	N/A	N/A	N/A	N/A
Dimethyl siloxane, reaction product with silica	67762-90-7	Data not available or insufficient for classification	N/A	N/A	N/A	N/A
Glycols,polyethylene,methyl 3-[1,3,3,3-tetramethyl-1-(trimethylsilyloxy)disiloxanyl]propyl ether	27306-78-1	Data not available or insufficient for classification	N/A	N/A	N/A	N/A
Aluminium Oxide	1344-28-1	Data not available or insufficient for classification	N/A	N/A	N/A	N/A
Titanium dioxide	13463-67-7	Data not available or insufficient for classification	N/A	N/A	N/A	N/A
Oils, mint, Mentha arvensis piperascensis, var. piperascens, Labiatae.	68917-18-0	Data not available or insufficient for classification	N/A	N/A	N/A	N/A
Silane, trimethyl-2-propenyl-	762-72-1	Data not available or insufficient for classification	N/A	N/A	N/A	N/A
Siloxanes and silicones, Di-Me, vinyl group-terminated	68083-19-2	Data not available or insufficient for classification	N/A	N/A	N/A	N/A

**12.3 : Bioaccumulative potential**

Material	CAS Nbr	Test type	Duration	Study Type	Test result	Protocol
Dimethyl methyl hydrogen silicone fluid	68037-59-2	Data not available or insufficient for classification	N/A	N/A	N/A	N/A
Dimethyl siloxane, reaction product with silica	67762-90-7	Data not available or insufficient for classification	N/A	N/A	N/A	N/A

**3M™ ESPE™IMPRINT™ 4 PENTA™ SUPER QUICK HEAVY BASE**

Glycols,polyethylene,methyl 3-[1,3,3,3-tetramethyl-1-(trimethylsiloxy)disiloxanyl]propyl ether	27306-78-1	Data not available or insufficient for classification	N/A	N/A	N/A	N/A
Aluminium Oxide	1344-28-1	Data not available or insufficient for classification	N/A	N/A	N/A	N/A
Titanium dioxide	13463-67-7	Experimental BCF-Carp	42 days	Bioaccumulation factor	9.6	Other methods
Oils, mint, Mentha arvensis piperascensis, var. piperascens, Labiatae.	68917-18-0	Data not available or insufficient for classification	N/A	N/A	N/A	N/A
Silane, trimethyl-2-propenyl-	762-72-1	Data not available or insufficient for classification	N/A	N/A	N/A	N/A
Siloxanes and silicones, Di-Me, vinyl group-terminated	68083-19-2	Data not available or insufficient for classification	N/A	N/A	N/A	N/A

**12.4. Mobility in soil**

Please contact manufacturer for more details

**12.5. Results of the PBT and vPvB assessment**

No information available at this time, contact manufacturer for more details

**12.6. Other adverse effects**

No information available.

**SECTION 13: Disposal considerations****13.1 Waste treatment methods**

See Section 11.1 Information on toxicological effects

Dispose of completely cured (or polymerised) material in a permitted industrial waste facility. As a disposal alternative, incinerate uncured product in a permitted waste incineration facility. Combustion products will include halogen acid (HCl/HF/HBr). Facility must be capable of handling halogenated materials. If no other disposal options are available, waste product that has been completely cured or polymerised may be placed in a landfill properly designed for industrial waste.

The coding of a waste stream is based on the application of the product by the consumer. Since this is out of the control of 3M, no waste code(s) for products after use will be provided. Please refer to the European Waste Code (EWC - 2000/532/EC and amendments) to assign the correct waste code to your waste stream. Ensure national and/or regional regulations are complied with and always use a licensed waste contractor.

## 3M™ ESPE™IMPRINT™ 4 PENTA™ SUPER QUICK HEAVY BASE

### EU waste code (product as sold)

180106\* Chemicals consisting of or containing dangerous substances.

### EU waste code (product container after use)

180107 Chemicals other than those mentioned in 18 01 06

## SECTION 14: Transportation information

ADR/IATA/IMDG: Not restricted for transport.

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### Carcinogenicity

<u>Ingredient</u>	<u>CAS Nbr</u>	<u>Classification</u>	<u>Regulation</u>
Titanium dioxide	13463-67-7	Grp. 2B: Possible human carc.	International Agency for Research on Cancer

#### Global inventory status

All applicable chemical ingredients in this material are listed on the European Inventory of Existing Chemical Substances (EINECS), or are exempt polymers whose monomers are listed on EINECS. Contact 3M for more information. The components of this material are in compliance with the China "Measures on Environmental Management of New Chemical Substance". Certain restrictions may apply. Contact the selling division for additional information. The components of this product are in compliance with the new substance notification requirements of CEPA.

### 15.2. Chemical Safety Assessment

Not applicable

## SECTION 16: Other information

### List of relevant H statements

H302	Harmful if swallowed.
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.
H332	Harmful if inhaled.

### List of relevant R-phrases

R20	Harmful by inhalation.
R22	Harmful if swallowed.
R36	Irritating to eyes.
R43	May cause sensitisation by skin contact.

### Revision information:

No revision information is available.

DISCLAIMER: The information on this Safety Data Sheet is based on our experience and is correct to the best of our knowledge at the date of publication, but we do not accept any liability for any loss, damage or injury resulting from its use (except as required by law). The information may not be valid for any use not referred to in this Data Sheet or use of the product in combination with other materials. For these reasons, it is important that customers carry out their own test to satisfy themselves as to the suitability of the product for their own intended applications.

3M United Kingdom MSDSs are available at [www.3M.com/uk](http://www.3M.com/uk)





## Säkerhetsdatablad

Upphovsrätt, 2013, 3M Company. Samtliga rättigheter förbehållna. Kopiering och/eller nedladdning av denna information i syfte att tillgodogöra sig 3M:s produkter på tillbörligt sätt är tillåten under förutsättning att: (1) informationen kopieras i sin helhet utan några ändringar om inte 3M, i förväg lämnar skriftligt godkännande därtill, och (2) vare sig kopian eller originalet säljs vidare eller på annat sätt distribueras i vinstsyfte.

**Dokumentnummer:** 31-6686-5 **Version:** 1.00  
**Datum (nytt eller omarbetat):** 2013-09-03 **Föregående datum:** Första upplagan  
**Version (avser transportinformation):** 1.00 (2013-09-03)

Säkerhetsdatabladet har sammanställts i enlighet med REACH (EG nr 1907/2006 med ändringar).

### Avsnitt 1: Namnet på ämnet/blandningen och bolaget/företaget

#### 1.1 Produktbeteckning

3M™ ESPE™ IMPRINT™ 4 PENTA™ SUPER QUICK Catalyst

#### 1.2 Relevanta identifierade användningar av ämnet eller blandningen och användningar som det avråds från

##### Identifierade användningar

Dentalprodukt.

#### 1.3 Närmare upplysningar om den som tillhandahåller säkerhetsdatablad

**Adress:** 3M Svenska AB, 191 89 Sollentuna  
**Telefon:** 08-92 21 00  
**e-post:** miljo.sv@mmm.com  
**Hemsida:** www.3M.se

#### 1.4 Telefonnummer för nödsituationer

Giftinformationscentralen: 08-33 12 31 eller akut 112

### Avsnitt 2: Farliga egenskaper

#### 2.1 Klassificering av ämnet eller blandningen

CLP-förordningen (EG) nr 1272/2008

##### Klassificering:

Denna produkt är undantagen krav på faroklassificering enligt CLP (EG nr 1272/2008 med ändringar).

Substans-(67/548/EEG) /Preparat-(1999/45/EG) direktivet

#### 2.2 Märkningsuppgifter

CLP-förordningen (EG) nr 1272/2008

Substans-(67/548/EEG) /Preparat-(1999/45/EG) direktivet

Farosymbol(er)

-

**Innehåller**

Inga ingående ämnen behöver anges på etiketten.

**Riskfraser -**

Skyddsfraser -

**Kommentarer angående märkning**

Denna produkt varningsmärks ej, eftersom den är sådan medicinteknisk produkt som undantas från märkningskrav genom artikel 1.5d i CLP (EG nr 1272/2008).

**2.3 Andra faror**

Inga kända

**Avsnitt 3: Sammansättning/information om beståndsdelar**

Beståndsdelar	CAS-nr	EG-nr	Vikt-%	Klassificering
Nefelinsyenit	37244-96-5		60 - 70	
Vinylsiloxaner och vinylsilikoner	68083-19-2		15 - 25	
Poly(dimetylsiloxan)	63148-62-9		5 - 15	
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	67762-90-7		1 - 5	

Se avsnitt 16 för lydelse av de R-fraser och faroangivelser (H) som nämns i avsnitt 3.

Se avsnitt 2 för information om de anmärkningar som har tillämpats på ovanstående ämnen.

Tabellen visar klassificeringar fastställda inom EU samt kompletterande egenklassificeringar respektive klassificeringar från råvaruleverantörer.

För information om beståndsdelars hygieniska gränsvärde eller PBT/vPvB-status, se avsnitt 8 och 12 av detta SDB.

**Avsnitt 4: Åtgärder vid första hjälpen****4.1 Beskrivning av åtgärder vid första hjälpen****Inandning**

Kräver ingen första hjälpen åtgärd.

**Hudkontakt**

Tvätta med tvål och vatten. Om symptom uppstår, kontakta läkare.

**Ögonkontakt**

Skölj med stora mängder vatten. Ta ur eventuella kontaktlinser om det går lätt. Fortsätt att skölja. Om symptom kvarstår, sök läkarhjälp.

**Vid förtäring**

Skölj munnen. Vid obehag, sök läkarhjälp.

**4.2 De viktigaste symptomen och effekterna, både akuta och fördröjda**

Se avsnitt 11.1 Information om de toxikologiska effekterna

**4.3 Angivande av omedelbar medicinsk behandling och särskild behandling som eventuellt krävs**

Ej tillämpligt

## Avsnitt 5: Brandbekämpningsåtgärder

### 5.1 Släckmedel

Vid brand: Släck branden med brandbekämpningsmedel lämpligt för vanligt brännbart material såsom vatten eller skum.

### 5.2 Särskilda faror som ämnet eller blandningen kan medföra

Ingen inneboende i denna produkt.

### Farliga sönderdelnings- eller biprodukter

#### Ämne

Formaldehyd

Kolmonoxid

Koldioxid

Irriterande gaser eller ångor

#### Betingelser

Vid förbränning

Vid förbränning

Vid förbränning

Vid förbränning

### 5.3 Råd till brandbekämpningspersonal

Inga ovanliga brand- eller explosionsrisker förväntas.

## Avsnitt 6: Åtgärder vid oavsiktliga utsläpp

### 6.1 Personliga skyddsåtgärder, skyddsutrustning och åtgärder vid nödsituationer

Observera försiktighetsåtgärder i andra avsnitt.

### 6.2 Miljöskyddsåtgärder

Undvik utsläpp till miljön.

### 6.3 Metoder och material för inneslutning och sanering

Samla upp så mycket som möjligt av spillet. Placera i en förslutbar behållare. Torka upp rester. Förslut behållaren. Kassera uppsamlat material så snart som möjligt.

### 6.4 Hänvisning till andra avsnitt

Se avsnitt 8 och avsnitt 13 för mer information

## Avsnitt 7: Hantering och lagring

### 7.1 Försiktighetsmått för säker hantering

Använd ej i begränsat utrymme med dålig eller ingen ventilation. Ät inte, drick inte och rök inte när du använder produkten. Tvätta grundligt efter användning.

### 7.2 Förhållanden för säker lagring, inklusive eventuell oförenlighet

Skyddas från solljus. Förvaras inte i stark värme.

### 7.3 Specifik slutanvändning

Se information i avsnitt 7.1 och 7.2 för rekommendationer kring hantering och förvaring. Se avsnitt 8 för rekommendationer avseende begränsning av exponering samt personlig skyddsutrustning.

## Avsnitt 8: Begränsning av exponeringen/personligt skydd

### 8.1 Kontrollparametrar

#### Hygieniska gränsvärden

Det finns inget gränsvärde för något av de ämnen som angivits under avsnitt 3 i detta säkerhetsdatablad.



## 8.2 Begränsning av exponeringen

### 8.2.1 Lämpliga tekniska kontrollåtgärder

Ej tillämbart.

### 8.2.2 Individuella skyddsåtgärder, t.ex. personlig skyddsutrustning

#### Ögon/ansiktsskydd

Skyddsglasögon av följande typ/typer rekommenderas: Skyddsglasögon med sidoskydd.

#### Hud/handskydd

Gör en exponeringsbedömning för att avgöra om det finns risk för hudkontakt. Välj vid behov ut och använd skyddshandskar och/eller skyddskläder för att förebygga hudkontakt. Konsultera tillverkare av skyddshandskar/skyddskläder för val av lämpligt material.

#### Andningsskydd

Vid normal användning förväntas inte luftföroreningar i sådan omfattning som kräver andningsskydd.

## Avsnitt 9: Fysikaliska och kemiska egenskaper

### 9.1 Information om grundläggande fysikaliska och kemiska egenskaper

Aggregationstillstånd	Fast ämne
Specifik fysikalisk form:	Pasta
Utseende/luft	Svag karaktäristisk lukt. Vitfärgad pasta.
Lukttröskel	<i>Inga data tillgängliga</i>
pH	<i>Inga data tillgängliga</i>
Kokpunkt/kokpunktsintervall	<i>Ej tillämpligt</i>
Smältpunkt	<i>Ej tillämpligt</i>
Brandfarlighet (fast form, gas)	Ej klassificerad
Explosiva egenskaper	Ej klassificerad
Oxiderande egenskaper	Ej klassificerad
Flampunkt	Ingen flampunkt
Självantändningstemperatur	<i>Inga data tillgängliga</i>
Undre brännbarhets-/explosionsgräns	<i>Inga data tillgängliga</i>
Övre brännbarhets-/explosionsgräns	<i>Inga data tillgängliga</i>
Ångtryck	<i>Inga data tillgängliga</i>
Relativ densitet	1,6 - 1,7 [Ref:vatten=1]
Löslighet i vatten	Försumbar
Löslighet, ej vatten	<i>Inga data tillgängliga</i>
Fördelningskoefficient: n-oktanol/vatten	<i>Inga data tillgängliga</i>
Avdunstningshastighet	<i>Ej tillämpligt</i>
Ångdensitet	<i>Inga data tillgängliga</i>
Sönderdelningstemperatur	<i>Inga data tillgängliga</i>
Viskositet	<i>Inga data tillgängliga</i>
Densitet	1,6 g/cm <sup>3</sup> - 1,7 g/cm <sup>3</sup>

### 9.2 Annan information

Flyktiga organiska föreningar	<i>Ej tillämpligt</i>
Flyktiga föreningar	<i>Ej tillämpligt</i>

Flyktiga organiska föreningar (utom undantagna lösningsmedel och vatten. US std)

*Ej tillämpligt*

## Avsnitt 10: Stabilitet och reaktivitet

### 10.1 Reaktivitet

Denna produkt kan vara reaktiv med vissa ämnen under vissa omständigheter - se övriga rubriker i detta avsnitt.

### 10.2 Kemisk stabilitet

Stabil.

### 10.3 Risken för farliga reaktioner

Farlig polymerisation sker ej

### 10.4 Förhållanden som ska undvikas

Värme

### 10.5 Oförenliga material

Inga kända

### 10.6 Farliga sönderdelningsprodukter

#### Ämne

Inga kända.

#### Betingelser

Se avsnitt 5.2 för farliga sönderdelningsprodukter vid förbränning.

## Avsnitt 11: Toxikologisk information

Nedanstående information överensstämmer inte alltid helt med produktens klassificering i avsnitt 2 i de fall då det finns en av myndighet fastställd ämnesklassificering. Dessutom avspeglas inte nödvändigtvis enskilda beståndsdelars toxikologiska data i produktens klassificering och/eller i symptom vid exponering, eftersom en beståndsdel kan ingå i halt som understiger gränsen för klassificering av blandningen, en beståndsdel är kanske inte tillgänglig i produkten på sådant sätt att exponering kan ske, eller så är viss tox.data inte relevant för produkten i sin helhet.

### 11.1 Information om de toxikologiska effekterna

#### Symptom och tecken på exponering

Baserat på testdata och/eller information om ingående beståndsdelar, så kan denna produkt ge följande hälsoeffekter:

#### **Inandning**

Inga skadliga hälsoeffekter förväntas vid inandning.

#### **Hudkontakt**

Kontakt med huden under produktens användning förväntas inte ge någon betydande irritation.

#### **Ögonkontakt**

Ögonkontakt vid användning av produkten förväntas ej orsaka nämnvärd irritation.

#### **Förtäring**

Kan vara skadligt vid förtäring. Irritation i mag/tarmkanalen: symptom kan vara magsmärtor, upprörd mage, illamående, kräkning och diarré.

## Toxikologiska data

### Akut toxicitet

Namn	Expo.väg	Art	Värde
Produkten	Förtäring		Data ej tillgänglig eller otillräcklig för klassificering; beräknad ATE3 752,6 mg/kg
Nefelinsyenit	Förtäring		LD50 beräknad att vara 2 000 - 5 000 mg/kg
Vinylsiloxaner och vinylsilikoner	Dermal	Kanin	LD50 > 15 440 mg/kg
Vinylsiloxaner och vinylsilikoner	Förtäring	Råtta	LD50 > 15 440 mg/kg
Poly(dimetylsiloxan)	Dermal	Kanin	LD50 > 19 400 mg/kg
Poly(dimetylsiloxan)	Förtäring	Råtta	LD50 > 17 000 mg/kg
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	Dermal	Kanin	LD50 > 5 000 mg/kg
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	Inandning-damm/dimma (4 h)	Råtta	LC50 > 0,691 mg/l
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	Förtäring	Råtta	LD50 > 5 110 mg/kg

ATE=uppskattad akut toxicitet (acute toxicity estimate)

### Frätande/irriterande på huden

Namn	Art	Värde
Nefelinsyenit		Data ej tillgänglig eller otillräcklig för klassificering.
Vinylsiloxaner och vinylsilikoner		Data ej tillgänglig eller otillräcklig för klassificering.
Poly(dimetylsiloxan)	Kanin	Ingen signifikant irritation
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	Kanin	Ingen signifikant irritation

### Allvarlig ögonskada/ögonirritation

Namn	Art	Värde
Nefelinsyenit		Data ej tillgänglig eller otillräcklig för klassificering.
Vinylsiloxaner och vinylsilikoner		Data ej tillgänglig eller otillräcklig för klassificering.
Poly(dimetylsiloxan)	Kanin	Ingen signifikant irritation
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	Kanin	Ingen signifikant irritation

### Hudsensibilisering

Namn	Art	Värde
Nefelinsyenit		Data ej tillgänglig eller otillräcklig för klassificering.
Vinylsiloxaner och vinylsilikoner		Data ej tillgänglig eller otillräcklig för klassificering.
Poly(dimetylsiloxan)		Data ej tillgänglig eller otillräcklig för klassificering.
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	Human och djur	Ej sensibiliserande

### Luftvägssensibilisering

Namn	Art	Värde
Nefelinsyenit		Data ej tillgänglig eller otillräcklig för klassificering.
Vinylsiloxaner och vinylsilikoner		Data ej tillgänglig eller otillräcklig för klassificering.
Poly(dimetylsiloxan)		Data ej tillgänglig eller otillräcklig för klassificering.
Dimetylsiloxan, reaktionsprodukt med kiseldioxid		Data ej tillgänglig eller otillräcklig för klassificering.

**Mutagenitet i könsceller**

Namn	Expo.väg	Värde
Nefelinsyenit		Data ej tillgänglig eller otillräcklig för klassificering.
Vinylsiloxaner och vinylsilikoner		Data ej tillgänglig eller otillräcklig för klassificering.
Poly(dimetylsiloxan)		Data ej tillgänglig eller otillräcklig för klassificering.
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	In vitro	Ej mutagen

**Cancerogenitet**

Namn	Expo.väg	Art	Värde
Nefelinsyenit			Data ej tillgänglig eller otillräcklig för klassificering.
Vinylsiloxaner och vinylsilikoner			Data ej tillgänglig eller otillräcklig för klassificering.
Poly(dimetylsiloxan)			Data ej tillgänglig eller otillräcklig för klassificering.
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	Ej specificerade	Mus	Data är ej tillräcklig för klassificering

**Reproduktionstoxicitet**

**Reproduktions- och/eller utvecklingseffekter**

Namn	Expo.väg	Värde	Art	Resultat	Expo.tid
Nefelinsyenit		Data ej tillgänglig eller otillräcklig för klassificering.			
Vinylsiloxaner och vinylsilikoner		Data ej tillgänglig eller otillräcklig för klassificering.			
Poly(dimetylsiloxan)		Data ej tillgänglig eller otillräcklig för klassificering.			
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	Förtäring	Ej reproduktionstoxisk (honlig)	Råtta	NOAEL 509 mg/kg/day	1 generation
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	Förtäring	Ej reproduktionstoxisk (hanlig)	Råtta	NOAEL 497 mg/kg/day	1 generation
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	Förtäring	Ej utvecklingstoxisk	Råtta	NOAEL 1 350 mg/kg/day	under organbildning

**Målorg.**

**Specifik organtoxicitet - enstaka exponering**

Namn	Expo.väg	Målorg.	Värde	Art	Resultat	Expo.tid
Nefelinsyenit			Data ej tillgänglig eller otillräcklig för klassificering.			
Vinylsiloxaner och vinylsilikoner			Data ej tillgänglig eller otillräcklig för klassificering.			
Poly(dimetylsiloxan)			Data ej tillgänglig eller otillräcklig för klassificering.			
Dimetylsiloxan,			Data ej tillgänglig eller otillräcklig			

**3M™ ESPE™ IMPRINT™ 4 PENTA™ SUPER QUICK Catalyst**

reaktionsprodukt med kiseldioxid			för klassificering.			
----------------------------------	--	--	---------------------	--	--	--

**Specifik organototoxicitet - upprepad exponering**

Namn	Expo.väg	Målorg.	Värde	Art	Resultat	Expo.tid
Nefelinsyenit			Data ej tillgänglig eller otillräcklig för klassificering.			
Vinylsiloxaner och vinylsilikoner			Data ej tillgänglig eller otillräcklig för klassificering.			
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	Inandning	andningsorgan   silikos	All data är negativ	Människa	NOAEL Ej tillgänglig	yrkesmässig exponering

**Fara vid aspiration**

Namn	Värde
Nefelinsyenit	Ingen fara vid aspiration
Vinylsiloxaner och vinylsilikoner	Ingen fara vid aspiration
Poly(dimetylsiloxan)	Ingen fara vid aspiration
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	Ingen fara vid aspiration

Vid frågor som gäller den toxikologiska informationen i detta SDB, vänligen se kontaktuppgifter på första sidan.

**Avsnitt 12: Ekologisk information**

Nedanstående information överensstämmer inte alltid helt med produktens klassificering i avsnitt 2 i de fall då det finns en av myndighet fastställd ämnesklassificering. Ytterligare information som ligger till grund för produktens klassificering i avsnitt 2 kan lämnas vid förfrågan. Information om en beståndsdelns uppträdande och effekt i miljön avspeglas dessutom kanske inte i detta avsnitt om ämnet ingår i en halt som är under gränsen för klassificering av blandningen, eller om ämnet inte förväntas vara tillgängligt för exponering eller om data inte bedöms som relevant för produkten i sin helhet.

**12.1 Toxicitet**

Inga testdata tillgängliga för produkten

Produkt/ämne	Cas-nr	Organism	Typ	Exponering	Slutpunkt för testet	Resultat
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	67762-90-7		Data ej tillgänglig eller otillräcklig för klassificering.			
Nefelinsyenit	37244-96-5		Data ej tillgänglig eller otillräcklig för klassificering.			
Poly(dimetylsiloxan)	63148-62-9		Data ej tillgänglig eller otillräcklig för klassificering.			
Vinylsiloxaner och	68083-19-2		Data ej tillgänglig eller			

**3MT™ ESPE™ IMPRINT™ 4 PENTA™ SUPER QUICK Catalyst**

vinylsilikoner			otillräcklig för klassificering.			
----------------	--	--	----------------------------------	--	--	--

**12.2 Persistens och nedbrytbarhet**

Inga testdata tillgängliga

Produkt/ämne	Cas-nr	Typ av test	Varaktighet	Typ av studie	Resultat	Protokoll
Poly(dimetylsiloxan)	63148-62-9	Data ej tillgänglig eller otillräcklig för klassificering.	N/A	N/A	N/A	N/A
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	67762-90-7	Data ej tillgänglig eller otillräcklig för klassificering.	N/A	N/A	N/A	N/A
Nefelinsyenit	37244-96-5	Data ej tillgänglig eller otillräcklig för klassificering.	N/A	N/A	N/A	N/A
Vinylsiloxaner och vinylsilikoner	68083-19-2	Data ej tillgänglig eller otillräcklig för klassificering.	N/A	N/A	N/A	N/A

**12.3 Bioackumuleringsförmåga**

Inga testdata tillgängliga

Produkt/ämne	Cas-nr	Typ av test	Varaktighet	Typ av studie	Resultat	Protokoll
Vinylsiloxaner och vinylsilikoner	68083-19-2	Data ej tillgänglig eller otillräcklig för klassificering.	N/A	N/A	N/A	N/A
Poly(dimetylsiloxan)	63148-62-9	Data ej tillgänglig eller otillräcklig för klassificering.	N/A	N/A	N/A	N/A
Nefelinsyenit	37244-96-5	Data ej tillgänglig eller otillräcklig för klassificering.	N/A	N/A	N/A	N/A
Dimetylsiloxan, reaktionsprodukt med kiseldioxid	67762-90-7	Data ej tillgänglig eller otillräcklig för klassificering.	N/A	N/A	N/A	N/A

**12.4 Rörligheten i jord**

Kontakta tillverkaren för mer information

**12.5 Resultat av PBT- och vPvB-bedömningen**

Ingen information tillgänglig, kontakta tillverkaren för mer detaljer.

**12.6 Andra skadliga effekter**

Ingen information tillgänglig

## **Avsnitt 13: Avfallshantering**

### **13.1 Avfallsbehandlingsmetoder**

Innehåll/förpackning avfallshandteras i enlighet med lokala/regionala/nationella/internationella regler.

Avfallskoden (EWC) baseras på vilken källa som givit upphov till avfallet. För bestämning av lämplig avfallskod i varje enskilt fall se Avfallsförordningen (SFS 2011:927 med ändringar) bilaga 4. Säkerställ även att eventuella ytterligare nationella och/eller regionala krav efterlevs. Samverka endast med godkända avfallshämtare.

#### **Avfallskod (produkt i överlåtet skick)**

18 01 06\* Kemikalier som består av eller som innehåller farliga ämnen

#### **Avfallskod (produktförpackning efter användning)**

18 01 07 Andra kemikalier än de som anges i 18 01 06

### **Förpackningsmaterial**

3M Svenska AB är anslutet till FTI (Förpacknings- och tidningsinsamlingen). Kunder kan därför lämna våra tomma förpackningar utan kostnad. För information om närmaste lämningsställe ring 0200-880310. Förpackningar som innehållit kemiska produkter ska vara väl tömda och dropporra. Undantag är förpackningar med symbol T+, T, N eller R52, vilka istället ska lämnas som farligt avfall.

## **Avsnitt 14: Transportinformation**

ADR/IMDG/IATA: Ej begränsat för transport

## **Avsnitt 15: Gällande föreskrifter**

### **15.1. Föreskrifter/lagstiftning om ämnet eller blandningen när det gäller säkerhet, hälsa och miljö**

#### **Status i globala kemikalieregister**

Kontakta 3M för mer information. Produktens beståndsdelar möter kraven i China "Measures on Environmental Management of New Chemical Substance". Vissa begränsningar kan förekomma. Kontakta 3M vid behov av närmare information. Produktens beståndsdelar möter CEPA:s krav på New Substance Notification.

### **15.2. Kemikaliesäkerhetsbedömning**

Ej tillämpligt

## **Avsnitt 16: Annan information**

### **Information om uppdateringar**

Ingen information om uppdateringar finns tillgänglig.

Informationen i detta säkerhetsdatablad är baserad på vår erfarenhet och är, så vitt vi känner till, korrekt vid tidpunkten för dess publicering, men vi åtar oss inget ansvar för någon ekonomisk, sak- eller personskada som uppstår till följd av användning av informationen (med förbehåll för vad som är föreskrivet i lag). Informationen skall inte tillämpas i fråga om sådan användning som inte anges i detta säkerhetsdatablad eller användning av produkten i kombination med andra material. Av dessa skäl är det viktigt att kunder genomför egna tester för att fastställa att produkten passar det tilltänkta användningsområdet.

Se [www.3M.se/sdb](http://www.3M.se/sdb) för 3M Svenska AB:s säkerhetsdatablad.

