

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 03.04.2024

Version number 2 (replaces version 1)

Revision: 03.04.2024


SECTION 1: Identification of the substance/mixture and of the company/undertaking

- **1.1 Product identifier**
- **Trade name: Ivotion Bond Polymer**
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**
No further relevant information available.
- **Application of the substance / the mixture** Bonding material for the cementation of denture teeth
- **1.3 Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan
PRINCIPALITY OF LIECHTENSTEIN

- **Tel:** +423 235 35 35
Fax: +423 235 33 60
- **Further information obtainable from:**
Regulatory Affairs
sds@ivoclar.com
- **1.4 Emergency telephone number:** +423 / 235 33 13 (Ivoclar Vivadent AG, 9494 Schaan, Liechtenstein)

SECTION 2: Hazards identification

- **2.1 Classification of the substance or mixture**
 - **Classification according to Regulation (EC) No 1272/2008**
Skin Sens. 1 H317 May cause an allergic skin reaction.

 - **2.2 Label elements**
 - **Labelling according to Regulation (EC) No 1272/2008**
The product is classified and labelled according to the GB CLP regulation.
 - **Hazard pictograms**
- 

GHS07
- **Signal word** Warning
 - **Hazard-determining components of labelling:**
dibenzoyl peroxide
 - **Hazard statements**
H317 May cause an allergic skin reaction.
 - **Precautionary statements**
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P302+P352 IF ON SKIN: Wash with plenty of soap and water.
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P362+P364 Take off contaminated clothing and wash it before reuse.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.
 - **2.3 Other hazards**
Dust generation
Particular danger of slipping on leaked/spilled product.
 - **Results of PBT and vPvB assessment**
 - **PBT:** Not applicable.

(Contd. on page 2)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 03.04.2024

Version number 2 (replaces version 1)

Revision: 03.04.2024

Trade name: Ivotion Bond Polymer

· **vPvB:** Not applicable.

(Contd. of page 1)

SECTION 3: Composition/information on ingredients

· **3.2 Mixtures**

· **Description:** > 95% polymethylmethacrylate

· **Dangerous components:**

CAS: 94-36-0	dibenzoyl peroxide	1-<2.5%
EINECS: 202-327-6	Org. Perox. B, H241; Eye Irrit. 2, H319; Skin Sens. 1, H317	
Reg.nr.: 01-2119511472-50-XXXX		

· **Additional information:** For the wording of the listed hazard phrases refer to section 16.

SECTION 4: First aid measures

· **4.1 Description of first aid measures**

· **General information:** No special measures required.

· **After inhalation:** Supply fresh air and to be sure call for a doctor.

· **After skin contact:**

Rinse with water.

If skin irritation continues, consult a doctor.

· **After eye contact:** Rinse opened eye for several minutes under running water. Then consult a doctor.

· **After swallowing:**

Rinse out mouth and then drink plenty of water.

If symptoms persist consult doctor.

· **4.2 Most important symptoms and effects, both acute and delayed** No further relevant information available.

· **4.3 Indication of any immediate medical attention and special treatment needed**

No further relevant information available.

SECTION 5: Firefighting measures

· **5.1 Extinguishing media**

· **Suitable extinguishing agents:** Use fire extinguishing methods suitable to surrounding conditions.

· **5.2 Special hazards arising from the substance or mixture** No further relevant information available.

· **5.3 Advice for firefighters**

· **Protective equipment:** No special measures required.

SECTION 6: Accidental release measures

· **6.1 Personal precautions, protective equipment and emergency procedures**

Wear protective equipment. Keep unprotected persons away.

· **6.2 Environmental precautions:**

No special measures required.

Do not allow to enter sewers/ surface or ground water.

· **6.3 Methods and material for containment and cleaning up:** Pick up mechanically.

· **6.4 Reference to other sections**

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

GB

(Contd. on page 3)

Safety data sheet according to 1907/2006/EC, Article 31

Printing date 03.04.2024

Version number 2 (replaces version 1)

Revision: 03.04.2024

Trade name: Ivotion Bond Polymer

(Contd. of page 2)

SECTION 7: Handling and storage

- **7.1 Precautions for safe handling**
Only adequately trained personnel should handle this product.
For use in dentistry only.
- **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** Store only in the original receptacle.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Protect from heat and direct sunlight.
- **7.3 Specific end use(s)** No further relevant information available.

SECTION 8: Exposure controls/personal protection

- **8.1 Control parameters**
 - **Ingredients with limit values that require monitoring at the workplace:**
- | | |
|-------------|--------------------------------------|
| CAS: | 94-36-0 dibenzoyl peroxide |
| WEL | Long-term value: 5 mg/m ³ |
- **Regulatory information** WEL: EH40/2011
 - **Additional information:** The lists valid during the making were used as basis.

- **8.2 Exposure controls**
- **Appropriate engineering controls** No further data; see section 7.
- **Individual protection measures, such as personal protective equipment**
- **General protective and hygienic measures:**
Usual hygienic measures for dental practice and dental laboratories.
Keep away from foodstuffs, beverages and feed.
Wash hands before breaks and at the end of work.
Avoid contact with the skin.
- **Respiratory protection:** Use respiratory protective device against the effects of dust.
- **Hand protection**



Protective gloves (EN 374)

After use of gloves apply skin-cleaning agents and skin cosmetics.

- **Material of gloves**
Butyl rubber, BR
Natural rubber, NR
Nitrile rubber, NBR
Fluorocarbon rubber (Viton)
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Penetration time of glove material**
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

(Contd. on page 4)

GB

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 03.04.2024

Version number 2 (replaces version 1)

Revision: 03.04.2024

Trade name: **Ivotion Bond Polymer**

(Contd. of page 3)

· Eye/face protection



Tightly sealed goggles (EN 166)

SECTION 9: Physical and chemical properties

· 9.1 Information on basic physical and chemical properties

· General Information

· Colour:	According to product specification
· Odour:	Characteristic
· Odour threshold:	Not determined.
· Melting point/freezing point:	150 °C
· Boiling point or initial boiling point and boiling range	Undetermined.
· Flammability	Not determined.
· Flash point:	250 °C
· pH	Not applicable.
· Viscosity:	
· Kinematic viscosity	Not applicable.
· Dynamic:	Not applicable.
· Solubility	
· water:	Nearly insoluble.
· Partition coefficient n-octanol/water (log value)	Not determined.
· Vapour pressure:	Not applicable.
· Density and/or relative density	
· Density at 20 °C:	1.2 g/cm ³
· Relative density	Not determined.
· Vapour density	Not applicable.

· 9.2 Other information

· Appearance:	
· Form:	Powder
· Important information on protection of health and environment, and on safety.	
· Ignition temperature:	Product is not selfigniting.
· Explosive properties:	Product does not present an explosion hazard.
· Change in condition	
· Evaporation rate	Not applicable.

· Information with regard to physical hazard classes

· Explosives	Void
· Flammable gases	Void
· Aerosols	Void
· Oxidising gases	Void
· Gases under pressure	Void
· Flammable liquids	Void
· Flammable solids	Void
· Self-reactive substances and mixtures	Void
· Pyrophoric liquids	Void
· Pyrophoric solids	Void
· Self-heating substances and mixtures	Void

(Contd. on page 5)

GB

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 03.04.2024

Version number 2 (replaces version 1)

Revision: 03.04.2024

Trade name: Ivotion Bond Polymer

(Contd. of page 4)

- | | |
|------------------------------------------------------------------------------------|------|
| · Substances and mixtures, which emit flammable gases in contact with water | Void |
| · Oxidising liquids | Void |
| · Oxidising solids | Void |
| · Organic peroxides | Void |
| · Corrosive to metals | Void |
| · Desensitised explosives | Void |

SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability** Stable under normal handling and storage conditions.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions** No dangerous reactions known.
- **10.4 Conditions to avoid** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** None under normal conditions of storage and use.

SECTION 11: Toxicological information

- **11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**
- **Acute toxicity** Based on available data, the classification criteria are not met.
- **Skin corrosion/irritation** Based on available data, the classification criteria are not met.
- **Serious eye damage/irritation** Based on available data, the classification criteria are not met.
- **Respiratory or skin sensitisation**
May cause an allergic skin reaction.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.
- **Additional toxicological information:** No further relevant information available.
- **11.2 Information on other hazards**

· **Endocrine disrupting properties**

None of the ingredients is listed.

SECTION 12: Ecological information

- **12.1 Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **12.2 Persistence and degradability** No further relevant information available.
- **12.3 Bioaccumulative potential** No further relevant information available.
- **12.4 Mobility in soil** No further relevant information available.
- **12.5 Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **12.6 Endocrine disrupting properties**
The product does not contain substances with endocrine disrupting properties.

(Contd. on page 6)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 03.04.2024

Version number 2 (replaces version 1)

Revision: 03.04.2024

Trade name: Ivotion Bond Polymer

(Contd. of page 5)

- **12.7 Other adverse effects**
- **Additional ecological information:**
- **General notes:**
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
- **Recommendation**
Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

SECTION 14: Transport information

- | | |
|-----------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------|
| · 14.1 UN number or ID number
· ADR/RID/ADN, ADN, IMDG, IATA | Void |
| · 14.2 UN proper shipping name
· ADR/RID/ADN, ADN, IMDG, IATA | Void |
| · 14.3 Transport hazard class(es)
· ADR/RID/ADN, ADN, IMDG, IATA
· Class | Void |
| · 14.4 Packing group
· ADR/RID/ADN, IMDG, IATA | Void |
| · 14.5 Environmental hazards:
· Marine pollutant: | No |
| · 14.6 Special precautions for user | Not applicable. |
| · 14.7 Maritime transport in bulk according to IMO instruments | Not applicable. |
| · Transport/Additional information: | Product is not classified as a dangerous good for transport (ADR, IMDG, IATA). |
| · UN "Model Regulation": | Void |

SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **National regulations:**
- **Other regulations, limitations and prohibitive regulations**
The product is a medical device according to Regulation 2017/745 on Medical Devices (MDR).

(Contd. on page 7)

GB

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 03.04.2024

Version number 2 (replaces version 1)

Revision: 03.04.2024

Trade name: Ivotion Bond Polymer

(Contd. of page 6)

· **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Relevant phrases**

H241 Heating may cause a fire or explosion.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

· **Abbreviations and acronyms:**

IMO: International Maritime Organization

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

Org. Perox. B: Organic peroxides – Type B

Eye Irrit. 2: Serious eye damage/eye irritation – Category 2

Skin Sens. 1: Skin sensitisation – Category 1

· *** Data compared to the previous version altered.**

GB

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 03.04.2024



Version number 2 (replaces version 1)

Revision: 03.04.2024

SECTION 1: Identification of the substance/mixture and of the company/undertaking

- **1.1 Product identifier**
- **Trade name: Ivotion Bond Monomer**
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**
No further relevant information available.
- **Application of the substance / the mixture** Bonding material for the cementation of denture teeth
- **1.3 Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan
PRINCIPALITY OF LIECHTENSTEIN
- **Tel:** +423 235 35 35
Fax: +423 235 33 60
- **Further information obtainable from:**
Regulatory Affairs
sds@ivoclar.com
- **1.4 Emergency telephone number:** +423 / 235 33 13 (Ivoclar Vivadent AG, 9494 Schaan, Liechtenstein)

SECTION 2: Hazards identification

- **2.1 Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008**
Flam. Liq. 2 H225 Highly flammable liquid and vapour.
Skin Irrit. 2 H315 Causes skin irritation.
Skin Sens. 1 H317 May cause an allergic skin reaction.
STOT SE 3 H335 May cause respiratory irritation.
- **2.2 Label elements**
- **Labelling according to Regulation (EC) No 1272/2008**
The product is classified and labelled according to the GB CLP regulation.
- **Hazard pictograms**

GHS02 GHS07
- **Signal word** Danger
- **Hazard-determining components of labelling:**
methyl methacrylate
1,4-butanediol dimethacrylate
- **Hazard statements**
H225 Highly flammable liquid and vapour.
H315 Causes skin irritation.
H317 May cause an allergic skin reaction.
H335 May cause respiratory irritation.
- **Precautionary statements**
P240 Ground and bond container and receiving equipment.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P302+P352 IF ON SKIN: Wash with plenty of water.
P312 Call a POISON CENTER/doctor if you feel unwell.

(Contd. on page 2)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 03.04.2024

Version number 2 (replaces version 1)

Revision: 03.04.2024

Trade name: Ivotion Bond Monomer

(Contd. of page 1)

P403+P235 Store in a well-ventilated place. Keep cool.

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Additional information:

Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

2.3 Other hazards

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

SECTION 3: Composition/information on ingredients

3.2 Mixtures

Description: Mixture of substances listed below with nonhazardous additions.

Dangerous components:

CAS: 80-62-6 EINECS: 201-297-1 Reg.nr.: 01-2119452498-28-XXXX	<i>methyl methacrylate</i> Flam. Liq. 2, H225; Skin Irrit. 2, H315; Skin Sens. 1, H317; STOT SE 3, H335	50-100%
CAS: 2082-81-7 EINECS: 218-218-1 Reg.nr.: 01-2119967415-30-XXXX	<i>1,4-butanediol dimethacrylate</i> Skin Sens. 1, H317	2.5-<10%

Additional information: For the wording of the listed hazard phrases refer to section 16.

SECTION 4: First aid measures

4.1 Description of first aid measures

General information: Immediately remove any clothing soiled by the product.

After inhalation:

Supply fresh air; consult doctor in case of complaints.

In case of unconsciousness place patient stably in side position for transportation.

After skin contact:

Immediately rinse with water.

If skin irritation continues, consult a doctor.

After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.

After swallowing:

Rinse out mouth and then drink plenty of water.

Do not induce vomiting; call for medical help immediately.

4.2 Most important symptoms and effects, both acute and delayed No further relevant information available.

4.3 Indication of any immediate medical attention and special treatment needed

No further relevant information available.

SECTION 5: Firefighting measures

5.1 Extinguishing media

Suitable extinguishing agents: CO₂, sand, extinguishing powder. Do not use water.

For safety reasons unsuitable extinguishing agents: Water with full jet

5.2 Special hazards arising from the substance or mixture No further relevant information available.

5.3 Advice for firefighters

Protective equipment: No special measures required.

(Contd. on page 3)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 03.04.2024

Version number 2 (replaces version 1)

Revision: 03.04.2024

Trade name: Ivotion Bond Monomer

(Contd. of page 2)

· **Additional information** Cool endangered receptacles with water spray.

SECTION 6: Accidental release measures

- **6.1 Personal precautions, protective equipment and emergency procedures**
Wear protective equipment. Keep unprotected persons away.
- **6.2 Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **6.3 Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Ensure adequate ventilation.
- **6.4 Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

SECTION 7: Handling and storage

- **7.1 Precautions for safe handling**
Only adequately trained personnel should handle this product.
Ensure good ventilation/exhaustion at the workplace.
For use in dentistry only.
- **Information about fire - and explosion protection:**
Keep ignition sources away - Do not smoke.
Protect against electrostatic charges.
- **7.2 Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:**
Store in a cool location.
Store only in the original receptacle.
- **Information about storage in one common storage facility:** Store away from oxidising agents.
- **Further information about storage conditions:**
Keep container tightly sealed.
Store receptacle in a well ventilated area.
Protect from heat and direct sunlight.
- **7.3 Specific end use(s)** No further relevant information available.

SECTION 8: Exposure controls/personal protection

· 8.1 Control parameters

· **Ingredients with limit values that require monitoring at the workplace:**

CAS: 80-62-6 methyl methacrylate

WEL	Short-term value: 416 mg/m ³ , 100 ppm
	Long-term value: 208 mg/m ³ , 50 ppm

- **Regulatory information** WEL: EH40/2011
- **Additional information:** The lists valid during the making were used as basis.

· 8.2 Exposure controls

- **Appropriate engineering controls** No further data; see section 7.
- **Individual protection measures, such as personal protective equipment**
- **General protective and hygienic measures:**
Usual hygienic measures for dental practice and dental laboratories.
Keep away from foodstuffs, beverages and feed.

(Contd. on page 4)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 03.04.2024

Version number 2 (replaces version 1)

Revision: 03.04.2024

Trade name: Ivotion Bond Monomer

(Contd. of page 3)

Immediately remove all soiled and contaminated clothing

Wash hands before breaks and at the end of work.

Avoid contact with the eyes and skin.

Do not inhale gases / fumes / aerosols.

· **Respiratory protection:**

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.

· **Recommended filter device for short term use:**

Filter A1

Filter A2

Filter A3

· **Hand protection**



Protective gloves (EN 374)

After use of gloves apply skin-cleaning agents and skin cosmetics.

· **Material of gloves**

Butyl rubber, BR

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Not suitable are gloves made of the following materials:**

Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

· **Eye/face protection**



Tightly sealed goggles (EN 166)

SECTION 9: Physical and chemical properties

· **9.1 Information on basic physical and chemical properties**

· **General Information**

· Colour:	<i>Colourless</i>
· Odour:	<i>Pungent</i>
· Odour threshold:	<i>Not determined.</i>
· Melting point/freezing point:	<i>-48 °C</i>
· Boiling point or initial boiling point and boiling range	<i>101 °C</i>
· Lower and upper explosion limit	
· Lower:	<i>2.1 Vol %</i>
· Upper:	<i>12.5 Vol %</i>
· Flash point:	<i>10 °C</i>
· Auto-ignition temperature:	<i>430 °C</i>
· pH	<i>Not determined.</i>
· Viscosity:	
· Kinematic viscosity	<i>Not determined.</i>
· Dynamic:	<i>Not determined.</i>
· Solubility	
· water at 20 °C:	<i>1.6 g/l</i>

(Contd. on page 5)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 03.04.2024

Version number 2 (replaces version 1)

Revision: 03.04.2024

Trade name: Ivotion Bond Monomer

(Contd. of page 4)

· Partition coefficient n-octanol/water (log value)	Not determined.
· Vapour pressure at 20 °C:	47 hPa
· Density and/or relative density	
· Density at 20 °C:	0.943 g/cm ³
· Relative density	Not determined.
· Vapour density	Not determined.
· 9.2 Other information	No further relevant information available.
· Appearance:	
· Form:	Fluid
· Important information on protection of health and environment, and on safety.	
· Ignition temperature:	Product is not selfigniting.
· Explosive properties:	Product is not explosive. However, formation of explosive air/vapour mixtures are possible.
· Change in condition	
· Evaporation rate	Not determined.
· Information with regard to physical hazard classes	
· Explosives	Void
· Flammable gases	Void
· Aerosols	Void
· Oxidising gases	Void
· Gases under pressure	Void
· Flammable liquids	Highly flammable liquid and vapour.
· Flammable solids	Void
· Self-reactive substances and mixtures	Void
· Pyrophoric liquids	Void
· Pyrophoric solids	Void
· Self-heating substances and mixtures	Void
· Substances and mixtures, which emit flammable gases in contact with water	Void
· Oxidising liquids	Void
· Oxidising solids	Void
· Organic peroxides	Void
· Corrosive to metals	Void
· Desensitised explosives	Void

SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability** Stable under normal handling and storage conditions.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions**
Forms explosive gas mixture with air.
Reacts with strong oxidising agents.
Exothermic polymerisation.
- **10.4 Conditions to avoid** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** None under normal conditions of storage and use.

GB

(Contd. on page 6)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 03.04.2024

Version number 2 (replaces version 1)

Revision: 03.04.2024

Trade name: **Ivotion Bond Monomer**

(Contd. of page 5)

SECTION 11: Toxicological information

· **11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**

· **Acute toxicity** Based on available data, the classification criteria are not met.

· **LD/LC50 values relevant for classification:**

CAS: 80-62-6 methyl methacrylate

Oral	LD50	7872 mg/kg (rat)
------	------	------------------

· **Skin corrosion/irritation**

Causes skin irritation.

· **Serious eye damage/irritation** Based on available data, the classification criteria are not met.

· **Respiratory or skin sensitisation**

May cause an allergic skin reaction.

· **Germ cell mutagenicity** Based on available data, the classification criteria are not met.

· **Carcinogenicity** Based on available data, the classification criteria are not met.

· **Reproductive toxicity** Based on available data, the classification criteria are not met.

· **STOT-single exposure**

May cause respiratory irritation.

· **STOT-repeated exposure** Based on available data, the classification criteria are not met.

· **Aspiration hazard** Based on available data, the classification criteria are not met.

· **Additional toxicological information:** No further relevant information available.

· **11.2 Information on other hazards**

· **Endocrine disrupting properties**

None of the ingredients is listed.

SECTION 12: Ecological information

· **12.1 Toxicity**

· **Aquatic toxicity:** No further relevant information available.

· **12.2 Persistence and degradability** No further relevant information available.

· **12.3 Bioaccumulative potential** No further relevant information available.

· **12.4 Mobility in soil** No further relevant information available.

· **12.5 Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

· **vPvB:** Not applicable.

· **12.6 Endocrine disrupting properties**

The product does not contain substances with endocrine disrupting properties.

· **12.7 Other adverse effects** No further relevant information available.

· **Additional ecological information:**

· **General notes:**

Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

SECTION 13: Disposal considerations

· **13.1 Waste treatment methods**

· **Recommendation**

Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.

(Contd. on page 7)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 03.04.2024

Version number 2 (replaces version 1)



Revision: 03.04.2024

Trade name: **Ivotion Bond Monomer**

(Contd. of page 6)

- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

SECTION 14: Transport information

· 14.1 UN number or ID number · ADR/RID/ADN, IMDG, IATA	UN1247
· 14.2 UN proper shipping name · ADR/RID/ADN · IMDG, IATA	1247 METHYL METHACRYLATE MONOMER, STABILIZED solution METHYL METHACRYLATE MONOMER, STABILIZED solution
· 14.3 Transport hazard class(es) · ADR/RID/ADN	
	
· Class · Label	3 (F1) Flammable liquids. 3
· IMDG, IATA	
	
· Class · Label	3 Flammable liquids. 3
· 14.4 Packing group · ADR/RID/ADN, IMDG, IATA	II
· 14.5 Environmental hazards: · Marine pollutant:	No
· 14.6 Special precautions for user · Hazard identification number (Kemler code): · EMS Number: · Stowage Category · Stowage Code	Warning: Flammable liquids. 339 F-E,S-D B SW2 Clear of living quarters.
· 14.7 Maritime transport in bulk according to IMO instruments	Not applicable.
· Transport/Additional information:	
· ADR/RID/ADN · Limited quantities (LQ) · Excepted quantities (EQ)	1L Code: E2 Maximum net quantity per inner packaging: 30 ml Maximum net quantity per outer packaging: 500 ml
· Transport category · Tunnel restriction code	2 D/E

(Contd. on page 8)

GB

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 03.04.2024

Version number 2 (replaces version 1)

Revision: 03.04.2024

Trade name: Ivotion Bond Monomer

(Contd. of page 7)

· IMDG	
· Limited quantities (LQ)	1L
· Excepted quantities (EQ)	Code: E2 Maximum net quantity per inner packaging: 30 ml Maximum net quantity per outer packaging: 500 ml
· UN "Model Regulation":	UN 1247 METHYL METHACRYLATE MONOMER, STABILIZED SOLUTION, 3, II

SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **Seveso category P5c FLAMMABLE LIQUIDS**
- **Qualifying quantity (tonnes) for the application of lower-tier requirements** 5000 t
- **Qualifying quantity (tonnes) for the application of upper-tier requirements** 50000 t
- **National regulations:**
- **Other regulations, limitations and prohibitive regulations**
The product is a medical device according to Regulation 2017/745 on Medical Devices (MDR).
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Relevant phrases**
H225 Highly flammable liquid and vapour.
H315 Causes skin irritation.
H317 May cause an allergic skin reaction.
H335 May cause respiratory irritation.
- **Classification according to Regulation (EC) No 1272/2008**
The classification of the mixture is generally based on the calculation method using substance data according to Regulation (EC) No 1272/2008.
- **Abbreviations and acronyms:**
IMO: International Maritime Organization
ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)
RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)
IMDG: International Maritime Code for Dangerous Goods
IATA: International Air Transport Association
GHS: Globally Harmonised System of Classification and Labelling of Chemicals
EINECS: European Inventory of Existing Commercial Chemical Substances
ELINCS: European List of Notified Chemical Substances
CAS: Chemical Abstracts Service (division of the American Chemical Society)
LC50: Lethal concentration, 50 percent
LD50: Lethal dose, 50 percent
PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative
Flam. Liq. 2: Flammable liquids – Category 2
Skin Irrit. 2: Skin corrosion/irritation – Category 2
Skin Sens. 1: Skin sensitisation – Category 1
STOT SE 3: Specific target organ toxicity (single exposure) – Category 3
- *** Data compared to the previous version altered.**

**Sicherheitsdatenblatt
gemäß 1907/2006/EG, Artikel 31**

Druckdatum: 03.04.2024

Versionsnummer 2 (ersetzt Version 1)

überarbeitet am: 03.04.2024

ABSCHNITT 1: Bezeichnung des Stoffs beziehungsweise des Gemischs und des Unternehmens

- **1.1 Produktidentifikator**
- **Handelsname: Ivotion Bond Modelling Liquid**
- **1.2 Relevante identifizierte Verwendungen des Stoffs oder Gemischs und Verwendungen, von denen abgeraten wird**
Keine weiteren relevanten Informationen verfügbar.
- **Verwendung des Stoffs / des Gemisches** Herstellung von Zahnersatz
- **1.3 Einzelheiten zum Lieferanten, der das Sicherheitsdatenblatt bereitstellt**
- **Hersteller/Lieferant:**
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan
FÜRSTENTUM LIECHTENSTEIN

Tel: +423 235 35 35
Fax: +423 235 33 60
- **Auskunftgebender Bereich:**
Regulatory Affairs
sds@ivoclar.com
- **1.4 Notrufnummer:**
145 (Schweiz: Tox Info Suisse, 24-h-Notfallnummer)
+423 / 235 33 13 (Ivoclar Vivadent AG, 9494 Schaan, Liechtenstein)

ABSCHNITT 2: Mögliche Gefahren

- **2.1 Einstufung des Stoffs oder Gemischs**
- **Einstufung gemäß Verordnung (EG) Nr. 1272/2008**
Flam. Liq. 2 H225 Flüssigkeit und Dampf leicht entzündbar.
Skin Irrit. 2 H315 Verursacht Hautreizungen.
Skin Sens. 1 H317 Kann allergische Hautreaktionen verursachen.
STOT SE 3 H335 Kann die Atemwege reizen.
- **2.2 Kennzeichnungselemente**
- **Kennzeichnung gemäß Verordnung (EG) Nr. 1272/2008**
Das Produkt ist gemäß CLP-Verordnung eingestuft und gekennzeichnet.
- **Gefahrenpiktogramme**
 
GHS02 GHS07
- **Signalwort** Gefahr
- **Gefahrbestimmende Komponenten zur Etikettierung:**
Methylmethacrylat
1,4-Butandioldimethacrylat
N,N-Bis(2-hydroxyethyl)-p-toluidin
- **Gefahrenhinweise**
H225 Flüssigkeit und Dampf leicht entzündbar.
H315 Verursacht Hautreizungen.
H317 Kann allergische Hautreaktionen verursachen.
H335 Kann die Atemwege reizen.

(Fortsetzung auf Seite 2)

Sicherheitsdatenblatt gemäß 1907/2006/EG, Artikel 31

Druckdatum: 03.04.2024

Versionsnummer 2 (ersetzt Version 1)

überarbeitet am: 03.04.2024

Handelsname: **Ivotion Bond Modelling Liquid**

(Fortsetzung von Seite 1)

- **Sicherheitshinweise**
 - P240 Behälter und zu befüllende Anlage erden.
 - P280 Schutzhandschuhe/ Schutzkleidung/ Augenschutz/ Gesichtsschutz/ Gehörschutz tragen.
 - P302+P352 BEI BERÜHRUNG MIT DER HAUT: Mit viel Wasser und Seife waschen.
 - P312 Bei Unwohlsein GIFTINFORMATIONSZENTRUM/Arzt anrufen.
 - P403+P235 An einem gut belüfteten Ort aufbewahren. Kühl halten.
 - P501 Entsorgung des Inhalts / des Behälters gemäß den örtlichen / regionalen / nationalen/ internationalen Vorschriften.
- **Zusätzliche Angaben:**
Handelsübliche medizinische Handschuhe bieten keinen Schutz gegen den sensibilisierenden Effekt von Methacrylaten.
- **2.3 Sonstige Gefahren**
- **Ergebnisse der PBT- und vPvB-Beurteilung**
- **PBT:** Nicht anwendbar.
- **vPvB:** Nicht anwendbar.

ABSCHNITT 3: Zusammensetzung/Angaben zu Bestandteilen

- **3.2 Gemische**
- **Beschreibung:** Gemisch aus nachfolgend angeführten Stoffen mit ungefährlichen Beimengungen.

- **Gefährliche Inhaltsstoffe:**

CAS: 80-62-6 EINECS: 201-297-1 Reg.nr.: 01-2119452498-28-XXXX	Methylmethacrylat Flam. Liq. 2, H225; Skin Irrit. 2, H315; Skin Sens. 1, H317; STOT SE 3, H335	50-100%
CAS: 2082-81-7 EINECS: 218-218-1 Reg.nr.: 01-2119967415-30-XXXX	1,4-Butandioldimethacrylat Skin Sens. 1, H317	2,5-10%
CAS: 3077-12-1 EINECS: 221-359-1 Reg.nr.: 01-2120791684-40-xxxx	N,N-Bis(2-hydroxyethyl)-p-toluidin Eye Dam. 1, H318; Acute Tox. 4, H302; Skin Sens. 1, H317; Aquatic Chronic 3, H412	0,1-<1%

- **Zusätzliche Hinweise:** Der Wortlaut der angeführten Gefahrenhinweise ist dem Abschnitt 16 zu entnehmen.

ABSCHNITT 4: Erste-Hilfe-Maßnahmen

- **4.1 Beschreibung der Erste-Hilfe-Maßnahmen**
- **Allgemeine Hinweise:** Mit Produkt verunreinigte Kleidungsstücke unverzüglich entfernen.
- **Nach Einatmen:**
Frischlufzufuhr, bei Beschwerden Arzt aufsuchen.
Bei Bewusstlosigkeit Lagerung und Transport in stabiler Seitenlage.
- **Nach Hautkontakt:**
Sofort mit Wasser abwaschen.
Bei andauernder Hautreizung Arzt aufsuchen.
- **Nach Augenkontakt:**
Augen bei geöffnetem Lidspalt mehrere Minuten unter fließendem Wasser abspülen und Arzt konsultieren.
- **Nach Verschlucken:**
Mund ausspülen und reichlich Wasser nachtrinken.
Kein Erbrechen herbeiführen, sofort Arzthilfe zuziehen.
- **4.2 Wichtigste akute und verzögert auftretende Symptome und Wirkungen**
Keine weiteren relevanten Informationen verfügbar.

(Fortsetzung auf Seite 3)

Sicherheitsdatenblatt gemäß 1907/2006/EG, Artikel 31

Druckdatum: 03.04.2024

Versionsnummer 2 (ersetzt Version 1)

überarbeitet am: 03.04.2024

Handelsname: Ivotion Bond Modelling Liquid

(Fortsetzung von Seite 2)

- **4.3 Hinweise auf ärztliche Soforthilfe oder Spezialbehandlung**
Keine weiteren relevanten Informationen verfügbar.

ABSCHNITT 5: Maßnahmen zur Brandbekämpfung

- **5.1 Löschmittel**
- **Geeignete Löschmittel:** CO₂, Sand, Löschpulver. Kein Wasser verwenden.
- **Aus Sicherheitsgründen ungeeignete Löschmittel:** Wasser im Vollstrahl
- **5.2 Besondere vom Stoff oder Gemisch ausgehende Gefahren**
Keine weiteren relevanten Informationen verfügbar.
- **5.3 Hinweise für die Brandbekämpfung**
- **Besondere Schutzausrüstung:** Keine besonderen Maßnahmen erforderlich.
- **Weitere Angaben** Gefährdete Behälter mit Wassersprühstrahl kühlen.

ABSCHNITT 6: Maßnahmen bei unbeabsichtigter Freisetzung

- **6.1 Personenbezogene Vorsichtsmaßnahmen, Schutzausrüstungen und in Notfällen anzuwendende Verfahren**
Schutzausrüstung tragen. Ungeschützte Personen fernhalten.
- **6.2 Umweltschutzmaßnahmen:** Nicht in die Kanalisation/Oberflächenwasser/Grundwasser gelangen lassen.
- **6.3 Methoden und Material für Rückhaltung und Reinigung:**
Mit flüssigkeitsbindendem Material (Sand, Kieselgur, Säurebinder, Universalbinder, Sägemehl) aufnehmen.
Für ausreichende Lüftung sorgen.
- **6.4 Verweis auf andere Abschnitte**
Informationen zur sicheren Handhabung siehe Abschnitt 7.
Informationen zur persönlichen Schutzausrüstung siehe Abschnitt 8.
Informationen zur Entsorgung siehe Abschnitt 13.

ABSCHNITT 7: Handhabung und Lagerung

- **7.1 Schutzmaßnahmen zur sicheren Handhabung**
Handhabung dieses Produkts nur durch ausreichend geschultes Personal.
Für gute Belüftung/Absaugung am Arbeitsplatz sorgen.
Nur für zahnärztlichen oder zahntechnischen Gebrauch.
- **Hinweise zum Brand- und Explosionsschutz:**
Zündquellen fernhalten - nicht rauchen.
Maßnahmen gegen elektrostatische Aufladung treffen.
- **7.2 Bedingungen zur sicheren Lagerung unter Berücksichtigung von Unverträglichkeiten**
- **Lagerung:**
- **Anforderung an Lagerräume und Behälter:**
An einem kühlen Ort lagern.
Nur im Originalgebinde aufbewahren.
- **Zusammenlagerungshinweise:** Getrennt von Oxidationsmitteln aufbewahren.
- **Weitere Angaben zu den Lagerbedingungen:**
Behälter dicht geschlossen halten.
Behälter an einem gut gelüfteten Ort aufbewahren.
Vor Hitze und direkter Sonnenbestrahlung schützen.
- **Klassifizierung nach Betriebssicherheitsverordnung (BetrSichV):** Leichtentzündlich
- **7.3 Spezifische Endanwendungen** Keine weiteren relevanten Informationen verfügbar.

DE

(Fortsetzung auf Seite 4)

Sicherheitsdatenblatt gemäß 1907/2006/EG, Artikel 31

Druckdatum: 03.04.2024

Versionsnummer 2 (ersetzt Version 1)

überarbeitet am: 03.04.2024

Handelsname: Ivotion Bond Modelling Liquid

(Fortsetzung von Seite 3)

* ABSCHNITT 8: Begrenzung und Überwachung der Exposition/Persönliche Schutzausrüstungen

· 8.1 Zu überwachende Parameter

· Bestandteile mit arbeitsplatzbezogenen, zu überwachenden Grenzwerten:

CAS: 80-62-6 Methylmethacrylat

AGW (Deutschland) Langzeitwert: 210 mg/m³, 50 ml/m³
2(I);DFG, EU, Y

MAK (Österreich) Kurzzeitwert: 420 mg/m³, 100 ml/m³
Langzeitwert: 210 mg/m³, 50 ml/m³

MAK (Schweiz) Kurzzeitwert: 420 mg/m³, 100 ml/m³
Langzeitwert: 210 mg/m³, 50 ml/m³
S SSc;

CAS: 2082-81-7 1,4-Butandioldimethacrylat

MAK (Deutschland) vgl.Abschn.IV

· Rechtsvorschriften

AGW (Deutschland): TRGS 900

MAK (Österreich): GKV 2020, 156. Verordnung, 09.04.2021, Teil II

MAK (Schweiz): MAK- und BAT-Liste

MAK (Deutschland): MAK- und BAT-Liste

· Zusätzliche Hinweise: Als Grundlage dienen die bei der Erstellung gültigen Listen.

· 8.2 Begrenzung und Überwachung der Exposition

· Geeignete technische Steuerungseinrichtungen Keine weiteren Angaben, siehe Abschnitt 7.

· Individuelle Schutzmaßnahmen, zum Beispiel persönliche Schutzausrüstung

· Allgemeine Schutz- und Hygienemaßnahmen:

Berufsübliche Hygienemaßnahmen einhalten.

Von Nahrungsmitteln, Getränken und Futtermitteln fernhalten.

Beschmutzte, getränkte Kleidung sofort ausziehen.

Vor den Pausen und bei Arbeitsende Hände waschen.

Berührung mit den Augen und der Haut vermeiden.

Gase/Dämpfe/Aerosole nicht einatmen.

· Atemschutz

Bei kurzzeitiger oder geringer Belastung Atemfiltergerät; bei intensiver bzw. längerer Exposition unluftunabhängiges Atemschutzgerät verwenden.

· Empfohlenes Filtergerät für kurzzeitigen Einsatz:

Filter A1

Filter A2

Filter A3

· Handschutz



Schutzhandschuhe (EN 374)

Nach der Verwendung von Handschuhen Hautreinigung- und Hautpflegemittel einsetzen.

· Handschuhmaterial

Butylkautschuk

Auswahl des Handschuhmaterials unter Beachtung der Durchbruchzeiten, Permeationsraten und der Degradation.

· Durchdringungszeit des Handschuhmaterials

Die genaue Durchbruchzeit ist beim Schutzhandschuhhersteller zu erfahren und zu beachten.

(Fortsetzung auf Seite 5)

Sicherheitsdatenblatt gemäß 1907/2006/EG, Artikel 31

Druckdatum: 03.04.2024

Versionsnummer 2 (ersetzt Version 1)

überarbeitet am: 03.04.2024

Handelsname: Ivotion Bond Modelling Liquid

(Fortsetzung von Seite 4)

- **Nicht geeignet sind Handschuhe aus folgenden Materialien:**
Handelsübliche medizinische Handschuhe bieten keinen Schutz gegen den sensibilisierenden Effekt von Methacrylaten.
- **Augen-/Gesichtsschutz**



Dichtschließende Schutzbrille (EN 166)

ABSCHNITT 9: Physikalische und chemische Eigenschaften

· 9.1 Angaben zu den grundlegenden physikalischen und chemischen Eigenschaften

- **Allgemeine Angaben**
- **Farbe** Farblos
- **Geruch:** Stechend
- **Geruchsschwelle:** Nicht bestimmt.
- **Schmelzpunkt/Gefrierpunkt:** -48 °C
- **Siedepunkt oder Siedebeginn und Siedebereich** 101 °C
- **Untere und obere Explosionsgrenze**
- **Untere:** 2,1 Vol %
- **Obere:** 12,5 Vol %
- **Flammpunkt:** 10 °C
- **Zündtemperatur** 430 °C
- **pH-Wert:** Nicht bestimmt.
- **Viskosität:**
- **Kinematische Viskosität** Nicht bestimmt.
- **Dynamisch:** Nicht bestimmt.
- **Löslichkeit**
- **Wasser bei 20 °C:** 1,6 g/l
- **Verteilungskoeffizient n-Oktanol/Wasser (log-Wert)** Nicht bestimmt.
- **Dampfdruck bei 20 °C:** 47 hPa
- **Dichte und/oder relative Dichte**
- **Dichte bei 20 °C:** 0,943 g/cm³
- **Relative Dichte** Nicht bestimmt.
- **Dampfdichte** Nicht bestimmt.

- **9.2 Sonstige Angaben** Keine weiteren relevanten Informationen verfügbar.
- **Aussehen:**
- **Form:** Flüssig
- **Wichtige Angaben zum Gesundheits- und Umweltschutz sowie zur Sicherheit**
- **Zündtemperatur:** Das Produkt ist nicht selbstentzündlich.
- **Explosive Eigenschaften:** Das Produkt ist nicht explosionsgefährlich, jedoch ist die Bildung explosionsgefährlicher Dampf-/Luftgemische möglich.
- **Zustandsänderung**
- **Verdampfungsgeschwindigkeit** Nicht bestimmt.

- **Angaben über physikalische Gefahrenklassen**
- **Explosive Stoffe/Gemische und Erzeugnisse mit Explosivstoff** entfällt
- **Entzündbare Gase** entfällt
- **Aerosole** entfällt
- **Oxidierende Gase** entfällt

(Fortsetzung auf Seite 6)

Sicherheitsdatenblatt gemäß 1907/2006/EG, Artikel 31

Druckdatum: 03.04.2024

Versionsnummer 2 (ersetzt Version 1)

überarbeitet am: 03.04.2024

Handelsname: Ivotion Bond Modelling Liquid

(Fortsetzung von Seite 5)

· Gase unter Druck	entfällt
· Entzündbare Flüssigkeiten	Flüssigkeit und Dampf leicht entzündbar.
· Entzündbare Feststoffe	entfällt
· Selbsterseztliche Stoffe und Gemische	entfällt
· Pyrophore Flüssigkeiten	entfällt
· Pyrophore Feststoffe	entfällt
· Selbsterhitzungsfähige Stoffe und Gemische	entfällt
· Stoffe und Gemische, die in Kontakt mit Wasser entzündbare Gase entwickeln	entfällt
· Oxidierende Flüssigkeiten	entfällt
· Oxidierende Feststoffe	entfällt
· Organische Peroxide	entfällt
· Gegenüber Metallen korrosiv wirkende Stoffe und Gemische	entfällt
· Desensibilisierte Stoffe/Gemische und Erzeugnisse mit Explosivstoff	entfällt

ABSCHNITT 10: Stabilität und Reaktivität

- **10.1 Reaktivität** Keine weiteren relevanten Informationen verfügbar.
- **10.2 Chemische Stabilität** Stabil bei normaler Handhabung und Lagerung.
- **Thermische Zersetzung / zu vermeidende Bedingungen:** Keine Zersetzung bei bestimmungsgemäßer Verwendung.
- **10.3 Möglichkeit gefährlicher Reaktionen**
Bildung explosiver Gasgemische mit Luft.
Reaktionen mit starken Oxidationsmitteln.
Polymerisation unter Wärmeentwicklung.
- **10.4 Zu vermeidende Bedingungen** Keine weiteren relevanten Informationen verfügbar.
- **10.5 Unverträgliche Materialien:** Keine weiteren relevanten Informationen verfügbar.
- **10.6 Gefährliche Zersetzungsprodukte:** Keine unter normalen Lager- und Handhabungsbedingungen.

ABSCHNITT 11: Toxikologische Angaben

- **11.1 Angaben zu den Gefahrenklassen im Sinne der Verordnung (EG) Nr. 1272/2008**
- **Akute Toxizität** Aufgrund der verfügbaren Daten sind die Einstufungskriterien nicht erfüllt.

· Einstufungsrelevante LD/LC50-Werte:

CAS: 80-62-6 Methylmethacrylat

Oral	LD50	7872 mg/kg (rat)
------	------	------------------

- **Ätz-/Reizwirkung auf die Haut**
Verursacht Hautreizungen.
- **Schwere Augenschädigung/-reizung**
Aufgrund der verfügbaren Daten sind die Einstufungskriterien nicht erfüllt.
- **Sensibilisierung der Atemwege/Haut**
Kann allergische Hautreaktionen verursachen.
- **Keimzellmutagenität** Aufgrund der verfügbaren Daten sind die Einstufungskriterien nicht erfüllt.
- **Karzinogenität** Aufgrund der verfügbaren Daten sind die Einstufungskriterien nicht erfüllt.
- **Reproduktionstoxizität** Aufgrund der verfügbaren Daten sind die Einstufungskriterien nicht erfüllt.
- **Spezifische Zielorgan-Toxizität bei einmaliger Exposition**
Kann die Atemwege reizen.
- **Spezifische Zielorgan-Toxizität bei wiederholter Exposition**
Aufgrund der verfügbaren Daten sind die Einstufungskriterien nicht erfüllt.
- **Aspirationsgefahr** Aufgrund der verfügbaren Daten sind die Einstufungskriterien nicht erfüllt.

(Fortsetzung auf Seite 7)

Sicherheitsdatenblatt
gemäß 1907/2006/EG, Artikel 31

Druckdatum: 03.04.2024

Versionsnummer 2 (ersetzt Version 1)

überarbeitet am: 03.04.2024

Handelsname: Ivotion Bond Modelling Liquid

(Fortsetzung von Seite 6)

- **Zusätzliche toxikologische Hinweise:** Keine weiteren relevanten Informationen verfügbar.
- **11.2 Angaben über sonstige Gefahren**

- **Endokrinschädliche Eigenschaften**

Keiner der Inhaltsstoffe ist enthalten.

ABSCHNITT 12: Umweltbezogene Angaben

- **12.1 Toxizität**
- **Aquatische Toxizität:** Keine weiteren relevanten Informationen verfügbar.
- **12.2 Persistenz und Abbaubarkeit** Keine weiteren relevanten Informationen verfügbar.
- **12.3 Bioakkumulationspotenzial** Keine weiteren relevanten Informationen verfügbar.
- **12.4 Mobilität im Boden** Keine weiteren relevanten Informationen verfügbar.
- **12.5 Ergebnisse der PBT- und vPvB-Beurteilung**
- **PBT:** Nicht anwendbar.
- **vPvB:** Nicht anwendbar.
- **12.6 Endokrinschädliche Eigenschaften**
Das Produkt enthält keine Stoffe mit endokrinschädlichen Eigenschaften.
- **12.7 Andere schädliche Wirkungen** Keine weiteren relevanten Informationen verfügbar.
- **Weitere ökologische Hinweise:**
- **Allgemeine Hinweise:**
Wassergefährdungsklasse 1 (Selbsteinstufung): schwach wassergefährdend
Nicht unverdünnt bzw. in größeren Mengen in das Grundwasser, in Gewässer oder in die Kanalisation gelangen lassen.

ABSCHNITT 13: Hinweise zur Entsorgung

- **13.1 Verfahren der Abfallbehandlung**
- **Empfehlung:**
Unter Beachtung der örtlichen, behördlichen Vorschriften einer geordneten Deponie oder Verbrennungsanlage zuführen.
- **Europäisches Abfallverzeichnis**
- 18 01 06* | Chemikalien, die aus gefährlichen Stoffen bestehen oder solche enthalten
- **Ungereinigte Verpackungen:**
- **Empfehlung:** Entsorgung gemäß den behördlichen Vorschriften.

ABSCHNITT 14: Angaben zum Transport

- **14.1 UN-Nummer oder ID-Nummer**
- **ADR/RID/ADN, IMDG, IATA** UN1247
- **14.2 Ordnungsgemäße UN-Versandbezeichnung**
- **ADR/RID/ADN** 1247 METHYLMETHACRYLAT, MONOMER, STABILISIERT, Lösung
- **IMDG, IATA** METHYL METHACRYLATE MONOMER, STABILIZED solution

(Fortsetzung auf Seite 8)

Sicherheitsdatenblatt
gemäß 1907/2006/EG, Artikel 31

Druckdatum: 03.04.2024

Versionsnummer 2 (ersetzt Version 1)

überarbeitet am: 03.04.2024

Handelsname: Ivotion Bond Modelling Liquid

(Fortsetzung von Seite 7)

· **14.3 Transportgefahrenklassen**

· **ADR/RID/ADN**



· **Klasse** 3 (F1) Entzündbare flüssige Stoffe
 · **Gefahrzettel** 3

· **IMDG, IATA**



· **Class** 3 Entzündbare flüssige Stoffe
 · **Label** 3

· **14.4 Verpackungsgruppe**

· **ADR/RID/ADN, IMDG, IATA** II

· **14.5 Umweltgefahren:**

· **Marine pollutant:** Nein

· **14.6 Besondere Vorsichtsmaßnahmen für den Verwender**

Achtung: Entzündbare flüssige Stoffe

· **Nummer zur Kennzeichnung der Gefahr (Kemler-Zahl):** 339

· **EMS-Nummer:** F-E,S-D

· **Stowage Category** B

· **Stowage Code** SW2 Clear of living quarters.

· **14.7 Massengutbeförderung auf dem Seeweg gemäß IMO-Instrumenten**

Nicht anwendbar.

· **Transport/weitere Angaben:**

· **ADR/RID/ADN**

· **Begrenzte Menge (LQ)** 1L

· **Freigestellte Mengen (EQ)** Code: E2

Höchste Nettomenge je Innenverpackung: 30 ml

Höchste Nettomenge je Außenverpackung: 500 ml

· **Beförderungskategorie** 2

· **Tunnelbeschränkungscode** D/E

· **IMDG**

· **Limited quantities (LQ)** 1L

· **Excepted quantities (EQ)** Code: E2

Maximum net quantity per inner packaging: 30 ml

Maximum net quantity per outer packaging: 500 ml

· **UN "Model Regulation":**

UN 1247 METHYLMETHACRYLAT, MONOMER, STABILISIERT, LÖSUNG, 3, II

DE

(Fortsetzung auf Seite 9)

Sicherheitsdatenblatt gemäß 1907/2006/EG, Artikel 31

Druckdatum: 03.04.2024

Versionsnummer 2 (ersetzt Version 1)

überarbeitet am: 03.04.2024

Handelsname: Ivotion Bond Modelling Liquid

(Fortsetzung von Seite 8)

ABSCHNITT 15: Rechtsvorschriften

· **15.1 Vorschriften zu Sicherheit, Gesundheits- und Umweltschutz/spezifische Rechtsvorschriften für den Stoff oder das Gemisch**

- Richtlinie 2012/18/EU
- Mengenschwelle (in Tonnen) für die Anwendung in Betrieben der unteren Klasse 5000 t
- Mengenschwelle (in Tonnen) für die Anwendung in Betrieben der oberen Klasse 50000 t
- VERORDNUNG (EG) Nr. 1907/2006 ANHANG XVII Beschränkungsbedingungen: 3

· Richtlinie 2011/65/EU zur Beschränkung der Verwendung bestimmter gefährlicher Stoffe in Elektro- und Elektronikgeräten – Anhang II

Keiner der Inhaltsstoffe ist enthalten.

· VERORDNUNG (EU) 2019/1148

· Anhang I - BESCHRÄNKTE AUSGANGSSTOFFE FÜR EXPLOSIVSTOFFE (Oberer Konzentrationsgrenzwert für eine Genehmigung nach Artikel 5 Absatz 3)

Keiner der Inhaltsstoffe ist enthalten.

· Anhang II - MELDEPFLICHTIGE AUSGANGSSTOFFE FÜR EXPLOSIVSTOFFE

Keiner der Inhaltsstoffe ist enthalten.

· Verordnung (EG) Nr. 273/2004 betreffend Drogenausgangsstoffe

Keiner der Inhaltsstoffe ist enthalten.

· Verordnung (EG) Nr. 111/2005 zur Festlegung von Vorschriften für die Überwachung des Handels mit Drogenaustauschstoffen zwischen der Gemeinschaft und Drittländern

Keiner der Inhaltsstoffe ist enthalten.

· **Nationale Vorschriften:**

· **Hinweise zur Beschäftigungsbeschränkung:**

822.115, Jugendarbeitsschutzverordnung - ArGV 5 und 822.115.2, Verordnung des WBF über gefährliche Arbeiten für Jugendliche sind zu beachten.

ArGV 1 und 822.111.52, Verordnung des WBF über gefährliche und beschwerliche Arbeiten bei Schwangerschaft und Mutterschaft sind nicht zutreffend.

· **Wassergefährdungsklasse: WGK 1 (Selbsteinstufung):** schwach wassergefährdend.

· **Sonstige Vorschriften, Beschränkungen und Verbotsverordnungen**

Das Produkt fällt unter die Verordnung (EU) 2017/745 über Medizinprodukte (MDR).

· **15.2 Stoffsicherheitsbeurteilung:** Eine Stoffsicherheitsbeurteilung wurde nicht durchgeführt.

ABSCHNITT 16: Sonstige Angaben

Die Angaben stützen sich auf den heutigen Stand unserer Kenntnisse, sie stellen jedoch keine Zusicherung von Produkteigenschaften dar und begründen kein vertragliches Rechtsverhältnis.

· **Relevante Sätze**

H225 Flüssigkeit und Dampf leicht entzündbar.

H302 Gesundheitsschädlich bei Verschlucken.

H315 Verursacht Hautreizungen.

H317 Kann allergische Hautreaktionen verursachen.

H318 Verursacht schwere Augenschäden.

H335 Kann die Atemwege reizen.

H412 Schädlich für Wasserorganismen, mit langfristiger Wirkung.

· **Einstufung gemäß Verordnung (EG) Nr. 1272/2008**

Die Einstufung der Mischung basiert generell auf der Berechnungsmethode unter Verwendung von Stoffdaten gemäß Verordnung (EG) Nr. 1272/2008.

· **Versionsnummer der Vorgängerversion: 1**

(Fortsetzung auf Seite 10)

Sicherheitsdatenblatt gemäß 1907/2006/EG, Artikel 31

Druckdatum: 03.04.2024

Versionsnummer 2 (ersetzt Version 1)

überarbeitet am: 03.04.2024

Handelsname: Ivotion Bond Modelling Liquid

(Fortsetzung von Seite 9)

· **Abkürzungen und Akronyme:**

MAK: Maximale Arbeitsplatz-Konzentration

AWG: Arbeitsplatzgrenzwert

DFG: Deutsche Forschungsgemeinschaft

TRGS: Technische Regeln für Gefahrstoffe

GKV: Grenzwerteverordnung

E/e: einatembare Fraktion

a: alveolengängige Fraktion

SSa: Eine Schädigung der Leibesfrucht kann auch bei Einhaltung des MAK-Wertes auftreten.

SSb: Eine Schädigung der Leibesfrucht kann auch bei Einhaltung des MAK-Wertes nicht ausgeschlossen werden.

SSc: Eine Schädigung der Leibesfrucht braucht bei Einhaltung des MAK-Wertes nicht befürchtet zu werden.

H: hautresorptiv

X: krebserzeugender Stoff der Kat. 1A oder 1B oder krebserzeugende Tätigkeit oder Verfahren nach § 2 Absatz 3 Nr. 4 der

Gefahrstoffverordnung – es ist zusätzlich § 10 GefStoffV zu beachten

Y: ein Risiko der Fruchtschädigung braucht bei Einhaltung des Arbeitsplatzgrenzwertes und des biologischen Grenzwertes (BGW) nicht befürchtet zu werden

Z: ein Risiko der Fruchtschädigung kann auch bei Einhaltung des AGW und des BGW nicht ausgeschlossen werden

ArGV: Verordnung zum Arbeitsgesetz (Schweiz)

WBF: Eidgenössische Departement für Wirtschaft, Bildung und Forschung (Schweiz)

IMO: International Maritime Organization

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

Flam. Liq. 2: Entzündbare Flüssigkeiten – Kategorie 2

Acute Tox. 4: Akute Toxizität – Kategorie 4

Skin Irrit. 2: Hautreizende/-ätzende Wirkung – Kategorie 2

Eye Dam. 1: Schwere Augenschädigung/Augenreizung – Kategorie 1

Skin Sens. 1: Sensibilisierung der Haut – Kategorie 1

STOT SE 3: Spezifische Zielorgan-Toxizität (einmalige Exposition) – Kategorie 3

Aquatic Chronic 3: Gewässergefährdend - langfristig gewässergefährdend – Kategorie 3

· *** Daten gegenüber der Vorversion geändert**