

Trade name: Kombiplast

Substance number: S0019 Version: 1 / GB Date revised: 09.06.2023

> Replaces Version: - / GB Print date: 09.06.2023

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Kombiplast

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/preparation

Dental splinting material long-term

1.3. Details of the supplier of the safety data sheet

Address/Manufacturer

Dreve Dentamid GmbH Max-Planck-Straße 31

59423 Unna

Telephone no. +49 2303 8807-0 Fax no. +49 2303 8807-29

Information provided Department Research & Development: Fax: +49 2303 8807-562

by / telephone E-mail address of

sicherheitsdatenblatt@dreve.com

person responsible

for this SDS

1.4. Emergency telephone number

Henkel Fire Department / 24h-Emergency-Contact-No.: +49 211 797-3350

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

This product is not classified hazardous in accordance with Regulation (EC) No 1272/2008.

2.2. Label elements

Labelling according to regulation (EC) No 1272/2008

Supplemental information

EUH210 Safety data sheet available on request.

2.3. Other hazards

No special hazards have to be mentioned.

The Substance does not meet PBT-criteria. This substance does not meet the vPvB-criteria. This substance does not have endocrine disrupting properties with respect to humans. This substance does not have endocrine disrupting properties with respect to non-target organisms.

SECTION 3: Composition/information on ingredients

3.2. Mixtures



Trade name: Kombiplast

Substance number: S0019 Version: 1 / GB Date revised: 09.06.2023

Replaces Version: - / GB Print date: 09.06.2023

Chemical characterization

Double coating plate of PETG and EVA

Hazardous ingredients

vinyl acetate

CAS No. 108-05-4 EINECS no. 203-545-4

Concentration \Rightarrow 0,1 < 1 %

Classification (Regulation (EC) No. 1272/2008)

Flam. Liq. 2 H225 Acute Tox. 4 H332 Carc. 2 H351 STOT SE 3 H335

Additional remarks:

CLP Regulation (EC) No 1272/2008, Annex VI, Note D

SECTION 4: First aid measures

4.1. Description of first aid measures

General information

No special measures required

After inhalation

Ensure supply of fresh air.

After skin contact

After contact with molten product cool quickly with cold water. Consult a doctor if skin irritation persists.

After eye contact

Separate eyelids, wash the eyes thoroughly with water (15 min.). In case of irritation consult an oculist.

After ingestion

Rinse out mouth and give plenty of water to drink.

Adhere to personal protective measures when giving first aid

First aider: Pay attention to self-protection!

4.2. Most important symptoms and effects, both acute and delayed

Until now no symptoms known so far.

4.3. Indication of any immediate medical attention and special treatment needed

Hints for the physician / hazards

In the case of swallowing with subsequent vomiting, aspiration of the lungs can occur which can lead to chemical pneumonia or asphyxiation.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media

Recommended: alcohol resistant foam, CO2, powders, water spray/mist, Extinguishing measures to suit surroundings



Trade name: Kombiplast

Substance number: S0019 Version: 1 / GB Date revised: 09.06.2023

Replaces Version: - / GB Print date: 09.06.2023

5.3. Advice for firefighters

Special protective equipment for fire-fighting

In case of combustion use a suitable breathing apparatus.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Avoid contact with skin, eyes and clothing.

6.2. Environmental precautions

Do not allow to enter drains or waterways.

6.3. Methods and material for containment and cleaning up

Pick up mechanically. Clean contaminated floors and objects thoroughly, observing environmental regulations.

6.4. Reference to other sections

Refer to protective measures listed in Sections 7 and 8.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on safe handling

Observe the usual precautions for handling chemicals.

Advice on protection against fire and explosion

No special measures required.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

Store product in closed containers.

Hints on storage assembly

Do not store together with foodstuffs.

Further information on storage conditions

Keep container tightly closed and dry.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Other information

Contains no substances with occupational exposure limit values.

8.2. Exposure controls

General protective and hygiene measures

Observe the usual precautions for handling chemicals.

Respiratory protection

Do not inhale vapours



Trade name: Kombiplast

Substance number: S0019 Version: 1 / GB Date revised: 09.06.2023

Replaces Version: - / GB Print date: 09.06.2023

Hand protection

Heat-resistant gloves when using hot material

Eye protection

Safety glasses

Body protection

Clothing as usual in the chemical industry.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state Thermoforming blank clear, transparent characteristic

Melting point

Value > 100 °C

Freezing point

Remarks not determined

Boiling point or initial boiling point and boiling range

Remarks not determined

Flammability

evaluation not determined

Upper and lower explosive limits

Remarks not determined

Flash point

Remarks Not applicable

Ignition temperature

Remarks not determined

Decomposition temperature

Remarks not determined

pH value

Remarks not determined

Viscosity

Remarks not determined

Solubility(ies)

Remarks not determined

Partition coefficient n-octanol/water (log value)

Remarks not determined

Vapour pressure

Remarks not determined

Density and/or relative density

Value 1,2 g/cm³

Temperature 20 °C

Relative vapour density

Remarks not determined



Trade name: Kombiplast

Substance number: S0019 Version: 1 / GB Date revised: 09.06.2023

Replaces Version: - / GB Print date: 09.06.2023

9.2. Other information

Odour threshold

Remarks not determined

Evaporation rate (ether = 1):

Remarks not determined

Solubility in water

Remarks virtually insoluble

Explosive properties

evaluation not determined

Oxidising properties

Remarks not determined

Other information None known

SECTION 10: Stability and reactivity

10.1. Reactivity

No hazardous reactions when stored and handled according to prescribed instructions.

10.2. Chemical stability

No hazardous reactions known.

10.3. Possibility of hazardous reactions

No hazardous reactions known.

10.4. Conditions to avoid

Protect from heat and direct sunlight

10.5. Incompatible materials

None known

SECTION 11: Toxicological information

11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008

Acute oral toxicity

Remarks Based on available data, the classification criteria are not met.

Acute oral toxicity (Components)

vinyl acetate

Species rat (male) LD50 3470

3470 mg/kg

Acute dermal toxicity

Remarks Based on available data, the classification criteria are not met.

Acute dermal toxicity (Components)

vinyl acetate

Species rabbit

LD50 7440 mg/kg



Trade name: Kombiplast

Substance number: S0019 Version: 1 / GB Date revised: 09.06.2023

Replaces Version: - / GB Print date: 09.06.2023

Acute inhalational toxicity

Remarks Based on available data, the classification criteria are not met.

Acute inhalative toxicity (Components)

vinyl acetate

Species rat

LC50 4000 ppm(V)

Duration of exposure 4 h

Administration/Form Vapors

Skin corrosion/irritation

Remarks Based on available data, the classification criteria are not met.

Serious eye damage/irritation

Remarks Based on available data, the classification criteria are not met.

Sensitization

Remarks Based on available data, the classification criteria are not met.

Subacute, subchronic, chronic toxicity

Remarks not determined

Mutagenicity

Remarks Based on available data, the classification criteria are not met.

Reproductive toxicity

Remarks Based on available data, the classification criteria are not met.

Carcinogenicity

Remarks Based on available data, the classification criteria are not met.

Carcinogenicity (Components)

vinyl acetate

evaluation Suspected of causing cancer.

Specific Target Organ Toxicity (STOT)

Single exposure

Remarks Based on available data, the classification criteria are not met.

Repeated exposure

Remarks Based on available data, the classification criteria are not met.

Specific Target Organ Toxicity (STOT) (Components)

vinyl acetate

Single exposure

evaluation May cause respiratory irritation.

Aspiration hazard

Based on available data, the classification criteria are not met.

11.2 Information on other hazards

Endocrine disrupting properties with respect to humans

This substance does not have endocrine disrupting properties with respect to humans.

Other information

No toxicological data are available.

SECTION 12: Ecological information



Trade name: Kombiplast

Substance number: S0019 Version: 1 / GB Date revised: 09.06.2023

Replaces Version: - / GB Print date: 09.06.2023

12.1. Toxicity

General information

not determined

Fish toxicity (Components)

vinyl acetate

Species Fathead minnow (Pimephales promelas)
NOEC 0,551 mg/l

Duration of exposure 34 d

Method OECD 210

Daphnia toxicity (Components)

vinyl acetate

Species Daphnia magna

EC50 12,6 mg/l

Duration of exposure 48 h

Method OECD 202

Algae toxicity (Components)

vinyl acetate

Species Pseudokirchneriella subcapitata

ErC50 12,7 mg/l

Duration of exposure 72 h

Method OECD 201

12.2. Persistence and degradability

General information

not determined

Biodegradability (Components)

vinyl acetate

Value 82 to 98 %

Duration of test 14 d

evaluation Readily biodegradable (according to OECD criteria)

12.3. Bioaccumulative potential

General information

not determined

Partition coefficient n-octanol/water (log value)

Remarks not determined

Octanol/water partition coefficient (log Pow) (Components)

vinyl acetate

log Pow 0,73

Temperature 25 °C

12.4. Mobility in soil

General information

not determined

12.5. Results of PBT and vPvB assessment

General information

not determined

Results of PBT and vPvB assessment



Trade name: Kombiplast

Substance number: S0019 Version: 1 / GB Date revised: 09.06.2023

Replaces Version: - / GB Print date: 09.06.2023

The Substance does not meet PBT-criteria. This substance does not meet the vPvB-criteria.

12.6 Endocrine disrupting properties

Endocrine disrupting properties with respect to the envrionment

This substance does not have endocrine disrupting properties with respect to non-target organisms.

12.7. Other adverse effects

General information

not determined

General information / ecology

Do not discharge product unmonitored into the environment.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal recommendations for the product

Must not be disposed together with household garbage.

Dispose of waste according to applicable legislation.

Disposal recommendations for packaging

Packaging that cannot be cleaned should be disposed off in agreement with the regional waste disposal company.

SECTION 14: Transport information

	Land transport ADR/RID	Marine transport IMDG/GGVSee	Air transport ICAO/IATA
14.1. UN number or ID number	The product does not constitute a hazardous substance in land transport.	The product does not constitute a hazardous substance in sea transport.	The product does not constitute a hazardous substance in air transport.
14.2. UN proper shipping name	-	-	-
14.3. Transport hazard class(es)		-	-
Label			
14.4. Packing group		-	-

SECTION 15: Regulatory information

15.2. Chemical safety assessment

For this preparation a chemical safety assessment has not been carried out.



Trade name: Kombiplast

Substance number: S0019 Version: 1 / GB Date revised: 09.06.2023

Replaces Version: - / GB Print date: 09.06.2023

SECTION 16: Other information

Hazard statements listed in Chapter 2/3

H225 Highly flammable liquid and vapour.

H332 Harmful if inhaled.

H335 May cause respiratory irritation. H351 Suspected of causing cancer.

CLP categories listed in Chapter 2/3

Acute Tox. 4 Acute toxicity, Category 4
Carc. 2 Carcinogenicity, Category 2
Flam. Liq. 2 Flammable liquid, Category 2

STOT SE 3 Specific target organ toxicity - single exposure, Category 3

Supplemental information

Relevant changes compared with the previous version of the safety data sheet are marked with: *** This information is based on our present state of knowledge. However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.