Ansell

MICRO-TOUCH[®] NITRA-TEX[™]

Performance Series

Nitrile, powder-free examination glove with strong barrier protection

This non-sterile glove is made from a softer synthetic formulation to offer a comfortable fit and feel. It also has an excellent resistance from a broad range of commonly used chemicals. High tear resistance ensures optimal protection, and a chlorinated inner coating facilitates ease of donning.

KEY FEATURES AND BENEFITS:

- Excellent resistance to a broad range of commonly used chemicals
- Protects from latex Type I allergy
- Tested for use with chemotherapy drugs*

* See product packaging or contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

WARNING: No glove completely protects against all chemicals. Users should test the suitability of this product against the specific chemicals and environment where used.

** This product is only available in Europe, the Middle East and Africa (EMEA).

Recommended For

- Tested for use with chemotherapy drugs
- Standard, moderate risk examination procedures
- Environmentally friendly facilities
- Protection from Type I latex allergy in HCW's or patients

FEATURED TECHNOLOGIES





MICRO-TOUCH[®] Everyday, Everywhere Protection

MICRO-TOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.



MICRO-TOUCH[®] NITRA-TEX[™]

PRODUCT INFORMATION

Material	Nitrile		Not Made From Natural Rubber Latex		l Yes	Yes	
Color	Blue		Cuff Length		Extended	Extended	
Powder Content	Powder-Free		External Glove Surface		Fully Textu	Fully Textured	
Freedom from Holes (Inspection level I)	1.5 AQL		Palm Thickness (mm/mil)		il) 0.10 / 3.9	0.10 / 3.9	
Finger Thickness (mm/mil)	0.12 / 4.7		Alle	Allergy Prevention Latex (Type I)			
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)		Tested For Use With Chemotherapy Drugs			Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)	
Sterile	No		Antistatic		Not Tested	Not Tested	
Glove Length (mm/inches)	240 / 9.5		Product Segmentation		Performanc	Performance Series	
ORDERING INFORMATIC	N						
Size	XS (5.5 - 6)	S (6.5 - 7))	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)	
Product Code	700111	700112		700113	700114	700115	
PACKAGING AND STOR	AGE						
Packaging	100 gloves per dis	100 gloves per dispenser 10 dispensers per case 1000 gloves per case					

Shelf Life **3** Years Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or **Storage Instructions**

ignition.

STANDARDS AND CERTIFICATIONS

CE 0493, Category III, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 421:2010, ISO 13485, ISO 10993-10:2002, ISO 9001, Food Contact

Vulcanization Chemical Accelerators

Zinc Diethyldithiocarbamate (ZDEC)

Zinc Dibutyldithiocarbamate (ZDBC)

Only a very small number of users may be sensitive to this ingredient(s) and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America US Tel: 800 952 9916 CA Tel: 1-844-494-7854 Email: insidesalesus@ansell.com Latin America & Caribbean Tel: +52(442) 296 20 50 Email: cslac@ansell.com Brazil Tel: +55-11-3356-3100 Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa Tel: +6 Tel: +32 (0) 2 528 74 00 Email: Email: info@ansell.eu

Asia Pacific Tel: +603 8310 6688 apac.medical@ansell.com

Australia & New Zealand Tel: +61 3 9270 7270 Email: protection@ap.ansell.com

NITRA-TEX

ansell.com Ansell, [®] and [™] are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking. © 2022 Ansell Limited. All Rights Reserved.