

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 21.11.2023

Version number 4 (replaces version 3)

Revision: 21.11.2023

- **1.1 Product identifier**
- **Trade name: *Omni Okklusionsspray grün***
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**  
No further relevant information available.
- **Application of the substance / the mixture** *Auxiliary for dental technology*
- **1.3 Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
*PRISMAN GmbH  
Otto Hahn Ring 6-18  
D-64653 Lorsch  
Germany*
- **Further information obtainable from:**  
*Abteilung Produktsicherheit  
Alexander.Metz@prisman.de*
- **1.4 Emergency telephone number:** ++49 (0)6251 866980-0, Mo - Fr 8-18 Uhr

## SECTION 2: Hazards identification

- **2.1 Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008**  
*Aerosol 1 H222 Extremely flammable aerosol.  
H229 Pressurised container: May burst if heated.*
- **2.2 Label elements**
- **Labelling according to Regulation (EC) No 1272/2008**  
*The product is classified and labelled according to the GB CLP regulation.*
- **Hazard pictograms**



GHS02

- **Signal word** *Danger*
- **Hazard statements**  
*H222 Extremely flammable aerosol.  
H229 Pressurised container: May burst if heated.*
- **Precautionary statements**  
*P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.  
P211 Do not spray on an open flame or other ignition source.  
P251 Do not pierce or burn, even after use.  
P410+P412 Protect from sunlight. Do not expose to temperatures exceeding 50 °C/122 °F.*
- **Additional information:**  
*Pressurized container: protect from sunlight and do not expose to temperatures exceeding 50°C. Do not pierce or burn, even after use.  
Do not spray on a naked flame or any incandescent material. Keep away from sources of ignition - No smoking.  
Keep out of the reach of children*
- **2.3 Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** *Not applicable.*
- **vPvB:** *Not applicable.*

GB

(Contd. on page 2)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 21.11.2023

Version number 4 (replaces version 3)

Revision: 21.11.2023

Trade name: **Omni Okklusionsspray grün**

(Contd. of page 1)

### SECTION 3: Composition/information on ingredients

#### 3.2 Mixtures

Description: Mixture of substances listed below with nonhazardous additions.

#### Dangerous components:

CAS: 75-28-5 EINECS: 200-857-2 Index number: 601-004-00-0 RTECS: TZ 4300000 Reg.nr.: 01-2119485395-27-xxxx	isobutane	 Flam. Gas 1, H220 Press. Gas (Comp.), H280	50-100%
------------------------------------------------------------------------------------------------------------------------	-----------	-----------------------------------------------------------------------------------------------------------------------------------	---------

Additional information: For the wording of the listed hazard phrases refer to section 16.

### SECTION 4: First aid measures

#### 4.1 Description of first aid measures

After inhalation: Supply fresh air; consult doctor in case of complaints.

#### After skin contact:

In cases of frost bites, rinse with plenty of water. Do not remove clothing.

Rinse with warm water.

After eye contact: Rinse opened eye for several minutes under running water.

After swallowing: Rinse out mouth and then drink plenty of water.

4.2 Most important symptoms and effects, both acute and delayed No further relevant information available.

4.3 Indication of any immediate medical attention and special treatment needed

No further relevant information available.

### SECTION 5: Firefighting measures

#### 5.1 Extinguishing media

#### Suitable extinguishing agents:

CO<sub>2</sub>, powder or water spray. Fight larger fires with water spray.

Use fire extinguishing methods suitable to surrounding conditions.

#### 5.2 Special hazards arising from the substance or mixture

In case of fire, the following can be released:

Carbon monoxide (CO)

#### 5.3 Advice for firefighters

Protective equipment: Do not inhale explosion gases or combustion gases.

### SECTION 6: Accidental release measures

#### 6.1 Personal precautions, protective equipment and emergency procedures

Keep away from ignition sources.

Ensure adequate ventilation

Wear protective equipment. Keep unprotected persons away.

#### 6.2 Environmental precautions: No special measures required.

#### 6.3 Methods and material for containment and cleaning up:

Allow to evaporate.

Ensure adequate ventilation.

#### 6.4 Reference to other sections

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

(Contd. on page 3)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 21.11.2023

Version number 4 (replaces version 3)

Revision: 21.11.2023

**Trade name: Omni Okklusionsspray grün**

(Contd. of page 2)

See Section 13 for disposal information.

### SECTION 7: Handling and storage

· **7.1 Precautions for safe handling** Keep away from heat and direct sunlight.

· **Information about fire - and explosion protection:**

Do not spray onto a naked flame or any incandescent material.

Highly volatile, flammable constituents are released during processing.

Fumes can combine with air to form an explosive mixture.

Emergency cooling must be available in case of nearby fire.

Keep ignition sources away - Do not smoke.

· **7.2 Conditions for safe storage, including any incompatibilities**

· **Storage:**

· **Requirements to be met by storerooms and receptacles:**

Store only in the original receptacle.

Observe official regulations on storing packagings with pressurised containers.

· **Information about storage in one common storage facility:** Not required.

· **Further information about storage conditions:** Keep container tightly sealed.

· **7.3 Specific end use(s)** No further relevant information available.

### SECTION 8: Exposure controls/personal protection

· **8.1 Control parameters**

· **Ingredients with limit values that require monitoring at the workplace:**

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

· **Additional information:** The lists valid during the making were used as basis.

· **8.2 Exposure controls**

· **Appropriate engineering controls** No further data; see section 7.

· **Individual protection measures, such as personal protective equipment**

· **General protective and hygienic measures:** Wash hands before breaks and at the end of work.

· **Respiratory protection:**

Not necessary if room is well-ventilated.

Not required.

Use suitable respiratory protective device in case of insufficient ventilation.

· **Hand protection**

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

-

· **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

-

(Contd. on page 4)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 21.11.2023

Version number 4 (replaces version 3)

Revision: 21.11.2023

Trade name: **Omni Okklusionsspray grün**

(Contd. of page 3)

- For the permanent contact in work areas without heightened risk of injury (e.g. Laboratory) gloves made of the following material are suitable:  
Rubber gloves
- For the permanent contact gloves made of the following materials are suitable: Neoprene gloves
- Eye/face protection Not required.

### SECTION 9: Physical and chemical properties

#### · 9.1 Information on basic physical and chemical properties

##### · General Information

- |                                                            |                                    |
|------------------------------------------------------------|------------------------------------|
| · Physical state                                           | Aerosol                            |
| · Colour:                                                  | Different according to colouring   |
| · Odour:                                                   | Pleasant                           |
| · Odour threshold:                                         | Not determined.                    |
| · Melting point/freezing point:                            | -159.4 °C                          |
| · Boiling point or initial boiling point and boiling range | -44 °C                             |
| · Flammability                                             | Extremely flammable liquefied gas. |
| · Lower and upper explosion limit                          |                                    |
| · Lower:                                                   | 1.8 Vol % (Propan/Butan)           |
| · Upper:                                                   | 8.5 Vol % (Propan/Butan)           |
| · Flash point:                                             | Not applicable, as aerosol.        |
| · Auto-ignition temperature:                               | 460 °C                             |
| · Decomposition temperature:                               | Not determined.                    |
| · pH                                                       | Not applicable                     |
| · Viscosity:                                               |                                    |
| · Kinematic viscosity                                      | Not applicable                     |
| · Dynamic:                                                 | Not applicable                     |
| · Solubility                                               |                                    |
| · water at 20 °C:                                          | 0.049 g/l (Isobutan)               |
| · Partition coefficient n-octanol/water (log value)        | Not determined.                    |
| · Vapour pressure at 20 °C:                                | 3,000 hPa                          |
| · Density and/or relative density                          |                                    |
| · Density at 20 °C:                                        | 0.56 g/cm <sup>3</sup>             |
| · Relative density                                         | Not determined.                    |
| · Vapour density                                           | Not determined.                    |

#### · 9.2 Other information

- |                                                                                 |                              |
|---------------------------------------------------------------------------------|------------------------------|
| · Appearance:                                                                   |                              |
| · Form:                                                                         | Aerosol                      |
| · Important information on protection of health and environment, and on safety. |                              |
| · Ignition temperature:                                                         | Product is not selfigniting. |
| · Explosive properties:                                                         | Not determined.              |
| · Solvent content:                                                              |                              |
| · VOC (EC)                                                                      | 95 %                         |
| · Solids content:                                                               | <1.0 %                       |
| · Change in condition                                                           |                              |
| · Evaporation rate                                                              | Not applicable.              |

#### · Information with regard to physical hazard classes

- |                   |                          |
|-------------------|--------------------------|
| · Explosives      | Void                     |
| · Flammable gases | Extremely flammable gas. |
| · Aerosols        | Void                     |

(Contd. on page 5)

**Safety data sheet**  
according to 1907/2006/EC, Article 31

Printing date 21.11.2023

Version number 4 (replaces version 3)

Revision: 21.11.2023

Trade name: **Omni Okklusionsspray grün**

(Contd. of page 4)

· <b>Oxidising gases</b>	Void
· <b>Gases under pressure</b>	Contains gas under pressure; may explode if heated.
· <b>Flammable liquids</b>	Void
· <b>Flammable solids</b>	Void
· <b>Self-reactive substances and mixtures</b>	Void
· <b>Pyrophoric liquids</b>	Void
· <b>Pyrophoric solids</b>	Void
· <b>Self-heating substances and mixtures</b>	Void
· <b>Substances and mixtures, which emit flammable gases in contact with water</b>	Void
· <b>Oxidising liquids</b>	Void
· <b>Oxidising solids</b>	Void
· <b>Organic peroxides</b>	Void
· <b>Corrosive to metals</b>	Void
· <b>Desensitised explosives</b>	Void

**SECTION 10: Stability and reactivity**

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions** No dangerous reactions known.
- **10.4 Conditions to avoid** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

**SECTION 11: Toxicological information**

- **11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**
- **Acute toxicity** Based on available data, the classification criteria are not met.

· **LD/LC50 values relevant for classification:**

75-28-5 isobutane

Inhalative	LC50/4 h	258,000 mg/l (rat)
------------	----------	--------------------

- **Skin corrosion/irritation** Based on available data, the classification criteria are not met.
- **Serious eye damage/irritation** Based on available data, the classification criteria are not met.
- **Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

· **11.2 Information on other hazards**

· **Endocrine disrupting properties**

None of the ingredients is listed.

GB

(Contd. on page 6)

**Safety data sheet**  
according to 1907/2006/EC, Article 31

Printing date 21.11.2023

Version number 4 (replaces version 3)

Revision: 21.11.2023

Trade name: **Omni Okklusionsspray grün**

(Contd. of page 5)

### SECTION 12: Ecological information

#### · 12.1 Toxicity

##### · Aquatic toxicity:

##### 75-28-5 isobutane


EC50	7.71 mg/l (A)
LC50/96h	27.98 mg/l (fish)

- **12.2 Persistence and degradability** No further relevant information available.
- **12.3 Bioaccumulative potential** No further relevant information available.
- **12.4 Mobility in soil** No further relevant information available.
- **12.5 Results of PBT and vPvB assessment**
  - **PBT:** Not applicable.
  - **vPvB:** Not applicable.
- **12.6 Endocrine disrupting properties**  
The product does not contain substances with endocrine disrupting properties.
- **12.7 Other adverse effects**
- **Remark:** The product is highly volatile.
- **Additional ecological information:**
- **General notes:** Not hazardous for water.

### SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
- **Recommendation**  
Must be specially treated adhering to official regulations.  
Must not be disposed together with household garbage. Do not allow product to reach sewage system.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

### SECTION 14: Transport information

- **14.1 UN number or ID number**
- **ADR, IMDG, IATA** UN1950
- **14.2 UN proper shipping name**
- **ADR** 1950 AEROSOLS
- **IMDG** AEROSOLS
- **IATA** AEROSOLS, flammable
- **14.3 Transport hazard class(es)**
- **ADR**
- 
- **Class** 2 5F Gases.

(Contd. on page 7)

GB

**Safety data sheet**  
according to 1907/2006/EC, Article 31

Printing date 21.11.2023

Version number 4 (replaces version 3)

Revision: 21.11.2023

Trade name: **Omni Okklusionsspray grün**

(Contd. of page 6)

· **Label** 2.1· **IMDG, IATA**

· **Class** 2.1 Gases.  
· **Label** 2.1

· **14.4 Packing group**  
· **ADR, IMDG, IATA** Void

· **14.5 Environmental hazards:**  
· **Marine pollutant:** No

· **14.6 Special precautions for user** Warning: Gases.  
· **Hazard identification number (Kemler code):** -  
· **EMS Number:** F-D,S-U

· **14.7 Maritime transport in bulk according to IMO instruments** Not applicable.

· **Transport/Additional information:**

· **ADR**  
· **Excepted quantities (EQ)** Code: E0  
Not permitted as Excepted Quantity

· **IMDG**  
· **Limited quantities (LQ)** 1L  
· **Excepted quantities (EQ)** Code: E0  
Not permitted as Excepted Quantity

### SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**  
· **Labelling according to Regulation (EC) No 1272/2008**  
GHS label elements  
The product is classified and labelled according to the GB CLP regulation.  
· **Hazard pictograms**



GHS02

- **Signal word** Danger  
· **Hazard statements**  
H222 Extremely flammable aerosol.  
H229 Pressurised container: May burst if heated.  
· **Precautionary statements**  
P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.  
P211 Do not spray on an open flame or other ignition source.  
P251 Do not pierce or burn, even after use.

(Contd. on page 8)

GB

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 21.11.2023

Version number 4 (replaces version 3)

Revision: 21.11.2023

**Trade name: Omni Okklusionsspray grün**

(Contd. of page 7)

P410+P412 Protect from sunlight. Do not expose to temperatures exceeding 50 °C/122 °F.

- **Directive 2012/18/EU**
- **Qualifying quantity (tonnes) for the application of lower-tier requirements** 150 t
- **Qualifying quantity (tonnes) for the application of upper-tier requirements** 500 t
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

### SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Relevant phrases**  
 H220 Extremely flammable gas.  
 H280 Contains gas under pressure; may explode if heated.
- **Recommended restriction of use** Product only for professional use
- **Department issuing SDS:** Abteilung Produktsicherheit
- **Contact:** Hr. Dr. Metz
- **Abbreviations and acronyms:**  
 RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)  
 IATA-DGR: Dangerous Goods Regulations by the "International Air Transport Association" (IATA)  
 ICAO: International Civil Aviation Organisation  
 ICAO-TI: Technical Instructions by the "International Civil Aviation Organisation" (ICAO)  
 ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)  
 IMDG: International Maritime Code for Dangerous Goods  
 IATA: International Air Transport Association  
 GHS: Globally Harmonised System of Classification and Labelling of Chemicals  
 EINECS: European Inventory of Existing Commercial Chemical Substances  
 ELINCS: European List of Notified Chemical Substances  
 CAS: Chemical Abstracts Service (division of the American Chemical Society)  
 VOC: Volatile Organic Compounds (USA, EU)  
 LC50: Lethal concentration, 50 percent  
 LD50: Lethal dose, 50 percent  
 PBT: Persistent, Bioaccumulative and Toxic  
 vPvB: very Persistent and very Bioaccumulative  
 Flam. Gas 1: Flammable gases – Category 1  
 Aerosol 1: Aerosols – Category 1  
 : Aerosols – Category 3  
 Press. Gas (Comp.): Gases under pressure – Compressed gas
- **\* Data compared to the previous version altered.**

GB