

MicroMega **GLIDEPATH**

One Glider
RECI Glider



Schaffung eines pfades durch jedes
hindernis

Produktvorteile

- Geringere spannung und pfadpräparation
- T.Wire : Flexibilität
- Flache spitze : Hervorragende führung und erhaltung der wurzelkanalanatomie
- Variabler querschnitt: Hohe widerstandsfähigkeit und verbesserung der debris-entfernung
- Erhältlich in kontinuierlicher rotation und reziproker bewegung

Ø
Durchmesser

16



Taper

3%



Verfügbare Längen

21mm

25mm

31mm

MicroMega One Glider

5 instrumente, Steriler Blister



Kontinuierliche rotation



Motordrehzahl

250 - 400 rpm



Maximales drehmoment

1.2 N.cm

Referenzen

B5 MicroMega One Glider n16 3% L21 Ste	60036825
B5 MicroMega One Glider n16 3% L25 Ste	60036827
B5 MicroMega One Glider n16 3% L31 Ste	60036829



Reziproke Winkel

170°/60°



Motordrehzahl

400rpm



Maximales drehmoment

4 N.cm

MicroMega RECI Glider

5 instrumente, Steriler Blister

Referenzen

B5 MicroMega RECI Glider n16 3% L21 Ste	60036826
B5 MicroMega RECI Glider n16 3% L25 Ste	60036828
B5 MicroMega RECI Glider n16 3% L31 Ste	60036830



80000516-A 02.25

Legal manufacturer :

Micro-Mega SA

12 rue du Tunnel
25000 BESANCON / France
T +33 3 81 54 42 42
mm.sales@coltene.com

MicroMega One Glider

MicroMega RECI Glider

Advertising | Medical device class IIa according to Directive 93/42/EEC | Notified body : GMED | Medical device for dental care, meant for professional dental use only. See product labelling and instructions for use on the website : micro-mega.com



COLTENEGROUP