

# Safety Data Sheet

according to UK REACH Regulation

## PERform Alginatisation

Revision date: 17.09.2024

Product code: 100520

Page 1 of 6

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

PERform Alginatisation

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

##### Use of the substance/mixture

Dental aid for prosthetics

#### 1.3. Details of the supplier of the safety data sheet

Company name:	Hedent GmbH	
Street:	Obere Zeil 6	
Place:	D 61440 Oberursel	
Telephone:	+49 (0) 6171 - 52036	Telefax: +49 (0) 6171 - 52090
E-mail:	info@retec-dent.de	
Contact person:	Research & Development Department	
E-mail:	info@retec-dent.de	
Internet:	www.retec-dent.de	

#### 1.4. Emergency telephone number:

+49 (0) 6171-52036

### SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture

##### GB CLP Regulation

This mixture is not classified as hazardous in accordance with GB CLP Regulation.

#### 2.2. Label elements

#### 2.3. Other hazards

No information available.

### SECTION 3: Composition/information on ingredients

#### 3.2. Mixtures

##### Chemical characterization

Aqueous solution of alginates and preservatives

##### Relevant ingredients

none (according to UK REACH Regulation)

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

##### After inhalation

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

##### After contact with skin

Wash with plenty of water.

##### After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

##### After ingestion

Rinse mouth immediately and drink plenty of water.

# Safety Data Sheet

according to UK REACH Regulation

## PERform Alginatisation

Revision date: 17.09.2024

Product code: 100520

Page 2 of 6

### 4.2. Most important symptoms and effects, both acute and delayed

No information available.

### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

#### Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings.

### 5.2. Special hazards arising from the substance or mixture

Non-flammable.

### 5.3. Advice for firefighters

In case of fire: Wear self-contained breathing apparatus.

#### Additional information

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

#### General advice

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

### 6.2. Environmental precautions

Do not allow to enter into surface water or drains.

### 6.3. Methods and material for containment and cleaning up

#### Other information

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

### 6.4. Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

#### Advice on safe handling

No special measures are necessary.

#### Advice on protection against fire and explosion

No special fire protection measures are necessary.

#### Advice on general occupational hygiene

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat or drink.

### 7.2. Conditions for safe storage, including any incompatibilities

#### Requirements for storage rooms and vessels

Keep container tightly closed.

#### Hints on joint storage

No special measures are necessary.

# Safety Data Sheet

according to UK REACH Regulation

## PERform Alginatisation

Revision date: 17.09.2024

Product code: 100520

Page 3 of 6

### 7.3. Specific end use(s)

Dental aid for prosthetics

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

#### Additional advice on limit values

none

### 8.2. Exposure controls

#### Individual protection measures, such as personal protective equipment

##### Eye/face protection

Wear eye/face protection.

##### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

##### Skin protection

Wear suitable protective clothing.

##### Respiratory protection

In case of inadequate ventilation wear respiratory protection.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state:	Liquid	
Colour:	see packing	
Odour:	neutral	
Melting point/freezing point:		not determined
Boiling point or initial boiling point and boiling range:		not determined
Flammability:		not applicable
Lower explosion limits:		not applicable
Upper explosion limits:		not determined
Flash point:		not determined
Decomposition temperature:		not determined
pH-Value:		not determined
Water solubility: (at 20 °C)		1000 g/L
Solubility in other solvents not determined		
Partition coefficient n-octanol/water:		not determined
Vapour pressure: (at 20 °C)		17 hPa
Density (at 20 °C):		1,1 g/cm <sup>3</sup>
Relative vapour density:		not determined

### 9.2. Other information

#### Information with regard to physical hazard classes

Self-ignition temperature

# Safety Data Sheet

according to UK REACH Regulation

## PERform Alginatisation

Revision date: 17.09.2024

Product code: 100520

Page 4 of 6

Solid:	not applicable
Gas:	not applicable
Oxidizing properties	
Not oxidising.	
<b>Other safety characteristics</b>	
Evaporation rate:	not determined
Solid content:	not determined

### SECTION 10: Stability and reactivity

#### 10.1. Reactivity

No hazardous reaction when handled and stored according to provisions.

#### 10.2. Chemical stability

The product is stable under storage at normal ambient temperatures.

#### 10.3. Possibility of hazardous reactions

No known hazardous reactions.

#### 10.4. Conditions to avoid

none

#### 10.5. Incompatible materials

No information available.

#### 10.6. Hazardous decomposition products

No known hazardous decomposition products.

### SECTION 11: Toxicological information

#### 11.1. Information on hazard classes as defined in GB CLP Regulation

##### Acute toxicity

Based on available data, the classification criteria are not met.

##### ATEmix calculated

ATE (oral) > 2000 mg/kg; ATE (dermal) > 2000 mg/kg; ATE (inhalation vapour) > 20 mg/l; ATE (inhalation dust/mist) > 5 mg/l

##### Irritation and corrosivity

Skin corrosion/irritation: Based on available data, the classification criteria are not met.

Serious eye damage/eye irritation: Based on available data, the classification criteria are not met.

##### Sensitising effects

Based on available data, the classification criteria are not met.

##### Carcinogenic/mutagenic/toxic effects for reproduction

Germ cell mutagenicity: Based on available data, the classification criteria are not met.

Carcinogenicity: Based on available data, the classification criteria are not met.

Reproductive toxicity: Based on available data, the classification criteria are not met.

##### STOT-single exposure

Based on available data, the classification criteria are not met.

##### STOT-repeated exposure

Based on available data, the classification criteria are not met.

##### Aspiration hazard

Based on available data, the classification criteria are not met.

### SECTION 12: Ecological information

#### 12.1. Toxicity

The product is not: Ecotoxic.

# Safety Data Sheet

according to UK REACH Regulation

## PERform Alginatisation

Revision date: 17.09.2024

Product code: 100520

Page 5 of 6

### **12.2. Persistence and degradability**

The product has not been tested.

### **12.3. Bioaccumulative potential**

The product has not been tested.

### **12.4. Mobility in soil**

The product has not been tested.

### **12.5. Results of PBT and vPvB assessment**

The substances in the mixture do not meet the PBT/vPvB criteria according to UK REACH.

The product has not been tested.

### **12.6. Endocrine disrupting properties**

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

### **12.7. Other adverse effects**

No information available.

### **Further information**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

## SECTION 13: Disposal considerations

### **13.1. Waste treatment methods**

#### **Disposal recommendations**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

#### **Contaminated packaging**

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

## SECTION 14: Transport information

### **Land transport (ADR/RID)**

#### **14.1. UN number or ID number:**

No dangerous good in sense of this transport regulation.

#### **14.2. UN proper shipping name:**

No dangerous good in sense of this transport regulation.

#### **14.3. Transport hazard class(es):**

No dangerous good in sense of this transport regulation.

#### **14.4. Packing group:**

No dangerous good in sense of this transport regulation.

### **Inland waterways transport (ADN)**

#### **14.1. UN number or ID number:**

No dangerous good in sense of this transport regulation.

#### **14.2. UN proper shipping name:**

No dangerous good in sense of this transport regulation.

#### **14.3. Transport hazard class(es):**

No dangerous good in sense of this transport regulation.

#### **14.4. Packing group:**

No dangerous good in sense of this transport regulation.

### **Marine transport (IMDG)**

#### **14.1. UN number or ID number:**

No dangerous good in sense of this transport regulation.

#### **14.2. UN proper shipping name:**

No dangerous good in sense of this transport regulation.

#### **14.3. Transport hazard class(es):**

No dangerous good in sense of this transport regulation.

#### **14.4. Packing group:**

No dangerous good in sense of this transport regulation.

### **Air transport (ICAO-TI/IATA-DGR)**

#### **14.1. UN number or ID number:**

No dangerous good in sense of this transport regulation.

#### **14.2. UN proper shipping name:**

No dangerous good in sense of this transport regulation.

#### **14.3. Transport hazard class(es):**

No dangerous good in sense of this transport regulation.

#### **14.4. Packing group:**

No dangerous good in sense of this transport regulation.

### **14.5. Environmental hazards**

# Safety Data Sheet

according to UK REACH Regulation

## PERform Alginatisation

Revision date: 17.09.2024

Product code: 100520

Page 6 of 6

ENVIRONMENTALLY HAZARDOUS: No

### 14.6. Special precautions for user

No information available.

### 14.7. Maritime transport in bulk according to IMO instruments

not applicable

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### EU regulatory information

Directive 2004/42/EC on VOC in paints and varnishes: 0,1 % (1,1 g/l)

Information according to Directive 2012/18/EU (SEVESO III): Not subject to 2012/18/EU (SEVESO III)

#### National regulatory information

Employment restrictions: Observe restrictions to employment for juveniles according to the 'juvenile work protection guideline' (94/33/EC).

Water hazard class (D): 1 - slightly hazardous to water

### 15.2. Chemical safety assessment

Chemical safety assessments for substances in this mixture were not carried out.

## SECTION 16: Other information

### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

### Further Information

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

*(The data for the relevant ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*