


Standard Gutta Percha Points



Standard Gutta Percha Points are faster and effective and are used for quick obturation after root canal preparation. These are true to size and excellent in workability with non-distorted crimp. It comes in three colours and shapes with specifications 2%, 4% and 6%.
RED Shape: Large tapered 0.06 Taper, PINK Shape: Large tapered 0.04 taper, WHITE Shape: ISO Standard 0.02 Taper.

Certification: Class IIa, CE, ISO

Specification	#15	#20	#25	#30	#35	#40	#15-#40	#45	#50	#55	#60	#65	#70	#75	#80	#45-#80
Taper	2%															
package	120pcs															

Specification	#15	#20	#25	#30	#35	#40	#15-#40	#45	#50	#55	#60	#65	#70	#75	#80	#45-#80
Taper	4%															
package	60pcs															

Specification	#15	#20	#25	#30	#35	#40	#15-#40	#45	#50	#55	#60	#65	#70	#75	#80	#45-#80
Taper	6%															
package	60pcs															

Paper Points



Paper Points are used for root canal filling and comes in three colours, shapes with specifications 2%, 4% and 6%.
RED Shape: Large tapered 0.06 Taper, PINK Shape: Large tapered 0.04 taper, WHITE Shape: ISO Standard 0.02 Taper. These points ensures long shelf life along with freshness to give flexible workability.

Certification: Class Ia,DOC,ISO

Specification	#15	#20	#25	#30	#35	#40	#15-#40	#45	#50	#55	#60	#65	#70	#75	#80	#45-#80
Taper	2%															
package	200pcs															

Specification	#15	#20	#25	#30	#35	#40	#15-#40	#45	#50	#55	#60	#65	#70	#75	#80	#45-#80
Taper	4%															
package	100pcs															

Specification	#15	#20	#25	#30	#35	#40	#15-#40	#45	#50	#55	#60	#65	#70	#75	#80	#45-#80
Taper	6%															
package	100pcs															