

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010

Dental Use



Date Issued : 11/17/2006

SDS No : 69-001.10

Date Revised : 01/23/2017

Revision No : 10

Peak® SE Primer , Part 1 of 2**SECTION 1 : Identification of the substance/preparation and of the company/undertaking****1.1. Product identifier**

Product code : BI/15076
 Product name : Peak® SE Primer , Part 1 of 2
 Product description : Professional Dental Etchant

1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses : Professional Dental Etchant

1.3. Details of the supplier of the safety data sheet**Manufacturer**

Ultradent Products, Inc.
 505 W. Ultradent Drive (10200 South)
 South Jordan, UT 84095

EC Responsible Person

Ultradent Products GmbH
 Am Westhover Berg 30
 51149 Cologne Germany
 Email: infoDE@ultradent.com
Emergency Phone : +49(0)2203-35-92-0

1.4. Emergency telephone number

CHEMTREC (NORTH AMERICA) :(800) 424 - 9300
 (INTERNATIONAL) :+1(703) 527 - 3887

SECTION 2: Hazards identification**2.1. Classification of the substance or mixture****Classification according to Regulation(EC) No 1272/2008 [CLP]**

: The mixture has been assessed and/or tested for its physical, health and environmental hazards and the following classification applies.

Health : Eye Irritation, Category 2
 Skin Irritation, Category 2
 Skin Sensitization, Category 1
 Acute Toxicity (Dermal), Category 3

Physical : Flammable Liquids, Category 2

2.2. Label elements**Classification according to Regulation(EC) No 1272/2008 [CLP]**

Hazard pictogram(s) :



Flame

Exclamation
mark

Signal Word : DANGER

Hazard statement(s) : H225: Highly flammable liquid and vapour.
 H315: Causes skin irritation.
 H317: May cause an allergic skin reaction.

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010

Dental Use



Date Issued : 11/17/2006

SDS No : 69-001.10

Date Revised : 01/23/2017

Revision No : 10

Peak® SE Primer , Part 1 of 2

H319: Causes serious eye irritation.

Precautionary statement(s)

Prevention

: P280: Wear protective gloves/protective clothing/eye protection/face protection.

Response

: P305: IF IN EYES:

P351: Rinse cautiously with water for several minutes.

P337+P313: If eye irritation persists: Get medical advice/attention.

P302+P352: IF ON SKIN: Wash with plenty of soap and water.

P332+P313: If skin irritation occurs: Get medical advice/attention.

Storage

: P273: Avoid release to the environment.

Disposal

: P501: Dispose of in compliance with governmental regulation. (EC 1975L0442-21/11/2003)

2.3. Other hazards

Conclusion: The substance does not fulfill the PBT and vPvB criteria for screening assessment; there are no indications of P or B properties.

SECTION 3: Composition / information on ingredients

3.1. Substances

Not Applicable

3.2. Mixtures

Chemical Name	CAS	EINECS No.	Wt.%	Classification according to Regulation(EC) No 1272/2008 [CLP]
Ethyl Alcohol	64-17-5	200-578-6	< 40	Flam. Liq.,Cat. 2; Eye Irr.,Cat. 2; Skin Irr.,Cat. 2; STOT SE,Cat. 3; H225; H319
2-hydroxyethyl Methacrylate	868-77-9	212-782-2	< 20	Eye Irr.,Cat. 2B; Skin Sens.,Cat. 1B; H317; H319

For full text of H-statements: see SECTION 16.

SECTION 4: First aid measures

4.1. Description of first aid measures

Following eyes

: Immediately flush eyes with plenty of water. Get medical attention, if irritation persists.

Following skin

: Wash with soap and water. Get medical attention if irritation develops or persists.

Following ingestion

: If swallowed Don't induce vomiting. Rinse mouth and seek medical advice if necessary.

Following inhalation

: If inhaled, remove to fresh air. If breathing becomes difficult, call a physician.

4.2. Most important symptoms and effects, both acute and delayed

Eyes

: Irritating to eyes.

Skin

: May cause skin irritation and sensitization.

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010
Dental Use



Date Issued : 11/17/2006
SDS No : 69-001.10
Date Revised : 01/23/2017
Revision No : 10

Peak® SE Primer , Part 1 of 2

Ingestion : May be harmful if swallowed.
Inhalation : May be harmful if inhaled.

4.3. Indication of any immediate medical attention and special treatment needed

Notes to physician : Not Applicable

SECTION 5: Fire fighting measures

5.1. Extinguishing media

Extinguishing media : Foam, dry chemical, carbon dioxide (CO₂).

5.2. Special hazards arising from the substance or mixture

General hazard : This product has been classified as a Flammable Liquids Category 2.

5.3. Advice for firefighters

Fire fighting procedures : General: Evacuate all personnel; use protective equipment for fire-fighting. Use self-contained breathing apparatus when the product is involved in fire.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General procedures : Refer to Section 8 for Personal Protective Equipment.

6.2. Environmental precautions

Water spill : Do not allow to enter sewers or drains that may lead to waterways.

6.3. Methods and material for containment and cleaning up

Small Spill : Absorb liquid and place in sealed container for disposal. Vapors can travel to an ignition source.

Large Spill : Absorb with inert, damp non-combustible material, then flush area with water.

6.4. Reference to other sections

Reference to other sections : Not Applicable

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Handling : Keep away from heat, sparks and flame.

Storage : See product labeling.

7.2. Conditions for safe storage, including any incompatibilities

Shelf life : See product labeling

7.3. Specific end use(s)

Specific end use(s) : Professional Dental Etchant

SECTION 8: Exposure controls / personal protection

8.1. Control parameters

Control parameters : Ethyl Alcohol CAS: 64-17-5

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010
Dental Use



Date Issued : 11/17/2006
SDS No : 69-001.10
Date Revised : 01/23/2017
Revision No : 10

Peak® SE Primer , Part 1 of 2

Control Parameters	Exposure Limits
OSHA PEL General Industry	1000 ppm (1900 mg/m3) TWA
OSHA PEL - Construction Industry	1000 ppm (900 mg/m3) TWA
OSHA PEL Shipyard Employment	1000 ppm(1900 mg/m3) TWA
NIOSH REL	1000 ppm(1900 mg/m3) TWA
ACGIH TLV	1000 ppm(1880 mg/m3) STEL
CAL/OSHA PEL	1000 ppm(1900 mg/m3) TWA

8.2. Exposure controls

Eye/face protection : Wear eye protection .
Skin protection : Wear suitable protective clothing and gloves.
Respiratory protection : Use only in well-ventilated areas.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Liquid
Appearance : Clear, low viscosity liquid
Colour : Yellow
Odour : Methacrylate odor
pH : 0.9 - 1.5 (This value is when mixed with part 2 of 2)
Flash point : < 27°C
Specific gravity : 1.015

9.2. Other information

Additional information : No further relevant information available.

SECTION 10: Stability and reactivity

10.1. Reactivity

Reactivity : Stable when stored and handled under recommended conditions.

10.2. Chemical stability

Chemical stability : Polymerization occurs when exposed to visible light, ultraviolet light or extreme heat.

10.3. Possibility of hazardous reactions

Hazardous polymerization : None

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010
Dental Use



Date Issued : 11/17/2006
SDS No : 69-001.10
Date Revised : 01/23/2017
Revision No : 10

Peak® SE Primer , Part 1 of 2

10.4. Conditions to avoid

Conditions to avoid : Heat, flames, ignition sources, direct light.

10.5. Incompatible materials

Incompatible materials : None known.

10.6. Hazardous decomposition products

Hazardous decomposition products : None known.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute

Dermal LD₅₀ : > 800 mg/kg rat/rabbit

Notes : ~91% of the mixture consists of components of unknown toxicity (Dermal).

Oral LD₅₀ : > 4100 mg/kg Mouse/Rat/Guinea Pig

Notes : ~32% of the mixture consists of components of unknown toxicity (Oral).

Inhalation LC₅₀ : > 490 mg/l/4H Vapors Rat/Guinea Pig

Notes : ~23% of the mixture consists of components of unknown toxicity (Inhalation).

Other information : Device is biocompatible when used as directed by dental professionals per ISO 10993-1.

SECTION 12: Ecological information

12.1. Toxicity

Aquatic toxicity (acute) : Do not allow to enter sewers or drains that may lead to waterways.

12.2. Persistence and degradability

Persistence and degradability : Not classified

12.3. Bioaccumulative potential

Bioaccumulative potential : Not classified

12.4. Mobility in soil

Mobility in soil : Not Available

12.5. Results of PBT and vPvB assessment

Results of PBT and vPvB assessment : Conclusion: The substance does not fulfill the PBT and vPvB criteria for screening assessment; there are no indications of P or B properties.

12.6. Other adverse effects

Environmental data : Do not flush to sewer.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal method : Dispose of in compliance with governmental regulation.(EC 1975L0442-20/11/2003)

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010

Dental Use



Date Issued : 11/17/2006

SDS No : 69-001.10

Date Revised : 01/23/2017

Revision No : 10

Peak® SE Primer , Part 1 of 2

SECTION 14: Transport information

14.1. UN number

UN number : 1987

14.2. UN proper shipping name

UN proper shipping name : Alcohols, n.o.s. (Ethyl alcohol mixture)

14.3. Transport hazard class(es)

Primary hazard class/division : 3

Hazard classification : 3

14.4. Packing group

Packing group : III

14.5. Environmental hazards

Marine pollutant #1 : NAP

14.6. Special precautions for user

ADR - road : NAP

RID - rail : NAP

IMDG - sea : NAP

IATA - air : NAP

14.7. Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code

Transport in bulk : NAP

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

General information : Substances of very high concern (SVHC) according to REACH, Article 57:
Contains no SVHC substances.

General comments : Please refer to Medical Devices Directive 93/42/EEC.

15.2. Chemical safety assessment

Chemical safety assessment : See Section 11

Additional information : Ethyl Alcohol is only considered carcinogenic by IARC and CA Prop 65 when it is found in beverages. None of the other ingredients in the formulation of this product are considered to be carcinogenic by IARC, NTP, or CA Prop 65.

SECTION 16: Other information

Relevant H-statements (number and full text) : Eye Irr., Cat. 2B: Eye Irritation, Category 2B
Eye Irr., Cat. 2: Eye Irritation, Category 2
Flam. Liq., Cat. 2: Flammable Liquids, Category 2
STOT SE, Cat. 3: Target Organ Toxicity (Single exposure), Respiratory Tract Category 3
Skin Irr., Cat. 2: Skin Irritation, Category 2
Skin Sens., Cat. 1B: Skin Sensitization, Category 1B
H225: Highly flammable liquid and vapour.

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010

Dental Use



Date Issued : 11/17/2006

SDS No : 69-001.10

Date Revised : 01/23/2017

Revision No : 10

Peak® SE Primer , Part 1 of 2

H317: May cause an allergic skin reaction.

H319: Causes serious eye irritation.

Prepared by

: Anu Kattoju

Revision summary

: This SDS replaces the 03/19/2015 SDS. Revised: **Section 1:** MSDS No. **Section 2:** EMERGENCY OVERVIEW - IMMEDIATE CONCERNS 2.3. Other hazards. **Section 4:** INGESTION, INHALATION, SIGNS AND SYMPTOMS OF OVEREXPOSURE (INHALATION). **Section 5:** EXPLOSION HAZARDS, EXTINGUISHING MEDIA, FIRE EXPLOSION, FIRE FIGHTING EQUIPMENT, GENERAL HAZARD, SENSITIVE TO STATIC DISCHARGE, SENSITIVITY TO IMPACT. **Section 8:** . **Section 9:** APPEARANCE, COLOR, COMMENTS, ODOR, PERCENT VOLATILE, pH, SPECIFIC GRAVITY. **Section 10:** CONDITIONS TO AVOID. **Section 11:** ACUTE (DERMAL LD₅₀ (rabbit), DERMAL LD₅₀ (rabbit), DERMAL LD₅₀ (rabbit), ORAL LD₅₀ (rat), INHALATION LC₅₀ (rat), ORAL LD₅₀ (rat), ORAL LD₅₀ (rat), INHALATION LC₅₀ (rat), INHALATION LC₅₀ (rat), INHALATION LC₅₀ (rat), ORAL LD₅₀ (rat) ACUTE, GENERAL COMMENTS. **Section 12:** CHEMICAL FATE INFORMATION, DISTRIBUTION, ECOTOXICOLOGICAL INFORMATION, ENVIRONMENTAL DATA, BIOACCUMULATION/ACCUMULATION, AQUATIC TOXICITY (ACUTE) (48-HOUR EC₅₀, 96-HOUR EC₅₀), SECTION 12: Ecological information. **Section 14:** VESSEL (IMO/IMDG) - MARINE POLLUTANT #1 ROAD AND RAIL (UK only) (CDG) - SPECIAL PROVISIONS SECTION 14: Transport information, ADR - road, IMDG - sea, IATA - air. **Section 15:** COMMENTS, GENERAL COMMENTS, RoHS. **Section 16:** GENERAL STATEMENTS.

General statements

: NAP = Not Applicable
NE = Not Established
TLV = Threshold Limit Value
PEL = Permittable exposure limits
MAK = Maximum Workplace Concentration
STEL = Short-term exposure limit
STEV = Short-term exposure value
TWA = Time Weighted Average
PPE = Personal Protective Equipment
AMU = Atomic Mass Unit

Manufacturer disclaimer

: FOR DENTAL USE ONLY: Use as directed. The information and recommendations are taken from sources (raw material SDS(s) and manufacturer's knowledge) believed to be accurate; however, the manufacturer, makes no warranty with respect to the accuracy of the information or the suitability of the recommendation and assumes no liability to any user thereof. Each user should review these recommendations in the specific context of the intended use and determine whether they are appropriate.

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010

Dental Use



Date Issued : 04/09/2014

SDS No : 343-001.3

Date Revised : 05/18/2017

Revision No : 2

Peak® SE Primer, Part 2 of 2

SECTION 1 : Identification of the substance/preparation and of the company/undertaking

1.1. Product identifier

Product code : BJ/15077
 Product name : Peak® SE Primer, Part 2 of 2
 Product description : Adhesive Primer, Part 2

1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses : Professional Dental Adhesive Primer, Part 2

1.3. Details of the supplier of the safety data sheet

Manufacturer

Ultradent Products, Inc.
 505 W. Ultradent Drive (10200 South)
 South Jordan, UT 84095

EC Responsible Person

Ultradent Products GmbH
 Am Westhover Berg 30
 51149 Cologne Germany
 Email: infoDE@ultradent.com
 Emergency Phone : +49(0)2203-35-92-0

1.4. Emergency telephone number

CHEMTREC (NORTH AMERICA) : (800) 424 - 9300
 (INTERNATIONAL) : +1(703) 527 - 3887

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation(EC) No 1272/2008 [CLP]

: The mixture has been assessed and/or tested for its physical, health and environmental hazards and the following classification applies.

Physical : Flammable Liquids, Category 2

2.2. Label elements

Classification according to Regulation(EC) No 1272/2008 [CLP]

Hazard pictogram(s) :



Flame

Signal Word : DANGER

Hazard statement(s) : H225: Highly flammable liquid and vapour.

Precautionary statement(s)

Prevention : P280: Wear protective gloves/protective clothing/eye protection/face protection.
 P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
 P233: Keep container tightly closed.
 P240: Ground and bond container and receiving equipment.
 P243: Take action to prevent static discharges.

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010
Dental Use



Date Issued : 04/09/2014
SDS No : 343-001.3
Date Revised : 05/18/2017
Revision No : 2

Peak® SE Primer, Part 2 of 2

Response : P241: Use explosion-proof [electrical/ventilating/lighting] equipment.
P242: Use non-sparking tools.

Storage : P370+P378: In case of fire: Use Water spray, alcohol-resistant foam, dry chemical or carbon dioxide to extinguish.
P303+P361+P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower].

Disposal : P403+P235: Store in a well-ventilated place. Keep cool.
P501: Dispose of in compliance with governmental regulation. (EC 1975L0442-21/11/2003)

2.3. Other hazards

Conclusion: The substance does not fulfill the PBT and vPvB criteria for screening assessment; there are no indications of P or B properties.

SECTION 3: Composition / information on ingredients

3.1. Substances

Not Applicable

3.2. Mixtures

Chemical Name	CAS	EINECS No.	Wt.%	Classification according to Regulation(EC) No 1272/2008 [CLP]
Ethyl Alcohol	64-17-5	200-578-6	< 20	Flam. Liq.,Cat. 2; Skin Irr.,Cat. 2; Eye Irr.,Cat. 2; Acute Tox. (O),Cat. 4; STOT SE,Cat. 1; H225; H319

For full text of H-statements: see SECTION 16.

SECTION 4: First aid measures

4.1. Description of first aid measures

Following eyes : Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, get medical advice/attention.

Following skin : Immediately wash skin with soap and plenty of water. Remove contaminated clothing. Get medical attention if symptoms occur. Wash clothing before reuse.

Following ingestion : If swallowed in small quantity the product is harmless, otherwise consult a doctor.

Following inhalation : No specific treatment is necessary since this material is not likely to be hazardous by inhalation.

4.2. Most important symptoms and effects, both acute and delayed

Eyes : May cause eye irritation.

Skin : May cause slight irritation.

Ingestion : May be harmful if swallowed in large quantities.

Inhalation : None expected for this product.

4.3. Indication of any immediate medical attention and special treatment needed

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010

Dental Use



Date Issued : 04/09/2014

SDS No : 343-001.3

Date Revised : 05/18/2017

Revision No : 2

Peak® SE Primer, Part 2 of 2

Notes to physician : NAP

SECTION 5: Fire fighting measures**5.1. Extinguishing media**

Extinguishing media : Water spray, alcohol-resistant foam, dry chemical or carbon dioxide

5.2. Special hazards arising from the substance or mixture

General hazard : This product has been classified as a Flammable Liquids Category 2.

5.3. Advice for firefighters

Fire fighting procedures : General: Evacuate all personnel; use protective equipment for fire-fighting. Use self-contained breathing apparatus when the product is involved in fire.

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures**

General procedures : Refer to Section 8 for Personal Protective Equipment.

6.2. Environmental precautions

Water spill : Do not allow to enter sewers or drains that may lead to waterways.

6.3. Methods and material for containment and cleaning up

Small Spill : Absorb liquid and place in sealed container for disposal. Vapors can travel to an ignition source.

Large Spill : Absorb with inert, damp non-combustible material, then flush area with water.

6.4. Reference to other sections

Reference to other sections : NAP

SECTION 7: Handling and storage**7.1. Precautions for safe handling**

General procedures : Avoid contact with eyes, skin and clothing.

Handling : Keep away from heat, sparks and flame.

Storage : See product labeling.

7.2. Conditions for safe storage, including any incompatibilities

Shelf life : See product labeling

7.3. Specific end use(s)

Specific end use(s) : Professional Dental Adhesive Primer, Part 2

SECTION 8: Exposure controls / personal protection**8.1. Control parameters**

Control parameters : Ethyl Alcohol CAS: 64-17-5

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010
Dental Use



Date Issued : 04/09/2014
SDS No : 343-001.3
Date Revised : 05/18/2017
Revision No : 2

Peak® SE Primer, Part 2 of 2

Control Parameters	Exposure Limits
OSHA PEL General Industry	1000 ppm (1900 mg/m3) TWA
OSHA PEL - Construction Industry	1000 ppm (900 mg/m3) TWA
OSHA PEL Shipyard Employment	1000 ppm(1900 mg/m3) TWA
NIOSH REL	1000 ppm(1900 mg/m3) TWA
ACGIH TLV	1000 ppm(1880 mg/m3) STEL
CAL/OSHA PEL	1000 ppm(1900 mg/m3) TWA

8.2. Exposure controls

Eye/face protection : Wear eye protection .
Skin protection : Wear suitable gloves.
Respiratory protection : Good general ventilation should be sufficient to control airborne levels.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Liquid
Appearance : Clear, low viscosity liquid
Odour : Ethanol
pH : NAP (non-aqueous)
Flash point : 43°C
Specific gravity : 0.971

9.2. Other information

Additional information : No further relevant information available

SECTION 10: Stability and reactivity

10.1. Reactivity

Reactivity : Stable

10.2. Chemical stability

Chemical stability : Stable when stored and handled under recommended conditions.

10.3. Possibility of hazardous reactions

Hazardous polymerization : None

10.4. Conditions to avoid

Conditions to avoid : Heat, flames and sparks. Extremes of temperature and direct sunlight.

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010

Dental Use



Date Issued : 04/09/2014

SDS No : 343-001.3

Date Revised : 05/18/2017

Revision No : 2

Peak® SE Primer, Part 2 of 2

10.5. Incompatible materials

Incompatible materials : None known.

10.6. Hazardous decomposition products

Hazardous decomposition products : None known.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute

Dermal LD₅₀ : Not able to calculate ATE mix (dermal) due to lack of data on the components.

Oral LD₅₀ : > 20000 mg/kg Mouse

Notes : There are no components of unknown toxicity (oral).

Inhalation LC₅₀ : > 240000 mg/L/4H vapors Mouse

Notes : There are no components of unknown toxicity (inhalation).

Germ cell mutagenicity : The components in this product are not reported to produce mutagenic effects in humans.

Carcinogenicity : None of the components of this material are listed by IARC, NTP, or CA Prop 65 as carcinogens.

Reproductive toxicity : The components in this product are not reported to cause reproductive effects in humans.

Other information : Device is biocompatible when used as directed by dental professionals per ISO 10993-1.

SECTION 12: Ecological information

12.1. Toxicity

Aquatic toxicity (acute) : Do not allow to enter sewers or drains that may lead to waterways.

12.2. Persistence and degradability

Persistence and degradability : Not classified

12.3. Bioaccumulative potential

Bioaccumulative potential : Not classified

12.4. Mobility in soil

Mobility in soil : Not Available

12.5. Results of PBT and vPvB assessment

Results of PBT and vPvB assessment : Conclusion: The substance does not fulfill the PBT and vPvB criteria for screening assessment; there are no indications of P or B properties.

12.6. Other adverse effects

Environmental data : No further relevant information available.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010
Dental Use



Date Issued : 04/09/2014
SDS No : 343-001.3
Date Revised : 05/18/2017
Revision No : 2

Peak® SE Primer, Part 2 of 2

Disposal method : Dispose of in compliance with governmental regulation.(EC 1975L0442-20/11/2003)

SECTION 14: Transport information

14.1. UN number

UN number : 1987

14.2. UN proper shipping name

UN proper shipping name : Alcohols, n.o.s. (Ethyl alcohol mixture)

14.3. Transport hazard class(es)

Hazard classification : 3

14.4. Packing group

Packing group : III

14.5. Environmental hazards

Marine pollutant #1 : NAP

14.6. Special precautions for user

ADR - road : NAP

RID - rail : NAP

IMDG - sea : NAP

IATA - air : NAP

14.7. Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code

Transport in bulk : NAP

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

General information : Substances of very high concern (SVHC) according to REACH, Article 57:
Contains no SVHC substances.

General comments : Please refer to Medical Devices Directive 93/42/EEC.

15.2. Chemical safety assessment

Chemical safety assessment : See Section 11

Additional information : None of the ingredients are listed as a carcinogen by IARC, NTP, or CA Prop 65.

SECTION 16: Other information

Relevant H-statements (number and full text) : Acute Tox. (O), Cat. 4: Acute Toxicity (Oral), Category 4
Eye Irr., Cat. 2: Eye Irritation, Category 2
Flam. Liq., Cat. 2: Flammable Liquids, Category 2
STOT SE, Cat. 1: Target Organ Toxicity (Single exposure), Category 1
Skin Irr., Cat. 2: Skin Irritation, Category 2
H225: Highly flammable liquid and vapour.
H319: Causes serious eye irritation.

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010

Dental Use



Date Issued : 04/09/2014

SDS No : 343-001.3

Date Revised : 05/18/2017

Revision No : 2

Peak® SE Primer, Part 2 of 2

Prepared by

: Anu Kattoju

Revision summary

: This SDS replaces the 09/03/2015 SDS. Revised: **Section 1:** GENERAL USE, MSDS No. **Section 2:** EMERGENCY OVERVIEW - IMMEDIATE CONCERNS 2.3. Other hazards. **Section 4:** EYES, INHALATION, NOTES TO PHYSICIAN, SIGNS AND SYMPTOMS OF OVEREXPOSURE - EYES SKIN. **Section 5:** EXPLOSION HAZARDS, EXTINGUISHING MEDIA, FIRE EXPLOSION, FIRE FIGHTING EQUIPMENT, GENERAL HAZARD, SENSITIVE TO STATIC DISCHARGE, SENSITIVITY TO IMPACT. **Section 6:** . **Section 7:** GENERAL PROCEDURES, . **Section 8:** PERSONAL PROTECTIVE EQUIPMENT - EYES AND FACE . **Section 9:** APPEARANCE, COMMENTS, PERCENT VOLATILE, pH, SPECIFIC GRAVITY. **Section 11:** ACUTE (ORAL LD₅₀ (rat), INHALATION LC₅₀ (rat), ORAL LD₅₀ (rat), ORAL LD₅₀ (rat), INHALATION LC₅₀ (rat), INHALATION LC₅₀ (rat), INHALATION LC₅₀ (rat), ORAL LD₅₀ (rat) ACUTE, GENERAL COMMENTS, CARCINOGENICITY, GERM CELL MUTAGENICITY, REPRODUCTIVE TOXICITY. **Section 12:** CHEMICAL FATE INFORMATION, DISTRIBUTION, ENVIRONMENTAL DATA, BIOACCUMULATION/ACCUMULATION, AQUATIC TOXICITY (ACUTE) (48-HOUR EC₅₀, 96-HOUR EC₅₀), SECTION 12: Ecological information. **Section 14:** VESSEL (IMO/IMDG) - MARINE POLLUTANT #1 ROAD AND RAIL (UK only) (CDG) - SPECIAL PROVISIONS SECTION 14: Transport information, Primary hazard class/division, Secondary hazard class/division, ADR - road, IMDG - sea, IATA - air. **Section 15:** COMMENTS, GENERAL COMMENTS, RoHS. **Section 16:** GENERAL STATEMENTS.

General statements

: NAP = Not Applicable
NE = Not Established
TLV = Threshold Limit Value
PEL = Permittable exposure limits
MAK = Maximum Workplace Concentration
STEL = Short-term exposure limit
STEV = Short-term exposure value
TWA = Time Weighted Average
PPE = Personal Protective Equipment

Manufacturer disclaimer

: FOR DENTAL USE ONLY: Use as directed. The information and recommendations are taken from sources (raw material SDS(s) and manufacturer's knowledge) believed to be accurate; however, the manufacturer, makes no warranty with respect to the accuracy of the information or the suitability of the recommendation and assumes no liability to any user thereof. Each user should review these recommendations in the specific context of the intended use and determine whether they are appropriate.